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When the «▼»-symbole appears, battery is exhausted and needs replacing. Please remove the compartment lid of the thermometers and replace the exhausted battery by positioning the «+» at the top.

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## IMPORTANT SAFETY INSTRUCTIONS

- The thermometer is suitable exclusively for measurement of the body temperature!
- The minimum measurement time until the signalling tone [beep] must be maintained without exception!
- The thermometer contains small parts (battery, etc.) that can be swallowed by children. Therefore never leave the thermometer unattended to children!
- Protect the thermometer from hitting and striking!
- Avoid bending the tip more than 45°!
- Avoid ambient temperatures above 60 °C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer.
- The manufacturer unattended to children!
- The thermometer must be intact when immersed in liquid disinfectant.
- The manufacturer recommends to verify the accuracy by an authorised laboratory every 2 years.
- Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

To turn on the thermometer, push the button on the top; a short signalling tone [beep] signals «thermometer ON». A display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C, «-» and a flashing «C» appear at the upper right of the display. The thermometer is now ready for use.

### TURNING ON THE THERMOMETER

Choose the preferred method. During a measurement, the current temperature is displayed continuously. The «C» symbol flashes during this period of time. If the signalling tone sounds (*beep-beep-deep-deep*), the «C» symbol ceases to flash, this means that the measured increase in temperature is now less than 0,1 °C in 8 seconds and that the thermometer is ready for reading.

### USING THE THERMOMETER

10 times) and the «C» is no longer flashing, this means that the measured increase in temperature is now less than 0,1 °C in 8 seconds and that the thermometer is ready for reading.

### STORAGE OF MEASURED VALUES

If you hold the button down for more than 3 seconds when turning the thermometer on, the maximum temperature stored automatically at the last measurement will be displayed. At the same time, a small «M» for memory will appear at the right of the display. About 2 seconds after release the button, this value will disappear and the thermometer will be turned off.

### TURNING ON THE THERMOMETER

Press the button on the top; a short signalling tone [beep] signals «thermometer ON». A display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C, «-» and a flashing «C» appear at the upper right of the display. The thermometer is now ready for use.

## TECHNICAL SPECIFICATIONS

Type / Type:	Maximum thermometer
Measurement range / Gamme de mesure:	32,0 °C to 42,9 °C Temp. < 32,0 °C: display «-» for low (too low) Temp. > 42,9 °C: display «H» for high (too high)
Measurement accuracy / Precision:	± 0,1 °C between 34 °C and 42 °C ± 0,1 °C entre 34 °C et 42 °C
Operating temperature / Température de fonctionnement:	18 °C to 28 °C
Self-test / Dispositif d'autocontrôle:	Automatic internal check at a test value of 37 °C. If there is a deviation of > 0,1 °C, «ERR» (error) is displayed.
Display / Affichage:	Liquid crystal display (LCD) with 3 digits. Smallest unit of display 0,1 °C
Signalling tone / Signal sonore:	For signalling that the thermometer is ready to use and that the temperature increase is less than 0,1 °C / 8 seconds. For indicating that the thermometer is ready to use and that the temperature increase is less than 0,1 °C / 8 seconds.
Memory / Mémoire:	For storing the last measured value.
Storage temperature / Température de stockage:	-10 °C to +60 °C
Battery / Pile:	1,5 / 1,55 V, LR41
Reference to standards / Référence aux normes:	EN 12470-3: 2000, clinical thermometers ASTM E-1112, electronic thermometers

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We grant you a lifetime guarantee after the date of purchase. Any damage caused by improper handling shall not be covered by the guarantee. Batteries and packaging are also excluded from the guarantee. All other damage claims excluded. A guarantee claim must be submitted with the purchase receipt. Please pack your defective thermometer well and send with sufficient postage to:

Date of purchase:

Date d'achat: