Europe / Middle-East / Africa

Microlife AG

Espenstrasse 139

9443 Widnau / Switzerland

Tel. +41 / 71 727 70 30

Fax +41 / 71 727 70 39

Email admin@microlife.ch

www.microlife.com

Asia

Microlife Corporation.

9F, 431, RuiGang Road, NeiHu

Taipei, 114, Taiwan, R.O.C.

Tel. 886 2 8797-1288

Fax 886 2 8797-1283

Email service@microlife.com.tw

www.microlife.com

North / Central / South America

Microlife USA, Inc.

424 Skinner Blvd., Suite C

Dunedin, FL 34698 / USA

Tel. +1 727 451 0484

Fax +1 727 451 0492

Email msa@microlifeusa.com

www.microlife.com

microlife BP 3AX1

Wrist Watch Blood Pressure Monitor

Instruction Manual (1-19)

Tensiomètre-bracelet

Mode d'emploi (20-39)

Tensiómetro de Muñeca Dot Matrix

Manual de instrucciones (40-59)

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microlife®

Wrist Watch Blood Pressure Monitor

Instruction Manual

Table of Contents

1. Introduction

- 1.1. Features
- 1.2. Important Information about Self-Measurement

2. Important information on the subject of blood-pressure and its measurement

- 2.1. How does high/low blood-pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular increased/low values are obtained?
- 2.4. MAM technology (Microlife Average Mode)

3. The various components of the blood-pressure monitor

4. Putting the blood-pressure monitor into operation

- 4.1. Inserting the batteries
- 4.2. Initial Setup
- 4.3. Setting Alarm Times
- 4.4. Sleep Mode

5. Carrying out a measurement

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the Cuff
- 5.4. Measurement Mode Selection
- 5.5. Measurement Procedure
- 5.6. Discontinuing a measurement

6. Memory handling and user setting

- 6.1. Viewing Measurements Saved in Memory
- 6.2. Setting the User
- 6.3. Clearing the Memory

7. Adjusting the Display Contrast

8. Error messages / malfunctions

9. Care and maintenance, recalibration

10. Guarantee

11. Technical Specifications

12. www.microlife.com

1. Introduction

1.1. Features

The blood-pressure monitor is a fully automatic, digital blood-pressure measuring device for use on the wrist. This enables very fast and reliable measurement of systolic and diastolic blood pressures as well as the pulse frequency by way of the oscillometric measuring method. The device offers a very high and clinically tested measurement accuracy and has been designed to provide a maximum of user-friendliness.

Before using, please read through this instruction manual carefully and then keep it in a safe place. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

Attention!

1.2. Important information about self-measurement

- Do not forget: self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.
- Clinical trials have shown that wrist measurements can reach similar accuracy to upper arm measurements. However, some individuals may get a different result between both methods. Therefore, we generally recommend to compare the readings of this device to data obtained by your doctor. If there is a difference, always consider the variance into your final reading.

Electromagnetic interference

The device contains sensitive electronic components. Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

2. Important information on the subject of blood-pressure and its measurement

2.1. How does high/low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory centre, and adapted to the respective situation by way of feedback via the nervous system. To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood-vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: During the «blood ejection» (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's «rest period» (Diastole) minimal (diastolic blood-pressure value).

The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood-pressure is over 140 mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

With blood-pressure values that are too low, i.e. systolic values under 100 mmHg and/or diastolic values under 60 mmHg, likewise, please consult your doctor.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.**

Table for classifying blood-pressure values (units mmHg) according to World Health Organization:

Range	Systolic Blood-pressure	Diastolic Blood-pressure	Recommendations
Hypotension	lower than 100	lower than 60	Consult your doctor
Normal range	between 100 and 140	between 60 and 90	Self-check
Mild hypertension	between 140 and 160	between 90 and 100	Consult your doctor
Moderate hypertension	between 160 and 180	between 100 and 110	Consult your doctor
Severe hypertension	higher than 180	higher than 110	Consult your doctor urgently

Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». Please consult your doctor if you suspect that this might be the case.
- Correctly measured diastolic blood-pressure values above 120mmHg require **immediate medical treatment**.

2.3. What can be done, if regular increased/low values are obtained?

- a) Please consult your doctor.
- b) Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result. Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.
- c) There are many different causes of the appearance of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take, not only for reducing a medically established high blood pressure, but also for prevention. These measures are part of your general way of life:

A) Eating habits

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

Follow consistently any medical instructions for treating previous illness such as:

- Diabetes (Diabetes mellitus)
- Fat metabolism disorder
- Gout

C) Habits

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)

D) Physical constitution

- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

2.4. MAM technology (Microlife Average Mode)

- «MAM technology» is a new type of concept for optimum reliability in self-measurement of blood pressure.
- An advanced measurement accuracy is achieved by the automatic analysis of three successive measurements.
- The new system provides reliable values for the doctor and can be used as the basis for reliable diagnostics and medication therapy for high blood pressure.

A) Why MAM?

- Human blood pressure is not stable
- Scattering of devices

B) Key advantages

Reduction of:

- Device scattering
- Insufficient rest prior to measurement
- Movement artefacts
- Cuff positioning influences

C) Medical benefits


- Improved accuracy
- Reliable patient self-measurement data for the doctor
- Safe hypertension diagnostic
- Reliable therapy control

D) Measurement sequence

- With the full measurement cycles, the total measurement time remains less than 3 minutes, compared to a single measurement time of 1 and a half minutes.
- Single results are not displayed.
- Due to the «Data Analysis» result, a 4th measurement may be applied.

3. The various components of the blood-pressure monitor

The illustration shows the blood-pressure monitor, consisting of:


Wrist cuff Type WC1 12–21 cm, for wrist circumference of 12 to 21 cm

4. Putting the blood-pressure monitor into operation

4.1. Inserting the batteries

After you have unpacked your device, first insert the batteries.

1. Remove the cover as illustrated (A).
2. Lay the plastic strip (B) inside the battery compartment, and insert 2 x size AAA 1.5V batteries (C), observing the indicated polarity according to the markings on the cuff (D).

Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- If the battery warning appears in the display, the batteries are empty and must be replaced by new ones.
- Use «AAA» Long-Life or Alkaline 1.5V Batteries. The use of rechargeable batteries is not recommended.
- If the blood pressure monitor is left unused for long periods, remove the batteries from the device.

4.2. Initial Setup

When the batteries are installed, the display shows the language selection screen.

Press ▲ or ▼ to select the display language. Press **OK** to confirm the selection. The default language is English.

The device prompts you to set the time and date. The year blinks highlighting the first field to change.

Press ▲ or ▼ to change the year. Press **OK** to confirm the date.

The month, day, hour, and minute are all consecutively highlighted. Press ▲ or ▼ to change these fields and press **OK** to confirm. You are then prompted to set the alarms.

4.3. Setting Alarm Times


After you have set the time, you can set one or two alarms to alert you to take your blood pressure measurements. An alarm icon appears in the display.

Press ▲ or ▼ to set the hour and minute for Alarm 1. Press **OK** to confirm each setting.

When the alarm icon is highlighted, press ▲ or ▼ to enable  or disable  Alarm 1. Press **OK** to confirm.

You are then prompted to set Alarm 2. Follow the instructions for Alarm 1 to set this alarm.

After setting the alarms, press **OK** to exit the Time and Date Setup. An alarm icon  appears in the display to indicate that the alarm is activated.

Attention!


The time display is set on a 24-hour clock.

When the final field DISPLAY OFF is highlighted, press ▲ or ▼ to select whether you want the time and date display to remain on when the device is switched off. Press **OK** to confirm the settings and exit the time setup screen.

4.4. Sleep Mode

The device is designed to stay in sleep mode as batteries are loaded. The illustration shows the display elements, when the BPM is in sleep mode, and the DISPLAY OFF feature is disabled (see Section 4.2).

5. Carrying out a measurement

5.1. Before the measurement

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about 5 minutes before the measurement.
- Measure always on the same wrist (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.


5.2. Common sources of error

Note: Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the wrist artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood-pressure will be measured! (Each 15cm difference in height results in a measurement error of 10mmHg!)
- A loose cuff causes false measurement values.

5.3. Fitting the cuff

- Remove all eventual objects and jewellery (e.g. wristwatch) from the wrist in question. Draw the cuff over the wrist.
- The distance between the cuff and the hand should be approx. 10 mm.
- Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight, whereby no space should remain between the cuff and the wrist.
- Lay the arm on a table, with the palm upwards. Support the arm a little with a rest (cushion), so that the cuff rests at about the same height as the heart. Take care, that the cuff lies free. Remain so for 5 minutes sitting quietly, before beginning with the measurement.


5.4. Measurement Mode Selection

Choose which measurement mode you want to use before taking a blood pressure reading. Press **OK** to enter the menus, press **▲** or **▼** to select Measuring mode, and press **OK** to enter the submenu. Select «1» or «3/MAM» by the selection buttons and confirm your selection by **OK**.

Average (MAM) Mode

In average mode, the unit takes 3 continuous measurements, and afterwards calculates the result based on an average using MAM technology (see Section 2.4).

There is a 15-second rest time between measurements; the unit beeps during the last 5 seconds to remind you of the upcoming measurement.

Single Mode


In single mode «1», only one measurement is taken.

5.5. Measurement Procedure

After the cuff and arm have been appropriately positioned and the measurement mode is set, the measurement can be taken as follows. Press the O/I button. The Microlife Welcome screen appears, and is followed by a screen which shows you how to properly hold your arms. If the display off function is enabled, the Sleep Mode screen appears for two seconds before measurements are taken.


The pump begins to inflate the cuff. And the increasing cuff pressure is displayed.

After the suitable inflation pressure has been reached, the pump stops, the measurement signals are displayed and the heart icon flashes.


This curve shows each heart beat activity, Movement artefacts would become visible by large deviations within the curve.

A long beep sounds when the measurement has been completed. The systolic and diastolic blood pressure and pulse rate appear in the display. The measurement reading remains on the display for approximately 30 seconds. If no button is pressed during the 30-second period the instrument switches itself off in order to preserve the batteries. When the unit is set to Average (MAM) Mode, three separate measurements take place as explained under 5.4.

	MAM
3/03	15:48
sys	124 mmHg
dia	73 mmHg
pul	77 /min.

If the data of each cycle differ too much from each other, a fourth measurement is performed before the result is displayed. In rare cases when the blood pressure is so unstable that even after four measurements the data vary too much, an error message is displayed. Refer to Error Messages and Malfunctions (Section 7) at the end of this manual.

5.6. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason, for example if the patient feels unwell, press the O/I button. The device immediately lowers the cuff-pressure and enters into sleep mode.

6. Memory handling and user setting

6.1 Viewing Measurements Saved in Memory

The blood pressure monitor automatically stores each of the last 30 measurement values. To view the stored measurements for a user, press **▲** or **▼** in sleep mode. The last measurement taken for the user is displayed.

The user, measurement sequence number, date, and time are displayed on the top. Data obtained in Average Mode are indicated by a MAM symbol in the upper right corner.

Press **▲** or **▼** to scroll through the list items. The first item in the list is the average of all stored data.

Press **OK** and the display shows table of all data.

The first entry is the average. Data obtained in Average Mode are indicated by an «A» symbol at the very right.


Press **▲** or **▼** to scroll through the list.

Press **OK** and the display shows graph of all data. The lines indicate the systolic (top end) and Diastolic (bottom end) values. Press **▲** or **▼** to move the cursor between measurements. The currently selected measurement is indicated by a blinking line. Press **O/I** to exit this mode.

Measurement
number

	9	MAM
3/03		15:48
sys	124	mmHg
dia	73	mmhg
pul	77	/min.

	3/03	15:48	
No.	sys	dia	pul
avg.	149	110	85
1	145	109	83
2	119	71	74 A
3	124	73	79
4	130	80	77 A

6.2. Setting the User

This monitor allows you to track readings for two individuals independently. Scroll to User and press **OK** to enter the submenu.

Press **▲** or **▼** to make your selection. Select  to set the device to not save readings.



Select  to select user 1. Any readings made when user 1 is selected will be saved separately from other readings.

Select  to select user 2. Any readings made when user 2 is selected will be saved separately from other readings.

After selecting user 1 or user 2 you are prompted to enter the user name. You can use up to six characters. To keep the default name «User_1», press **O/I**. To enter a different name, follow these instructions:

Press **▲** or **▼** to scroll through the alphanumeric characters A ~ Z and 0 ~ 9. There is also one underscore and one space character. Press **OK** to accept the character and move to the next space.


After you have finished entering the name, press **OK** to exit the user setup menu. The user name is displayed on the screen.

6.3. Clearing the Memory

This item enables you to clear readings for user 1 or user 2. You can erase single readings, or erase all readings at once.

Select user 1 or user 2 as explained under 6.2. Scroll to Clear Memory and press **OK**. The user you selected is shown in the upper left corner.


Clear Single Readings

Select «SINGLE» by the selection buttons and press **OK** to erase single readings. Press **▲** or **▼** to select a reading from the list. An arrow **▶** highlights the selected reading. Press **OK** to erase the selected reading.

CLEAR MEMORY				
		3/03	15:48	
No.	SYS	DIA	PUL	
1	124	84	76	
2	139	95	74	
3	139	98	80	
4	146	106	72	


Clear All Readings

Select «All» by the selection buttons and press **OK** to erase all readings of the selected user. You are prompted to confirm the action. Press **OK** to erase all readings of the selected user. Press the 0/I button to cancel deletion of all readings.

7. Adjusting the Display Contrast

This item enables you to adjust the contrast. Scroll to Contrast and press **OK** to enter the submenu. Press **▲** to increase the display contrast. Press **▼** to decrease the display contrast. After you have adjusted the contrast, press **OK** to confirm the setting.

8. Error messages / malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 1).

Error No.	Possible cause(s)
ERR 1	No pulse has been detected.
ERR 2	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the Measurement (Artefact).
ERR 3	The inflation of the cuff takes too long. The cuff is not correctly seated.
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.
ERR 6	Single data differ too much during Average Mode even after 4 cycles. No average result can be displayed.
HI	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
LO	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

* Please consult your doctor, if this or any other problem occurs repeatedly.

Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched on although the batteries are in place.	<ol style="list-style-type: none"> 1. Check batteries for correct polarity and if necessary insert correctly. 2. If the display is unusual, re-insert batteries or exchange them.
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none"> 1. Check the positioning of the cuff. 2. Measure the blood-pressure again in peace and quiet under observance of the details made under point 5.
Every measurement produces a different value although the instrument functions normally and the values displayed are normal.	<ul style="list-style-type: none"> • Please read the following information and the points listed under «Common sources of error». Repeat the measurement.
Blood pressure measured differ from those values measured by the doctor.	<ul style="list-style-type: none"> • Record the daily development of the values and consult your doctor.

Further Information

The level of blood-pressure is subject to fluctuations even with healthy people. Therefore it is important to ensure that to compare **measurements always requires the same conditions (Quiet conditions)!**

You must consult your specialist dealer or chemist if there are technical problems with the blood-pressure instrument. **Never attempt to repair the instrument yourself!**

Any unauthorised opening of the instrument invalidates all guarantee claims!

9. Care and maintenance, recalibration

- Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight compartment. Handle this carefully and avoid all types of straining through twisting or buckling.
- Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!

- d) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- e) Never open the device! Otherwise the manufacturer calibration becomes invalid!

Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display every 2 years.

Your specialist dealer would be pleased to provide more extensive information about this.

10. Guarantee

The blood-pressure monitor is guaranteed for **3 years** from date of purchase. This guarantee includes the instrument and the cuff. The guarantee does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument by third parties.

The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

Name and company address of the responsible dealer:

11. Technical specifications

Weight:	140g (with batteries)
Size:	67 x 70 x 80 mm
Storage temperature:	-5 to +50°C 15 to 90% relative humidity maximum
Operation temperature:	10 to 40°C
Display:	dot-matrix LCD (L iquid C rystal D isplay)
Measuring method:	oscillometric
Pressure sensor:	solid-state
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 per minute
Measuring resolution:	1 mmHg
Accuracy:	Pressure within ± 3 mmHg Pulse ± 5 % of the reading
Memory:	Storing the last 2x30 measurements automatically
Power source:	2 x 1.5 V Batteries; size AAA
Accessories:	Storage case
Reference to standards:	EU Directives 93/42/EEC; NIBP requirements: EN 1060-1 /-3 /-4, ANSI / AAMI SP10

Technical and colour alterations reserved!

12. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com

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