Europe / Middle-East / Africa

Microlife AG

Espenstrasse 139

9443 Widnau / Switzerland

Tel. +41 / 71 727 70 30

Fax +41 / 71 727 70 39

Email admin@microlife.ch

www.microlife.com

Asia

Microlife Corporation.

9F, 431, RuiGang Road, NeiHu

Taipei, 114, Taiwan, R.O.C.

Tel. 886 2 8797-1288

Fax 886 2 8797-1283

Email service@microlife.com.tw

www.microlife.com

North / Central / South America

Microlife USA, Inc.

424 Skinner Blvd., Suite C

Dunedin, FL 34698 / USA

Tel. +1 727 451 0484

Fax +1 727 451 0492

Email msa@microlifeusa.com

www.microlife.com

microlife BP 3BTO-AP

Automatic Blood Pressure Monitor

Instruction Manual (1-11)

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Table of contents

1. Introduction

- 1.1. Features
- 1.2. Important information about self-measurement

2. Important information on the subject of blood pressure and its measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?
- 2.4. PAD - Pulse Arrhythmia Detection

3. Components of your blood pressure monitor

4. Using your blood pressure monitor

- 4.1. Inserting the batteries
- 4.2. Using an AC adapter (special accessory)
- 4.3. Cuff connection
- 4.4. Setting the time and date

5. Carrying out a measurement

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Discontinuing a measurement
- 5.6. Memory – recall of the measurements
- 5.7. Memory – cancellation of all measurements

6. Error messages/malfunctions

7. Care and maintenance, recalibration

8. Guarantee

9. Technical specifications

10. www.microlife.com

1. INTRODUCTION

1.1. Features

Your blood pressure monitor is a fully automatic, digital blood pressure measuring device with integrated PAD technology for use on the upper arm. This monitor enables very fast and reliable measurement of the systolic and diastolic blood pressure as well your pulse by use of the oscillometric method.

This device is easy to use and is proven in clinical studies to provide excellent accuracy. The large display makes it easy to visualize operational status and cuff pressure during a measurement.

Please read through this instruction manual carefully before use and then keep it in a safe place. For further questions on the subject of blood pressure and its measurement, please contact your doctor.

Attention!

1.2. Important information about self-measurement

- Self-measurement means **control**, not diagnosis or treatment. Your values must always be discussed with your doctor. You should never alter the dosages of any medication without direction from your doctor.
- The pulse reading is **not** suitable for checking heart pacemakers!
- If you have been diagnosed with a severe Arrhythmia or irregular heartbeat, measurements made with this instrument should only be evaluated after consultation with the doctor.

Electromagnetic interference:

This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens. Electrical interference can lead to temporary impairment of the measuring accuracy).

2. IMPORTANT INFORMATION OF BLOOD PRESSURE AND ITS MEASUREMENT

2.1. How does high/low blood pressure arise?

Your blood pressure level is determined in the circulatory center of your brain. Your nervous system allows your body to adapt or alter blood pressure in response to different situations. Your body alters your pulse and the width of blood vessels through changes in muscles in the walls of blood vessels.

Your blood pressure reading is highest when your heart pumps or ejects blood. This stage is called your Systolic Blood Pressure. Your blood pressure is lowest when the heart rests. (in-between beats) This is called your diastolic blood pressure.

It is critical to maintain blood pressure values within a «normal» range(s) in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90mmHg and/or the systolic blood pressure is over 140mmHg. If you obtain readings in this range, consult your doctor immediately. High blood pressure values over time damage blood vessels, vital organs such as the kidney and even your heart.

When blood pressure values are too low, i.e. systolic values under 100mmHg and/or diastolic values under 60mmHg, consult your doctor.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to independently alter drug doses prescribed by your doctor.**

Table for classifying blood pressure values (units mmHg): World Health Organization:

Range	Systolic blood pressure	Diastolic blood pressure	Recommendations
Hypotension	lower than 100	lower than 60	Check with your doctor
Normal range	between 100 and 140	between 60 and 90	Self-check
Mild hypertension	between 140 and 160	between 90 and 100	Consult your doctor
Moderately serious hypertension	between 160 and 180	between 100 and 110	Consult your doctor
Serious hypertension	higher than 180	higher than 110	Urgent! See your doctor.
Specific systolic hypertension	higher than 140	lower than 90	Consult your doctor

Further information

- If your values are mostly «normal» under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». In any case, please discuss the values with your doctor.
- Correctly measured diastolic blood pressure values above 120mmHg require **immediate medical treatment**.

2.3. What can be done, if regular increased/low values are obtained?

- a) Please consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) over time are associated with considerable risks to health. Blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). This can result in a deficient supply of blood to important organs (heart, brain, muscles). Additionally, the heart will become structurally damaged.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. Secondary hypertension can cause organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are lifestyle changes you can make to prevent and reduce high blood pressure. These measures must be part of a healthy lifestyle and include:

A) Eating habits

- Maintain a normal weight as prescribed by your doctor.
- Avoid excessive consumption of common salt. Please note many «packaged foods» contain high levels of salt.
- Avoid fatty foods. (Packaged foods are frequently high in fats)

B) Previous illnesses

Consistently follow any medical instructions for treating previous illness such as:

- Diabetes (Diabetes mellitus or sugar diabetes)
- Fat metabolism disorder
- Gout

C) Habits

- Eliminate smoking
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee, tea, chocolate, etc.)

D) Physical conditioning

- After a preliminary medical examination, exercise regularly.
- Choose sports which require endurance and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. Your physician will help you develop an exercise routine that is appropriate for you.

2.4. PAD - Pulse Arrhythmia Detection

Appearance of the Arrhythmia indicator

The appearance of this symbol  signifies that a certain pulse irregularity was detected during the measurement. The result can vary from your normal blood pressure. As a rule this is not a cause for concern; however, if the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend you inform your doctor. Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Arrhythmia indicator

The Microlife device BP 3BT0-AP is an oscillometric blood pressure measuring device that also analyses the pulse frequency during measurement. The accuracy of this device has been clinically validated.


If pulse irregularities occur during measurement, the arrhythmia symbol is displayed after the measurement. If the symbol appears **more frequently** (e.g. several times per week on measurements performed daily) **or if it suddenly appears more often than usual**, we recommend the patient to seek medical advice.

The device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

3. COMPONENTS OF YOUR BLOOD PRESSURE MONITOR

The illustration shows the blood pressure monitor, consisting of:


a) Measuring unit

b) Upper arm cuff

M-size-cuff (22 - 32 cm)

L-size-cuff (32 - 42 cm) – available as special accessory


4. USING YOUR BLOOD PRESSURE MONITOR

4.1. Inserting the batteries


Insert batteries immediately after unpacking the device. The battery compartment is located on the bottom of the device (see illustration).

- Remove cover as illustrated
- Insert the batteries (4 x size AA 1,5V), thereby observing the indicated polarity.
- If the battery warning appears in the display, the batteries are empty and must be replaced.

Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AA» Long-Life or Alkaline 1.5 V Batteries. The use of 1.2 V Accumulators is not recommended.
- If the blood pressure monitor is left unused for long periods, please remove the batteries from the device.


Functional check: Hold on the START button down to test all the display symbols. When functioning correctly all symbols must appear.

4.2. Using an AC adapter (special accessory)

This blood pressure monitor can be operated with the Microlife mains adapter (output 6 V DC / 600 mA with plug).

- a) Plug the plug into the socket at the left side of the instrument.
- b) Plug the mains adapter into a 230V or 110V power socket. Test that power is available by pressing the start button.


Note:

- No power is taken from the batteries while the mains adapter is connected to the instrument.
- If power is interrupted during the measurement (e.g. by accidental removal of the mains adapter from the socket, the instrument must be reset by removing the plug from its socket and re-inserting the connections.
- Please consult your dealer if you have questions relating to the mains adapter.

4.3. Cuff connection

Insert the cuff tube into the opening provided on the left side of the instrument, as shown in the diagram.

4.4. Setting the time and date

- 1) After the new batteries are fitted, the year number flashes in the display. You can set the year by pressing the M-button. To confirm and then set the month, press the Time button.
- 2) You can now set the month using the M-button. Press the time button to confirm and then set the day.
- 3) Please follow the instructions above to set the day, hour and minutes.
- 4) Once you have set the minutes and pressed the Time button, the date and time are set and the time is displayed.
- 5) If you want to change the date and time, press and hold the Time button down for approx. 3 seconds until the year number starts to flash. Now you can enter the new values as described above.

5. CARRYING OUT A MEASUREMENT

5.1. Before the measurement

- Avoid eating, smoking as well as all forms of exertion directly before measurement. All these factors influence the measurement result. Try to find time to relax by sitting in an armchair in a quiet atmosphere for about 5 minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Take measurements regularly at the same time of day, since blood pressure changes during the course of the day.




5.2. Common sources of error

Note: Comparable blood pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support their arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an false reading will be obtained. Each 15cm difference in height results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is of extraordinary importance. The cuff size is dependent upon the circumference of the arm (measured in the centre). The permissible range is printed on the cuff. If this is not suitable for your use, please contact your dealer. **Note: Only use clinically approved Microlife Cuffs!**
- A loose cuff or a sideways protruding air pocket causes false measurement values.

5.3. Fitting the cuff

- Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and that the rubber tube leaves the cuff on the inner side of the arm.
- Tighten the free end of the cuff and close the cuff by affixing the Velcro.
- The cuff should be snug on your upper arm so that you can fit 2 fingers between the cuff and your upper arm. Any piece of clothing which restricts the arm (e.g. a pullover) must be taken off.
- Secure the cuff with the Velcro closer in such a way that it lies comfortably and is not too tight. Lay your arm on a table (palm upwards) so that the cuff is at the same height as the heart. Do not bend the tube!


Comment:

If it is not possible to fit the cuff to your left arm, it can also be placed on the right. However, all measurements should be made using the same arm.

5.4. Measuring procedure

After the cuff has been appropriately positioned, the measurement can begin:

- Press the START button. The pump begins to inflate the cuff. The rising pressure in the cuff is shown in the display.
- After the suitable inflation pressure has been reached, the pump stops and the pressure gradually falls. The cuff pressures are displayed. In case that the inflation pressure is not sufficient, the monitor automatically re-inflates to a higher level.


- c) When the instrument detects a pulse, the heart symbol in the display starts to flash and a beep is heard for every heartbeat.
- d) A longer beep is sounded when the measurement has been completed. The systolic, and diastolic blood pressures and pulse rate now appear in the display.
- e) The appearance of the PAD Symbol  means that a certain pulse irregularity was detected during the measurement. Please refer to item 2.4. of this manual and consult your doctor with the related information.
- f) The measurement readings remain on the display until you switch off the instrument. If no button is pressed for a period of 5 minutes the instrument switches itself off in order to preserve the batteries.

5.5. Discontinuing a measurement


If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the START button can be pressed at any time.

The device immediately lowers the cuff pressure automatically.

5.6. Memory - recall of measurements

This blood pressure monitor automatically stores the last 99 measurement values. By pressing the M-button, the last Measurement, as well as previous measurements can be displayed one after the other.

(MR30: Value of the last measurement)


(MR29: Value of the measurement before MR 30)

5.7. Memory– cancellation of all measurements


Attention!

Before you delete all readings stored in the memory, make sure you will not need refer to the readings at a later date.

Keeping a written record is prudent and may provide additional information for your doctor's visit.


In order to delete all stored readings, depress the M-button for at least 7 seconds, the display will show the symbol «CL» and 3 short beep sounds will be heard to indicate deletion of stored readings.

6. ERROR MESSAGES/MALFUNCTIONS

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.
(Example Error no. 1)

Error No.	Possible cause(s)
ERR 1	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. The tube may have become unplugged after the systolic blood pressure was measured. Further possible cause: No pulse could be detected
ERR 2	Unnatural pressure impulses impair the measurement result. Possible cause: The arm was moved during the measurement (artefact).
ERR 3	Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.
HI	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
LO	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

* Please consult your doctor, if this or any other problem occurs repeatedly.

Other possible malfunctions and possible solutions – If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched. Batteries have been installed.	<ol style="list-style-type: none">1. Check the polarity of the batteries (+/-)2. If the display is erratic or unusual, remove the batteries and re-install new batteries.
The pressure does not rise although the pump is running.	<ul style="list-style-type: none">• Check the connection of the cuff tube and connect properly if necessary.
The instrument frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none">1. Ensure the cuff fits correctly.2. Make sure that the cuff is not fitted too tightly. Ensure tight clothing such as a rolled-up sleeve isn't exerting pressure on the arm above the measuring position. Take off articles of clothing if necessary.3. Measure blood pressure again in complete peace and quiet.
Every measurement produce different value although the instrument functions normally and the values displayed are normal	<ul style="list-style-type: none">• Please read the following information and the points listed under «Common sources of error». Repeat the measurement.
Blood pressure values measured differ from those measured by the doctor.	<ul style="list-style-type: none">• Record the daily development of the values and consult your doctor.



Further information

The level of blood pressure is subject to fluctuations even in healthy people. It is important to compare measurements taken under the same conditions and at the same time of day. (Quiet conditions!)

For licensing, this device has been subjected to strict clinical tests and supervised by experienced Cardiac Specialists. If you have any questions regarding the use of this blood pressure monitor, please ask your dealer or pharmacist for the Microlife Service representative in your country. The Microlife Service Team will be happy to help you. **Never attempt to repair the instrument yourself!** Any unauthorized opening of the instrument invalidates all guarantee claims!

7. CARE AND MAINTENANCE, RECALIBRATION

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle the cuff carefully and avoid all types of straining through twisting or buckling.
- c) Handle the tube carefully. Avoid stress from pulling, pinching and sharp edges.
- d) Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. Carefully remove spots on the cuff with a damp cloth and soapsuds. The cuff **should not be washed in a clothes or dishwasher!**
- e) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the device!** This can negatively impact the calibration! (accuracy)

Periodical re-calibration

Sensitive measuring devices must be checked for accuracy from time to time. We recommend a periodical inspection of your unit by an authorized microlife dealer **every 2 years**.

Your factory authorized microlife dealer would be pleased to provide more information regarding calibration.

8. GUARANTEE

This blood pressure monitor is **guaranteed for 3 years** from date of purchase. This guarantee includes the instrument and the cuff. The guarantee does not apply to damage caused by improper handling, damage from leaking batteries, accidents, not following the operating instructions or alterations made to the instrument by third parties.

The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

Name and company address of the responsible dealer:

9. TECHNICAL SPECIFICATIONS

Weight:	460 g (with batteries)
Size:	115 (W) x 182 (L) x 76 (H) mm
Storage temperature:	-20 to +50°C
Humidity:	15 to 90% relative humidity maximum
Operation temperature:	10 to 40°C
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 beats per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Automatically stores the last 99 measurements
Measuring resolution:	1 mmHg
Accuracy:	Pressure within ± 3 mmHg
Pulse	± 5 % of the reading
Power source:	a) 4 x 1.5 V Batteries; size AA b) Mains adapter DC 6, V 600 mA (optional)
Cuffs:	M-size for arm circumference 22-32 cm or L-size for arm circumference 32-42 cm (optional)
Reference to standards:	EU directive 93/42/EEC; NIBP - requirements: EN 1060-1 /-3 /-4 & ANSI / AAMI SP10

Technical alterations reserved!

10. www.microlife.com

Detailed user information about our products as well as services can be found at www.microlife.com.

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