
Guarantee Card

IR 100


Hypertension

Human Fever

Asthma

Flexible Heating

Europe / Middle-East / Africa

Microlife AG
Espenstrasse 139
9443 Widnau / Switzerland
Tel. +41 / 71 727 70 30
Fax +41 / 71 727 70 39
Email admin@microlife.ch
www.microlife.com

Asia

Microlife Corporation.
9F, 431, RuiGang Road, NeiHu
Taipei, 114, Taiwan, R.O.C.
Tel. 886 2 8797-1288
Fax.886 2 8797-1283
Email service@microlife.com.tw
www.microlife.com

North / Central / South America

Microlife USA, Inc.
424 Skinner Blvd., Suite C
Dunedin, FL 34698 / USA
Tel. +1 727 451 0484
Fax +1 727 451 0492
Email msa@microlifeusa.com
www.microlife.com

CE0044

IM 04

IB IR 100-S EN-ES 1909


Microlife IR 100

EN → 2
ES → 8







microlife®

Microlife IR 100

Microlife IR 100

Microlife Ear Thermometer IR 100

EN

- ① Measuring Sensor
- ② START Button
- ③ Display
- ④ ON/OFF Button
- ⑤ Battery Compartment Cover
- ⑥ All Segments displayed
- ⑦ Memory
- ⑧ Ready for Measurement
- ⑨ Measurement complete
- ⑩ Out-of-ear Temperature Indication
- ⑪ Low Battery Indication

- ⑫ Changing between Celsius and Fahrenheit
- ⑬ Recall Mode
- ⑭ Recall the last 12 Readings
- ⑮ Measured Temperature too high
- ⑯ Measured Temperature too low
- ⑰ Ambient Temperature too high
- ⑱ Ambient Temperature too low
- ⑲ Error Function Display
- ⑳ Blank Display
- ㉑ Flat Battery
- ㉒ Replacing the Battery

The Microlife Ear Thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, this thermometer can provide a stable, heat-interference-free reading with each measurement. The instrument performs a self-test every time it is switched on to always guarantee the specified accuracy of measurements. The Microlife Ear Thermometer is intended for the periodic measurement and monitoring of human body temperature in the home. It is intended for use on people of all ages. **This thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.** Please read through these instructions carefully so that you understand all functions and safety information.

- Read the instructions carefully before using this device.
- Type BF applied part

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Specialist Dealer
Vendedor especializado

Table of Contents

1. The Advantages of this Ear Thermometer
2. Important Safety Instructions
3. How this Ear Thermometer measures Temperature
4. Control Displays and Symbols
5. Directions for Use
6. Changing between Celsius and Fahrenheit
7. How to recall 12 readings in Memory Mode
8. Error Messages
9. Cleaning and Disinfecting
10. Battery Replacement
11. Guarantee
12. Technical Specifications
13. www.microlife.com
Guarantee Card (see Back Cover)

1. The Advantages of this Ear Thermometer

Multiple Uses (Wide Range Measurement)

This thermometer offers a wide measurement range feature from 0 °C to 100.0 °C (32.0 °F to 212.0 °F), meaning the unit can be used as an ear thermometer to measure

body temperature or it can be used to measure surface temperature of the following:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

Probe Cover Free

This thermometer is more user-friendly and more cost-effective since a probe cover is not required.

Measurement in 1 second

The innovative infrared technology allows measurement of ear temperature in only 1 second.

Accurate and reliable

The unique probe assembly construction, incorporating an advanced infrared sensor, ensures that each measurement is accurate and reliable.

Gentle and Easy to Use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can even be used on a sleeping child, thereby causing no disruption.
- This thermometer is quick, therefore child-friendly.

Auto-Display Memory

The last reading is automatically displayed for 2 seconds when the unit is switched ON.

Multiple Reading Recall

Users will be able to recall the last 12 readings when entering the recall mode, enabling efficient tracking of temperature variations.

Safe and Hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Cleaning of the probe can be done with an alcohol-moistened cotton tissue, making this thermometer completely hygienic for use by the whole family.

Fever Alarm

10 short beeps alert the patient that he/she may have fever.

2. Important Safety Instructions

- This instrument may only be used for the purposes described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.

- **Never immerse this instrument in water or other liquids. For cleaning please follow the instructions in the «Cleaning and Disinfecting» section.**
- Do not use the instrument if you think it is damaged or notice anything unusual.
- Never open the instrument.
- Earwax in the ear canal may cause a lower temperature reading. Therefore it is important to ensure the subject's ear canal is clean.
- This instrument comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Protect it from:
 - extreme temperatures
 - impact and dropping
 - contamination and dust
 - direct sunlight
 - heat and cold
- If the instrument is not going to be used for a prolonged period the battery should be removed.
- Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.

WARNING: Use of this instrument is not intended as a substitute for consultation with your physician. This instrument is NOT waterproof so do not immerse in liquids.

3. How this Ear Thermometer measures Temperature

This thermometer measures infrared energy radiated from the eardrum and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. The measured reading obtained directly from the eardrum (Tympanic Membrane) ensures the most accurate ear temperature.

Measurements taken from the surrounding tissue of the ear canal generate lower readings and may result in misdiagnosis of fever.

To avoid an inaccurate measurement:

1. Switch on the thermometer by pressing the ON/OFF button ④.
2. After one beep is heard (and the temperature scale icon is flashing), straighten the ear canal by gently pulling the middle of the ear back and up.

3. Place the probe ① firmly into the ear canal, press the START button ② and keep the probe in the ear until the thermometer beeps to signal the completion of the measurement.

4. Control Displays and Symbols

- **All segments displayed** ⑥: Press the ON/OFF button ④ to turn on the unit; all segments will be shown for 2 seconds.
- **Memory** ⑦: The last reading will automatically be shown on the display for 2 seconds.
- **Ready for measurement** ⑧: The unit is ready for measurement, the «°C» or «°F» icon will keep flashing.
- **Measurement complete** ⑨: The reading will be shown on the display ③ with the «°C» or «°F» icon; the unit is ready for the next measurement, when the «°C» or «°F» icon is flashing again.
- **Out-of-ear temperature indication** ⑩: A crossed-ear-icon will appear on the display ③ if the reading falls outside the range 32.0 ~ 42.2 °C (89.6 ~ 108.0 °F).
- **Low battery indication** ⑪: When the unit is turned on, the battery icon will keep flashing to remind the user to replace the battery.

5. Directions for Use

1. Press the ON/OFF button (4). The display (3) is activated to show all segments for 2 seconds.
2. The last measurement reading will automatically be shown on the display for 2 seconds with the «M» icon (7).
3. When the «°C» or «°F» icon is flashing, a beep sound is heard and the thermometer is ready for measurement (8).
4. Straighten the ear canal by pulling the ear up and back to give a clear view of the eardrum.
 - For children under 1 year; pull the ear straight back.
 - For children of 1 year to adult; pull the ear up and back.Also refer to the short instructions at the front.
5. While gently pulling the ear, insert the probe snugly into the ear canal (**max. 1 sec.**) and **immediately** press the START button (2). Release the button and wait for the beep sound. This is the indication that confirms the end of measurement.
6. Remove the thermometer from the ear canal. The display shows the measured temperature (9).

NOTE:

- **To ensure accurate readings, wait at least 30 sec. after 3-5 continuous measurements.**

- Accumulation of ear wax on the probe can result in less accurate temperature readings or cross infection between users. **Therefore, it is essential that the probe is clean before each measurement.** For cleaning, follow the instructions in the «**Cleaning and Disinfecting**» section.
- **After cleaning the measuring sensor (1) with alcohol, wait 5 minutes before taking the next measurement,** in order to allow the thermometer to reach its operating reference temperature.
- 10 short beeps will sound when the temperature is higher than 37.5 °C (99.5 °F) in order to alert the patient that he/she may have fever.
- For an infant, it is best to have the child lying flat with his/her head sideways so the ear is facing upwards. For an older child or adult, it is best to stand behind and slightly to the side of the patient.
- Always take the temperature in the same ear, since the temperature readings may be different from ear to ear.
- Wait for a few minutes to take the ear temperature after sleeping.

- In the following situations it is recommended that three temperatures in the same ear be taken and the highest one taken as the reading:
 1. New born infants in the first 100 days.
 2. Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 3. When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.
 4. If the measurement is surprisingly low.

6. Changing between Celsius and Fahrenheit

This thermometer can display temperature measurements in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn OFF the unit, **press and hold** the START button (2) for 5 seconds; after 5 seconds, the current measurement scale («°C» or «°F» icon) will be flashing on the display (12). Change the measurement scale between °C and °F by pressing the START button (2). When the measurement scale has been chosen, wait for 5 seconds and the unit will automatically enter the «ready for measuring» mode.

7. How to recall 12 readings in Memory Mode

This thermometer can recall the last 12 readings.

- **Recall mode** ⑬: Press the START button ② to enter Recall mode when power is off. The memory icon «M» will flash.
- **Reading 1 - the last reading** ⑭: Press and release the START button ② to recall the last reading. Display 1 alone with memory icon.
- **Reading 12 - readings in succession:** Press and release the START button ② consecutively to recall the last 12 readings in succession.

Pressing and releasing the START button ② after the last 12 readings have been recalled will resume the above sequence from reading 1.

8. Error Messages

- **Measured temperature too high** ⑮: Displays «H» when measured temperature is higher than 100.0 °C or 212.0 °F.
- **Measured temperature too low** ⑯: Displays «L» when measured temperature is lower than 0 °C or 32.0 °F.

- **Ambient temperature too high** ⑰: Displays «H» in conjunction with the «▲» when ambient temperature is higher than 40.0 °C or 104.0 °F.
- **Ambient temperature too low** ⑱: Display «L» in conjunction with the «▼» when ambient temperature is lower than 5.0 °C or 41.0 °F.
- **Error function display** ⑲: The system has a malfunction.
- **Blank display** ⑳: Check if the battery has been loaded correctly. Also check polarity (<+> and <->) of the battery.
- **Dead battery indication** ㉑: If the steady battery icon is the only symbol shown on the display, the battery should be replaced immediately.

9. Cleaning and Disinfecting

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the probe lens and the display.

10. Battery Replacement

This instrument is supplied with one lithium battery, type CR2032. Replace with a new CR2032 battery when the flashing battery symbol appears on the display ㉑.

Remove the battery cover by sliding it in the direction shown. Remove the battery and replace with a new one ㉒.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

11. Guarantee

This instrument is covered by a **2 year guarantee** from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) confirming date of purchase or the receipt.

- The guarantee covers the instrument. Batteries and packaging are not included.
- Opening or altering the instrument invalidates the guarantee.
- The guarantee does not cover damage caused by improper handling, discharged batteries, accidents or non-compliance with the operating instructions.

Please contact Microlife-service.

12. Technical Specifications

Type: Ear Thermometer IR 100

Measurement range: 0 °C to 100.0 °C (32.0 °F to 212.0 °F)

Resolution: 0.1 °C / °F

Measurement accuracy: Laboratory:
±0.2 °C, 32.0 ~ 42.2 °C
(±0.4 °F, 89.6 ~ 108.0 °F)

Display: Liquid Crystal Display, 4 digits plus special icons

Acoustic:

- The unit is turned ON and ready for the measurement: 1 short beep
- Complete the measurement: 1 long beep
- System error or malfunction: 3 short beeps
- Fever alarm: 10 short beeps

Memory:

- Auto-Display the last measured temperature
- 12 readings recall in the Memory Mode

Operating temperature: 5 °C to 40 °C (41.0 °F to 104 °F)

Storage temperature: -25 °C to +55 °C (-13 °F to 131 °F)

Automatic Switch-off: Approx. 1 minute after last measurement has been taken.

Battery: CR2032 BATTERY (X1) 3V - at least 1000 measurements

Dimensions: 120 x 35 x 60 mm

Weight: 53 g (with battery), 50 g (w/o battery)

Reference to standards: Complies with ASTM E1965 requirements

Technical alterations reserved.

According to the Medical Product User Act a biennial technical inspection is recommended for professional users. Please observe the applicable disposal regulations.

13. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com.

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