microlife®

Instruction Manual

GoldTemp Thermometer – MT 1622

IB MT 1622 GT VarA 2804

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IMPORTANT SAFETY INSTRUCTIONS

- The thermometer is suitable exclusively for measurement of the body temperature!
- The minimum measurement time until the signalling tone [beep] must be maintained without exception!
- The thermometer contains small parts (battery, etc.) that can be swallowed by children. Therefore never leave the thermometer unattended to children.
- Protect the thermometer from hitting and striking!
- Avoid ambient temperatures above 60°C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer. The thermometer must be intact when immersed in liquid disinfectant.
- Warning:** The thermometer contains a long-life battery. For disposal, please comply with local regulations for special waste.
- The manufacturer recommends to verify the accuracy by an authorised laboratory every 2 years.

EN

TURNING ON THE THERMOMETER

To turn on the thermometer, push the button over the display window; a short signalling tone [beep] signals «thermometer ON.» A display test is performed. Then at an ambient temperature of less than 32°C, an «L» and a flashing «C» appear at the upper right of the display. The thermometer is now ready for use.

FUNCTION TEST

Proper functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccuracy), this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

USING THE THERMOMETER

During a measurement, the current temperature is displayed continuously. The «C» symbol flashes during this period of time. If the signalling tone sounds (beep-beep-beep 10 times) and the «C» is no longer flashing, this means that the measured increase in temperature is now less than 0,01 °C in 8 seconds and that the thermometer is ready for reading.

To prolong the life of the battery, turn it off after use by pressing the operating button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

STORAGE OF MEASURED VALUES

If you hold the button down for more than 3 seconds when turning the thermometer on, the maximum temperature stored automatically at the last measurement will be displayed. At the same time, a small «M» for memory will appear at the right of the display. After 3 more seconds, this value will disappear and the thermometer will switch to the normal measurement mode. Then the stored value is erased.

CLEANING AND DISINFECTING

Name: Isopropyl alcohol 70%; dilution with water: 95,0 %
Use a cotton swab moistened with alcohol (70% Isopropyl) to clean and disinfect the thermometer casing and the measuring probe.

BATTERY REPLACEMENT

When the «▼» symbol (upside-down triangle) appears in the lower right of the display, the battery is exhausted and needs replacing. For the battery replacement please remove the compartment lid of the thermometer. Replace the exhausted battery by positioning the + at the top. Make sure a battery of the same type is on hand. Batteries can be purchased at any electrical goods shops. For battery disposal please refer to «safety instructions».

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GARANTEE

We grant you a lifetime warranty after the date of purchase. Any damage caused by improper handling shall not be covered by the warranty. Batteries and packaging are also excluded from the warranty. All other damage claims excluded. A warranty claim must be submitted with the purchase receipt. Please pack your defective thermometer well and send with sufficient postage to:

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