SPRAGUE RAPPAORT ST-77

microlife®

Features

- Functions 5-way convertible purposes in one unit
- Furnish 3 different sizes of bell chest-piece
- Large/Small diaphragms for heart sounds of low/high frequency
- Adult/Medium/Infant size bells for L/M pitched sounds, murmurs and ribs detection
- Name card included

Design for all purposes, you will find these combinations are astonishingly well designed and their high sensitivities will enable you to clearly discern low and high frequencies and to hear the faintest murmurs with ease. Also spare eartips of 3 different sizes are furnished enabling the user to choose a pair that best fits him for comfortable use. They are packed in compact small case put with 3 bell type chest-pieces, and 2 spare transparent plastic diaphragms, large and small.

Sprague Rappaport Type Improved Stethoscope comes as a complete set with following components:

1. Chest-piece with large and small diaphragm..... 1 pc
2. Binaural, chrome plated with white plastic medium size eartips 1 pair
3. Heavy wall tubing with metal clip..... 1 pair
4. Interchangeable bell type chest-piece, adult size 1 pc
5. Interchangeable bell type chest-piece, medium size 1 pc
6. Interchangeable bell type chest-piece, infant size..... 1 pc
7. Spare eartips, large and small size 2 pairs
8. Spare transparent plastic diaphragms, large and small 2 pairs
9. Plastic container for accessories 1 pc

The Structure of The Chestpiece

Large Diaphragm: This is to detect heart sounds of low frequency, diastolic and the third heart sounds, as well as first and second heart sounds. The high pitched murmurs emitting from the heart is also detectable.

Small Diaphragm: This design is to hear infant heart sounds of high pitched nature.

Adult Size Bell: To pinpoint low and medium pitched sounds and murmurs this bell is very useful.

Medium Size Bell: Concentrates on narrow spots such as between ribs or any other small places to detect low and medium pitched sounds or murmurs.


Infant size Bell: Best suits the examination of infants and small children. Low and medium pitched sound and murmurs emitting from infant heart are best sensed by the bell.

How To Use

The Sprague Rappaport is dual headed and ready furnished with large and small size diaphragms fitted on revolvable twin head drum which contains leakage free stem valve made in hair thin precision.

The Operation Procedure Are as Below.

1. First of all, put the binaural tips into ears.
2. A soft tapping of diaphragm by your finger will let the sounds come through and you can tell if it is in working position.
3. If you don't hear the tapping sound, rotate the head 180° where it clicks settle in the ready position, facing opposite side.
4. Then, tap again and you will hear the tapping sounds which means the device is correctly set for use.
5. Now you can diagnose the examinee with the stethoscope.
6. Go through the same procedures when you switch from one chest-piece to another.

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