microlife®

BP A6 BT

Blood Pressure Monitor

EN →	1	SL →	76
RU →	13	SR →	88
BG →	27	HU →	100
RO →	40	HR →	112
CZ →	52	PL →	123
SK →	64		

Microlife Corporation
9F, No. 431, RuiGuang Road, Nei-Hu
Taipei, 114, Taiwan, China
www.microlife.com

CE REP Microlife UAB
P. Lukšio g. 32
08222 Vilnius
Lithuania

CE 1639

IB BP A6 BT 4G E-V11 2725
Revision Date: 2025-06-04
Issue Date: 2024-11-20


Preparation

Sit on a back-supported chair and keep your legs uncrossed. / Сядьте, не скрещивая ноги, на стул со спинкой. / Сядьте на стол с облегалка и не кръстосвайте краката си. / Așezați-vă într-un scaun cu spatar și nu încrucișați picioarele / Sediňte si na židli s opěrkou a nepřekřížujte nohy. / Sadnite si na stoličku s opierkou a neprekrížujte si nohy. / Usedite se na stol s podprtím hrbtom in nogami na tleh (ne prekrížajte jih). / Sedite na stolicu koja podupire leđa i nemojte prekrštati noge. / Üljön egy háttámlás székre és lábait ne rakja keresztbe! / Sjednite na stolicu koja podupire leđa i neka Vam noge ne budu prekrížene. / Ustądz na krzesle wspieranym plecami i nie rozstawiaj nóg.

Avoid thick or close-fitting garments on the upper arm. / Освободите плечо от плотной или плотно облегающей одежды. / Избягвайте дебели или плътно прилепващи дрехи върху горната част на ръката. / Evitați îmbracaminte groasă sau stâmta pe braț / Vyhněte se těsnému nebo přiléhavému oděvu na horní části paže. / Vyhněte sa tesnému alebo priliehavému odevu na hornej časti ramena. / Alzogibajte se tesnim oblačilom. / Izbegavajte usku odeću na nadlaktici. / Kerülje a vastag vagy szoros ruhákat a felkaron! / Izbjegavajte usku odjeću na nadlaktici. / Unikaj grubej lub dopasowanej odzieży na ramieniu.

Place the artery-mark on the cuff over your artery. / Поместите манжету так, чтобы значок артерии на манжете находился над артерией руки в сгибе локтя. / Поставете маркера на маншета върху артерията. / Plasați semnul artera al manșetei peste artera dumneavoastră / Umístěte manžetu se značkou tepny nad Vaši tepnou. / Namestite znak za arterijo, ki je na manšeti, na vašo arterijo. / Postavite oznaku za arteriju na manžetini preko Vaše arterije. / Helyezze a mandzsetta artériajelzőjét az artériája fölé! / Postavite oznaku arterije na manžeti preko Vaše arterije. / Umieść mankiecie na tętnicy.


Fit the cuff closely, but not too tight. / Поместите манжету так, чтобы она плотно, но не туго прилегалла. / Стегнете добре маншета, но не твърде много. / Fixați manșeta ferm dar nu prea stransă / Nasadite manžetu těsně, ale ne příliš těsně. / Nasadite manžetu tesne, ale nie príliš tesne. / Manšeto dobro namestite, vendar ne pretesno. / Dobro zategnite manžetnu, ali ne previse stegnuto. / Rögzítse a mandzsettát szorosan, de ne feszítse meg túl erősen! / Postavite manžetu blizu, ali ne pretjensno. / Zacišnij markiet dokladnie, jednak niezbyt silnie.

Position the cuff 1-2 cm above your elbow. / Расположите манжету на 1-2 см выше локтя. / Поставете маншета 1-2 см над лакътя. / Plasați manșeta cu 1-2 cm deasupra articulației / Manžetu umístěte 1-2 cm nad loket. / Manžetu umiestnite 1-2 cm nad laket. / Namestite manšeto 1-2 cm nad komolcem. / Postavite manžetnu 1-2 cm iznad Vašeg lakta. / A mandzsetta 1-2 centiméterrel legyen a könyökétől feljebb! / Postavite manžetu 1-2 cm iznad Vašeg lakta. / Ustaw mankiet na 1-2 cm powyżej łokcia.

Keep your arm still and do not speak during the measurement. / Держите руку неподвижно и не разговаривайте во время измерения. / Дръжте ръката си неподвижна и не говорете по време на измерването. / Mějte ruku uvolněnou a během měření nemlve. / Majte ruku uvolnenú a počas merania nerozprávajte. / Držite ruku uvolnioną a počas merania nerozprávajte. / Med merjenjem ne premikajte roke in ne govorite. / Držite ruku mimo i nemojte pričati tokom merenja. / Ne mozgassa a karját és ne beszéljen mérés közben! / Držite ruku mimo i nemojte pričati tijekom mjerenja. / Trzymaj rękę nieruchomo i nie mów podczas pomiaru.

1.

Avoid eating, bathing, smoking or caffeine (approx. 30 min).
Избегайте приема пищи, купания, курения или приема кофеиносодержащих продуктов (прибл. за 30 минут до измерения).
Избягавайте ядене, къпане, пушене или кофеин (около 30 минути).
Evitați sa mâncați, sa faceți baie, sa fumați, sau sa beți cafea (aproximativ 30 min).
Vyhněte se jídlu, koupání, kouření nebo kofeinu (přibližně 30 minut).
Vyhňte sa jedlu, kúpaniu, fajčeniu alebo kofeínu (približne 30 minút).
Izogibajte se hrani, kopanju in kajenju (približno 30 minut).
Izbegavajte hranu, kupanje, pušenje i kofein (približno 30 minuta).
Kerülje az evést, fürdést, dohányzást vagy koffeinfogyasztást (a megelőző kb. 30 percben)!
Izbjegavajte jesti, kupati se, pušiti i kofein (oko 30 min.)
Unikaj jedenia, kapieli, palenia lub kofeiny (około 30 minut).

2.

Avoid activity and relax for 5-10 min.
В течение 5-10 мин. оставайтесь в расслабленном состоянии и избегайте физической нагрузки.
Избягвайте активности и се отпуснете в продължение на 5-10 минути.
Evitați orice activitate și relaxați-vă 5-10 min.
Vyhněte se fyzické aktivitě a odpočívajte 5-10 minut.
Vyhňte sa fyzickej aktivity a odpočívajte 5-10 minút.
Izogibajte se fizični aktivnosti in počivajte 5-10 minut.
Izbegavajte fizičku aktivnost i opustite se tokom 5-10 minuta.
Kerülje a fizikai tevékenységet, pihenjen 5-10 percig!
Izbjegavajte aktivnost i opustite se 5-10 min.
Unikaj aktywności i zrelaksuj się przez 5-10 minut.

3.

Measure before medication intake.
Измеряйте до приема лекарства.
Измерете преди приема на лекарството ви.
Měření proveďte před přijetím léku.
Meranie uskutočnite pred užitím liekov.
Merjenje izvedite pred zaužitjem zdravil.
Obavite merenje pre uzimanja leka.
Mérjen a gyógyszer bevétele előtt!
Izvedite mjerjenje prije nego što uzmete lijek.
Zmierz przed przyjęciem leku.

Name of Purchaser / Ф.И.О. покупателя /
Име на купувача / Numele cumpărătorului / Imię
i nazwisko nabywcy / Jméno kupujícího / Meno
zákazníka / Ime in priimek kupca / Ime i prezime
kupca / Vásártó neve / Ime i prezime kupca

Serial Number / Серийный номер / Серийен
номер / Număr de serie / Numer seryjny /
Výrobní číslo / Výrobné číslo / Serijska številka /
Serijski broj / Sorozatszám / Serijski broj

Date of Purchase / Дата покупки / Дата на
закупуване / Data cumpărării / Data zakupu /
Datum nákupu / Dátum kúpy / Datum nakupa /
Datum kupovine / Vásárlás dátuma / Datum
kupovine

Specialist Dealer / Специализированный дилер /
Специалист дистрибутор / Distributor de spe-
cialitate / Przedstawiciel / Specializovaný dealer /
Špecializovaný predajca / Szpecializirani trgovec /
Ovlašteni diler / Forgalmazó / Ovlašteni prodavač

- ① START/STOP Button
- ② Display
- ③ M-button (memory)
- ④ Time button
- ⑤ Cuff socket
- ⑥ AFIB/MAM Switch
- ⑦ Traffic light indicator
- ⑧ USB Port
- ⑨ Mains Adapter Socket
- ⑩ Battery compartment
- ⑪ Cuff
- ⑫ Cuff connector
- ⑬ - «Backward» Button
- ⑭ + «Forward» Button
- ⑮ Lock switch

Display

- ⑯ Systolic value
- ⑰ Diastolic value
- ⑱ Pulse rate
- ⑲ Battery display
- ⑳ Cuff check indicator
- ㉑ Arm movement indicator
- ㉒ AFIB/MAM Interval time
- ㉓ Atrial Fibrillation Indicator (AFIB)
- ㉔ AFIB/MAM Mode
- ㉕ Stored value
- ㉖ Date/Time
- ㉗ Pulse indicator

Dear Customer,

This device was developed in collaboration with physicians and clinical tests carried out prove its measurement accuracy to be of a very high standard.*

Microlife AFIBsens is the world's leading digital blood pressure measurement technology for the detection of atrial fibrillation (AF) and arterial hypertension. These are the two top risk factors of getting a stroke or heart disease. It is important to detect AF and hypertension at an early stage, even though you may not experi-

ence any symptoms. AF screening in general and thus also with the Microlife AFIB algorithm, is recommended for people of 65 years and older. The AFIB algorithm indicates that atrial fibrillation may be present. For this reason, it is recommended that you visit your doctor when the device gives an AFIB signal during your blood pressure measurement. The AFIB algorithm of Microlife has been clinically investigated by several prominent clinical investigators and showed that the device detects patients with AFIB at a certainty of 97-100%.^{1,2}

If you have any questions, problems or want to order spare parts please contact your local Microlife-Customer Service. Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit the internet at www.microlife.com where you will find a wealth of invaluable information on our products.

Stay healthy – Microlife Corporation!

** This device uses the same measuring technology as the award winning «BP 3BTO-A» model tested according to the British and Irish Hypertension Society (BIHS) protocol.*

¹ Kearley K, Selwood M, Van den Briel A, Thompson M, Mant D, Hobbs FR et al.: *Triage tests for identifying atrial fibrillation in primary care: a diagnostic accuracy study comparing single-lead ECG and modified BP monitors. BMJ Open 2014; 4:e004565.*

² Wiesel J, Arbesfeld B, Schechter D: *Comparison of the Microlife blood pressure monitor with the Omron blood pressure monitor for detecting atrial fibrillation. Am J Cardiol 2014; 114:1046-1048.*

Table of contents

1. **Introduction**
 - Document scope
 - Disclaimers
2. **Important information**
 - Device description
 - Intended use
 - Intended user
 - Intended patient
 - Intended use environment and conditions
 - Indications
 - Contra-indications
 - Side effects

- Warning
- Caution
- Electromagnetic compatibility
- Adverse events and reporting
- 3. **Device information**
 - Device accessories
 - Low battery
 - Flat battery - replacement
- 4. **Device installation and setup**
 - Inserting the batteries
 - Setting the date and time
 - Selecting the correct cuff
 - Connecting the cuff to the device
 - Selecting standard or AFIB/MAM mode
 - AFIB/MAM mode
- 5. **Measurement preparation**
 - Before taking a measurement
 - Correct cuff fitting and posture for taking a measurement
- 6. **Measurement operation**
 - Starting measurement
 - How not to store a reading
 - Manual inflation
- 7. **Measurement interpretation**
 - How do I evaluate my blood pressure
 - Appearance of the Atrial Fibrillation Indicator for early Detection (Active only in MAM mode)
- 8. **Data memory function**
 - Viewing the stored values
 - Memory full
 - Clearing all values
- 9. **Bluetooth® Function**
- 10. **PC-Link functions**
- 11. **Device error and troubleshooting**
- 12. **Device maintenance and disposal**
 - Storage
 - Calibration and support
 - Disposal
- 13. **Specifications and compliance**
 - Technical specifications
 - Compliance information
- 14. **Supplement information for users and patients**
 - Guarantee
 - Symbols and definitions

1. Introduction

Document scope

Read the instructions carefully before using this device.

This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.

Disclaimers

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Microlife Corp. is under license. Other trademarks and trade names are those of their respective owners.

Apple, the Apple logo, iPad, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Android and Google Play are both trademarks of Google Inc. The Windows word mark and logos are registered trademarks of Microsoft Corporation and its affiliated companies.

Microlife® is a registered trademark of Microlife Corporation. Trademarks and trade names are those of their respective owners.

2. Important information

Device description

A digital home-use blood pressure monitor is a medical device that utilizes the principles of cuff-based oscillometric method and digital signal process to compute and provide a blood pressure measurement.

Intended use

This device is intended to measure brachial blood pressures (systole and diastole) and pulse rate with an optional function for detecting the presence of atrial fibrillation (AFIB) during measurement.

Intended user

The device is intended to be operated by adults and adolescents with adequate vision, motor functions, and education, capable of understanding the instructions for use and operating general household electrical appliances.

Intended patient

The intended patients are normotensive and hypertensive adults and adolescents (aged 12 years or older) of the general population.

The intended patients also include those with conditions of diabetes, pregnancy, pre-eclampsia, atherosclerosis.

Intended use environment and conditions

The device is intended for use in a home healthcare environment (e.g. general household without medically trained personnel) by patients (e.g. for self-measurement) or by a care giver.

Indications

This device measures blood pressures for indications of:

- Diagnosis of white-coat hypertension and masked hypertension and identifying white-coat effect and masked uncontrolled hypertension.
- Evaluate blood pressure in response to treatment.
- Confirming the diagnosis of resistant hypertension.
- Detecting morning hypertension.

The device is also indicated for screening the presence of irregular pulses suggestive of Atrial Fibrillation (AF).

Contra-indications

- The device is not intended for measuring blood pressure in pediatric patients of age younger than 12 years old (children, infant, or neonates).
- The device measures blood pressure using a pressured cuff. If the measuring limb suffers from injuries (for example open wounds) or under conditions or treatments (for example intravenous drip) making it unsuitable for surface contact or pressurization, do not use the device, to avoid worsening of the injuries or conditions.
- Avoid taking measurements of patients with conditions, diseases, and susceptible to environment conditions that lead to uncontrollable motions (e.g. trembling or shivering) and inability to communicate clearly (for example children and unconscious patients).
- The device uses oscillometric method to determine blood pressure and requires the measured limb with normal perfusion. The device is not intended to be used on a limb with restricted or impaired blood circulation. Consult with your doctor if you have severe perfusion or blood disorders before using the device.

Side effects

In rare cases, slight bruising may result after measurement due to pressurization of the arm.

Warning

NOTE: Warning items indicate potentially hazardous situations, if not avoided, may result in death, critical or serious injury to the user or patient.

- Avoid taking measurement on the arm on the side of a mastectomy or lymph node clearance.
- Avoid taking measurements on the arm with intravascular access or therapy or an arterio-venous (A-V) shunt. Cuff and pressurization may temporarily interfere with blood flow and could result in injury.
- Presence of significant cardiac arrhythmia during measurement may interfere with blood pressure measurement and affect the reliability of blood pressure readings. Consult with your doctor about whether the device is suitable for use in this case.
- DO NOT use this device in a moving vehicle (for example in a car or on an aircraft).
- DO NOT use this device for purposes beyond described in this instructions for use. The manufacturer cannot be held liable for damage caused by incorrect application.
- The measurement result of this device is not a medical diagnosis and not intended to substitute consultation and diagnosis by a qualified professional healthcare provider (e.g., physician, pharmacist, or other licensed health-care professionals).
- DO NOT use this device for self-diagnosis or for self-treatment of a medical condition. Seek advice from a health-care professional immediately if the patient is clearly unwell and/or having physiological or medical symptoms.
- Inspect the device, cuff and other parts for damage. DO NOT USE the device, cuff or parts if they appear damaged or operating abnormally.
- Blood flow of the arm is temporarily interrupted during measurement from cuff pressurization. Extended periods of cuff pressurization reduces peripheral circulation. Beware of signs (e.g tissue discoloration) of impeded peripheral circulation when taking prolonged or multiple measurements. It is recommended to rest between measurements. Abort measurement, loosen

the cuff (or disconnect the cuff and device) and rest to restore perfusion.

- DO NOT use this device in oxygen rich environment or near flammable gas.
- DO NOT use this device with other medical electrical (ME) equipment simultaneously. This may cause device malfunction or measurement inaccuracies.
- Use and store the device, cuff and parts in temperature and humidity conditions specified in the «Specifications and compliance». Usage and storage of the device, cuff and parts in conditions outside ranges given in the «Specifications and compliance» may results in device malfunction and the safety of usage.
- The atrial fibrillation (AF) screening function of this device is not meant to diagnose atrial fibrillation. Diagnosis of atrial fibrillation must be made by a cardiologist based on interpretation of electrocardiograms (ECG).
- Keep the device away from children and people incapable of operating the device. Beware of the risks of accidental ingestion of small parts and of strangulation with the cables and tubes of this device and accessories.
- DO NOT let children operate the device alone.

Caution

NOTE: Caution items indicate potentially hazardous situations, if not avoided, may result in minor or negligible injury to the user or patient, or damage to the property or environment.

- The device is not intended to measure pulse rate to check the frequency of a pacemaker.
- This device may incorrectly detect atrial fibrillation (AF) in people with pacemakers or defibrillators.
- DO NOT disassemble or attempt to service the device, accessory, and parts, during use or in storage. Access to the device internal hardware and software is prohibited. Unauthorized access and servicing of the device, during use or in storage, may compromise the safety and performance of the device.
- The device is intended only for measuring blood pressure on your upper arm. DO NOT measure other sites because the reading does not reflect your blood pressure accurately.
- When measuring patients of arm circumference of 50 cm or above, please ensure the cuff is fitted and secured tightly on the

patient's arm. Measurement errors may occur more frequently if the cuff is fitted loosely; it's recommended to re-fit and tighten the cuff, then re-attempt measurement in such case.

- After a measurement is completed, loosen the cuff, and rest the arm to restore limb perfusion, before taking another measurement.
- Avoid kinking, pressing, and moving of the cuff tube during device operation, as this affects reading reliability and may cause injury if the cuff pressurization is prolonged, and deflation interrupted.
- Use this device only with compatible accessories and parts from Microlife, including cuffs, connectors, and AC adapters. Using non-compatible accessories may compromise the safety and performance of the device.
- Protect the device and accessories from the following to avoid damaging the device:
 - water, other liquids, and moisture
 - extreme temperatures
 - impacts and vibrations
 - direct sunlight
 - contamination and dust
- This device is reusable. It is recommended to clean the device and the accessory before and after use if the device is dirty from use or after storage.
- Always use the arm cuff of range appropriate for the mid arm circumference of the patient (upper arm only).
- Stop using this device and cuff and consult with your doctor if you experience skin irritation or discomfort.
- DO NOT use this device, cuff, or parts after the expiration of its stated service life.

Electromagnetic compatibility

- This device is compliant with electromagnetic disturbances standard.

Further documentation in compliance with EN 60601-1-2 EMC standard is available from Microlife on www.microlife.com/electro-magnetic-compatibility.

- DO NOT use this device in proximity of equipment that may cause electromagnetic disturbance (EMD), such as high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, and computerized tomography (CT) scanners. This device is not certified for operation near these

equipments, which may cause device malfunction and measurement inaccuracies.

- DO NOT use this device close to strong electromagnetic fields and portable radio frequency communication devices (for example microwave oven and mobile devices). Keep a minimum distance of 0.3 m from such devices when using this device.

 **Caution:** The use of non-Microlife or non-compatible accessories may result in increased emissions or decreased immunity of the equipment or system.

Adverse events and reporting

Please report any serious incident, injury or adverse event that has occurred in relation to the device to the manufacturer/ European authorized representative (EC REP), and to the competent authority.

3. Device information

Package Contents

- 1 x Microlife BP A6 BT
- 1 x Microlife Soft Cuff M-L
- 1 x Instruction manual
- 1 x BP diary
- 1 x Cuff size marker
- 1 x Cuff holder
- 1 x Micro USB cable
- 1 x Storage bag
- 1 x Microlife AC-1024D (DC 6V, 600 mA) optional
- 4 x 1.5 V alkaline batteries; type AAA

 **CAUTION:** Inspect the device, cuff, and other parts for damage. DO NOT USE the device, cuff or parts if they appear damaged or operating abnormally.

Device accessories

Blood pressure cuffs

Microlife offers cuffs, covering a wide range of arm sizes.

Microlife Soft Cuff S	Range 17-22 cm
Microlife Soft Cuff M	Range 22-32 cm
Microlife Soft Cuff M-L	Range 22-42 cm
Microlife Soft Cuff L	Range 32-42 cm
Microlife Soft Cuff L-XL	Range 32-52 cm

Contact your local authorized Microlife distributor if the standard cuff of the device is not the correct size for your arm.

AC adapter

You can operate this device using the Microlife AC adapter model AC-1024D (DC 6V, 600 mA).

 **Warning:** Do not use the AC adapter if the adapter or the cable is damaged. If the device, adapter, or cable is damaged, turn off the power and unplug the AC adapter immediately.

 **Warning:** Only use the AC adapter with outlets of compatible voltage rating.

 **Warning:** Do not plug or unplug the AC adapter from the outlet with wet hands.

 **Warning:** Do not damage the AC Adapter. Handle the AC adapter with care. Avoid pulling, bending, and tempering of the adapter cable.

 **Warning:** Unplug the AC adapter before cleaning this device.

 **Warning:** The mains adapter is not waterproof. DO NOT pour or spray liquid on the mains adapter.

 **Note:** When device is powered by the AC adapter the batteries are not drained.

1. Plug the adapter jack into a suitable adapter socket (⑨). Check to ensure the adapter or cable are not damaged.
2. Plug the adapter plug into the mains socket.

Batteries

Use 4 new 1.5 V, size AAA alkaline batteries.

 **Caution:** Do not use expired batteries or mix new and used batteries together.

 **Caution:** Remove batteries if the device is not going to be used for a prolonged period.

You can also operate this device using rechargeable batteries.

 Only use «NiMH» type reusable batteries.

 Batteries must be removed and recharged when the flat battery symbol appears. They should not remain inside the device as they may become damaged (total discharge as a result of low use of the device, even when switched off).

 Batteries cannot be charged in the blood pressure monitor. Recharge batteries in an external charger and observe the information regarding charging, care and durability.

Low battery

When the batteries are approximately $\frac{3}{4}$ empty the battery symbol  will flash as soon as the device is switched on (partly filled battery displayed). Although the device will continue to measure reliably, you should obtain replacement batteries.

Flat battery - replacement

When the batteries are flat, the battery symbol  will flash as soon as the device is switched on (flat battery displayed). You cannot take any further measurements and must replace the batteries.

1. Open the battery compartment  on the bottom of the device.
2. Replace the batteries – ensure correct polarity as shown by the symbols in the compartment.
3. To set date and time, follow the procedure described in Section «Setting the date and time».

 The memory retains all values although date and time must be reset – the year number therefore flashes automatically after the batteries are replaced.

4. Device installation and setup

Inserting the batteries

Switch the lock switch  to «unlock» position. The battery compartment  is on the bottom of the device. Insert the batteries (4 x 1.5 V, size AAA), thereby observing the indicated polarity.

Setting the date and time

1. After the new batteries are fitted, the year number flashes in the display. You can set the year by pressing either the «+»   or the «-»   button. To confirm and then set the month, press the time button .
2. Press the «+»   or the «-»   button to set the month. Press the time button  to confirm and then set the day.
3. Follow the instructions above to set the day, hour and minutes.
4. Once you have set the minutes and pressed the time button, the date and time are set and the time is displayed.
5. If you want to change the date and time, press and hold the time button for approx. 3 seconds until the year number starts to flash. Now you can enter the new values as described above.

Caution: Make sure date and time settings are correct on the device. Incorrect settings results in misleading data and time records of the measurements.

Selecting the correct cuff

Check if the cuff size is suitable for the circumference of your upper arms. The upper arm circumference can be measured using a tape measure around the mid-point of the upper arm.

Please see cuff range in chapter «Device accessories».

Caution: Using an undersized or oversized cuff for measurement can result in inaccurate blood pressure values. Use the correctly sized cuff for measurement to ensure the readings are reliable.

Connecting the cuff to the device

Connect the cuff to the device by inserting the cuff connector  into the cuff socket  as far as it will go.

Note: A loose connection will result in inaccurate readings, and an error message («Err 3»).

Selecting standard or AFIB/MAM mode

This device enables you to select either standard (standard single measurement) or AFIB/MAM mode (automatic triple measurement). To select standard mode, slide the AFIB/MAM switch  on the side of the device downwards to position «1» and to select AFIB/MAM mode, slide this switch upwards to position «3».

AFIB/MAM mode

In AFIB/MAM mode, 3 measurements are automatically taken in succession and the result is then automatically analysed and displayed. Because blood pressure constantly fluctuates, a result determined in this way is more reliable than one produced by a single measurement. AF detection is only activated in AFIB/MAM mode.

- When you select the 3 measurements, the AFIB/MAM-symbol  appears in the display.
- The bottom, right hand section of the display shows a 1, 2 or 3 to indicate which of the 3 measurements is currently being taken.
- There is a break of 15 seconds between the measurements. A count down indicates the remaining time.
- The individual results are not displayed. Your blood pressure will only be displayed after all 3 measurements are taken.
- Do not remove the cuff between measurements.

- If one of the individual measurements was questionable, a fourth one is automatically taken.

5. Measurement preparation

Before taking a measurement

- ▶ Avoid heavy activity, eating or smoking immediately before the measurement.
- ▶ Empty your bladder prior to measurement.
- ▶ Sit down on a back-supported chair and relax for 5 minutes. Keep your feet flat on the floor and do not cross your legs.
- ▶ **Always measure on the same arm** (normally left). It is recommended that doctors perform double arm measurements on a patient's first visit in order to determine which arm to measure in the future. The arm with the higher blood pressure should be measured.

Correct cuff fitting and posture for taking a measurement

- ▶ Always ensure that the correct cuff size is used (marking on the cuff).
- ▶ Remove close-fitting garments from the upper arm. To avoid constriction, shirt sleeves should not be rolled up - they do not interfere with the cuff if they are laid flat.
- ▶ Fit the cuff closely, but not too tight.
- ▶ Make sure that the cuff is positioned 1-2 cm above the elbow.
- ▶ The **artery mark** on the cuff (ca. 3 cm long bar) must lie over the artery which runs down the inner side of the arm.
- ▶ Support your arm so it is relaxed.
- ▶ Ensure that the cuff is at the same height as your heart.

6. Measurement operation

Starting measurement

1. Slide the lock switch (15) down to the «unlock» position. Press the START/STOP button (1) to start measuring.
2. The cuff will now pump up automatically. Relax, do not move and do not tense your arm muscles until the measurement result is displayed. Breathe normally and do not talk.
3. When the correct pressure is reached, the pumping stops and the pressure falls gradually. If the required pressure was not reached, the device will automatically pump some more air into the cuff.
4. During the measurement, the pulse indicator (27) flashes in the display.

5. The result, comprising the systolic (16) and the diastolic (17) blood pressure and the pulse rate (18) are displayed. Note also the explanations on further display symbols in this booklet.
6. When the device has finished measuring, remove the cuff.
7. Switch off the device. (The monitor does switch off automatically after approx. 1 min.).

Caution: Remain still and do not move or talk during measurement. Motions caused by talking, moving, trembling and other vibrations may interfere with the measurement and affect the measurement accuracy!

Caution: You can stop the measurement at any time by pressing the ON/OFF button or open the cuff (e.g. if you feel uneasy or an unpleasant pressure sensation).

This monitor is specially tested for use in pregnancy and pre-eclampsia. When you detect unusual high readings in pregnancy, you should measure after a short while again (eg. 1 hour). If the reading is still too high, consult your doctor or gynecologist.

In pregnancy the AFIB symbol can be ignored.

How not to store a reading

As soon as the reading is displayed press and hold the ON/OFF button (1) until «M» (25) is flashing. Confirm to delete the reading by pressing the time button (4).

«CL» is displayed when the reading is deleted from the memory successfully.

Manual inflation

In case of high systolic blood pressure, it can be an advantage to set the pressure individually. Press the ON/OFF button after the monitor has been pumped up to a level of approx. 30 mmHg (shown on the display). Keep the button pressed until the pressure is about 40 mmHg above the expected systolic value – then release the button.

7. Measurement interpretation

How do I evaluate my blood pressure

The bars on the left-hand edge of the display (7) show you the range within which the indicated blood pressure value lies. Depending on the height of the bar, the readout value is either within the optimum (green), elevated (yellow), too high (orange) or dangerously high (red) range. The classification corresponds to

the 4 ranges in the table as defined by the international guidelines (ESH, ESC, JSH).

Range	Systolic	Diastolic	Recommendation
1. blood pressure normal	< 120	< 74	Self-check
2. blood pressure optimum	120 - 129	74 - 79	Self-check
3. blood pressure elevated	130 - 134	80 - 84	Self-check
4. blood pressure too high	135 - 159	85 - 99	Seek medical advice
5. blood pressure dangerously high	≥ 160	≥ 100	Urgently seek medical advice!

The higher value is the one that determines the evaluation.

Example: a blood pressure value of **140/80** mmHg or a value of **130/90** mmHg indicates «blood pressure too high».

NOTE: The blood pressure classification is a general guideline of blood pressure level at home, but diagnosis of hypertension should be made by a healthcare professional based on specific conditions of the patient. Consult with your doctor for questions about the interpretation and classification of your blood pressure values.

Appearance of the Atrial Fibrillation Indicator for early Detection (Active only in MAM mode)

This device is able to detect atrial fibrillation (AF). This symbol ②③ indicates that atrial fibrillation was detected during the measurement. Please refer to the next paragraph for information regarding the consultation with your doctor.

Information for the doctor on frequent appearance of the atrial fibrillation indicator

This device is an oscillometric blood pressure monitor that also analyses pulse irregularity during measurement. The device is clinically tested.

The AFIB Symbol ②③ is indicated if atrial fibrillation is detected during measurement.

If the AFIB symbol appears after having performed a full blood pressure measurement episode (triplicate measurements), the patient is advised to perform another measurement episode (triplicate measurements). If the AFIB symbol appears again, we recommend the patient to seek medical advice.

If the AFIB-symbol appears on the screen of the blood pressure monitor, it indicates the possible presence of atrial fibrillation.

The atrial fibrillation diagnosis however, **must** be made by a **cardiologist** based on ECG interpretation.

Caution: In the presence of atrial fibrillation the diastolic blood pressure value may not be accurate.

Caution: In the presence of atrial fibrillation using AFIB/MAM-mode is recommended for more reliable blood pressure measurement.

Information available on: www.microlife.com/afib

8. Data memory function

This device automatically stores the last **99** measurement values.

Viewing the stored values

Switch the lock switch ①⑤ to «unlock» position. Press the M-button ③ briefly. The display first shows «M» ②⑤ and «28A», which stands for the average of all stored values.

Pressing the «+» ①④ or the «-» ①③ button repeatedly enables you to move from one stored value to another. Press the M-button again to exit the memory mode.

Memory full

☞ Pay attention that the maximum memory capacity of 99 memories per user is not exceeded. **When the 99 memory is full, the oldest value is automatically overwritten with the 100th value.** Values should be evaluated by a doctor before the memory capacity is reached – otherwise data will be lost.

Clearing all values

If you are sure that you want to permanently remove all stored values, hold down the M-button (the device must have been switched off beforehand) until «CL» appears and then release the button. To permanently clear the memory, press the M-button while «CL» is flashing. **Individual values cannot be cleared.**

☞ **Cancel deletion:** press START/STOP button ① while «CL» is flashing.

9. Bluetooth® Function

This device can be used in conjunction with a smartphone running the «Microlife Connected Health+» App. The Bluetooth® connection is automatically active after the measurement has been completed.

1. Activate Bluetooth® on your smartphone.
2. Activate the Bluetooth® function on the device. Press the «+» button ⑭ while the device is in stand-by mode. The icon «bt» starts flashing on the display.

For more detailed information visit www.microlife.com/connect.

10. PC-Link functions

This device can be used in conjunction with a personal computer (PC) running the Microlife Blood Pressure Analyzer+ (BPA+) software. The memory data can be transferred to the PC by connecting the monitor via a cable.

If no download-voucher and cable is included download the BPA+ software from www.microlife.com/software and use a USB cable with a Mini-B 5 pin connector.

☞ During the connection, the device is completely controlled by the computer.

11. Device error and troubleshooting

If an error occurs during the measurement, the measurement is interrupted and an error message, e.g. «Err 3», is displayed.

Error	Description	Potential cause and remedy
«Err 1»	Signal too weak	The pulse signals on the cuff are too weak. Re-position the cuff and repeat the measurement.*
«Err 2» ②①	Error signal	During the measurement, error signals were detected by the cuff, caused for instance by movement or muscle tension. Repeat the measurement, keeping your arm still.
«Err 3» ②②	Abnormal cuff pressure	An adequate pressure cannot be generated in the cuff. A leak may have occurred. Check that the cuff is correctly connected and is not too loose. Replace the batteries if necessary. Repeat the measurement.
«Err 5»	Abnormal result	The measuring signals are inaccurate and no result can therefore be displayed. Read through the checklist for taking a reliable measurement and then repeat the measurement.*
«Err 6»	AFIB/MAM Mode	There were too many errors during the measurement in AFIB/MAM mode, making it impossible to obtain a final result. Read through the checklist for taking a reliable measurement and then repeat the measurement.*
«HI»	Pulse or cuff pressure too high	The pressure in the cuff is too high (over 299 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
«LO»	Pulse too low	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*
«FL»	Problem with Bluetooth® connection	If «FL» flashes on the display, the device failed to establish a connection with your smartphone. Make sure that Bluetooth® is activated and repeat the procedure.

** Please immediately consult your doctor, if this or any other problem occurs repeatedly.*

12. Device maintenance and disposal

Cleaning the device

The device can be cleaned when necessary (e.g., between uses by different patients).

Use a soft cloth, dry or wet with detergent, to gently wipe the exterior of the device remove dust or stains.

Cleaning the cuff

Use a soft cloth, dry or wet with mild detergent, to carefully wipe the cuff to remove dust or stains.

Caution: Do not wash the cuff in a washing machine or dishwasher!

Caution: Do not use fabric softener.

Caution: Do not dry tumble dry or iron the cuff cover

Storage

When not in use:

- Disconnect the cuff and parts from the device.
- Keep the device and accessories in a dry, cool place away from sunlight, with ambient conditions within the temperature and humidity ranges described in the «Specifications and compliance» section.
- Remove the batteries from the device if the device will not be used for an extended period.

Warning: Storing the device **unused** for an extended period without removing batteries increases the chance of battery fluid leakage, which may lead to device damage and skin irritation when in contact. If your eye or skin is exposed to battery fluid, wash the exposed part immediately with ample clean water. Consult a doctor if irritation or discomfort persists.

Calibration and support

The device is calibrated during manufacture. In general, it is recommended to have the device verified by your local designated Microlife device distributor every two years, or after mechanical impact, liquid ingress, and/or device malfunctions. For questions related to device measurement accuracy, please contact your local designated Microlife device distributor.

Caution: Do not attempt to service or calibrate the device and accessories yourself.

Disposal

This device is medical electrical equipment. Dispose this device and batteries in accordance with the Waste Electrical and Electronic Equipment (WEEE) directive and applicable local regulations. DO NOT dispose of the device and batteries with domestic or commercial waste.

13. Specifications and compliance

Technical specifications

NOTE: Technical specifications subjected to change without notice.

Device Type:	Digital non-invasive blood pressure monitor
Model number:	BP3GU1-7B
Reference number	BP A6 BT
Operating conditions:	10 - 40 °C / 50 - 104 °F 15 - 90 % relative maximum humidity 700 hPa – 1060 hPa
Storage and transport conditions:	-20 - +55 °C 15 - 90 % relative maximum humidity
Weight:	354 g (including batteries)
Dimensions:	160 x 80 x 32 mm
Measuring procedure:	oscillometric, corresponding to Korotkoff method: Phase I systolic, Phase V diastolic
Pressure resolution:	1 mmHg
Cuff pressure display range:	0 - 299 mmHg
Measurement range:	SYS: 60 - 255 mmHg DIA: 40 - 200 mmHg Pulse: 40 - 199 beats per minute
Static accuracy:	± 3 mmHg
Pulse accuracy:	± 5 % of the readout value
Wireless Communication:	Bluetooth® 4.0
Power source - internal:	4 x 1.5 V AAA batteries

Power source - external (optional):	AC Adapter model: Microlife AC-1024D Input: 100-240 V Output: 6.0 V, 0.6 A, 3.6 W
Ingress protection (IP) rating:	IP20: Protected against solid objects with a diameter of ≥ 12.5 mm. No protection against ingress of water.
Applied part type reference:	 Type BF
Service life - device:	5 years or 10000 measurements, whichever comes first
Service life - cuff:	2 years or 5000 measurements, whichever comes first
Battery lifetime:	approx. 400 measurements (1.5 V alkaline batteries; size AAA)

Compliance information

This device complies with the requirements of the Medical Device Regulation (EU)2017/745.

Compliant standards:

EN 60601-1
EN 60601-1-2
EN 60601-1-11
EN IEC 80601-2-30
EN ISO 81060-2

14. Supplement information for users and patients

Guarantee

This device is covered by a **5 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by using non-Microlife specified accessories or parts, incorrect application or non-compliance with the instruction for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).

- Accessories and wearing parts: Batteries, power adapter (optional).

The cuff is covered by a functional guarantee (bladder tightness) for 2 years.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website:

www.microlife.com/support

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

Symbols and definitions

	Medical device
	CE Marking of Conformity
	Importer
	Authorized representative in the European Community
	Manufacturer
	Country of manufacture (Date of manufacture if date printed next to symbol)
	Model number
	Reference number
	Serial number (YYYY-MM-DD-SSSSS; year-month-day-serial number)
	Lot number (YYYY-MM-DD; year-month-day)
	Unique Device Identifier
	Caution

General warning sign

Type BF applied part

Direct current

IP20

Protected against solid objects with a diameter of ≥ 12.5 mm. No protection against ingress of water.

Keep dry

Temperature limitation for operating **or** storage

Humidity limitation for operating **and** storage

Atmospheric pressure limitation

Read instructions for use before operating the device.

Dispose in accordance with waste electrical and electronic equipment (WEEE) directive.

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¹ Kearley K, Selwood M, Van den Bruel A, Thompson M, Mant D, Hobbs FR et al.: Triage tests for identifying atrial fibrillation in primary care: a diagnostic accuracy study comparing single-lead ECG and modified BP monitors. *BMJ Open* 2014; 4:e004565.

² Wiesel J, Arbesfeld B, Schechter D: Comparison of the Microlife blood pressure monitor with the Omron blood pressure monitor for detecting atrial fibrillation. *Am J Cardiol* 2014; 114:1046-1048.

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¹ *Kearley K, Selwood M, Van den Bruel A, Thompson M, Mant D, Hobbs FR et al.: Triage tests for identifying atrial fibrillation in primary care: a diagnostic accuracy study comparing single-lead ECG and modified BP monitors. BMJ Open 2014; 4:e004565.*

² *Wiesel J, Arbesfeld B, Schechter D: Comparison of the Microlife blood pressure monitor with the Omron blood pressure monitor for detecting atrial fibrillation. Am J Cardiol 2014; 114:1046-1048.*

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