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C E 0044

WHAT IS A BASAL TEMPERATURE?

A basal temperature is the normal body temperature of a healthy person immediately upon awakening after a restful nights sleep. More precisely, it is the body temperature measured under so called basal conditions (12 hours after eating, after a restful sleep, no exercise, no emotional excitement, normal room temperature). For optimum accuracy, your basal temperature should be taken at the same time every morning upon waking. Changes in basal temperature are related to changes in the female reproductive cycle. By taking your basal temperature daily and charting it, you can notice slight changes and patterns. These charts of your basal temperature, in addition to other information, can help you understand your menstrual cycle and ovulation. Recording and interpreting basal temperature changes is called the «**Basal Body Temperature Method**». The Basal Temperature Method can be used to help in family planning. Charting of temperature and prediction of ovulation can be used to time sexual intercourse during fertile days to aid in increasing the likelihood of pregnancy.

CHARTING YOUR BASAL TEMPERATURE

By charting your basal temperature daily, you can see a pattern and notice the changes that occur in every ovulation cycle. You can record your daily basal temperature on the enclosed charts. Keeping a charted record for at least three months will usually provide the information you need to determine your approximate day of ovulation during a regular menstrual cycle. Most women keep both the Microlife Basal Thermometer and Temperature Tracking Chart at their bedside, for convenience.

HOW TO TAKE ACCURATE BASAL TEMPERATURES

To obtain an accurate basal temperature, you must take your temperature when you first awake in the morning. For best results, this waking temperature should be taken at the same time each morning. Since any activity may tend to raise your temperature, do not get out of bed until you have taken your basal temperature. Postpone going to the bathroom, eating, smoking or drinking until after taking your basal temperature. Cold drinks lower mouth temperature; hot drinks, smoking and exercise cause higher readings. Using an electric blanket or heating pad can affect your basal temperature. If you use one, it should be kept at the same setting each night throughout the time you are taking your basal temperature.

Should I use the oral, rectal or vaginal method?

Since temperatures taken at different sites can vary, it is important to select one method and then take your temperature the same way every day. Use either the oral, vaginal or rectal method. Using more than one method will provide inaccurate results.

TAKING YOUR BASAL TEMPERATURE

Select oral, rectal or vaginal as your measurement based on professional recommendation.

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Oral Method

1. Push the on/off button located besides the display window. Hold the button down to recall the previous measured temperature. After releasing the button the thermometer will beep and the display will read: «88.88 °C» This is a function check indicating the unit is working properly.
2. The °C symbol will flash continuously when the unit is ready to measure temperature. Please note that the «▼» indicates a low battery. See instructions for replacing the battery before measuring temperature. Do not expect an accurate temperature reading if this low battery indicator is present. An «ERR» message in the window indicates your unit is not functioning properly.
 - Do not drink hot or cold fluids, exercise, smoke, or perform other activities that will raise or lower temperature readings when compared to your normal, average temperature. The mouth should remain closed up to 5 minutes before attempting a reading.
 - Place the sensor tip well under the tongue as indicated by the «✓» (fig. 1).
 - The peak temperature should be reached in approximately 60 seconds. Opening the mouth or improper placement of the probe tip could result in a longer time for a reading.
3. When the peak temperature is reached, the °C symbol will stop flashing, and the thermometer will beep repeatedly. **NOTE: It is not necessary to hear the beep to ensure the peak temperature has been reached. When the °C symbol stops flashing and the display reads a consistent temperature, the peak temperature has been reached.**
4. Read and record the temperature and time for reference. The reading will automatically be stored in the memory.
5. Push the on/off button to turn the unit off. If you do not turn the unit off, it will shut off automatically in approximately 10 minutes.


fig. 1

Rectal Method

Lie on your side with knees slightly bent. With one hand, gently slide the tip of the thermometer no more than 1/2 inch into the rectum. Once used rectally, the thermometer should not be used orally for sanitary reasons.

Vaginal Method

Please seek guidance from your health care professional. Taking vaginal temperatures is beyond the scope of this manual.

IMPORTANT SAFETY INSTRUCTIONS

- The thermometer is suitable exclusively for measurement of body temperature!
- The minimum measurement time until the signal tone (beep) must be maintained without exception.
- Protect the thermometer from hitting and striking!
- The thermometer contains small parts (battery, etc.) that can be swallowed by children. Therefore never allow children to use the thermometer unattended!
- The manufacturer recommends verifying the accuracy by an authorized laboratory every 2 years.

WARNING: Battery should not be charged or placed into extreme heat as it may explode.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

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
your ovary releases an egg. This helps identify the best time to increase your likelihood of becoming pregnant. Note Jane timed intercourse with the rise in temperature on day #6 and #7 as indicated by a «•». Jane's temperature dropped back to 36.38 °C the day following ovulation and remained consistent through day 14. On day 15 through 28, Jane's temperature rose and remained elevated. This second rise during her cycle is likely related to our body's release of the hormone progesterone. This hormone prepares your uterus for a fertilized egg. If you don't become pregnant, the level of progesterone falls so your temperature does too and your cycle begins again with the start of your period as Jane's did on day 30 as indicated by the «X».

CHANGING THE BATTERY

When the «▼» appears in the lower right of the display, the battery is exhausted and needs replacing. Make sure a battery of the same type is on hand and proceed as follows: Remove battery compartment lid. Use a small screwdriver to remove battery. Insert new equivalent battery type (1.55V, LR41 or SR41). Ensure battery is correctly positioned with the «+» facing up or towards the back of the unit. Carefully replace the battery compartment lid. Ensure «O» Ring (washer) is in place to maintain water resistance.

IMPORTANT INFORMATION TO HELP YOU UNDERSTAND THE BASAL TEMPERATURE METHOD

Hormonal changes associated with ovulation produce a slight rise in a woman's basal body temperature, (approximately 0.3°C/0.5°F), as well as other recognizable signs such as changes in cervical mucus. Learning to record and interpret basal temperature changes is known as the Basal Temperature Method. By monitoring the basal temperature changes, couples, in conjunction with their physician or family planning counselor, may be able to determine when ovulation takes place. Monitoring other symptoms, such as cervical mucus changes, can assist in determining the fertile days prior to ovulation. When temperature changes are used in conjunction with other fertility symptoms, such as the cervical mucus, this is known as the Sympto-Thermal Method. Determination of ovulation can be used to time sexual intercourse during fertile days to aid in the likelihood of achieving pregnancy.

THE FEMALE REPRODUCTIVE CYCLE

In order to understand the Basal Temperature Method it is helpful to review the basics of human female reproductive activity.

A. The role hormones play in ovulation

At the beginning of each menstrual cycle, the pituitary gland secretes a hormone (FSH or Follicle Stimulating Hormone), that acts upon the ovaries to prepare for ovulation. As one or more follicles in an ovary prepare for ovulation, they secrete another hormone called estrogen, that prepares the uterus for pregnancy by causing the inner lining of the uterus to thicken. Estrogen also causes some glands in the cervix to secrete a mucus discharge that is necessary for fertility. This mucus changes the environment in the vagina to make it more favorable to sperm life and provides a medium in which the sperm move upward into the uterus. The mucus typically starts as a sticky substance. During the time of fertility, it becomes more liquid and plentiful usually attaining a consistency very similar to raw egg whites. Sometimes the mucus becomes so watery

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that it will no longer stretch like raw egg whites, but continues to produce sensations of lubrication or wetness at the outer lips of the vagina.

B. What is ovulation?

At ovulation, the ovarian follicle releases an egg which is picked up by one of the Fallopian tubes and transported toward the uterus. In order for pregnancy to occur, the egg must be fertilized in the fallopian tube by sperm within 24 hours after ovulation. After 24 hours the egg disintegrates and cannot be fertilized.

C. After ovulation?


After ovulation, the ovarian follicle that released the egg secretes another hormone called progesterone. This hormone suppresses further ovulation. The progesterone continues to build the vaginal lining and maintains it after ovulation. It also causes the basal temperature to rise and the cervical mucus to thicken so much that it appears to dry up. Ovarian progesterone secretion continues for approximately 14 days. When it stops, the inner lining of the uterus can no longer be maintained, so it is discharged off in a process called menstruation (your period), and the menstrual cycle starts again.

D. How to determine the date of ovulation?

By carefully taking basal temperatures upon waking and recording them on a chart, the characteristic rise in temperature after ovulation can be observed. The Basal Digital Thermometer makes it easy to read and record the slight temperature variations which occur when ovulation has taken place.

E. Temperature ranges (oral) before and after ovulation

Your basal temperature before ovulation is usually in the range of 35.8 - 36.7 °C, (96.5 - 98 °F). The most frequent range is between 36.3 - 36.5 °C (97.4 - 97.7 °F). After ovulation the basal temperature rises about 0.3°C (0.5°F) and stays in that higher range until the beginning of the next menstrual period. The highest probability of conception is during the days immediately prior to ovulation as well as on the day of ovulation itself. Refer to the chart below for Fahrenheit/Celsius conversions.

°F	95.0	95.5	96.0	96.5	97.0	97.5
°C	35.0	35.3	35.6	35.8	36.1	36.4

F. Other important information

Remember that the thermometer is designed only to help determine when ovulation has occurred. Mucus observations, as noted in «The female productive cycle», can also be helpful when recorded as shown on the sample chart. Notation of these changes, in addition to your temperature pattern, will assist your physician or counselor in providing you with proper guidance. It is commonly suggested that couples who have had difficulty in achieving pregnancy should not have intercourse more than once a day nor on consecutive days, in order to maintain satisfactory quantities of sperm. The most fertile time continues up through the first day of drying up of the mucus or the first day of temperature elevation, whichever comes later. The temperature record should be maintained on a daily basis. If pregnancy is achieved, the basal temperature will remain elevated for several months. A temperature pattern that remains elevated for 21 days (a week or more than the usual elevated pattern) is a good indicator that pregnancy has been achieved. If you do not become pregnant after 6 months, see your doctor. If you believe

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you are pregnant, see your doctor immediately.

AVOIDING PREGNANCY WITH NATURAL FAMILY PLANNING

In natural family planning, there are several ways of determining the end of pre-ovulation fertility and the beginning of post-ovulation infertility. Professionals teaching natural family planning claim most unplanned pregnancies occur during pre-ovulation infertility. This brochure does not include guidelines for natural family planning, since it is beyond the scope of this publication. Your basal thermometer only helps you determine when ovulation has occurred. Calculating infertile periods involves other bodily signs in addition to temperature change and requires professional guidance when avoidance of pregnancy is your objective. **NOTE: This product is not for contraceptive use!** This device is to be used to aid in conception and is not to be used for contraceptive purposes. If you do not become pregnant after 6 months, see your doctor.

CAUTION

The effectiveness of the methods outlined for seeking and avoiding pregnancy depend on factors such as a woman's individual reactions, the accuracy with which daily temperatures and other symptoms are recorded, and the absence of outside factors which may affect body temperature. Consequently, no representation of any kind can be made as to the reliability of the results obtained in any individual case, and the patient is urged to consult with and be guided by the advice of her physician or counselor in the use of this thermometer.

TECHNICAL SPECIFICATIONS

Measuring range: 32.00 - 42.99 °C (90.00 - 109.99 °F)
Resolution: 0.01 °C (0.01 °F)
Accuracy: ± 0.1 °C (± 0.2 °F)
Operating temp.: 10 °C to 40 °C
Storage temp.: 15 - 95 % relative maximum humidity
-25 °C to +60 °C
Power Source: LR 41x1 / 1.5V or SR41x1 / 1.55V
Beep Function: a. 1 Beep when power on or off
b. 10 beeps after measurement completed
Memory function: Press and hold the button for more than 2 sec., the unit will recall the last reading.
Auto-off function: ± 10 min. after measurement is completed
Low battery indication: Displays the symbol «▼» to remind user to replace the battery.
No measurement is available before the battery is replaced.
Self-test function: Displays «ERR» when system malfunctions
Response time: 60 sec. for oral measurements (reference)
Reference to standards: EN 12470-3; ASTM E1112; IEC 60601-1; IEC 60601-1-2 (EMC)
This device complies with the requirements of the Medical Device Directive 93/42/EEC. Technical alterations reserved.

GUARANTEE

We grant you a lifetime guarantee after the date of purchase. Any damage caused by improper handling shall not be covered by the guarantee. Batteries and packaging are also excluded from the guarantee. All other damage claims excluded. A guarantee claim must be submitted with the purchase receipt. Please pack your defective instrument well and send with sufficient postage to the Microlife distributor.

**Europe / Middle-East / Africa
Europa / Medio Oriente / África**
Microlife AG, Espenstrasse 139,
9443 Widnau / Switzerland
Email admin@microlife.ch
www.microlife.com

Microlife USA, Inc.
424 Skinner Blvd.
Suite B, Dunedin
FL 34698 / USA
Email msa@microlifeusa.com
www.microlife.com

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