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Measuring sensor  
Embout thermosensible

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### IMPORTANT SAFETY INSTRUCTIONS

EN

- The thermometer is suitable exclusively for measurement of the body temperature!
- The battery cover is sonic-welded with the main housing of the thermometer; the battery is well sealed inside the thermometer to reduce the risk to be swallowed by children.
- Protect the thermometer from hitting and striking!
- Avoid ambient temperatures above 60 °C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer. The thermometer must be intact when immersed in liquid disinfectant.
- The manufacturer recommends to verify the accuracy by an authorised laboratory every 2 years.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

### TURNING ON THE THERMOMETER

To turn on the thermometer, push the button on top; the display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C, an «L» and a flashing «°C» appear at the upper right of the display. The thermometer is now ready for use.

### USING THE THERMOMETER

Choose the preferred measuring method. During a measurement, the current temperature is displayed continuously. The «°C» symbol flashes during this period of time. If the «°C» is no longer flashing, this means that the measured increase in temperature is now less than 0,1 °C in 16 seconds and that the thermometer is ready for reading.

To prolong the life of the battery, turn it off after use by pressing the operating button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

### STORAGE OF MEASURED VALUES

If you hold the button down for more than 3 seconds when turning the thermometer on, the maximum temperature stored automatically at the last measurement will be displayed. At the same time, a small «M» for memory will appear at the right of the display. After 3 more seconds, this value will disappear and the thermometer will switch to the normal measurement mode. Then the stored value is erased.

### MEASURING METHODS

#### In the mouth (oral)

Position the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue. The measuring sensor must be in good contact with the tissue. Close your mouth and breath evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air.

Approx. measuring time: 1 minute.

#### In the anus (rectal)

This is the most reliable measuring method, and is especially appropriate for infants and small children. Carefully insert the tip of the thermometer 2 to 3 cm into the anal aperture.

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## In the armpit (axillary)

From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required. Approx. measuring time: 5 to 10 minutes!

## CLEANING AND DISINFECTING

The thermometer is 100 % waterproof and may be immersed in water or disinfectant solution for cleaning with no problem. Then wipe it off with a moist towel. Please use one of the following disinfectants for disinfecting, following the manufacturer's instructions in each case.

Name of disinfectant	Manufacturer	Dilution with water	Max. duration Method
Isopropyl alcohol 70%			Immerse for 24 hours
Gigasept FF	Schilke & Mayr GmbH	5.0 %	Immerse for 24 hours
Lysoformin 3000	Rosemann GmbH	6.0 %	Immerse for 24 hours
Aseptisol	Bode Chemie	4.0 %	Immerse for 24 hours

**Note:** Any other cleaning solution, dilution or cleaning method can cause a malfunction of the thermometer!

## BATTERY EXHAUSTION

When the «▼» symbol (upside-down triangle) appears in the lower right of the display, the battery is exhausted and thermometer needs to be replaced. For battery disposal please refer to «safety instructions».

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## GUARANTEE

We grant you a lifetime guarantee after the date of purchase. Any damage caused by improper handling shall not be covered by the guarantee. Batteries and packaging are also excluded from the guarantee. All other damage claims excluded. A guarantee claim must be submitted with the purchase receipt. Please pack your defective thermometer well and send with sufficient postage to:

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