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WHAT IS A BASAL TEMPERATURE?

A basal temperature is the normal body temperature of a healthy person immediately upon awakening after a restful nights sleep. More precisely, it is the body temperature measured under so called basal conditions (12 hours after eating, after a restful sleep, no exercise, no emotional excitement, normal room temperature). For optimum accuracy, your basal temperature should be taken at the same time every morning upon waking. Changes in basal temperature are related to changes in the female reproductive cycle. By taking your basal temperature daily and charting it, you can notice slight changes and patterns. These charts of your basal temperature, in addition to other information, can help you understand your menstrual cycle and ovulation. Recording and interpreting basal temperature changes is called the «**Basal Body Temperature Method**». The Basal Temperature Method can be used to help in family planning. Charting of temperature and prediction of ovulation can be used to time sexual intercourse during fertile days to aid in increasing the likelihood of pregnancy.

CHARTING YOUR BASAL TEMPERATURE

By charting your basal temperature daily, you can see a pattern and notice the changes that occur in every ovulation cycle. You can record your daily basal temperature on the enclosed charts. Keeping a charted record for at least three months will usually provide the information you need to determine your approximate day of ovulation during a regular menstrual cycle. Most women keep both the Microlife Basal Thermometer and Temperature Tracking Chart at their bedside, for convenience.

HOW TO TAKE ACCURATE BASAL TEMPERATURES

To obtain an accurate basal temperature, you must take your temperature when you first awake in the morning. For best results, this waking temperature should be taken at the same time each morning. Since any activity may tend to raise your temperature, do not get out of bed until you have taken your basal temperature. Postpone going to the bathroom, eating, smoking or drinking until after taking your basal temperature. Cold drinks lower mouth temperature; hot drinks, smoking and exercise cause higher readings. Using an electric blanket or heating pad can affect your basal temperature. If you use one, it should be kept at the same setting each night throughout the time you are taking your basal temperature.

Should I use the oral, rectal or vaginal method?

Since temperatures taken at different sites can vary, it is important to select one method and then take your temperature the same way every day. Use either the oral, vaginal or rectal method. Using more than one method will provide inaccurate results.

TAKING YOUR BASAL TEMPERATURE

Select oral, rectal or vaginal as your measurement based on professional recommendation.

Oral Method

1. Push the on/off button located besides the display window. Hold the button down to recall the previous measured temperature. After releasing the button the thermometer will beep and the display will read: «88.88 °C». This is a function check indicating the unit is working properly.
2. The °C symbol will flash continuously when the unit is ready to measure. Please note that the «▼» indicates a low battery. See

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PL

instructions for replacing the battery before measuring temperature. Do not expect an accurate temperature reading if this low battery indicator is present. An «ERR» message in the window indicates your unit is not functioning properly.

- Do not eat or drink anything hot or cold 10 minutes before the measurement. The mouth should remain closed up to 2 minutes before starting a reading.
 - Place the sensor tip well under the tongue as indicated by the «✓» (fig. 1).
 - The peak temperature should be reached in approximately 60 seconds. Opening the mouth or improper placement of the probe tip could result in a longer time for a reading.
3. When the peak temperature is reached, the °C symbol will stop flashing, and the thermometer will beep repeatedly. **NOTE: It is not necessary to hear the beep to ensure the peak temperature has been reached. When the °C symbol stops flashing and the display reads a consistent temperature, the peak temperature has been reached.**
 4. Read and record the temperature and time for reference. The reading will automatically be stored in the memory.
 5. Push the ON/OFF button to turn off the unit. The unit will shut off automatically in approx. 10 minutes if you do not turn it off.


fig. 1 / Rys.1

Rectal Method

Lie on your side with knees slightly bent. With one hand, gently slide the tip of the thermometer 2 to 3 cm into the rectum. Once used rectally, the thermometer should not be used orally for sanitary reasons.

Vaginal Method

Please seek guidance from your health care professional. Taking vaginal temperatures is beyond the scope of this manual.

IMPORTANT SAFETY INSTRUCTIONS

- Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.
- The instrument may be used only for measuring body temperature!
- Do not use this device if you think it is damaged or notice anything unusual.
- We recommend cleaning this device according to the cleaning instructions before first use for personal hygiene.
- Do not attempt rectal measurements on persons with rectal disorders. Doing so may aggravate or worsen the disorder.
- The minimum measurement time until the beep is heard must be maintained without exception!
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.
- The function of this device may be compromised when used close to strong electromagnetic fields such as mobile phones or radio installations and we recommend a distance of at least 1 m. In cases where you suspect this to be unavoidable, please verify if the device is working properly before use.
- Protect the instrument from impact and dropping!
- Avoid ambient temperatures above 60 °C. NEVER boil the instrument!
- Use only the commercial disinfectants listed in the section «Cleaning and Disinfecting» to clean the device to avoid damage to the device.
- We recommend this device is tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact your local MicroLife-Service to arrange the test.

WARNING: Battery should not be charged or placed into extreme heat as it may explode.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

RECORDING BASAL TEMPERATURE

Please refer to the sample recording and blank temperature

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tracking charts enclosed. Prepare your chart following the example. Day 1 for each cycle is the first day of menstruation / «Period». Above the pre-printed days of the cycle enter corresponding calendar dates. During your menstrual flow, it is not necessary to record temperature (unless you have short cycles). Instead, mark an X in each space provided to indicate the number of days your period lasts. Begin recording your basal temperature the first morning following your period.

For each daily temperature, trace down the column below the cycle day until it intersects with your recorded temperature. Place a dot at the intersection. A notation should be made for any rise in temperature due to illness, emotional stress, or sleeplessness. Intercourse can be indicated by placing a circle around the dot of the appropriate day. Chart your readings for at least three months in order to be sure that you have determined your approximate time of ovulation during any regular menstrual cycle. Careful temperature recording will help establish your particular pattern. Completed charts should then be presented to your physician or family planning counselor, who will assist you in accurately interpreting your results.

HOW TO INTERPRET THE BASAL BODY TEMPERATURE (BBT) CURVE

Let's follow the sample BBT curve on this page to help better interpret typical results. «Jane Smith» has noted the last 4 days of her period with «X». Day #4 was Jane's last day of menstrual flow so on day #5 she recorded her temperature as 36.33 °C and made a «*» mark on the chart. On day #6 Jane's temperature rose only 0.04 °C to 36.37 °C. On day #7 her temperature rose to 36.67 °C, a full 0.3 °C increase in one day. This increase in temperature can be caused by hormonal changes that occur in your body 12-24 hours before your ovary releases an egg. This helps identify the best time to increase your likelihood of becoming pregnant. Note Jane timed intercourse with the rise in temperature on day #6 and #7 as indicated by a «•». Jane's temperature dropped back to 36.38 °C the day following ovulation and remained consistent through day 14. On day 15 through 28, Jane's temperature rose and remained elevated. This second rise during her cycle is likely related to our body's release of the hormone progesterone. This hormone prepares your uterus for a fertilized egg. If you don't become pregnant, the level of progesterone falls so your temperature does too and your cycle begins again with the start of

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

your period as Jane's did on day 30 as indicated by the «X».

CHANGING THE BATTERY

When the «▼» appears in the lower right of the display, the battery is exhausted and needs replacing. Make sure a battery of the same type is on hand and proceed as follows: Remove battery compartment lid. Use a small screwdriver to remove battery. Insert a new equivalent battery type LR41 (1.5V) / SR41 (1.55V). Ensure the battery is correctly positioned with the «+» facing up or towards the back of the unit. Carefully replace the battery compartment lid. Ensure the «O» Ring (washer) is in place to maintain water resistance.

IMPORTANT INFORMATION TO HELP YOU UNDERSTAND THE BASAL TEMPERATURE METHOD

Hormonal changes associated with ovulation produce a slight rise in a woman's basal body temperature, (approximately 0.3°C/0.5°F), as well as other recognizable signs such as changes in cervical mucus. Learning to record and interpret basal temperature changes is known as the Basal Temperature Method. By monitoring the basal temperature changes, couples, in conjunction with their physician or family planning counselor, may be able to determine when ovulation takes place. Monitoring other symptoms, such as cervical mucus changes, can assist in determining the fertile days prior to ovulation. When temperature changes are used in conjunction with other fertility symptoms, such as the cervical mucus, this is known as the Sympto-Thermal Method. Determination of ovulation can be used to time sexual intercourse during fertile days to aid in the likelihood of achieving pregnancy.

THE FEMALE REPRODUCTIVE CYCLE

In order to understand the Basal Temperature Method it is helpful to review the basics of human female reproductive activity.

A. The role hormones play in ovulation

At the beginning of each menstrual cycle, the pituitary gland secretes a hormone (FSH or Follicle Stimulating Hormone), that acts upon the ovaries to prepare for ovulation. As one or more follicles in an ovary prepare for ovulation, they secrete another hormone called estrogen, that prepares the uterus for pregnancy by causing the inner lining of the uterus to thicken. Estrogen also causes some glands in the cervix to secrete a mucus discharge that is necessary for fertility. This mucus changes the environment in the vagina to make it more favorable to sperm life and provides a medium in which the sperm move upward into the uterus. The mucus typically starts as a sticky substance. During the time of fertility, it becomes more liquid and plentiful usually attaining a consistency very similar to raw egg whites. Sometimes the mucus becomes so watery that it will no longer stretch like raw egg whites, but continues to produce sensations of lubrication or wetness at the outer lips of the vagina.

B. What is ovulation?

At ovulation, the ovarian follicle releases an egg which is picked up by one of the Fallopian tubes and transported toward the uterus. In order for pregnancy to occur, the egg must be fertilized in the fallopian tube by sperm within 24 hours after ovulation. After 24 hours the egg disintegrates and cannot be fertilized.

C. After ovulation?

After ovulation, the ovarian follicle that released the egg secretes another hormone called progesterone. This hormone suppresses further ovulation. The progesterone continues to build the vaginal lining and maintains it after ovulation. It also causes the basal temperature to rise and the cervical mucus to thicken so much that it appears to dry up. Ovarian progesterone secretion continues for approximately 14 days. When it stops, the inner lining of the uterus can no longer be maintained, so it is discharged off in a process called menstruation (your period), and the menstrual cycle starts again.

D. How to determine the date of ovulation?

By carefully taking basal temperatures upon waking and re-cording them on a chart, the characteristic rise in temperature after ovulation can be observed. The Basal Digital Thermometer makes it easy to read and record the slight temperature variations which occur when

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ovulation has taken place.


E. Temperature ranges (oral) before and after ovulation

Your basal temperature before ovulation is usually in the range of 35.8 - 36.7 °C / 96.5 - 98 °F. The most frequent range is between 36.3 - 36.5 °C / 97.4 - 97.7 °F. After ovulation the basal temperature rises about 0.3°C / 0.5°F and stays in that higher range until the beginning of the next menstrual period. The highest probability of conception is during the days immediately prior to ovulation as well as on the day of ovulation itself.

Refer to the chart below for Fahrenheit/Celsius conversions.

F. Other important information

Remember that the thermometer is designed only to help determi-

°F	95.0	95.5	96.0	96.5	97.0	97.5	98.0	98.5	98.6	99.0	99.5	100
°C	35.0	35.3	35.6	35.8	36.1	36.4	36.7	36.9	37.0	37.2	37.5	37.8

ne when ovulation has occurred. Mucus observations, as noted in «The female productive cycle», can also be helpful when recorded as shown on the sample chart. Notation of these changes, in addition to your temperature pattern, will assist your physician or counselor in providing you with proper guidance. It is commonly suggested that couples who have had difficulty in achieving pregnancy should not have intercourse more than once a day nor on consecutive days, in order to maintain satisfactory quantities of sperm. The most fertile time continues up through the first day of drying up of the mucus or the first day of temperature elevation, whichever comes later. The temperature record should be maintained on a daily basis. If pregnancy is achieved, the basal temperature will remain elevated for several months. A temperature pattern that remains elevated for 21 days (a week or more than the usual elevated pattern) is a good indicator that pregnancy has been achieved. If you do not become pregnant after 6 months, see your doctor. If you believe you are pregnant, see your doctor immediately.

AVOIDING PREGNANCY WITH NATURAL FAMILY PLANNING

In natural family planning, there are several ways of determining the end of pre-ovulation fertility and the beginning of post-ovulation infertility. Professionals teaching natural family planning claim most unplanned pregnancies occur during pre-ovulation infertility. This brochure does not include guidelines for natural family planning, since it is beyond the scope of this publication. Your basal thermometer only helps you determine when ovulation has occurred. Calculating infertile periods involves other bodily signs in addition to temperature change and requires professional guidance when avoidance of pregnancy is your objective. **NOTE: This product is not for contraceptive use!** This thermometer is to be used to aid in conception and is not to be used for contractive purposes.

CAUTION

The effectiveness of the methods outlined for seeking and avoiding pregnancy depends on factors such as a woman's individual reactions, the accuracy with which daily temperatures and other symptoms are recorded, and the absence of outside factors which may affect body temperature. Consequently, no representation of any kind can be made as to the reliability of the results obtained in any individual case. Please consult with and be guided by the advice of your physician or counselor in the use of this thermometer.

CLEANING AND DISINFECTING

For disinfection in home use environment, use a 70% Isopropyl alcohol swab or a cotton tissue moistened with 70% Isopropyl alcohol to wipe surface pollutants off the thermometer probe. Always start wiping from the end of the thermometer probe (approx. at the middle of the thermometer) towards the thermometer tip. Afterwards the entire thermometer should be immersed in 70% Isopropyl alcohol for at

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least 5 minutes (max. 24 hours). After immersion, let the disinfectant dry off for 1 minute before next use.

For professional use: You may apply probe covers (available as accessory). Contact Microlife for further details concerning suitable disinfectants for your thermometer model.

TECHNICAL SPECIFICATIONS

Measuring range: 32.00 - 42.99 °C / 90.00 - 109.99 °F
Resolution: 0.01 °C / 0.01 °F
Accuracy: ± 0.1 °C / ± 0.2 °F
Operating conditions: 10 - 40 °C; 15-95 % relative maximum humidity
Storage conditions: -25 - +60 °C; 15-95 % relative maximum humidity
Battery: LR41 (1.5V) / SR41 (1.55V)
Battery lifetime: approx. 4500 measurements (using a new battery)

Beeper Function: a. 1 Beep when power on or off
b. 10 beeps after measurement completed

Memory function: Press and hold the button for more than 2 sec., the unit will recall the last reading.

Auto-off function: ± 10 min. after measurement is completed

Low battery indication: Displays the symbol «▼» to remind user to replace the battery.
No measurement is available before the battery is replaced.

Self-test function: Displays «ERR» when system malfunctions
Response time: 60 sec. for oral measurements (reference)

IP Class: IP67
Reference to standards: EN 12470-3; ASTM E1112; IEC 60601-1; IEC 60601-1-11; IEC 60601-1-2 (EMC)

Expected service life: 5 years or 10000 measurements
This device complies with the requirements of the Medical Device Directive 93/42/EEC. Technical alterations reserved.

GUARANTEE

This device is covered by a 2 year guarantee from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.
Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Battery.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website: www.microlife.com/support

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

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