IB MT 3010 E-V11 4923, Revision Date: 2023-11-13

Description of this Thermometer

- ON/OFF button
- Display
- Battery compartment cover
- Measuring sensor / measuring tip
- Cleaning and disinfecting area (thermometer probe only)

Antimicrobial Copper Properties

Surfaces made or covered by special copper alloys, have strong anti-bioicrobial properties against a wide variety of microorganisms\*, thus destroying these microorganisms. This activity reduces the microbial flora on the coated area and – due to the «halo phenomenon» – simultaneously causes a drastic reduction in pathogens on the remaining body of the thermometer. Thermometers containing parts with antimicrobial copper alloy drastically reduce microbial flora, minimizing the dispersion of contagious microorganisms and thus providing high safety to the end user.

\*Establishou A. Panos «The Role of Antimicrobial Copper Surfaces in Reducing Healthcare-associated Infections», European Infectious Disease, Volume 5, Issue 2, Autumn 2011.

Important Safety Instructions

- Follow instructions for use. This document provides important product and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.
- This device is only to be used for measuring human body temperature through oral, rectal or axillary. Do not attempt to take measurements at other sites, such as in the ear, as it may result in false readings and may lead to injury.
- Do not use this device if you think it is damaged or notice anything unusual.
- We recommend cleaning this device according to the cleaning instructions before first use for personal hygiene.
- The minimum measurement time until the beep is heard must be maintained with accuracy.
- Consider that different measurement locations may require continued measuring even after the beep, see section «Measuring methods / Normal body temperature».
- For safety reasons (risk of rectal perforation) rectal measurement in children younger than 3 years must be performed only by trained health professionals (professional users). Use another measurement method instead. For rectal fever measurement in small children younger than 3 years, thermometers with a flexible tip are available.
- Do not attempt rectal measurements on persons with rectal disorders. Doing so may aggravate or worsen the disorder.
- Ensure that the device is not used unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.
- Protect the device from impact and dropping!
- Avoid ambient temperatures above 60 °C. NEVER boil this device!
- Use only the commercial disinfectants listed in the section «Cleaning and disinfecting» to clean the device to avoid damage to the device.
- We recommend this device is tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact your local MicroLife-Service to arrange the test.

- WARNING:** The measurement result given by this device is not a diagnosis! Do not rely on the measurement result only.

- Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

- Read the instructions carefully before using this device.

- Type BF applied part

- Keep dry

- Turning on the Thermometer

To turn on the thermometer, press the ON/OFF button (1); a short beep signals «thermometer ON». A display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C, an «L» and a flashing «C» appear at the display field (2). The thermometer is now ready for use.

Function Test

Correct functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement is inaccurate), this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

Using the Thermometer

Before use, keep this thermometer out of any physical contact, for at least 2 hours. This utilizes the antimicrobial copper properties, as described in chapter «Antimicrobial Copper Properties». The use of antimicrobial copper is a supplement to and not a substitute for standard infectious control practices. Users must continue to follow all current infection control and cleaning practices. We recommend cleaning the thermometer as described in chapter «Cleaning and disinfecting».

Choose the preferred measuring method. When taking a measurement, the current temperature is continuously displayed and the «C» symbol flashes. If the beep is heard 10 times and the «C» is no longer flashing, this means that the measured increase in temperature is less than 0.1 °C in 16 seconds.

To prolong the battery life, turn off the thermometer by briefly pressing the ON/OFF button (1). Otherwise the thermometer will automatically turn off after about 10 minutes.

**Storage of Measured Values**

If the ON/OFF button (1) is pressed for more than 3 seconds when turning on the thermometer, the automatically stored maximum temperature during the last measurement will be displayed. At the same time, a «M» for memory will appear on the display. About 3 seconds after the button is released, the temperature value disappears and the thermometer is ready for measurement.

**Measuring methods / Normal body temperature**

**In the armpit (axillary) / 34.7 - 37.3 °C**

Wipe the underarm with a dry towel. Place the measuring sensor (4) under the arm into the center of the armpit so that it is touching the skin and position the patient’s arm next to the patient’s body. This ensures that the room air does not affect the reading. Because the axillary takes more time to reach its stable temperature, wait **at least 5 minutes**, regardless of the beep sound.

**In the mouth (oral) / 35.5 - 37.5 °C**

Do not eat or drink anything hot or cold 10 minutes before the measurement. The mouth should remain closed up to 2 minutes before starting a reading.

Position the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue. The measuring sensor (4) must be in good contact with the tissue. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. If this is not possible due to a blocked airways, another method for measuring should be used.

- Approx. measuring time: 1 minute!**
- In the anus (rectal) / 36.6 - 38.0 °C**
- Attention:** For the prevention of rectal perforation in children (younger than 3 years), we recommend using another measuring method, or the use of a thermometer with a flexible tip. Carefully insert the measuring sensor (4) of the thermometer 2 to 3 cm into the anal aperture.

The use of a probe cover and the use of a lubricant is recommended. If you are unsure of this measurement method, you should consult a professional for guidance.

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