*microlife*<sup>®</sup>

### Microlife IR 150

---





<b>EN</b>	➔	1
<b>ES</b>	➔	6
<b>FR</b>	➔	12
<b>IT</b>	➔	18
<b>DE</b>	➔	24
<b>TR</b>	➔	31
<b>PT</b>	➔	38
<b>NL</b>	➔	44
<b>GR</b>	➔	50
<b>AR</b>	➔	56
<b>FA</b>	➔	62

Microlife AG  
Esenstrasse 139  
9443 Widnau / Switzerland  
[www.microlife.com](http://www.microlife.com)

**CE0044**

IB IR 150 S-V11 5020  
Revision Date: 2020-11-24

*microlife*<sup>®</sup>

Name of Purchaser / Nom de l'acheteur / Nome del rivenditore / Name des Käufers / Alicinin adı / Nombre del comprador / Nome do comprador / Naam koper / Ονοματεπώνυμο αγοραστή / نام خریدار / اسم المشتري

Serial Number / Numéro de série / Numero di serie / Serien-Nr. / Seri numarasi / Número de serie / Número de série / Seriennummer / Αριθμός σειράς / رقم التسلسل / مدل

Date of Purchase / Date d'achat / Data d'acquisto / Kaufdatum / Satın alma tarihi / Fecha de compra / Data da compra / Datum van aankoop / Ημερομηνία αγοράς / تاريخ الشراء / شماره سریال

Specialist Dealer / Revendeur / Categoria rivenditore / Fachhändler / Uzman satıcı / Distribuidor especializado / Revendedor autorizado / Specialist Dealer / Εξειδικευμένος αντιπρόσωπος / تاريخ خرید / التاجر المختص

- ① Measuring sensor
- ② Display
- ③ START button
- ④ ON/OFF button
- ⑤ Battery compartment cover
- ⑥ M-button (memory)
- ⑦ All segments displayed
- ⑧ Memory
- ⑨ Ready for measuring
- ⑩ Measurement complete
- ⑪ Out-of-ear temperature indication
- ⑫ Low battery indicator
- ⑬ Changing between Celsius and Fahrenheit
- ⑭ Recall mode
- ⑮ Recall the last 30 readings
- ⑯ Measured temperature too high
- ⑰ Measured temperature too low
- ⑱ Ambient temperature too high
- ⑲ Ambient temperature too low
- ⑳ Error function display
- ㉑ «CLEAN ME» display
- ㉒ Blank display
- ㉓ Flat battery
- ㉔ Date/Time
- ㉕ Beeper function setting
- ㉖ Replacing the battery

Read the instructions carefully before using this device.

Type BF applied part

This Microlife thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, this device can provide a stable, heat-interference-free reading with each measurement. The device performs a self-test every time it is switched on to always guarantee the specified accuracy of any measurement. The Microlife Ear Thermometer is intended for the periodic measurement and monitoring of human body temperature. It is intended for use on people of all ages.

**This thermometer has been clinically tested and proven to be safe and accurate when used in accordance to the operating instruction manual.**

Please read through these instructions carefully in order for you to understand all functions and safety information.

## Table of Contents

### 1. The Advantages of this Thermometer

- Measurement in 1 second
- Multiple uses (wide range of measurement)
- Probe cover free
- Probe LED
- Accurate and reliable
- Gentle and easy to use
- Multiple readings recall
- Safe and hygienic
- Fever alarm

### 2. Important Safety Instructions

### 3. How this Thermometer measures Temperature

### 4. Control Displays and Symbols

### 5. Setting Date, Time and Beeper Functions

### 6. Directions for Use

### 7. Changing between Celsius and Fahrenheit

### 8. How to recall 30 readings in Memory Mode

### 9. Error Messages

### 10. Cleaning and Disinfecting

### 11. Battery Replacement

### 12. Guarantee

### 13. Technical Specifications

### 14. [www.microlife.com](http://www.microlife.com)

Guarantee Card (see Back Cover)

## 1. The Advantages of this Thermometer

### Measurement in 1 second

The innovative infrared technology allows measurement of ear temperature in only 1 second.

### Multiple uses (wide range of measurement)

This thermometer offers a wide range of measurement from 0 - 100.0 °C / 32 - 212.0 °F, meaning the unit can be used to measure body temperature or it also has a feature allowing it to be used to measure surface temperature of the following examples:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

### Probe cover free

This thermometer is more user-friendly and more cost-effective since a probe cover is not required.

### Probe LED

This thermometer includes a probe LED light which enables the user to find the correct ear position in the dark and is used as a reminder (flashing LED and probe light) to clean the thermometer after each temperature measurement.

### Accurate and reliable

The unique probe assembly construction incorporates an advanced infrared sensor, ensuring that each measurement is accurate and reliable.

### Gentle and easy to use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can even be used on a sleeping child without causing any interruption.
- This thermometer is quick, therefore child-friendly.

### Multiple readings recall

Users will be able to recall the last 30 readings with a record of both time and date when entering the recall mode, enabling efficient tracking of temperature variations.

### Safe and hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Cleaning of the probe can be done with an alcohol-moistened cotton tissue, making this thermometer completely hygienic for use by the whole family.

### Fever alarm

10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature equal to or higher than 37.5 °C.

## 2. Important Safety Instructions

- Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.
- This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.
- **Never immerse this device in water or other liquids. For cleaning please follow the instructions in the «Cleaning and Disinfecting» section.**
- Do not use this device if you think it is damaged or notice anything unusual.
- Never open this device.
- Earwax in the ear canal may cause a lower temperature reading. Therefore it is important to ensure the subject's ear canal is clean.
- If the measurement result is not consistent with the patient's finding or unusually low, repeat the measurement every 15 minutes or double check the result by another core body temperature measurement.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.
- Protect it from:
  - extreme temperatures
  - impact and dropping
  - contamination and dust
  - direct sunlight
  - heat and cold
- If the device is not going to be used for a prolonged period the battery should be removed.

**WARNING:** The measurement results given by this device is not a diagnosis. It is not replacing the need for the consultation of a physician, especially if not matching the patient's symptoms. Do not rely on the measurement result only, always consider other potentially occurring symptoms and the patient's feedback. Calling a doctor or an ambulance is advised if needed.

### 3. How this Thermometer measures Temperature

This thermometer measures infrared energy radiated from the eardrum and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. The measured reading obtained directly from the eardrum (Tympanic Membrane) ensures the most accurate ear temperature. Measurements taken from the surrounding tissue of the ear canal generate lower readings and may result in misdiagnosis of fever.

#### To avoid an inaccurate measurement

1. Switch on the thermometer by pressing the ON/OFF button ④.
2. After one beep is heard (and the temperature scale icon is flashing), straighten the ear canal by gently pulling the middle of the ear back and up.
3. Place the probe ① firmly into the ear canal, press the START button ③ and keep the probe in the ear until the thermometer beeps to signal the completion of the measurement.

### 4. Control Displays and Symbols

- **All segments displayed** ⑦: Press the ON/OFF button ④ to turn on the unit; all segments will be shown for 1 second.
- **Ready for measurement** ⑨: When the unit is ready for measurement, the «°C» or «°F» icon will keep flashing. The probe LED light is activated for 10 seconds.
- **Measurement complete** ⑩: The reading will be shown on the display ② with the «°C» or «°F» icon; the unit is ready for the next measurement, when the «°C» or «°F» icon is flashing again.
- **Out-of-ear temperature indication** ⑪: A crossed-ear-icon will appear on the display ② if the reading falls outside the range 32.0 ~ 42.2 °C / 89.6 ~ 108.0 °F.
- **Low battery indicator** ⑫: When the unit is turned on, the «▼» icon will keep flashing to remind the user to replace the battery.

### 5. Setting Date, Time and Beeper Functions

#### Setting the date and time

1. After the new battery is fitted, the year number flashes in the display ②. You can set the year by pressing the START button ③. To confirm and then set the month, press the M-button ⑥.
2. Press the START button ③ to set the month. Press the M-button ⑥ to confirm and then set the day.
3. Follow the previously mentioned instructions to set the day, 12 or 24 hour mode, hours and minutes.
4. Once you have set the minutes and pressed the M-button ⑥, the date and time are set and the time is displayed.

☞ If no button is pressed for 20 seconds, the device automatically switches to ready for measuring ⑨.

☞ **Cancel time setup:** Press the ON/OFF button ④ during time setup. The LCD will show Date/Time icons with «--:--». After that press the ON/OFF button ④ to start the measurement. If no further action is taken within 60 seconds, the device will automatically turn off.

☞ **Change current date and time:** Press and hold the M-button ⑥ for approx. 3 seconds until the year number starts to flash ②. Now you can enter the new values as described above.

#### Setting the beeper

1. When the device is switched off, press and hold the ON/OFF button ④ for 5 seconds to set the beeper ⑫.
2. Press the ON/OFF button ④ again to either turn the beeper on or off. The beeper is activated when the beeper icon ⑫ is not crossed-out.

☞ If no button is pressed for 5 seconds, the device automatically switches to ready for measuring ⑨.

### 6. Directions for Use

1. Press the ON/OFF button ④. The display ② is activated to show all segments for 1 second.
2. When the «°C» or «°F» icon is flashing, a beep sound is heard and the thermometer is ready for measuring ⑨.
3. The probe LED light is activated and will automatically turn off after 10 seconds.
4. Straighten the ear canal by pulling the ear up and back to give a clear view of the eardrum.
  - For children under 1 year; pull the ear straight back.
  - For children of 1 year to adult; pull the ear up and back.Also refer to the short instructions at the front.
5. While gently pulling the ear, insert the probe snugly into the ear canal (**about 1 sec.**) and **immediately** press the START button ③. Release the button and wait for the beep sound. This indicates the end of the measurement.
6. Remove the thermometer from the ear canal. The display shows the measured temperature ⑩.

#### ☞ NOTE:

- **Patients and thermometer should stay in similar room condition for at least 30 minutes.**

- To ensure accurate readings, wait at least 30 sec. after 3-5 continuous measurements.
  - Accumulation of ear wax on the probe can result in less accurate temperature readings or cross infection between users.
  - **It is essential that the probe is cleaned after every measurement.** Therefore, this device reminds the user to clean the probe when turning off the device. «CLEAN ME» (2) is displayed and the probe LED light will flash for 3 seconds. For cleaning, follow the instructions in the «Cleaning and Disinfecting» section.
  - **After cleaning the measuring sensor (1) with alcohol, wait 5 minutes before taking the next measurement,** in order to allow the thermometer to reach its operating reference temperature.
  - 10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature equal to or higher than 37.5 °C.
  - For an infant, it is best to have the child lying flat with his/her head sideways so the ear is facing upwards. For an older child or adult, it is best to stand behind and slightly to the side of the patient.
  - Always take the temperature in the same ear, since the temperature readings may be different from ear to ear.
  - In the following situations it is recommended that three temperatures in the same ear be taken and the highest one taken as the reading:
    1. New born infants in the first 100 days.
    2. Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
    3. When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the device and obtains consistent readings.
    4. If the measurement is surprisingly low.
  - Don't take a measurement while or immediately after nursing a baby.
  - Don't use the thermometer in high humidity environments.
  - Patients should not drink, eat or exercise before/while taking the measurement.
  - Doctors recommend rectal measurement for newborn infants within the first 6 months, as all other measuring methods might lead to ambiguous results.
  - **Readings from different measuring sites should not be compared as the normal body temperature varies by measuring site and time of day,** being highest in the evening and lowest about one hour before waking up.
- Normal body temperature ranges:
- Axillar: 34.7 - 37.3 °C / 94.5 - 99.1 °F
  - Oral: 35.5 - 37.5 °C / 95.9 - 99.5 °F
  - Rectal: 36.6 - 38.0 °C / 97.9 - 100.4 °F
  - Microlife IR 150: 35.4 - 37.4 °C / 95.7 - 99.3 °F

## 7. Changing between Celsius and Fahrenheit

This thermometer can display temperature readings in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn OFF the unit, **press and hold** the START button (3) for 5 seconds; after 5 seconds, the current measurement scale («°C» or «°F» icon) will flash on the display (13). Change the measurement scale between °C and °F by pressing the START button (3). When the measurement scale has been chosen, wait for 5 seconds and the unit will automatically enter the «ready for measuring» mode.

## 8. How to recall 30 readings in Memory Mode

This thermometer can recall the last 30 readings with a record of both time and date.

- **Recall mode (14):** Press the M-button (6) to enter recall mode when the power is off. The memory icon «M» will flash.
- **Reading 1 - the last reading (15):** Press and release the M-button (6) to recall the last reading. Number «1» and a flashing «M» are displayed.
- **Reading 30 - readings in succession:** Press and release the M-button (6) consecutively to recall the last 30 readings in succession.

Pressing and releasing the M-button (6) after the last 30 readings have been recalled will resume the above sequence from reading 1.

## 9. Error Messages

- **Measured temperature too high (16):** Displays «H» when measured temperature is higher than 100.0 °C / 212.0 °F.
- **Measured temperature too low (17):** Displays «L» when measured temperature is lower than 0 °C / 32.0 °F.
- **Ambient temperature too high (18):** Displays «H» when ambient temperature is higher than 40 °C / 104 °F.
- **Ambient temperature too low (19):** Display «L» when ambient temperature is lower than 10.0 °C / 50.0 °F.
- **Error function display (20):** The system has a malfunction.
- **Blank display (22):** Check if the battery has been loaded correctly. Also check polarity (<+> and <->) of the battery.
- **Flat battery indication (23):** If only the «▼» icon is shown on the display, the battery should be replaced immediately.

## 10. Cleaning and Disinfecting

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring sensor. Ensure that no liquid enters the interior of the device. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the device in water or other cleaning liquids. Take care not to scratch the surface of the sensor lens and the display.

## 11. Battery Replacement

This device is supplied with one lithium battery, type CR2032. The battery needs replacing when only the «▼» icon (23) is shown on the display.

To remove the battery cover you must first undo the screw (26), then remove the battery cover by sliding it in the direction shown. Insert the new battery with the + at the top. To close, replace the battery cover and tighten the screw.

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

## 12. Guarantee

This device is covered by a **5 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Batteries, Probe cover (optional).

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website:

[www.microlife.com/support](http://www.microlife.com/support)

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

## 13. Technical Specifications

<b>Type:</b>	Ear Thermometer IR 150
<b>Measurement range:</b>	0 - 100.0 °C / 32 - 212.0 °F
<b>Resolution:</b>	0.1 °C / °F
<b>Measurement accuracy:</b>	Laboratory: ±0.2 °C, 32.0 ~ 42.2 °C ±0.4 °F, 89.6 ~ 108.0 °F

<b>Display:</b>	Liquid Crystal Display, 4 digits plus special icons
<b>Acoustic:</b>	The unit is turned ON and ready for the measurement: 1 short beep. Complete the measurement: 1 long beep. System error or malfunction: 3 short beeps. Fever alarm: 10 short beeps.
<b>Memory:</b>	30 readings recall in the memory mode with a record of both time and date.
<b>Backlight:</b>	The display light will be GREEN for 1 second, when the unit is turned ON. The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5 °C / 99.5 °F. The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5 °C / 99.5 °F.
<b>Operating conditions:</b>	10 - 40 °C / 50.0 - 104 °F
<b>Storage conditions:</b>	15 - 95 % relative maximum humidity -25 - +55 °C / -13 - +131 °F
<b>Automatic Switch-off:</b>	15 - 95 % relative maximum humidity Approx. 1 minute after last measurement has been taken.
<b>Battery:</b>	1 x CR2032 battery 3V
<b>Battery lifetime:</b>	approx. 1000 measurements (using a new battery)
<b>Dimensions:</b>	139 x 39 x 42 mm
<b>Weight:</b>	54.5 g (with battery), 51 g (w/o battery)
<b>IP Class:</b>	IP22
<b>Reference to standards:</b>	EN 12470-5; ASTM E 1965; IEC 60601-1; IEC 60601-1-2 (EMC); IEC 60601-1-11
<b>Expected service life:</b>	5 years or 12000 measurements

This device complies with the requirements of the Medical Device Directive 93/42/EEC.

Technical alterations reserved.

According to the Medical Product User Act a biennial technical inspection is recommended for professional users. Please observe the applicable disposal regulations.

## 14. [www.microlife.com](http://www.microlife.com)

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at [www.microlife.com](http://www.microlife.com).

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