



Compressor Nebuliser

EN

Sterilizing with steam
All nebuliser components (except air hose and masks) can be heat steam sterilized up to 121 °C (20 min.) or 134 °C (7 min.). EN554/ISO11134.

Nebuliser «NEB PRO»
After sterilization always let all components cool down to ambient temperature before further use. Do not repeat sterilization cycle when components are still warm.

5. Maintenance, Care, and Service

Order all spare parts from your dealer or pharmacist, or contact Microlife-Service (see forward).

We recommend replacing the nebuliser, mouthpiece, and the face masks after 3 months use.

Check the filter after each use for clogging and replace the filter or filter + mouthpiece after 3 months use. Spare filters are provided with the device.

To replace the filter, open the air filter compartment (on the top of the device marked "filter") and exchange the filter (1).

Never unscrew the plastic lid.

6. Malfunctions and Actions to take

The device cannot be switched on:

- Ensure the power lead (2) is correctly plugged into the socket.
- The ON/OFF switch (3) is in the position «On».
- The device has automatically switched off due to overheating. Wait until it has cooled down and try again.

The nebuliser functions poorly or not at all:

- Ensure the air hose (6) is correctly connected at both ends.
- Ensure the air hose is not snatched, bent, dirty or blocked. If necessary, replace with a new one.
- Ensure the nebuliser (12) is fully assembled and the colored impactor (12-B) is placed correctly.
- Ensure required medication has been added.

7. Guarantee

This device is covered by a 3 year guarantee from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) confirming date of purchase or the receipt.

8. Technical Specifications

Nebulisation rate: 0.5 ml/min. (NaCl 0.9%)
Particle size: 75% < 5 μm NaCl 0.9%
2.86 μm (MMD) with 0.9% NaCl using Cascade Impactor
Compressor flow: 13 l/min.
Acoustic noise level: 53 dB
Power source: 230V - 50Hz / 150VA
Power lead length: 1.5 m
Nebuliser capacity: min. 2 ml; max. 12 ml
Residual volume: 0.7 ml
Operating limits: 30 min. On / 30 min. Off
Operating temperature: 10 - 40 °C / 50 - 104 °F
10-95% relative maximum humidity
Storage temperature: -25 - +70 °C / -13 - +58 °F
-55 - +55% relative maximum humidity
Weight: 1300 g
Dimensions: 209 x 140 x 124 mm
Reference to standards: EN 13544-1; EN 60601-1; EN 60601-1-2; CE0434
Technical alterations reserved:

Dear Customer,
This nebulizer is a high-quality device for inhalation therapy for asthma, chronic bronchitis, and other respiratory illnesses.

The operation and usage of the device is very straightforward. All common liquid medication can be used for inhalation therapy.

Please read through these instructions carefully so that you understand all functions and safety information. We want you to be happy with your Microlife product. If you have any questions, problems or want to order spare parts please contact Microlife-Customer Service. Your dealer or pharmacist will also give you the address and the phone number of our service center. You can visit the internet at www.microlife.com where you will find a wealth of available information on our products. Retain instructions in a safe place for future reference. Stay healthy - Microlife AG!

Read the instructions carefully before using this device.

Type B applied part

1. Important Safety Instructions

This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.

This device is not suitable for anaesthesia and lung ventilation.

This device should only be used with original accessories as shown in these instructions.

Do not use this device if you think it is damaged or not working correctly.

Never open this device.

This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.

Protect the device from:

- water and moisture

- extreme temperatures

- impact and dropping

- direct sunlight

- heat and cold

Only use the medication prescribed for you by your doctor and follow your doctor's instructions with regard to dosage, duration and frequency of the therapy.

Never heat the nebuliser over 60°.

Never leave children alone with this device unless supervised; some parts are small enough to be swallowed.

Electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Use of this device is not intended as a substitute for a consultation with your doctor.

2. Preparation and Usage of this Device

Prior to using the device for the first time, we recommend cleaning it as described in the section «Cleaning and Disinfecting».

1. Assemble the nebuliser kit (12). Ensure that all parts are complete.

2. Fill the nebuliser (13) with the inhalation solution as per your doctor's instructions. Ensure that you do not exceed the maximum level.

3. Connect the nebuliser with the air hose (6) to the compressor (1) and plug the power lead (2) into the socket (230V / 50Hz).

4. Switch the ON/OFF switch (3) to position «On» to turn on the device and place the mouthpiece (7) to the mouth or fit one of the face masks over mouth and nose.

- Take the mask (8) and give you a better drug delivery to the lungs.

- Choose between adult (9) or child face mask (9) and make sure that it encloses the mouth and nose area completely.

- Use the nebuliser (12) by inserting the nose piece (10) as prescribed by your doctor.

5. Breathe in and out calmly during therapy. Sit in a relaxed position with the upper body upright. Do not lie down while inhaling. Stop treatment if you feel unwell.

6. After completing the inhalation period recommended by your doctor, switch the ON/OFF switch (3) to position «Off» to turn off the device and unplug it from the socket.

7. Empty the remaining medication from the nebuliser and clean the device as described in the section «Cleaning and Disinfecting».

8. Nebuliser «NEB PRO» treatment options

This device is equipped with the innovative nebuliser «NEB PRO» (5), which includes the breath-enhanced valve system (14). All common medication for inhalation therapy can be nebulised effectively.

Advanced treatment

This treatment allows the maximum amount of medication to be delivered thus reaching the respiratory capacity of each patient. The system optimally adapts the flow of medication during inhalation (14-A), reducing medication loss during exhalation (14-B).

Use the breath-enhanced valve system (synchronized with breathing) to minimize the therapeutic effects of the medication and get the best results.

► Mount the membrane on the valve-holder and use the mouthpiece with respiratory valve (12-A).

► Follow the usage steps as described in Section 2.4.

► Only using the mouthpiece with respiratory valve (7) enables the synchronized functioning of the device.

Short treatment

The therapeutic session will be reduced if the membrane is not mounted on the valve-holder.

► A synchronized nebulisation is not possible

► Please use this treatment option.

4. Cleaning and Disinfecting

Thoroughly clean all components to remove medication residues and possible impurities after each treatment.

The compressor (1) and the air hose (6) should be cleaned with a clean, moist cloth.

► Always wash your hands well before cleaning and disinfecting the accessories.

► Replace the air hose for each treatment with a new patient or in case of impurities.

► Unplug the device before cleaning.

Cleaning with water

► Wash all nebuliser components (except air hose)

under warm tap water (max 60 °C) for about 5 minutes adding if necessary a small quantity of detergent. Following dosage and use limitations as provided by the manufacturer.

► Rinse thoroughly making sure that all deposits are washed away and leave to dry.

Disinfecting

All nebuliser components (except air hose) can be disinfected with chemical disinfectants following dosage and use limitations as provided by disinfectant manufacturer.

► Disinfectants are usually available at pharmacies.

5. Maintenance, Care, and Service

Observe all spare parts from your dealer or pharmacist, or contact Microlife-Service (see forward).

We recommend replacing the nebuliser, mouthpiece and the face masks after 3 months use.

Check the filter after each use for clogging and replace the filter or filter + mouthpiece after 3 months use. Spare filters are provided with the device.

To replace the filter, open the air filter compartment (on the top of the device marked "filter") and exchange the filter (1).

Never unscrew the plastic lid.

6. Malfunctions and Actions to take

The device cannot be switched on:

► Ensure the power lead (2) is correctly plugged into the socket.

► The ON/OFF switch (3) is in the position «On».

► The device has automatically switched off due to overheating. Wait until it has cooled down and try again.

The nebuliser functions poorly or not at all:

► Ensure the air hose (6) is correctly connected at both ends.

► Ensure the air hose is not snatched, bent, dirty or blocked. If necessary, replace with a new one.

► Ensure the nebuliser (12) is fully assembled and the colored impactor (12-B) is placed correctly.

► Ensure required medication has been added.

7. Guarantee

This device is covered by a 3 year guarantee from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) confirming date of purchase or the receipt.

8. Technical Specifications

Nebulisation rate: 0.5 ml/min. (NaCl 0.9%)

Particle size: 75% < 5 μm NaCl 0.9%

2.86 μm (MMD) with 0.9% NaCl using Cascade Impactor

Compressor flow: 13 l/min.

Acoustic noise level: 53 dB

Power source: 230V - 50Hz / 150VA

Power lead length: 1.5 m

Nebuliser capacity: min. 2 ml; max. 12 ml

Residual volume: 0.7 ml

Operating limits: 30 min. On / 30 min. Off

Operating temperature: 10 - 40 °C / 50 - 104 °F

Humidity range: 10-95% relative maximum humidity

Storage temperature: -25 - +70 °C / -13 - +58 °F

Weight: 1300 g

Dimensions: 209 x 140 x 124 mm

Reference to standards: EN 13544-1; EN 60601-1; EN 60601-1-2; CE0434

Technical alterations reserved:

9. Characteristics of the nebuliser «NEB PRO»

This device features an innovative nebuliser system (5), which is based on a cascade system of valves.

The compressor (1) and the air hose (6) are connected to the cascade system.

The nebuliser «NEB PRO» (5) includes the breath-enhanced valve system (14).

All common medication for inhalation therapy can be nebulised effectively.

Advanced treatment

This treatment allows the maximum amount of medication to be delivered thus reaching the respiratory capacity of each patient.

The system optimally adapts the flow of medication during inhalation (14-A),

reducing medication loss during exhalation (14-B).

Use the breath-enhanced valve system (synchronized with breathing) to minimize the therapeutic effects of the medication and get the best results.

► Mount the membrane on the valve-holder and use the mouthpiece with respiratory valve (12-A).

► Follow the usage steps as described in Section 2.4.

► Only using the mouthpiece with respiratory valve (7) enables the synchronized functioning of the device.

Short treatment

The therapeutic session will be reduced if the membrane is not mounted on the valve-holder.

► A synchronized nebulisation is not possible

► Please use this treatment option.

4. Cleaning and Disinfecting

Thoroughly clean all components to remove medication residues and possible impurities after each treatment.

The compressor (1) and the air hose (6) should be cleaned with a clean, moist cloth.

► Always wash your hands well before cleaning and disinfecting the accessories.

► Replace the air hose for each treatment with a new patient or in case of impurities.

► Unplug the device before cleaning.

Cleaning with water

► Wash all nebuliser components (except air hose)

under warm tap water (max 60 °C) for about

5 minutes adding if necessary a small quantity of

detergent. Following dosage and use limitations as provided by the manufacturer.

► Rinse thoroughly making sure that all deposits are

washed away and leave to dry.

Disinfecting

All nebuliser components

Kompresorový inhalátor

SK

- Piestový kompresor
- Napájaci kábel
- Prepnac ON/OFF
- Priestor pre vzduchový filter
- Rozprášovač «NEB PRO»
- Vzduchová hadica
- Náustok s ventilom
- Tvárová maska pre dospelých
- Tvárová maska pre deti
- Nosový nadstavec
- Výmena vzduchového filtra
- Zostavanie jednotlivých častí inhalátora
- Naplnenie inhalátora
- Inhalácia synchronizovaná s dýchaním

Väčší základník. Väčší komprimor Microlife je vysokovýkonný zariadenie na inhaláciu terapeut. astmy, chronickej bronchitídy a ostatných respiračných ochorení. Obsluha a používanie je veľmi jednoduchá. Na inhaláciu je možné použiť väčší jednoduchý. Prečítajte si prosím tento návod pozorné tak, aby ste porozumeli všetkým funkciám a funkciám ktorých sa jučia. Chcete aj počítať s výrobcom? Akožte akékoľvek otázky, voprosy alebo o chote objednajte našejednej diely, kontaktujte prosím zákaznického servisu Microlife. Väčší predajca alebo lekár. Väčší predajca vás poskytne adresu distribútoru Microlife vo Vašej krajinie. Pripádejte internetom stránku www.microlife.com, kde môžete nájsť množstvo neoceniteľných informácií o výrobku.

Pred použitím prístroja si pozorne prečítajte návod a uschovajte ho pre budúcu upotrebu.

Zostavte zdraví - Microlife AG!

Pred použitím zariadenia si pozorne prečítajte návod.

Použité súčasti typu B

1. Dôležité bezpečnostné pokyny

- Tento prístroj sa nášme používať iba na výdechu popisané v tomto návode. Výrobca nemôže niesť zodpovednosť za poškodenie spôsobené nespravnou používaniom.
- Prístroj je určený len pre jednotlivé používania.
- Prístroj by mal byť používaný len s originálnym príslušenstvom ake je uvedené v tomto návode.
- Nepoužívajte rúriky, ak sú myšliené, že je poškodený, alebo ak spoznajete niečo nezvyčajné.
- Nikdy rúriky nerobiť.
- Tentoraz prístroj obsahuje cítilce súšasťky, preto je potrebné s nimi zohľadniť opatrenie. Dodržiavať podmienky skladovania a prepravy diely je uvedené v kapitole «Technické údaje»!
- Pri používaní prístroja:
 - voda a výdech
 - extrémneho teplotami
 - nárazom a pádom
 - znečisteniu a prachom
 - priamym snehovým svetlom
 - tepmom a chladičom
- Na inhaláciu používajte len, ktoré sú výrobcom predpisane.
- Prístroj by mal byť používaný len s originálnym príslušenstvom ake je uvedené v tomto návode.
- Nepoužívajte rúriky, ak sú myšliené, že je poškodený, alebo ak spoznajete niečo nezvyčajné.
- Nikdy rúriky nerobiť.
- Tentoraz prístroj obsahuje cítilce súšasťky, preto je potrebné s nimi zohľadniť opatrenie. Dodržiavať podmienky skladovania a prepravy diely je uvedené v kapitole «Technické údaje»!

- Pri používaní prístroja:
 - voda a výdech
 - extrémneho teplotami
 - nárazom a pádom
 - znečisteniu a prachom
 - priamym snehovým svetlom
 - tepmom a chladičom

Na inhaláciu používajte len, ktoré sú výrobcom predpisane.

Prístroj je určený len pre jednotlivé používania.

Na inhaláciu používajte len, ktoré sú výrobcom predpisane.

Na inhaláciu použív