
Microlife NEB 50B

NEB 50B

IB NEB 50B S 9/2716

Guarantee Card

Name of Purchaser / Nom de l'acheteur / Nome del comprador / Nombre del comprador / Naam koper / نام خریدار

Serial Number / Numéro de série / Numero de serie / Serie-Nr. / Número de serie / Número de serie / Serienummer / Αριθμός σειράς / رقم التسلسلي

Date of Purchase / Date d'achat / Data d'acquisto / Kaufdatum / Fecha de compra / Data da compra / Datum van aankoop / حیثیت میو ایکوپ / تاریخ الشراء

Specialist Dealer / Revendeur / Categorie venditore / Fachhändler / Distributor especializado / Revendedor autorizado / Specialized Dealer / Εφορητής αποκύπων / الموزع المختص

3A Health Care S.r.l.
Via Marche Cattolica Loc. San Tommaso
25010 Novara (BS) - Italy
Tel. 0339 030 313177
Fax 0339 030 391914

CE 0434

Advanced Compressor Nebuliser

EN

- ① Piston compressor
- ② Power lead
- ③ ON/OFF switch
- ④ Air filter compartment
- ⑤ Nebuliser with switchable nebulisation rate
- ⑥ Air hose
- ⑦ Mouthpiece
- ⑧ Adult face mask
- ⑨ Child face mask
- ⑩ Nose piece
- ⑪ Replacing air filter
- ⑫ Assembling nebuliser kit

Dear Customer,

This professional nebuliser is a high-quality device for inhalation therapy for the treatment and care of upper and medium/lower respiratory tract.

The operation and usage of the device is very straightforward. All common liquid medication can be used for inhalation therapy.

Please follow these instructions carefully so that you understand all functions and safety information. We want to give you a great experience with your Microlife products. If you have any questions, problems or want to order spare parts please contact your local Microlife-Customer Service.

Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit our internet at www.microlife.com where you will find a wealth of invaluable information on our products. Retain instructions in a safe place for future reference.

Stay healthy – Microlife AG!

Read the instructions carefully before using this device.

Type BF applied part

Keep dry

1. Important Safety Instructions

• This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.

• This device is not suitable for anaesthesia and lung ventilation.

• This device should only be used with original accessories as shown in these instructions.

• Do not use this device if you think it is damaged or notice anything unusual.

• Never use this device near water.

• The device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.

• Protect it from:

- water and moisture
- extreme temperatures
- impact and dropping
- contact with dust
- direct sunlight
- heat and cold

• Only use the medication prescribed by you by your doctor and follow your doctor's instructions with regard to dosage, duration and frequency of the therapy.

• Never bend the nebuliser over 60°.

• The function of this device may be compromised when used close to strong electromagnetic fields such as mobile phones or radio installations and we recommend a distance of at least 1 m. In cases where there is no choice in the use of this device, please verify if the device is working properly before use.

• Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

• Electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

• Use of this device is not intended as a substitute for a consultation with your doctor.

2. Preparation and Usage of this Device

Prior to using the device for the first time, we recommend cleaning it as described in the section «Cleaning and Disinfecting».

1. Assemble the nebuliser kit (2). Ensure that all parts are complete. Use the wheel on the upper part of the nebuliser to select the position most suitable for the therapy to be performed (see chapter 3).

2. Fill the nebuliser with the inhalation solution as per your doctor's instructions. Ensure that you do not exceed the maximum level.

3. Connect the nebuliser with the air hose (6) to the compressor (1) and plug the power lead (2) into the socket (230V - 50Hz).

4. Switch the ON/OFF switch (3) to position «On» to turn on the device and place the mouthpiece (7) in the mouth or fit one of the face masks over mouth and nose.

- The mouthpiece gives you a better drug delivery to the lungs.

- Choose between adult (8) or child face mask (9) and make sure that it encloses the mouth and nose completely.

- Use all accessories including the nose piece (10) as prescribed by your doctor.

5. Breathe in and out calmly during the therapy. Sit in a relaxed position with the upper body upright. Do not lie down while inhaling. Stop inhalation if you feel unwell.

6. After completing the inhalation period recommended by your doctor, switch the ON/OFF switch (3) to position «Off» to turn off the device and unplug it from the socket.

7. Empty the remaining medication from the nebuliser and clean the device as described in the section «Cleaning and Disinfecting».

• This device was designed for intermittent use of 30 min. On / 30 min. Off. Switch off the device after 30 min. use and wait for another 30 min. before you resume treatment.

3. Choosing between 3 different Nebulisation Rates

The nebuliser generates an aerosol with variable characteristics to adapt to the patient's respiratory pathology. It can produce aerosols with aerodynamic characteristics that will deposit the maximum amount of the drug in the respiratory tract being treated.

Choose one of the 3 following therapeutic treatment settings by regulating the wheel on the nebuliser's (5) upper body.

• **Upper air passages (nose, throat, etc.):** Position the wheel towards the bottom at number 3; the aerodynamic diameter of the aerosol particles will be from 10.33 - 14.56 microm.

• **Traeths / bronchial tract:** Position the wheel in the center at number 2; the aerodynamic diameter of the aerosol particles will be from 6.21 - 9.45 micron.

• **Deep air passages:** Position the wheel towards the top, at number 1; the aerodynamic diameter of the aerosol particles will be from 1.14 - 5.32 micron.

4. Cleaning and Disinfecting

Thoroughly clean all components to remove medication residuals and possible impurities after each treatment.

The compressor (1) and the air hose (6) should be cleaned with a moist cloth.

• Always wash your hands well before cleaning and disinfecting the accessories.

• Do not expose the compressor to water or heat.

• Replace the air hose for each treatment with a new patient or in case of impurities.

• Unplug the device before cleaning.

Cleaning with water:

• Wash all nebuliser components (except air hose) under warm tap water (max. 60 °C) for about 5 minutes adding if necessary a small quantity of detergent following dosage and use limitations as provided by detergent manufacturer.

• Rinse thoroughly making sure that all deposits are washed away and leave to dry.

Disinfecting:

• All nebuliser components (except air hose) can be disinfected with chemical disinfectants following

dosage and use limitations as provided by disinfectant manufacturer.

• Disinfectants are usually available at pharmacies.

• Sterilizing with steam:

– All nebuliser components (except air hose and masks) can be heat sterilized up to 121 °C (20 min.) or 134 °C (7 min.). EN554/ISO11134.

– The sterilization packaging must conform to EN868/ ISO11607 and be suitable for steam sterilisation.

• After sterilization always let all components cool down to ambient temperature before further use.

• Do not repeat sterilization cycle when components are still warm.

5. Maintenance, Care, and Service

Order all spare parts from your dealer or pharmacist, or contact Microlife Service (see forward).

• We recommend replacing the nebuliser, mouthpiece, and the face masks after 3 months use.

• Check the filter continually for cleanliness and replace it if dirty, or after a maximum of 3 months use. Spare filters (11) are provided with the device.

• To replace the filter, open the air filter compartment (4).

• Do not repeat sterilization cycle when components are still warm.

6. Malfunctions and Actions to take

The device cannot be switched on:

• Ensure the air hose (6) is correctly plugged into the socket.

• Ensure the ON/OFF switch (3) is in the position «». The device has automatically switched off due to overheating. Wait until it has cooled down and try again.

The nebuliser functions poorly or not at all:

• Ensure the air hose is not crushed, bent, dirty or blocked. If necessary, replace with a new one.

• Ensure the nebuliser (2) is fully assembled and the colored vapouriser head (12-A) is placed correctly.

• Ensure the required medication has been added.

7. Guarantee

This device is covered by a 3 year warranty from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer confirming date of purchase or receipt.

• This device is not suitable for anaesthesia and lung ventilation.

• This device should only be used with original accessories as shown in these instructions.

• Do not use this device if you think it is damaged or notice anything unusual.

• Never use this device near water.

• The device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.

• Protect it from:

- water and moisture

- extreme temperatures

- impact and dropping

- contact with dust

- direct sunlight

- heat and cold

• Only use the medication prescribed by you by your doctor and follow your doctor's instructions with regard to dosage, duration and frequency of the therapy.

• Never bend the nebuliser over 60°.

• The function of this device may be compromised when used close to strong electromagnetic fields such as mobile phones or radio installations and we recommend a distance of at least 1 m. In cases where there is no choice in the use of this device, please verify if the device is working properly before use.

• Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

• Electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

• Use of this device is not intended as a substitute for a consultation with your doctor.

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