Shenzhen Jumper Medical Equipment Co., Ltd  
D Building, No. 71, Xintian Road, Fuyong Street,  
Baoan, Shenzhen, Guangdong 518103, China

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



10

V2.12	Limit Setup	*	a)
Alert	Alert	on	
Beep	Beep	off	
Demo	Demo	off	
Restore	Restore	ok	
Brightness	Brightness	4	
	Exit		

V2.12	Limit Setup	*	b)
SpO2	Hi	100	
SpO2	Lo	94	
PR Hi		130	
PR Lo		50	
+/-		+	
	Exit		

Name of Purchaser / Nom de l'acheteur / Nome del rivenditore /  
 Name des Käufers / Alicinin adı / Nombre del comprador / Nome do  
 comprador / Naam koper / Ονοματεπώνυμο αγοραστή / اسم المشتري

Serial Number / Numéro de série / Numero di serie / Serien-Nr. /  
 Seri numarası / Número de serie / Número de série / Seriennummer /  
 Αριθμός σειράς / رقم التسلسل

Date of Purchase / Date d'achat / Data d'acquisto / Kaufdatum /  
 Satın alma tarihi / Fecha de compra / Data da compra / Datum van  
 aankoop / Ημερομηνία αγοράς / تاريخ الشراء

Specialist Dealer / Revendeur / Categoria riven-ditore / Fachhändler /  
 Uzman satıcı / Distribuidor especializado / Revendedor autorizado /  
 Specialist Dealer / Εξειδικευμένος αντιπρόσωπος / التاجر المختص

- ① ON/OFF button
- ② Oxygen saturation (value as percentage)
- ③ Pulse rate (value in beats per minute)
- ④ Pulse bar
- ⑤ Low battery indicator
- ⑥ Inserting the batteries
- ⑦ Attaching the lanyard
- ⑧ Operation principle
- ⑨ Display modes (6 different)
- ⑩ Settings menu
  - a Interface 1
  - b Interface 2

This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products. Correct disposal helps to protect the environment and human health.

Read the important information in these instructions for use before using this device. Follow the instructions for use for your safety and keep it for future reference.

Type BF applied part

Class II equipment

**IP22**

Manufacturer

Low battery indicator

Serial number

Protected against dripping water

Authorized representative in the European Community

Medical device

Importer


Caution

% SpO<sub>2</sub>


Oxygen saturation (value as percentage)

♥ /Min

Pulse rate (value in beats per minute)

Operating conditions:  
5 – 40 °C / 41 – 104 °F

Storage conditions:  
-10 – +50°C / 14 – 122 °F

Humidity limitation for operating  
and storage

Temperature limitation for operating  
or storage

Ambient pressure limitation

Unique Device Identifier

Type number

CE 0598

CE Marking of Conformity

Dear Customer,

This Microlife fingertip pulse oximeter is a portable non-invasive device intended for spot-checking of the oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate of adults and pediatric patients. It is suitable for private use (at home, or on the go) as well as for use in the medical sector (hospitals, hospital-type facilities). It has been clinically proven to be of high precision during repeatability.

If you have any questions, problems or want to order spare parts please contact your local Microlife-Customer Service. Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit the internet at [www.microlife.com](http://www.microlife.com) where you will find a wealth of invaluable information on our products.

Retain instructions in a safe place for future reference.  
Stay healthy – Microlife Corporation!

#### Intended use:

The Fingertip Pulse Oximeter is a non-invasive device intended for spot-checking of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate. The portable fingertip device is indicated for adult and pediatric patients, for use at home and hospital environments.

#### Intended user:

Professional or lay person.

#### Table of contents

---

1. Important safety instructions
2. General description
3. Measurement principles
4. Directions for use
5. Inserting the batteries
6. ON/OFF button / Function button
7. Display Mode
8. Using the lanyard
9. Malfunctions and actions to take
10. Cleaning and disinfecting
11. Guarantee
12. Technical specifications

## 1. Important safety instructions

- Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.
- This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.

**Precaution:** Do not self diagnose or treat based on measurement results, please always consult a doctor. The product has no contraindications. The device has no side-effects if used correctly and residual risk is acceptable.

- **Never immerse this device in water or other liquids. For cleaning please follow the instructions in the « Cleaning and disinfecting» section.**
- Do not use this device if you think it is damaged or notice anything unusual.
- Never open this device.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical specifications» section.
- Protect it from:
  - water and moisture

- extreme temperatures
- impact and dropping
- contamination and dust
- direct sunlight
- heat and cold

- The function of this device may be compromised when used close to strong electromagnetic fields such as mobile phones or radio installations and we recommend a distance of at least 1 m (according to 60601-1-2 table 5). In cases where you suspect this to be unavoidable, please verify if the device is working properly before use.
- Do not use the device in an MRI or CT environment.
- This device is not intended for continuous monitoring.
- Do not sterilize this device using autoclaving or ethylene oxide sterilizing. This device is not intended for sterilization.
- If the device is not going to be used for a prolonged period the batteries should be removed.

Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

Use of this device is not intended as a substitute for a consultation with your doctor.

## 2. General description

Oxygen saturation indicates the percentage of hemoglobin in arterial blood that is loaded with oxygen. This is a very important parameter for the respiratory circulation system. Many respiratory diseases can result in lower oxygen saturation within human blood.

**Following factors can reduce oxygen saturation:** Automatic regulation of organ dysfunction caused by anesthesia, intensive postoperative trauma, injuries caused by some medical examinations. These situations may result in light-headedness, asthenia and vomiting. Therefore, it is very important to know the oxygen saturation of a patient so that doctors can detect problems in a timely manner.

## 3. Measurement principles

**Principle of this fingertip pulse oximeter:** A mathematical formula is established making use of Lambert Beer Law according to spectrum absorption characteristics of deoxygenated hemoglobin (Hb) and oxyhemoglobin (HbO<sub>2</sub>) in red and near-infrared zones.

**Operation principle of this device:** Based on full digital technology, the fingertip pulse oximeter non-invasively measures the actual content (oxygen saturation) of oxyhemoglobin (HbO<sub>2</sub>) in arterial blood using the optical transmittance method. The fingertip pulse oximeter measures the blood oxygen saturation and pulse rate of a human

body via finger artery. The fingertip pulse oximeter works by applying a sensor to a pulsating arteriolar vascular bed. The sensor contains a dual light source and photo detector. The one wavelength of light source is 660 nm, which is red light; the other is 905 nm, which is infrared-red light. Skin, bone, tissue and venous vessels normally absorb a constant amount of light over time. The photo detector in finger sensor collects and converts the light into an electronic signal which is proportional to the light intensity. The arteriolar bed normally pulsates and absorbs variable amounts of light during systole and diastole, as blood volume increases and decreases.

#### 4. Directions for use

---

1. Insert the batteries as described in the «Inserting the batteries» section.
2. Insert one finger (nail side up; index or middle finger is recommended) into the finger opening of the device. Be sure to fully insert the finger so that the sensors are completely covered by the finger.
3. Release the device allowing it to clamp down on the finger.
4. Press the ON/OFF button ① to turn the device on.
5. **Do not shake your finger during the test.** It is recommended that you do not move your body whilst taking a reading.

6. Your measurement values will appear on the screen after a few seconds.
7. Remove your finger from the device.
8. The device will automatically switch off after approx.  $10 \pm 2$  seconds after the finger is removed from the device.

The height of the bar graph ④ is an indication of the pulse and signal strength. The bar should be greater than 30 % for a proper reading.

The device must be able to measure the pulse properly to obtain an accurate SpO<sub>2</sub> measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO<sub>2</sub> measurement.

**The maximum application time at a single site should be less than 30 minutes**, in order to ensure correct sensor alignment and skin integrity.

#### Inaccurate measurements may occur if:

- The patient suffers from significant levels of dysfunctional hemoglobin (such as carboxyhemoglobin or methemoglobin).
- Intravascular dyes such as indocyanine green or methylene blue have been injected into the patient.

- Used in the presence of high ambient light (e.g. direct sunlight). Shield the sensor area with a surgical towel if necessary.
- There is excessive patient movement.
- The patient experiences venous pulsations.
- The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia.
- The patient is in cardiac arrest or is in shock.
- Fingernail polish or false fingernails are applied.

#### 5. Inserting the batteries

---

After you have unpacked your device, first insert the batteries. The battery compartment is on the bottom of the device. Remove the battery cover by sliding it in the direction shown. Insert the batteries (2 x 1.5 V, size AAA), thereby observing the indicated polarity.

Replace the batteries when the low power indicator ⑤ appears on the display.

Always replace both batteries at the same time.

#### 6. ON/OFF button / Function button

---

Press and release the ON/OFF button ① to turn on, hold the button for about one second. The device shows the

settings menu ⑩. Press or hold the ON/OFF button to perform corresponding operations. Hold it to set an item, or press it to switch an option or switch the display mode. Press means no more than 0.5 seconds, while hold means more than 0.5 seconds.

### Alert Sound Setting

Hold ON/OFF button ① while the device is switched on. Settings menu (interface 1) is displayed ⑩-a. Move «\*» to the corresponding option, and hold the function button to set **Alert to on** and set **Beep to off**. When **Alert** is set to **on** and the measured values of the blood oxygen saturation and pulse rate go beyond the upper limit or lower limit, the device gives off an alert sound. When **Beep** is set to **on**, a tick will be heard along with pulse beats during pulse rate measurement. While the «\*» symbol stays on the **Restore** option, hold the functional button to restore default settings.

### Brightness Setting

Press ON/OFF button ① to select the **Brightness** option and then hold the ON/OFF button to set the brightness to a value ranging from 1 to 5. The greater the value, the greater the brightness of the screen.

### Alert Range Setting

On settings menu (interface 2) ⑩-b, press the ON/OFF button ① to switch between options. On this interface, you can set the upper limit and lower limit of **SpO2 Alert** and

**PR Alert**. While the «\*» symbol stays on the +/- option, hold the functional button to set the option to + or -.

In + mode, select the corresponding option and hold the ON/OFF button to increase the upper or lower limit; in - mode, hold the ON/OFF button to decrease the upper or lower limit. Move «\*» to the **Exit** option, and hold the ON/OFF button to return to the monitoring interface.

## 7. Display Mode

When the device is switched on, shortly press the ON/OFF button ① to switch to another display mode to select your desired display mode ⑨. There are 6 different display modes. The default setting is mode 1.

## 8. Using the lanyard

1. Thread the thinner end of the lanyard through the hanging hole at the rear end of the device.
2. Thread the thicker end of the lanyard through the threaded end before pulling it tightly.

## 9. Malfunctions and actions to take

Description	Symptom/ Possible causes	Solutions
SpO <sub>2</sub> or pulse rate do not display normally.	<ol style="list-style-type: none"> <li>1. Finger is not inserted correctly.</li> <li>2. Patient SpO<sub>2</sub> value is too low to be measured.</li> <li>3. There is excessive illumination.</li> </ol>	<ol style="list-style-type: none"> <li>1. Retry inserting the finger.</li> <li>2. &amp; 3. Measure more times. If you determine the product is working correctly, consult your doctor.</li> </ol>
SpO <sub>2</sub> or pulse rate is shown unstable.	<ol style="list-style-type: none"> <li>1. Finger might not be inserted deep enough.</li> <li>2. Excessive patient movement.</li> </ol>	<ol style="list-style-type: none"> <li>1. Retry inserting the finger.</li> <li>2. Sit calmly and retry.</li> </ol>

Description	Symptom/ Possible causes	Solutions
The device cannot be powered on.	<ol style="list-style-type: none"> <li>1. No batteries or low battery power.</li> <li>2. Batteries are not installed correctly.</li> <li>3. The device may be damaged.</li> </ol>	<ol style="list-style-type: none"> <li>1. Replace the batteries.</li> <li>2. Remove and re-install the batteries.</li> <li>3. Contact your local Microlife-Customer Service.</li> </ol>
The display suddenly switches off.	<ol style="list-style-type: none"> <li>1. The device is automatically powered off, when no signal was detected after 10±2 seconds.</li> <li>2. The battery power is too low to operate.</li> </ol>	<ol style="list-style-type: none"> <li>1. Normal.</li> <li>2. Replace the batteries.</li> </ol>

## 10. Cleaning and disinfecting

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the silicone that touches the finger inside of the device. Also clean the finger being tested using alcohol before and after each test. Allow the device to dry thoroughly before use.

Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the device in water or other cleaning liquids.

## 11. Guarantee

This device is covered by a **2 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee. The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Battery.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website: [www.microlife.com/support](http://www.microlife.com/support) Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

## 12. Technical specifications

<b>Model:</b>	Pulse Oximeter JPD-500E
<b>Type:</b>	OXY 210
<b>Display:</b>	OLED display
<b>Display range:</b>	Oxygen saturation: 35 ~ 100% Pulse rate: 25 ~ 250 bpm
<b>Resolution:</b>	Oxygen saturation: 1% Pulse rate: 1 bpm
<b>Measurement precision:</b>	Oxygen saturation: ±3% (70 – 100%) No requirement: (≤ 69%) Pulse rate: ± 2 bpm

<b>Alert range:</b>	Upper limit: 50 – 100%
<b>Oxygen saturation:</b>	Lower limit: 50 – 100%
<b>Pulse rate:</b>	Upper limit: 25 – 250 bpm Lower limit: 25 – 250 bpm
<b>Alert error:</b>	Oxygen saturation: $\pm 1\%$ of the preset value Pulse rate: the greater of $\pm 10\%$ of the preset value and $\pm 5$ bpm PI (Perfusion Index) Weak PI Min. 0.2 %
<b>Operating conditions:</b>	5 – 40 °C / 41 – 104 °F 15 – 80% relative maximum humidity
<b>Storage conditions:</b>	-10 – +50°C / 14 – 122 °F 10 – 93% relative maximum humidity
<b>Automatic switch-off:</b>	Automatically shut down in 10 $\pm$ 2 seconds, when no or low signal is detected.
<b>Battery:</b>	2 x 1.5 V alkaline batteries; size AAA
<b>Battery lifetime:</b>	approx. 30 hours (using new batteries)

<b>Weight:</b>	42.5 g (including batteries)
<b>Dimensions:</b>	62 x 37 x 32 mm
<b>IP Class:</b>	IP22
<b>Reference to standards:</b>	EN ISO10993-1/-5/-10; IEC 60601-1; EN 60601-1-2; ISO 80601-2-61; EN 62304; EN 60601-1-6
<b>Expected service life:</b>	5 years (when used 15 times/day; 20 minutes for each measurement)

Technical alterations reserved.

**Note 1:** The fingertip pulse oximeter measurements are statistically distributed, only about two-thirds of pulse oximeter equipment measurements can be expected to fall within  $\pm 2\% A_{\text{rmse}}$  (Average-root-mean-square-error) of the value measured by a co-oximeter.

**Note 2:** The statistic conclusion of a controlled desaturation study which is guided by «ISO 80601-2-61, Annex EE, guideline for evaluating and documenting SpO<sub>2</sub> accuracy in human subjects». The statistic result displayed the accuracy distribution between the range of 70% – 100%, which is 2.83%.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

The device is a Class IIa medical device. **Device complies with European Regulation on medical devices EU MDR 2017/745.**

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## EMC Information-Guidance and Manufacturer's Declaration

### 1\* WARNING:

Use of this equipment adjacent to or stacked with other equipment should be avoided because it may result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

### 2\* WARNING:

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and may result in improper operation.

### 3\* WARNING:

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ME equipment, including cables specified by the manufacturer. Otherwise, it may result in the degradation of the performance of this equipment could result.

Declaration – Electromagnetic Immunity	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable

Declaration – Electromagnetic Immunity		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air
Electrical fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/ output lines	Not applicable
Surge IEC 61000-4-5	$\pm 0.5$ kV, $\pm 1$ kV line(s) to lines $\pm 0.5$ kV, $\pm 1$ kV, $\pm 2$ kV line(s) to earth	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycles	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m
NOTE: UT is the a.c. mains voltage prior to application of the test level.		

Declaration – Electromagnetic Immunity		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz	Not applicable
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m

Declaration – IMMUNITY to proximity fields from RF wireless communications equipment					
Immunity test	IEC60601 test level				Compliance level
	Test frequency	Modulation	Maximum power	Immunity level	
Radiated RF IEC 61000- 4-3	385 MHz	**Pulse Modulation: 18Hz	1.8 W	27 V/m	27 V/m
	450 MHz	*FM+ 5Hz deviation: 1kHz sine	2 W	28 V/m	28 V/m
	710 MHz, 745 MHz, 780 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m
	810 MHz, 870 MHz, 930 MHz	**Pulse Modulation: 18Hz	2 W	28 V/m	28 V/m
	1720 MHz, 1845 MHz, 1970 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m
	2450 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m
	5240 MHz, 5500 MHz, 5785 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m

**Note\* – As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.**

**Note\*\* – The carrier shall be modulated using a 50% duty cycle square wave signal.**