

WatchBP Analyzer

User Manual



Compatible models:

WatchBP Office series: BP3SK1-3B, TWIN200 ABI (with pulse wave presentation), WatchBP Office Vascular (TWIN200 VSR)

WatchBP O3 series: BP3SZ1-1, BP3MZ1-1, BP3MZ1-1A

microlife®

Content

1	Install the program.....	3
2	User Interface.....	4
	(1) Device Information	4
	(2) Patient information area	5
	(3) Measurement Folder area	5
	(4) Measurement Tab and Report Tab.....	5
3	Patient management.....	7
	(1) Add a new patient.....	7
	(2) Edit a patient.....	7
	(3) Delete a patient	7
4	Connecting the WatchBP device with the WatchBP analyzer.....	8
5	How to program device	10
	(1) Program ABPM.....	10
	(2) Program Office.....	12
6	How to perform a measurement controlled by WatchBP Analyzer	16
7	How to download measurement data from WatchBP devices	17
8	View the measurement data	18
9	How to generate PDF report and Excel report.....	20
	(1) Settings of 24-h ABPM report (for WatchBP O3 models).....	20
	(2) Settings of WatchBP Office report.....	20
	(3) Customization and report generation	21
10	How to view measurement history and delete a Measurement Folder	28
	(1) View measurement history.....	28
	(2) Delete a Measurement Folder	28
11	How to activate AFIB or Central BP feature	29
12	Code, note and potential cause and remedy.....	31

1 Install the program

Download the WatchBP Analyzer software from the Microlife website:

<http://www.microlife.com/support/software-professional-products>. Double click on downloaded installation program, the installation wizard dialog box will appear as shown in Fig. 1




Fig. 1 The installation wizard dialog box.

Click **“Install”** to start the installation of the WatchBP Analyzer. After successful installation, the dialog box will appear. (Fig 2)




Fig. 2

Click **“Finish”** to complete the installation.

2 User Interface




Fig. 3

(1) Device Information

Device information will show up when a WatchBP device is connected.

Device name : The name of the WatchBP model

Device ID : A unique device ID of the unit

Batteries : The area displays the voltage and status of batteries of the WatchBP O3. If the voltage is low and thus may not last for 24-hour session, the **“Renew battery”** message will be highlighted

Type : shows the version of the device

Options : shows the extra options of the device. AFIB – device with AF detector; Central BP – device with the possibility of taking central blood pressure measurement

[Functions]

Program device : For WatchBP O3 and Office, program the settings to the device.

The button is enabled when a compatible WatchBP device is connected. Click the button to program the device

*For WatchBP Vascular, this button cannot be pressed.

Measure : Start or Stop measurement (WatchBP Office only)

Set the measurement parameters, and then start a measurement

Download : Download BP data from the device memory to the software

(2) Patient information area

You can manage the patient information in the Patient information area. First select a patient from the patient list before starting to program the device, initiate a measurement or download measurement data from device to PC.

[Functions]

- New** : Create a new patient account
- Edit** : Edit the patient account
- Delete** : Delete patient data in the patient list
- Search** : Key in and search the name, ID or the date of birthday of the patients in the software database

(3) Measurement Folder area

The displayed Measurement Folder is for the selected patient. The new folders created according to the dates of the recorded measurement data in the device. Measurements on the same dates will be merged into the same folders.

[data type]

- WatchBP O3 : Measurements of WatchBP O3 series
- WatchBP Office : Measurements of WatchBP Office series
- WatchBP Office ABI : Measurements of WatchBP Office ABI (W/ waveform presentation)
- WatchBP Office Vascular : Measurements of WatchBP Office Vascular (W/ waveform presentation)
- 7 day : Measurements of HOME mode of WatchBP O3 Ambulatory (BP3MZ1-1)

(4) Measurement Tab and Report Tab

Measurement Tab contains measurements of a selected Measurement Folder.

[Caption]

- Date** : the date and time of the measurement
- SYS** : systolic blood pressure
- DIA** : diastolic blood pressure
- MAP** : mean arterial pressure
- PP** : pulse pressure
- cSYS** : central systolic blood pressure
- cDIA** : central diastolic blood pressure
- cPP** : central pulse pressure
- ABI** : Ankle brachial index

PWV : pulse wave velocity

Index of cycle: the number of measurements in a set

Limb : the specific limb of the measurement (only for simultaneous double cuffs measurements)

Posture : patient's position

HR : heart rate (beats per minute)

AFIB : atrial fibrillation

Exclude : option to exclude measurements from either PDF or Excel report

PVP wave: brachial pulse waveform of the measurement

CODE : events or errors

NOTE : brief description of the CODE, refer "Code, Note and potential cause and remedy" for detail

[Functions]

Setting Office Blood Pressure Threshold: Choose the threshold of systolic and diastolic from the drop-down menus. The values over the chosen threshold will be shown in red.

Report option : Click to hide AFIB result

Customization : You can type the name of centre/hospital, and upload an image file as your customized logo, the information you provide above will be displayed in the PDF report, folder, and XLSX report. You can choose the storage path here.

Generate PDF Report : Click to generate a PDF report for the selected Measurement Folder

Open Folder : Open the folder that contains the reports of the selected patient

Generate XLSX report : Click to generate an Excel report of the selected Measurement Folder.

3 Patient management

You can manage the patients' information through the patient information area (Fig 4).

Patient						Search
Name	ID	Gender	Age	DOB	Physician	
Green Lantern	GL	Male	30	1988.07.27	Dr. 1	New
Harley Quinn	HQ	Female	68	1950.07.02	Dr. Quinzel	Edit
Super Man	SM	Male	30	1988.09.05	doctor2	Delete

Fig. 4

(1) Add a new patient

1. Click **New** button, a dialog box appears (Fig 5).
2. Enter patient ID, Name, Sex and Date of Birth; Enter Physician ID, Email (optional) and Phone number (optional) then click **Okay** button to add a new patient account.

Patient Information

Patient

Sex Male Female

Physician

Fig. 5

(2) Edit a patient

Select a patient and click **Edit** button, the patient information dialog box with the record appears. Edit the information and click **Okay** button to save the changes.

(3) Delete a patient

If you would like to remove a patient from the list, select the patient's account and click **Delete** button. A confirmation screen pops up, and choose "**yes**" to delete the selected patient account from the list.

4 Connecting the WatchBP device with the WatchBP analyzer

(1) Connecting the device with WatchBP software cable:

Follow the instruction manual of the device to connect the device to the PC.

(2) Connecting the device using Bluetooth connection:

Prepare the device:

Model	Pairing – device	Ready to pair	Connected
WatchBP Office (BP3SK1-3B)	Press and hold the MODE button for 7 seconds	Bluetooth indicator flashes	Shows Bluetooth indicator on device. Displays device information on the Analyzer
WatchBP Office Vascular (TWIN200 VSR)	Press and hold the MODE button for 7 seconds	Bluetooth indicator flashes	Shows Bluetooth indicator on device. Displays device information on the Analyzer
WatchBP O3 (BP3SZ1-1)	Press and hold the Start/Stop button for 7 seconds	Bluetooth indicator flashes	Shows Bluetooth indicator on device. Displays device information on the Analyzer

**Bluetooth connection of the WatchBP Analyzer supports Microsoft Windows 10.*

Connecting the WatchBP device from the WatchBP Analyzer using Bluetooth:

Click **Bluetooth** at the upper right on the WatchBP Analyzer to view the Bluetooth device. Select the WatchBP device and click the **Connect** button to connect it with Analyzer. (Fig. 6)




Fig. 6

Device information will be shown up when the connection is successful.

*If the Bluetooth connection keeps failing, please reset the Bluetooth bonding by pressing and holding the Start/Stop button of the device for 7 seconds and start the connection again. Refer to user manual of the device for clearing bonding.

5 How to program device

** Program device will automatically clear all measurement data from the device. Make sure the measurement data on the device are downloaded to the PC before **program** the device!*

(1) Program ABPM

Click **Program device** button while WatchBP O3 is connected, shows ‘**Program ABPM**’ dialogue window.

Fig. 7a

Fig. 7b

Fig. 7c

Fig. 7d

<Setting measurement periods>

Click the number of periods you want to have.

Set the start and stop hours to the period by the drop-down menu.

Select the interval between each measurement.

Click the Central BP if you want to measure the central blood pressure.

<Setting Highest Inflation Pressure>

The “Highest Inflation Pressure” can be programmed to the device. The ‘Auto’ option is recommended because then the device will automatically inflate the cuff to the optimal cuff pressure. With the ‘Manual’ option you can select 180, 200, 220 or 240mmHg. The suggested Inflation Pressure is 30 to 40 mmHg above the expected systolic value of the patient. If the selected Highest Inflation Pressure selected is too low to measure a patient’s blood pressure it may result in re-pumping or erroneous blood pressure measurement values.

<Setting Ambulatory Option>

Select the ‘Hide BP reading’ to hide the readings of the measurement results during monitoring. This is the recommended option (default).

Select ‘Silent mode’ to disable the warning beep before starting the measurement. The patient will still get a warning 1 minute before the measurement by means of a short inflation of the cuff. The patient should keep the arm still during blood pressure measurement.

(2) Program Office (for programming the WatchBP Office device)

Fig. 8 presents the dialog window that is shown after clicking **Program device** button while the WatchBP Office (model No. BP3SK1-3B) is connected to the PC.




Fig. 8

The **number of measurements** can be set from 1 to 6 measurements. Both the **Count-down time before first measurement** and **Interval time between measurements** can be set at 15, 30, 60, 120, 180, 240 and 300 seconds. **Average without first measurement** can be selected if **Number of measurements** is 3 measurements or more. When selected, the first measurement will be discarded from the average blood pressure value.

<Child Mode (safeguard)>

On default this mode is on **Disable** – the device automatically inflates the cuff to an optimal cuff pressure for blood pressure measurement. While **Enable** is selected and confirmed by **Program**, the device automatically inflates the cuff to an optimal cuff pressure not higher than the pressure selected for **Child Mode**, however, the device may re- inflate the cuff to a pressure higher than the selected pressure if needed.

<Measurement options>

Hide BP readings and **Enable central BP measurement** can be programmed to the device. Please refer to the instruction manual of the device for detail of the features.

(3) Program Office ABI

Fig. 9 presents the dialog window that is shown after clicking **Measure** button while the WatchBP Office ABI (model No. TWIN200 ABI) is connected to the PC.




Fig. 9 Example of the selection of an arm (left), both arms (middle) and Ankle-brachial Index measurement (right).

Choose the intended limb(s) for the measurement at the top area.

Note: only allow the selection of an arm, both arms, or an arm with a leg (ABI measurement).

The **Number of measurements** can be set from 1 to 6 measurements when you want to measure the blood pressure of an arm or both arms. Both the **Count-down time before first measurement** and **Interval time between measurements** can be set at 15, 30, 60, 120, 180, 240 and 300 seconds.

Click **Measure** to start the measurement.

(4) Program Office Vascular

Fig. 10 presents the dialog window that is shown after clicking **Measure** button while the WatchBP Office Vascular (model No. TWIN200 VSR) is connected to the PC.




Fig. 10

Choose the intended limb(s) for the measurement at the top area.

The **Number of measurements** can be set from 1 to 6 measurements when you want to measure the blood pressure of an arm or both arms. Both the **Count-down time before first measurement** and **Interval time between measurements** can be set at 15, 30, 60, 120, 180, 240 and 300 seconds.

The **patient's height** can be set by typing or by adjusting the bar of height.

*The **patient's height** is crucial for baPWV. Please set the height correctly when taking the baPWV measurement.

<Child Mode (safeguard)>

On default this mode is on **Enable** – the device automatically inflates the cuff to an optimal cuff pressure for blood pressure measurement. While **Disable** is selected, the device inflates the cuff to an optimal cuff pressure not higher than the selected

pressure in the **Child Mode**, however, the device may re-inflate the cuff to a pressure higher than the selected pressure if needed.

<Measurement options>

Hide BP readings and **Enable Central BP Measurement** can be programmed to the device. Please refer to the instruction manual of the device for detail of the features. Click **Measure** to start the measurement.

6 How to perform a measurement controlled by WatchBP Analyzer

Connect the **WatchBP Office** blood pressure monitor to the WatchBP Analyzer using Microlife software cable or Bluetooth. **Program Device** first if you want to use a measurement program that differs from the one that was selected before. Click **Measure** button and then click **Yes** to confirm and to start taking measurements. The status of the measurement procedure is shown in the device information area (Fig. 11), measurement data are automatically downloaded after each measurement.

For **WatchBP Office Vascular**, Select parameters in **Measure** and click **Measure** at the bottom of the **Measure** screen. The status of the measurement procedure is shown in the device information area (Fig. 11), measurement data are automatically downloaded after each measurement.



Fig. 11

During the measurement process **Measure** button is switched to **Stop** button; Click **Stop** button if you would like to stop the measurement. (Fig. 12)




Fig. 12

7 How to download measurement data from WatchBP devices

Connect a WatchBP device to Analyzer using the WatchBP software cable or through a Bluetooth connection.* Select a patient by clicking anywhere of that patient bar and click the Download button to transfer the BP data from the device to the **Folder** of that patient in WatchBP Analyzer.

**Bluetooth connection of the WatchBP Analyzer supports Windows 10 operating system.*

If the ID of the patient selected on the software differs from the patient ID that has been programmed to the device, a message pops up so that the correct ID can be selected to assign the data. (Fig. 13, WatchBP O3 only). Select the patient and click **Continue** button.




Fig. 13

(HQ and SM are the patient ID examples)

If the device was programmed to also perform central blood pressure measurement, another dialogue shows up. Check the tick box before **transfer waveform data** and click **Continue** button to download the waveform. (Fig 14)




Fig. 14

Click the **Continue** button to start the download process. A new Measurement Folder with the date and time of the moment of downloading will be added to the assigned patient.

8 View the measurement data

Select a Folder named by the date of the measurement of the patient. Click **Measurement** Tab to view the measured records of the Folder. Fig. 15a is an example of a WatchBP O3 measurement list. Fig. 15b is an example of WatchBP Office measurement list. Fig. 15c is an example of measurements simultaneously taken by the WatchBP Office series with both cuffs. Fig. 15d is an example of measurements simultaneously taken by the WatchBP Office Vascular

The Limb indicates the specific limb of the measurement.

Measurement		Report											
	Date	SYS	DIA	HR	MAP	PP	cSYS	cDIA	cPP	AFIB	Exclude	PVP wave	CODE
3	07.30.2018 10:21	119	76	68	82	43	119	73	46		<input type="checkbox"/>		4,8
4	07.30.2018 10:22	119	78	68	86	41	120	75	45		<input type="checkbox"/>		4,8
5	07.30.2018 10:23	117	82	76	93	35	116	77	39	*	<input type="checkbox"/>		4,8
6	07.30.2018 10:25	115	78	72	87	37	116	75	41		<input type="checkbox"/>		4,8

Fig. 15a, an example of WatchBP O3 measurement list

Measurement		Report											
<input type="checkbox"/> Average without first measurement													
	Date	SYS	DIA	HR	MAP	PP	cSYS	cDIA	cPP	AFIB	PVP wave	CODE	NOTE
Average		120	76	69	86	44	121	72	49				
1	04.17.2019 10:44	128	77	66	85	51	130	71	59				
2	04.17.2019 10:45	114	73	70	79	41	115	71	44				
3	04.17.2019 10:47	118	78	70	93	40	119	75	44				

Fig. 15b, an example of WatchBP Office measurement list

Measurement		Report											
	Date	SYS	DIA	MAP	PP	ABI	Limb	HR	AFIB	Exclude	PVP wave	CODE	NOTE
1-1	05.20.2020 15:14	131	79	87	52		RightArm	93		<input type="checkbox"/>			
1-2	05.20.2020 15:14	140	92	98	48	1.07	RightLeg	93		<input type="checkbox"/>			
2-1	05.20.2020 15:15	129	72	82	57		RightArm	88		<input type="checkbox"/>			
2-2	05.20.2020 15:15	136	88	99	48	1.05	RightLeg	88		<input type="checkbox"/>			
3-1	05.20.2020 15:18	128	73	82	55		LeftArm	90		<input type="checkbox"/>			
3-2	05.20.2020 15:18	142	88	98	54	1.11	LeftLeg	90		<input type="checkbox"/>			

Fig. 15c, an example of the measurement list for simultaneous double arm measurement or ABI assessment

	Date	SYS	DIA	MAP	PP	cSYS	cDIA	cPP	ABI	PWV	Index of cycle	Limb	Posture	HR	AFIB	Exclude	PVP wave	CODE	NOTE
1	2020.12.02 15:29										1	LeftArm	---	124		<input type="checkbox"/>		ER 5,	Abnormal result,
2	2020.12.02 15:31	104	61	70	43	107	61	46			2	LeftArm	---	65		<input type="checkbox"/>			
3-1	2020.12.03 13:22	93	52	71	41				1.58	123	1	LeftArm	Supine	56		<input type="checkbox"/>			
3-2	2020.12.03 13:22	147	92	115	55				1.58	123	1	LeftLeg	Supine	56		<input type="checkbox"/>			
4-1	2020.12.03 13:24	92	51	63	41				1.59	123	2	LeftArm	Supine	58		<input type="checkbox"/>			
4-2	2020.12.03 13:24	146	87	106	59				1.59	123	2	LeftLeg	Supine	58		<input type="checkbox"/>			
5	2020.12.03 13:37	152	98	110	54	135	88	47			1	LeftArm	Supine	58		<input type="checkbox"/>			
6	2020.12.03 13:38	146	94	112	52						2	LeftArm	Supine	61		<input type="checkbox"/>		ER 15,	Abnormal central blood pressure reading,

Fig. 15d, an example of the WatchBP Office Vascular measurement list

*The first number indicates the order of recording. The dash and number after the

first number indicate that the blood pressure of different limbs was measured at the same time. (The examples are framed by red blocks)

Exclude:

If you select the tick box under **Exclude** the corresponding measurement value will exclude data from the analysis and the report. (Fig. 15a)

Average without first measurement:

You may check the box to average the data excluding first measurement (Fig. 15b)

Pulse Volume Plethysmographic (PVP) wave:

When double clicking on the data showing a small wave icon, the PVP wave of the corresponding measurement will be displayed (only for the central BP option). (Fig. 16a and 16b)




Fig. 16a PVP wave for an arm or both arms




Fig. 16b PVP wave for an arm with a leg

Code and Note:


Displays the codes of errors or events and brief description of the codes. Refer “Code, note and potential cause and remedy” for detail

9 How to generate PDF report and Excel report

Select a Measurement Folder of the patient. Click **Report** Tab to view the settings and generate report of the Measurement Folder

(1) Settings of 24-h ABPM report (for WatchBP O3 models)

When a WatchBP O3 Measurement Folder is selected, the **Report Tab** shows the report settings for WatchBP O3 24-hour ABPM measurement (Fig 17).



Setting Blood Pressure Threshold					
24-hour SYS	130	mmHg	24-hour DIA	80	mmHg
Day SYS	135	mmHg	Day DIA	85	mmHg
Night SYS	120	mmHg	Night DIA	70	mmHg
Setting Actual Awake / Asleep Time					
Awake Time	6		Asleep Time	22	

Fig. 17

(Sys indicates systolic blood pressure, Dia, diastolic blood pressure)

Setting Blood Pressure Threshold:

You may use the default threshold values (as recommended by the guidelines) or select another threshold value for 24h-hour, day and night time blood pressures.

Setting Actual Awake/ Asleep Time:

You can modify the awake and asleep time according to patient's actual lifestyle. The awake time can be selected from 00:00 to 23:00, asleep time can be selected from 00:00 to 23:00, too.

Report option:

Hide error and event messages to keep the report free of error messages (not recommended because the disadvantage is that you cannot see the cause of erroneous measurements from the report). **Hide AFIB result** can be used for e.g. children, young adults and pregnant women as is not recommended for these patient groups (Fig. 18).



Report Option	
<input checked="" type="checkbox"/>	Hide error and event messages
<input type="checkbox"/>	Hide AFIB result

Fig. 18

(2) Settings of WatchBP Office/Office Vascular report

When a WatchBP Office/Office Vascular Measurement Folder is selected, the **Report Tab** shows the report settings for WatchBP Office/Office Vascular.

Setting Office Blood Pressure Threshold:

Use the default threshold for defining hypertension or select another blood pressure threshold (Fig. 19).

Setting Office Blood Pressure Threshold			
SYS	140	▼	mmHg
DIA	90	▼	mmHg

Fig. 19

Report option

Hide error and event messages to keep the report free of error messages (This is not recommended because the disadvantage is that you cannot see the cause of erroneous measurements in the report). **Hide AFIB result** can be used for e.g. children, young adults and pregnant women as is not recommended for these patient groups. Check **Average without first measurement** to exclude the first measurement for averaging the blood pressure value. (Fig 20a, b).

Report Option	
<input checked="" type="checkbox"/>	Hide error and event messages
<input type="checkbox"/>	Hide AFIB result
<input type="checkbox"/>	Average without first measurement

Fig. 20a Options for WatchBP Office

Report Option	
<input type="checkbox"/>	Hide AFIB result

Fig. 20b Options for WatchBP Office Vascular

Customization and report generation


Enter the **physician, centre/hospital** and place a **Customized logo** to show on the PDF report, and choose the storage path (Fig 21).

Customization	
Physician	microlife
Centre/Hospital	Hospital
Customized logo	D:\Users\Morris.Huang\Desktop\squirtle.png Select Image File
Report Folder	C:\Microlife\WatchBP Analyzer\Report Path

Fig. 21

Generate PDF Report:

(1) For WatchBP O3 models: Click the **Generate PDF Report** button, an **Auto comment** dialog appears with the comments automatically generated. (Fig. 22) Comments can be modified or added manually.



Auto comment


Daytime Normotension, Isolated 24-h Diastolic Hypertension, Isolated NightTime Diastolic Hypertension, White Coat Hypertension, Non-Dipper

Save Cancel

Fig. 22

Click **Save** button to generate the report. A dialog appears after the report is generated successfully. Click **Yes** button to open the report.

(2) For WatchBP Office models/Office Vascular: Click the **Generate PDF Report** button, then choose the **Yes** button to open the report. (Fig. 23)



Generating Report

The report is ready. Do you want to open it?

Yes No

Open Folder:

Click **Open Folder** button to open the folder that contains the reports of the patient selected. By default, the reports are stored in the folder:

C:\Microlife\WatchBP Analyzer\Report\Patient ID

The file name of the report is built with the patient ID, device name, date and time the report is created for example: PatientID_WatchBPO3_2018_11_05_10_54.pdf.

Generate XLSX report:

Click the **Generate XLSX report** button to generate an Excel report of the selected Measurement Folder.

Example of the WatchBP O3, 24-hour ABPM report (Fig. 24).




Fig. 24

Example of the WatchBP Office report (Fig. 25a, b)




Fig. 25a

Example of the WatchBP Office report (Fig. 25b)





Fig. 25b

Example of the WatchBP Office Vascular report (Fig. 26a)

WatchBP Office											
Hypertension management				Patient ID	eee			Health facility			Hospital
				Name	www			Physician			microlife
				Sex	Female			Email			haha
				Age	30			Phone			777
				Date of birth	1990.12.01			Study date			2020.12.01
List of individual readings											
Right Arm											
DATE	TIME	SYS	DIA	MAP	PP	cSYS	cDIA	cPP	HR	Posture	AFIB
Average		-	-	-	-	-	-	-	-	---	-
Left Arm											
DATE	TIME	SYS	DIA	MAP	PP	cSYS	cDIA	cPP	HR	Posture	AFIB
2020.12.01	14:19	126	55	70	71	124	56	68	57	---	-
Average		126	55	70	71	124	56	68	57	---	-

Example of the WatchBP Office Vascular report (Fig. 26b)



Page 1 of 1

10 How to view measurement history and delete a Measurement

Folder

(1) View measurement history

Select a patient in the patient list. The Measurement Folder of the patient appears in the **Measurement Folder area**

Click Measurement Folder, then the measurements of the selected Folder are displayed on the **Measurement Tab**

(2) Delete a Measurement Folder

Right click the mouse on a Folder, appears the Delete option (Fig. 27). Click **Delete** to remove the Folder from the database.




Fig. 27

11 How to activate the device

Click **About** button on the upper right on the WatchBP Analyzer to view the detail of the device. (Fig. 28)




Fig. 28

AFIB or Central BP feature

The Atrial Fibrillation Detector (AFIB) and Central Blood Pressure measurement of the WatchBP O3 (model BP3SZ1-1) and WatchBP Office (model BP3SK1-3B) can be activated through the WatchBP Analyzer. (Fig. 29) There are three different versions of the device:

- **Standard:** standard ABPM
- **AFIB:** standard ABPM with Microlife Atrial Fibrillation Detector
- **AFIB & Central BP:** standard ABPM with Microlife Atrial Fibrillation Detector and Central Blood Pressure measurement




Fig. 29

ABI or PWV & Central BP feature

The ankle-brachial index (ABI), brachial-ankle pulse wave velocity (PWV), and Central Blood Pressure measurement of the WatchBP Office Vascular (model TWIN200 VSR) can be activated through the WatchBP Analyzer. (Fig. 30) There are three different versions of the device:




Fig. 30

Activation key is required specifically to match the device ID for the AFIB or Central BP feature activation. Please contact Microlife or the local distributor for additional information. Please click on the question mark (?) button at the right side of the **Activate** button to copy the device information. Contact Microlife or the local distributor and pass the information for the activation key. (Fig. 31)




Fig. 31

Click the **Activate** button specific for **AFIB, AFIB & Central BP, ABI, and PWV & Central BP**. then a dialog appears. Insert the received activation key and click **Confirm** button to activate the feature on the device. (Fig. 32)



Fig. 32

12 Code, note and potential cause and remedy

The CODE is shown at measurement list on Measurement Tab, PDF and Excel report.

The NOTE (the message provided by the software) is shown on Analyzer.

CODE	NOTE	Potential cause and remedy
1	Pill record	
2	device power on	
3	device power off	
4	Low battery	
5	Empty battery	
6	Re-try measurement	
7	Measurement manually terminated	
8	Manual measurement	
Er 1	Signal too weak	The pulse signals on the cuff are too weak. Reposition the cuff and repeat the measurement.
Er 2	Error signal	During the measurement, error signals were detected by the cuff, caused by e.g. movement or muscle tension. Repeat the measurement, keeping your arm still.
Er 3	No pressure in the cuff	An adequate pressure cannot be generated in the cuff. A leak may have occurred. Replace the cuff if necessary. Repeat the measurement.
Er 5	Abnormal result	The measuring signals are inaccurate, and no result can therefore be displayed. Read through the checklist for performing reliable measurements and then repeat the measurement.
Er 11	Signal too weak during central blood pressure measurement	The pulse signals on the cuff are too weak. Re-position the cuff and repeat the measurement.
Er 12	Error signal during central blood pressure measurement	During the measurement, error signals were detected by the cuff, caused by e.g. movement or muscle tension. Repeat the measurement, keeping your arm still
Er 13	Cuff pressure errors during central blood pressure measurement	An adequate pressure cannot be generated in the cuff. A leak may have occurred. Check if the cuff is correctly connected and is not too loose. Replace the batteries if necessary. Repeat the measurement
Er 15	Abnormal central blood pressure reading	The measuring signals are inaccurate and no result can therefore be displayed. Read through the checklist for performing reliable measurements and then repeat the measurement.
Er 21	Error signal during collecting pulse wave signals	Check if the cuff is correctly connected and is not too loose. Replace the blood pressure cuff if necessary. Repeat the measurement.
Er 23	Cuff pressure errors during collecting pulse wave signals	Check if the cuff is correctly connected and is not too loose. Replace the blood pressure cuff if necessary. Repeat the measurement.

Er 25	Abnormal result of baPWV reading	Check if the cuff is correctly connected and is not too loose. Replace the blood pressure cuff if necessary. Repeat the measurement.
Er F	The device has gone into "single fault condition"	Re-position the cuff and repeat the measurement. Replace the batteries if necessary. If the error persists, contact Microlife or local distributor
Er A	Flash memory error	Possible hardware fault. Try again. If the error persists, contact Microlife or local distributor
HI	Pulse or cuff pressure too high	The pressure in the cuff is too high (over 299 mmHg) OR the pulse is too high (over 239 beats per minute). Relax for 5 minutes and repeat the measurement.
LO	Pulse too low	The pulse is too low (less than 30 beats per minute). Repeat the measurement.