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---

- ① Piston compressor
- ② Power lead
- ③ ON/OFF Switch
- ④ Air filter compartment
- ⑤ Nebuliser
  - A: Advanced treatment (position A)
  - B: Short treatment (position B)
- ⑥ Nasal washer
- ⑦ Air tube
- ⑧ Mouthpiece with valve
- ⑨ Nose piece
- ⑩ Adult face mask
- ⑪ Child face mask
- ⑫ Air filter
- ⑬ Filling the nebuliser
- ⑭ Assembling and usage of the nasal washer

This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products. Correct disposal helps to protect the environment and human health.

Read the instructions carefully before using this device.

Type BF applied part

Class II equipment

Serial number

Reference number

Manufacturer

Date of manufacture

ON

OFF

Protection against solid foreign objects and harmful effects due to the ingress of water

Medical device

Caution

Single patient multiple use (for applied parts only)

Humidity limitation for operating **and** storage

Temperature limitation for operating **or** storage

Ambient pressure limitation

Unique Device Identifier

Type number

CE Marking of Conformity

### Intended use:

The device is an aerosoltherapy system suitable for domestic use. The device is designed for the production of compressed air to operate a small volume nebuliser for the production of medical aerosol for lung respiratory disorders.

### Patient population:

The device is intended for use with children from 2 years old, adolescent and adult patients.

## Intended user:

The use of the device does not require a specific knowledge or professional ability. The patient is the intended operator except in case of child and patient that required special assistance.

## Indications:

Lung acute or chronic diseases of respiratory disorders organs, or inflammation of the upper respiratory system.

## Contra-indications:

The device is not indicated to be used with quick-relief medications during life-threatening asthma attacks. No contraindications exist to the administration of aerosols by inhalation. Contraindications related to the medicine used must be checked on the medicine package leaflet. Consult your physician in case of doubts.

Dear Customer,

This nebuliser is an aerosol therapy system suitable for domestic use. This device is used for the nebulisation of liquids and liquid medication (aerosols) and for the treatment of the upper and lower respiratory tract.

If you have any questions, problems or want to order spare parts please contact your local Microlife-Customer Service. Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit the internet at [www.microlife.com](http://www.microlife.com) where you will find a wealth of invaluable information on our products.

Stay healthy – Microlife AG!

## Table of contents

---

<b>1. Important safety instructions</b>
<b>2. Preparation and usage of this device</b>
<b>3. NEB 2 in 1 - Aerosol therapy</b>
Advanced treatment
Short treatment
<b>4. Nasal washer</b>
<b>5. Cleaning and disinfecting</b>
<b>6. Maintenance, Care and Service</b>
Replacement of the nebuliser
Replacement of the air filter
<b>7. Malfunctions and actions to take</b>
The device cannot be switched on
The nebuliser functions poorly or not at all
<b>8. Guarantee</b>
<b>9. Disposal</b>
<b>10. Technical specifications</b>

## 1. Important safety instructions

---

- Use the device only as described in this manual and therefore as an aerosoltherapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
- Keep the user manual for future reference.
- Do not operate the unit in presence of any anesthetic mixture inflammable with oxygen or nitrogen protoxide.
- The correct functioning of the equipment can be affected by electromagnetic interferences which exceed the limits indicated by the European standards in force. In case this device interferes with other electrical devices, move it and plug it into a different power socket.
- In case of failure and/or malfunction, read the «Malfunctions and actions to take» section in the user manual. Do not handle or open the compressor housing.
- For repair operations address only to a technical service centre authorized by the manufacturer and require the use of original spare parts. The non-observation of the above mentioned indications can compromise the device safety.
- Comply with the safety regulations concerning the electrical devices and in particular:
  - use only original applied parts and components;
  - never submerge the unit in water;
  - never wet the device, it is not protected against water penetration;
  - never touch the unit with wet or moist hands;
  - do not leave the unit exposed to the weather elements;
  - place the unit on a stable and horizontal surface during its operation;
  - the use of this device by children and people with disabilities requires always the close supervision by an adult with full mental faculties;
  - do not pull the power cord or the device itself to unplug it from the power socket;
  - the power plug is the separation element from the grid power; keep the plug accessible when the device is in use.

- Before plugging in the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
- In case the power plug provided with the device does not fit your wall socket, contact qualified personnel for a replacement plug with that of a suitable one. In general, the use of adapters, simple or multiple, and/or extension cables is not recommended. If their use is indispensable, it is necessary to use types complying with safety regulations, paying attention that they do not exceed the maximum power limits, indicated on adapters and extension cables.
- Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not in operation.
- The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause damage to persons, animals or objects, for which the manufacturer cannot be held responsible.
- The power cord of this device cannot be replaced by the user. In case of a power cord damage, address to a technical service center authorized by the manufacturer for its replacement.
- The power supply cord should always be fully unwound in order to prevent dangerous overheating.
- Before performing any maintenance or cleaning operation, turn off the device and disconnect the plug from the main supply.
- Some parts of the unit are so small that they may be swallowed by children; keep the equipment out of children's reach.
- If you decide not to use the device any longer, it is recommended to dispose of it according to the current regulations.
- Make sure to:
  - use this device only with medicines prescribed by your doctor;
  - make the treatment using only the applied parts recommended by your doctor depending on the pathology;
  - use the nosepiece applied part only if expressly indicated by your doctor and paying attention to NEVER introduce the bifurcations in the nose, but only bringing them as close as possible.
- Check in the medicine instruction leaflet for possible contraindications for use with common aerosol therapy systems.
- To avoid strangulation and entanglement, keep cable and air tubes out of reach of young children.
- Do not position the equipment so that it is difficult to operate the disconnection device.
- Nebuliser and applied parts are single patient use. Device is multi-patient use.
- This nebulising system is not suitable for use in anaesthetic breathing system or a ventilator breathing system. The patient cannot undergo MRI scan while using this device.
- This device is designed to nebulise solution and suspension liquids.
- Do not use this device if you think it is damaged or notice anything unusual.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical specifications» section.
- The use of this device is not a replacement for visiting the doctor. In case of no health improvement after treatment consult your doctor again.

## 2. Preparation and usage of this device

- The device must be checked before each use in order to identify any malfunctions and/or damage caused during transportation and/or storage. During inhalation, sit upright and relaxed to avoid compressing the airways and impairing treatment effectiveness. The applied parts must be used only with a single patient, it is not recommended to use them with several patients.
  - After unpacking the device, check it for visible damages or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check if the applied parts are intact.
  - Before using the device, proceed with the cleaning operations as described in the chapter «Cleaning and disinfecting».
1. Assemble the nebuliser kit (14). Ensure that all parts are complete.
  2. Fill the nebuliser with the inhalation solution as per your doctor's instructions. Ensure that you do not exceed the maximum level.
  3. Connect the nebuliser with the air tube (7) to the compressor (1) and plug the power lead (2) into the socket (230V 50 Hz AC).
  4. To start the treatment, set ON/OFF switch (3) into the «I» position.
    - The mouthpiece (8) gives you a better drug delivery to the lungs.
    - Choose between adult (10) or child face mask (11) and make sure that it encloses the mouth and nose area completely.

- Use all applied parts including the nose piece ⑨ as prescribed by your doctor.
- 5. During inhalation, sit upright and relaxed at a table and not in an armchair, in order to avoid compressing your respiratory airways and impairing the treatment effectiveness. **Do not lie down while inhaling.** Stop inhalation if you feel unwell.
- 6. Inhale the aerosol solution using the prescribed applied part.
- 7. Once you have finished the treatment, switch the appliance off by pressing the ON/OFF switch ③.
- 8. Empty the remaining medication from the nebuliser and clean the device as described in the section «Cleaning and disinfecting»

 The device requires no calibration.

 Tampering with the device is absolutely forbidden.

 No modification to the device is permitted.

### 3. NEB 2 in 1 - Aerosol therapy

The advanced «NEB 2 in 1» aerosol therapy allows effective treatments of the medium and lower respiratory tract (asthma, chronic bronchitis, cystic fibrosis etc.).

#### Advanced treatment

This treatment option allows the maximum amount of medication to be delivered thus reaching the respiratory capacity of each patient. The valve system optimally adapts the flow of medication during inhalation, reducing medication loss during exhalation. Use the breath-enhanced valve system (synchronized with breathing) to optimize the therapeutic effects of the medication and to get the maximal benefit.

- ▶ Close the valve on the nebuliser top ⑤-A (position A) and use the mouthpiece with expiratory valve ⑧.
- ▶ Fill with medication ⑬ and connect the nebuliser top to the nebuliser bottom by rotating it clockwise.
- ▶ Follow the usage steps as described in «Section 2.».

 To enable the synchronized with breathing function of the nebuliser, use only the mouthpiece with the expiration valve ⑧. Do not use with masks or nose piece.

#### Short treatment

The therapeutic session will be reduced if the valve is open (position B) ⑤-B.

 A synchronized nebulisation is not possible using this treatment option.

### 4. Nasal washer

The nasal washer ⑥ is a medical device suitable for cleaning nasal cavities. The nasal washer helps to moisturize the nasal mucosa and creates a micronized jet for the treatment of respiratory problems of the upper airways (allergic and nonallergic rhinitis, rhinosinusitis, nasal polyps etc.).

1. Remove the nasal washer cover ⑭-1 by rotating counter clockwise.
2. Fill the container of the nasal washer with isotonic solution ⑭-2. Ensure that you do not exceed the maximum level.
3. Close the cover tightly by rotating clockwise ⑭-3.
4. Connect the nasal washer with the air hose ⑭-4 to the compressor ① and plug the power lead ② into the socket (230V 50 Hz AC).
5. Close the non-treated nostril with your finger.
6. Position the nasal washer cover on nostril and inhale through the nose ⑭-5.
7. Switch the ON/OFF switch ③ to position «I» to turn on the device and place your finger over the aperture on the nasal washer to start the treatment ⑭-6.

 We recommend changing from one nostril to the other every 10 – 20 seconds.

 Children and people in need of care must be assisted for usage.


 Only use with isotonic saline solution.

## 5. Cleaning and disinfecting

<b>Device(s): Aerosoltherapy nebuliser applied parts</b>	
	<p><b>Before the first use and after each treatment</b> Follow carefully the cleaning and disinfection instructions of the applied parts as they are very important for the device performances and the therapy success.</p> <ul style="list-style-type: none"> <li>• <b>Use original applied parts only.</b></li> <li>• <b>Do not clean or disinfect the air tube.</b></li> <li>• <b>The applied parts cannot be cleaned and disinfected by automated method.</b></li> </ul>
<b>Limitations on processing</b>	<p>The nebuliser must be replaced after a long period of inactivity, in case it shows deformations or breakings, or when the nebuliser nozzle is obstructed by dry medicine, dust, etc. We recommend to replace the nebuliser after a period between 6 months and 1 year, depending on the usage. The maximum times of cleaning and disinfection of the nebulizer are 360 times.</p> <p>The maximum times of cleaning and disinfection of the mask, the nosepiece and the mouthpiece are 360 times.</p>

<b>Instructions</b>	
<b>Preparation before cleaning</b>	<ul style="list-style-type: none"> <li>• Detach the air tube from the nebuliser.</li> <li>• Detach the mouthpiece or nosepiece or mask from the nebuliser.</li> <li>• Make sure that all residual volume is removed from the nebuliser.</li> <li>• Disassemble the nebuliser  by turning the top counterclockwise and remove the medicine conduction cone.</li> </ul> <p><b>Note:</b> To avoid microbial growth and residual drug drying, clean and disinfect the applied parts immediately after each use.</p>
	

<b>Cleaning</b>	<p>Briefly rinse all parts in running tap water for at least 10 seconds beforehand.</p> <p>Mix little dish detergent and warm tap water (i.e. FAIRY brand dish washing liquid for hand washing in a ratio of 2 ml : 1 L) in a clean container.</p> <p>Immerse the components of disassembled nebuliser, the mask, the mouthpiece and the nosepiece into the mixed water for about 5 min. And then scrub all the surface of all the components with clean and small brush for at least 8 times.</p> <p>Afterwards rinse all parts thoroughly in running tap water at least 30 seconds, to fully remove any possible dish washing liquid residue.</p>
<b>Disinfection</b>	<ul style="list-style-type: none"> <li>• After cleaning, disinfect all of the disassembled parts (only parts that have been cleaned can be disinfected effectively).</li> <li>• Boil the disassembled nebuliser, the masks, the mouthpiece and the nosepiece for 5 minutes in boiling tap water.</li> </ul>
<b>Drying</b>	<ul style="list-style-type: none"> <li>• Reassemble the nebuliser components and connect it to the air-outlet, switch the device on and let it work for 10 – 15 minutes.</li> <li>• Let all parts dry completely before reassembling and reuse to avoid risk of germ grow.</li> </ul> <p>Care needs to be taken not to contaminate the parts after they have been cleaned and disinfected. Contamination can be avoided by good hand washing and not touching the inside sections of the device when laying them out to dry or when reassembling.</p>
<b>Inspection</b>	<p>Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discolored parts.</p>
<b>Packaging</b>	<p>Pack dry parts in a clean and sealed container when not in use. Do <b>NOT</b> pack wet or damp parts.</p>

<b>Storage</b>  	Storage conditions refer to «Technical specifications». <b>Note:</b> re-clean and re-disinfect the parts if they are stored more than one day.
<b>Transportation</b>  	After cleaning and disinfection, always transport the parts with clean and sealed container. Contamination can be avoided by good hand washing and not touching the inside sections of the parts when taking out and re-assembling the parts for use.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

## 6. Maintenance, Care and Service

Order all spare parts from your dealer or pharmacist, or contact Microlife-Service (see foreword).

### Replacement of the nebuliser

Replace the nebuliser ⑤ after a long period of inactivity, in cases where it shows deformities, breakage, or when the vaporiser head ⑤-a is obstructed by dry medicine, dust, etc. We recommend to replace the nebuliser after a period between 6 months and 1 year depending on the usage.

Only use original nebulisers!

### Replacement of the air filter

In normal conditions of use, the air filter ⑫ must be replaced approximately after 200 working hours or after each year. We recommend to periodically check the air filter (10 – 12 treatments) and if the filter shows a grey or brown colour or is wet, replace it. Extract the filter and replace it with a new one.

Do not try to clean the filter for reusing it.

The air filter shall not be serviced or maintained while in use with a patient.

Only use original filters! Do not use the device without filter!

## 7. Malfunctions and actions to take

### The device cannot be switched on

- Ensure the power lead ② is correctly plugged into the socket.
- Ensure the ON/OFF switch ③ is in the position «I».

### The nebuliser functions poorly or not at all

- Ensure the air tube ⑦ is correctly connected at both ends.
- Ensure the air tube is not squashed, bent, dirty or blocked. If necessary, replace with a new one.
- Ensure the nebuliser ⑤ is fully assembled and the vaporiser head ⑤-a is placed correctly and not obstructed.
- Ensure the required medication has been added.

## 8. Guarantee

This device is covered by a **5 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Applied and wearing parts: Nebuliser, masks, mouthpiece, nose piece, tube, filters, nasal washer (optional).

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website:

[www.microlife.com/support](http://www.microlife.com/support)

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

## 9. Disposal

Item	Component	Disposal Description
1	Device ①	The component mainly includes plastic and electronic components. All comply with RoHS and REACH, and all could be safely disposed. This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products.
2	Applied parts: air tube ⑦, nebuliser ⑤, nose piece ⑨, mouth-piece ⑧, masks ⑩ / ⑪	The components are plastic. All comply with ROHS and REACH. Before disposal, all the applied parts should be cleaned according to manual instructions and then disinfected by boiling for 5 minutes.

## 10. Technical specifications

**Model:** NEB PRO  
**Type:** GCE841

**AEROSOL PERFORMANCES ACCORDING TO EN ISO 27427:2019 based on adult ventilatory pattern with Sodium Fluoride (NaF):**

**Aerosol output:** 0.25 ml  
**Aerosol output rate:** 0.089 ml/min.  
**Percentage of fill volume emitted per min:** 4.5 %  
**Residual Volume:** 1.74 ml  
**Particle size (MMAD):** 3.5 µm  
**GSD (geometric standard deviation):** 2.9 µm  
**RF (respirable fraction < 62 % 5 µm):**

**Large particle range (> 5 µm):** 38 %  
**Middle particle range (2-5 µm):** 34 %  
**Small particle range (< 2 µm):** 28 %  
**Operating air flow:** 5.31 l/min.  
**Acoustic noise level:** 52 dBA  
**Power source:** 230V 50 Hz AC  
**Current:** ≤ 1000 mA  
**Power lead length:** 1.4 m  
**Nebuliser capacity:** min. 2 ml; max. 8 ml  
**Operating limits:** Continuous use  
**Operating conditions:** 10 - 40 °C / 50 - 104 °F  
 10 - 95 % relative maximum humidity  
 700 - 1060 hPa Atmospheric pressure  
**Storage and shipping conditions:** -25 - +70 °C / -4 - +140 °F  
 10 - 95 % relative maximum humidity  
 700 - 1060 hPa Atmospheric pressure  
**Weight:** approx. 1765 g  
**Dimensions:** 289 x 180 x 127 mm  
**IP Class:** IP21  
**Reference to standards:** EN ISO 27427; EN 60601-1; EN 60601-1-2; EN 60601-1-6; IEC 60601-1-11  
**Expected service life:** 1000 hours

Class II device as regards protection against electric shocks. Nebuliser, mouthpiece and masks are type BF applied parts.

The technical specifications may change without prior notice.

Please report any serious incident that has occurred in relation to the device, injury or adverse event to the local competent authority and to the manufacturer or to the European authorised representative (EC REP).

Vigilance contact point:  
<https://ec.europa.eu/growth/sectors/medical-devices/contacts>  
**AEROSOL CHARACTERISTICS IN ACCORDANCE WITH REGULATION EN ISO27427**

Standards applied:

Electric safety standards CEI EN 60601-1 electromagnetic compatibility according to CEI EN 60601-1-2.

The device is a Class IIa medical device. **Device complies with European Regulation on medical devices EU MDR 2017/745.**

**Important information regarding Electromagnetic Compatibility (EMC):** This product, manufactured by Globalcare Medical Technology Co., Ltd., conforms to the Electromagnetic Compatibility (EMC) standard EN 60601-1-2:2015/A1:2021. Further documentation in compliance with this EMC standard is available from Microlife on [www.microlife.com/electro-magnetic-compatibility](http://www.microlife.com/electro-magnetic-compatibility).

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