 Microlife Corporation
 9F, No. 431, RuiGuang Road,
 NeiHu, Taipei, 114, Taiwan, China
www.microlife.com

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


microlife®

- ① ON/OFF button
- ② Display
- ③ Data compartment
- ④ Measuring tube
- ⑤ Mouthpiece
- ⑥ Memory scroll buttons
- ⑦ Battery compartment

Display

- ⑧ Date/Time
- ⑨ Traffic light display
- ⑩ AM/PM
- ⑪ PEF/FEV1 result
- ⑫ Active Bluetooth®
- ⑬ Connection
- ⑭ Data memory
- ⑮ PEF/FEV1 value
- ⑯ PEF/FEV1 units
- ⑰ Ready for measuring
- ⑱ Battery display

Data compartment

- ⑲ Time button
- ⑳ USB Port

Table of contents

1. **Introduction**
 - Document scope
 - Disclaimers
2. **Important information**
 - Device description
 - Intended purposes
 - Intended user
 - Intended patient
 - Intended use environment and conditions
 - Indications
 - Contra-indications
 - Warning
 - Caution
 - Electromagnetic compatibility
 - Adverse events and reporting
3. **Device information**
 - Package contents
 - Device accessories
 - Device optional parts
4. **Device installation and setup**
 - Mouthpiece installation
 - Battery installation
 - Battery replacement
 - Device date and time setup
 - Computer software installation
 - Mobile software installation
 - Settings function operation
 - Bluetooth® automatic activation
 - PEF traffic light reference selection
 - PEF traffic light reference value setup
 - Device date and time setup

5. **Measurement preparation**
6. **Measurement operation**
 - Standalone measurement mode operation
7. **Measurement interpretation**
 - Self-assessment using the PEF traffic light indicator
 - Baseline reference of the PEF traffic light indicator
8. **Memory function operation**
 - Review measurement readings in memory
 - Memory capacity low
 - Memory full
 - Delete the most recent reading in memory
 - Delete all readings in memory
9. **Bluetooth® function operation**
 - Turn on/off the Bluetooth® datalink function
 - Software data transfer via Bluetooth® datalink
10. **USB datalink function operation**
 - Turn on/off the USB datalink function
 - Software data transfer via USB datalink
11. **Device error and troubleshooting**
 - Device error and troubleshooting
 - Common issues and troubleshooting
12. **Device maintenance and disposal**
 - Cleaning
 - Cleaning the mouthpiece and the measuring tube
 - Cleaning the device
 - Storage
 - Calibration and support
 - Disposal
13. **Specifications and compliance**
 - Technical specifications
 - Compliance information
 - Compliant standards
14. **Supplement information for users and patients**
 - Information about measurement and monitoring of peak flow values
 - Which PEF values are normal?
 - Guarantee
 - Symbols and definitions

1. Introduction

Document scope

Please read this instructions for use, the device label, and all other information provided with this product, before using the device.

This instruction of use contains important information about the device safety, performance, and operation.

Disclaimers

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Microlife Corp. is under license. Other trademarks and trade names are those of their respective owners.

Apple, the Apple logo, iPad, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android and Google Play are both trademarks of Google Inc.

Microlife® is a registered trademark of Microlife Corporation. Trademarks and trade names are those of their respective owners.

2. Important information

Device description

MicroLife digital peak flow meter is a medical device that utilizes the principles of flow turbine, with optical interruption detection and digital signal processing, to compute and provide a measurement of peak expiratory flow and other pulmonary functional parameters. Measuring your pulmonary functions, as instructed by your physician or other qualified healthcare professionals, is useful for monitoring your airway and lung conditions or diseases.

Intended purposes

This device is intended to measure the following pulmonary function parameters:

- Peak Expiratory Flow (PEF): The maximum flow rate a person can exhale forcibly after taking a deep as possible breath.
- Forced Expiratory Volume 1-Second (FEV1): The air volume of forcible exhalation a person can exhale in 1 second after taking a deep as possible breath.

Intended user

The device is intended to be operated by adult users with the vision and motor functions as well as the literacy and basic education, capable of understanding the content of this instructions for use and operating general household electrical appliances, for use of self-monitoring or monitoring of other subjects.

Intended patient

The device is intended for measurements of pediatric and adult subjects.

Intended use environment and conditions

The device is intended for use in a home healthcare environment (for example general household) and professional healthcare environment by patients (for self-measurement) or by a care giver.

Indications

This device offers pulmonary function measurements for use in screening, diagnosis, and monitoring of the following indications under healthcare professional supervision:

- Asthma
- Chronic Obstructive Pulmonary Disease (COPD)
- Other respiratory diseases or conditions affecting airway flow, as instructed by healthcare professionals

Contra-indications

The device is not suitable for subjects who are unable to conduct the measurement correctly, for example:

- Patients unconscious or in disabling conditions that does not allow correct self-measurement operating conditions.
- Bedridden patients unable to stand or sit upright as required for correct measurement.
- Patients whose mouth cannot cover the mouthpiece as required for correct measurement.

Warning

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient, or cause damage to the device or other property.

- DO NOT use this device for purposes beyond described in this instructions for use. The manufacturer cannot be held liable for damage caused by incorrect application.
- Inspect the device and other parts for damage. DO NOT USE the device or parts if they appear damaged or operating abnormally.
- DO NOT use this device in oxygen rich environment or near flammable gas.
- DO NOT use this device in a moving vehicle (for example in a car or on an aircraft).
- Keep the device away from children and people incapable of operating the device. Beware of the risks of accidental ingestion of small parts and of strangulation with the cables and tubes of this device and accessories.
- DO NOT let children operate the device alone.
- The measurement result of this device is not a medical diagnosis and not intended to substitute consultation and diagnosis by a qualified professional healthcare provider (e.g., physician, pharmacist, or other licensed health-care professionals).

Do not use this device for self-diagnosis or for self-treatment of a medical condition. Seek advice from a health-care professional immediately if the patient is clearly unwell and/or having physiological or medical symptoms.

- DO NOT change the patient medication and treatment based on the result of one or multiple measurements. Treatment and medication changes should be prescribed only by a medical professional.
- Use and store the device and parts in temperature and humidity conditions specified in the «Technical specifications». Usage and storage of the device and parts in conditions outside ranges given in the «Technical specifications» may result in device malfunction and/or impact the safety of usage.

Caution

Indicates a potentially hazardous situation, which if not avoided, could result in death or serious injury.

- Protect the device and accessories from the following to avoid damaging the device:
 - water, other liquids, and moisture
 - extreme temperatures
 - impacts and vibrations
 - direct sunlight
 - contamination and dust
- Stop using this device and consult with your doctor if you experience skin irritation or discomfort.
- DO NOT use this device or parts after the expiration of its stated service life.
- Do not disassemble or attempt to service the device, accessory and parts, during use or in storage. Access to the device internal hardware and software is prohibited. Unauthorized access and servicing of the device, during use or in storage, may compromise the safety and performance of the device.
- Using non-compatible accessories may compromise the safety and performance of the device.
- This device is reusable. It is recommended to clean and disinfect the device and the accessory prior to use if contamination or cross contamination are possible.
- Ensure the device has acclimated to the environment before taking measurement. It is recommended to wait approximately 4 hours after device storage at the maximum or minimum storage and transport temperature for device to cool down or warm up.

Electromagnetic compatibility

- This device is compliant with electromagnetic disturbances standard.

Further documentation in compliance with EN 60601-1-2 EMC standard is available from MicroLife on www.microlife.com/electro-magnetic-compatibility.

- Do not use this device close to strong electromagnetic fields and portable radio frequency communication devices (for example microwave oven and mobile devices). Keep a minimum distance of 0.3 m from such devices when using this device.
- DO NOT use this device in proximity of equipment that may cause electromagnetic disturbance (EMD), such as high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, and computerized tomography (CT) scanners. This device is not certified for operation near these equipments, which may cause device malfunction and measurement inaccuracies.
- This device features Bluetooth® that emits radio frequency (RF) in the 2.4GHz band. Do not use this device in locations where RF is restricted (for example, on an aircraft). Turn off the device and remove the power source if necessary when in RF restricted locations.
- This device operates in an unlicensed ISM band at 2.4GHz. In case this device is used near other wireless devices (for example wireless LAN) which operates on the same frequency band as this device, there is a possibility that interference may occur. If interference occurs, stop the operation of other devices or relocate this product away from other wireless devices before using it.

Adverse events and reporting

Please report any serious incident, injury or adverse event that has occurred in relation to the device to the manufacturer/ European authorized representative (EU REP), and to the competent authority.

3. Device information

Package contents

- 1 x Microlife PF 200 BT
- 2 x 1.5 V alkaline batteries; type AAA
- 2 x Spare mouthpiece
- 1 x Data cable USB-C to USB-A
- 1 x Instruction manual
- 1 x Medication card
- 1 x Storage bag

Device optional parts

- Mouthpiece
- 2 x 1.5 V alkaline batteries; size AAA

4. Device installation and setup

Mouthpiece installation

The mouthpiece (5) is fitted into the measurement tube (4) and can be inserted or removed with a gentle push or pull. Check if the mouthpiece is inserted all the way into the measurement tube prior to use.

- ☞ Prior to the initial use, it's recommended to clean the mouthpiece and flow tube.

Battery installation

1. Open the battery compartment (7) and remove the cover.
2. Insert 2 batteries (1.5V, size AAA) following polarity symbols in the compartment.
3. Reinstall the battery compartment cover (7).

Battery replacement

When the battery low symbol (18) appears in the display (2), the device is inoperable until the batteries have been replaced with 2 new batteries (1.5 V size AAA).

- ☞ During a battery exchange the data stored in the device memory is saved and will NOT be lost. After battery replacement, the device date / time setting needs to be re-set.

- ☞ If the device is left unused for long periods, remove the batteries.

Device date and time setup

The device needs to be setup with the correct date and time to take measurements and record results in memory. After installing new batteries in the device, the device will enter date and time setting operation automatically, with the date and time display symbols blinking.

Follow the steps to set the device date/time:

1. Open the data compartment (3) and remove the cover.
2. If the year is not blinking in the display, press the TIME button next to the clock symbol with a paper clip and the year will start blinking; release the switch
3. Using the two arrow keys (6) on the front of the device you can decrease (left arrow button) or increase (right arrow button) the number. Release the arrow button when the correct number is reached.
4. To confirm the year and then set the month, press the TIME button.
5. You can now set the month using the arrow buttons (6). (Example: pressing 2 times the right arrow button advances to 06 for June). To confirm and then set the day, press the TIME button.
6. Please follow the previous sequence to set the day, hour and minutes.
7. Once you have set the last minute and pressed the TIME button, the date and time are set, and the time is displayed.
8. Close the data compartment (3).

The device date and time setting can be changed using the settings mode function manually or updated using compatible mobile or computer software.

If you hold down the arrow key for more than 2 seconds, it speeds up.

Date/time can also be set easily from computer when you run the Microlife Analyzer Software and Connected Health+ App.

Computer software installation

The device can transfer data to your personal computer (PC) using the USB data link function with the included USB data cable and compatible Microlife software. The free downloadable Microlife Asthma Analyzer Software allows you to transfer, store, and review measurement results on your Windows OS or macOS computer.

Software download and information are available on: www.microlife.com/support/software-downloads

Mobile software installation

The device can connect to compatible Microlife mobile application software on an Apple iOS or Google Android OS smartphones via Bluetooth® to transfer data to the software via Bluetooth® datalink mode.

The free downloadable Microlife Connected Health+ App allows you to transfer, store, and review measurement results on your Apple iOS or Google Android OS smartphone.

Software download and information are available on: www.microlife.com/technologies/connect

Settings function operation

1. Turn on device and wait for device to become ready for measurement.
2. Press TIME button (19) to enter settings mode.
3. In settings mode, the following settings operations are available in the settings menu in the listed sequences:
 1. Bluetooth® automatic activation on/off selection
 2. PEF Traffic light reference selection
 3. PEF Traffic light reference value setup
 4. Device date and time setup

Bluetooth® automatic activation

The first setting in settings menu is the on and off of automatic activation of Bluetooth® function after completion of measurement, for the ease of use with data transfer from device to mobile software. In this setting, device display shows Bluetooth® symbol, «Auto», «SEt», and «On» or «OFF» blinking.

1. When this function is set to «On», the device automatically turns on Bluetooth® after measurement completion; when this function is set to «OFF», the device will return to stand-by after measurement completion. By default, this is set to «On». Press Memory Scroll buttons (6) to change the selection.
2. Press time button (19) to confirm selection and proceed to the next settings. Press O/I button (1) to keep current setting and exit settings menu.

PEF traffic light reference selection

1. The second setting in settings menu is the selection of reference maximum used for PEF traffic light. Device display will show «PEF», «rEF», «SEt», memory symbol (19), with «On» or «OFF» blinking.
2. There are 2 reference methods available for selection, press memory scroll buttons (6) to change selection:
 - Select «On»: Device automatically use the highest PEF measurement result in memory (you personal best PEF) as the reference maximum PEF for the Traffic light feature
 - Select «OFF»: Device uses customized PEF valve input by the user as the reference maximum PEF for the Traffic light feature
3. Press time button (19) to confirm selection and proceed to the next settings menu. Press O/I button (1) to keep the current setting and exit settings menu.

PEF traffic light reference value setup

If customized PEF traffic light reference is selected, the device will enable adjustment of customized maximum PEF value. Device display will show «aPEF» and «SEt» with the PEF value blinking.

Press memory scroll buttons (6) to set the PEF value to be used as the reference maximum of PEF traffic light.

Press time button to confirm selection and proceed to the next setting menu. Press O/I button (1) to keep the current setting and exit settings menu.

Device date and time setup

The last setting in settings menu is device date and time, the date and time display symbols will blink to allow adjustments, starting with year, then month, day, hour, and minute.

Press Memory Scroll buttons (6) to adjust the value. Press TIME button (19) to confirm value and proceed to the next adjustment. Press O/I button (1) to keep the current setting and exit settings menu.

 Device date and time can also be set using companion computer or mobile software.

5. Measurement preparation

Check the following items in preparation of measurement:

- Take the measurement in positions of standing or sitting upright to obtain accurate readings.
- Check if the device has been cleaned or disinfected in case the device was previously used by another patient.

 Adult supervision is required for taking measurements of children, elderly person, and patients with disabilities.

6. Measurement operation

Standalone measurement mode operation

1. Press O/I button (1) to turn the device on.
 2. When device turns on, the previous measurement result (11) will appear on the display (0 if there is no data). The device will then sound two short beeps and show blinking arrows (17) on display. The display will blink to indicate it is ready for measurement.
 3. Follow these steps to perform a measurement correctly.
 - You can perform the measurement while standing or sitting upright. For better comparison of your data, you should always perform the measurement in a similar position.
 - Hold the device with both hands on the two textured compartments.
 - Cover the mouthpiece tightly with your lips. Ensure the lips form an airtight seal around the mouthpiece (5).**NOTE:** Do not bend your neck nor obstruct the mouthpiece with your tongue.
 - Inhale completely and hold your breath for a moment.
 - Blow into the measuring tube (4) as hard and as fast as you can.
 - NOTE:** Avoid slow and sustained exhalation to obtain accurate reading.
 - Device will sound a long beep to indicate the exhalation is measured; the PEF value is displayed for 3 seconds followed by FEV1 value (13). Then, the device will sound two short beeps to indicate that it is ready for a new measurement.
 - Repeat measurement as needed. It is recommended to perform three or more measurements sequentially.
4. After each measurement, the actual measurement result is shown followed by the highest reading of your current measurement session.
 5. Press the O/I button (1) to finish the measurement session and record the highest measurement value in the device memory. The device will display the saved PEF and FEV1 values with the memory slot number (14).
 6. After saving the measurement result, the device will activate Bluetooth® datalink if automatic Bluetooth® activation is set to «On». To Turn device off without Bluetooth® data transfer, press O/I button (1) again.
 7. Clean the measuring tube and mouthpiece after use.

7. Measurement interpretation

Self-assessment using the PEF traffic light indicator

The colored bars on the left of the display is the PEF Traffic Light Indicator (9), which provides a quick visual reference comparing your PEF result versus a reference baseline value. As your PEF result raises or lowers, the arrow on the display will show the measurement as either within the normal (green), borderline (yellow) or danger (red) range.

- Green Zone – OK: If your measurement is in the green zone, your condition appears to be under control. Continue your treatment plan.
- Yellow Zone – Caution: If your measurement is in the yellow zone, measure more frequently and follow your treatment plan.
- Red Zone – Alert: If your measurement is in the red zone, your condition is serious. Act as discussed with your physician or seek emergency medical treatment.

Baseline reference of the PEF traffic light indicator

There are two different settings of baseline reference value available for the PEF traffic light indicator:

- **Memory:** Using your personal highest PEF result in the device memory as the baseline for the traffic light indicator zones.
- **Customized:** Set a customized PEF baseline (e.g., predicted reference) value manually on the device or use the compatible software for the traffic light indicator zones.

 This device is preset with a customized baseline of 500 L/min.

- The green zone ranges between baseline or predicted value (MAX) and 80% of that reading. To determine the yellow/green border zone, multiply the baseline MAX by 0.8 [Example: 500 L/min x 0.8 = 400 L/min].
- The yellow zone ranges between baseline or predicted value (MAX) and 60% of that reading. To determine the yellow/red border zone, multiply the Baseline MAX by 0.6. [Example: 500 L/min x 0.6 = 300 L/min].
- The red zone ranges below 50% of baseline or predicted value.

 Setting a predicted reference PEF value, creating a treatment plan and modifying your traffic light zone settings should only be done under direction of a physician (or other licensed health care practitioner). If you would like to use this feature, determine your customized zones with your doctor and record them on the traffic light card (enclosed).

8. Memory function operation

During a measurement session, the device will measure multiple exhalation readings and automatically record the highest reading when the session ends (for example: 1 session is 1 on-off cycle). The device memory stores 400 readings.

Review measurement readings in memory

To review the measurement readings in device memory, switch the device on and press Memory Scroll buttons (6). Press «<» of the Memory Scroll button to review the second latest readings; press «>» works in the opposite direction.

Memory capacity low

When the number of readings stored in the device memory reaches 390 sets or more, the device will display (14) blinking after switched on to indicate that the device memory capacity is low.

Memory full

When the number of readings stored in the device memory reaches 400 sets, the device will beep and display (14) «400» blinking when switched on to indicate that the device memory is full. From this point on, any new measurement readings will replace and overwrite the oldest reading in the device memory.

Delete the most recent reading in memory

To delete the most recent data in the device memory, while the device is in memory function, press Memory Scroll (both «<» and «>») buttons (6) simultaneously for 5 seconds until the display shows «CLR» blinking. Press the Memory Scroll buttons (6) again to delete the most recent reading in device memory. Wait until the «CLR» message on the display elapses after 3 seconds to exit without change.

Delete all readings in memory

To delete all data in the device memory, while the device is in memory function, press Memory Scroll (both «<» and «>») buttons (6) simultaneously for 5 seconds until the display shows «CLR» blinking. Press O/I button (1) to delete all readings in the device memory. Wait until the «CLR» message on the display elapses after 3 seconds exit without change.

9. Bluetooth® function operation

Bluetooth® pairing between the device and app is required when using the Bluetooth® datalink function for the first time. Please refer to the FAQ of the app for information about using the app and Bluetooth® pairing with the device.

Turn on/off the Bluetooth® datalink function

1. Press O/I button (1) for 4 seconds while device is off to activate the Bluetooth® datalink function; device display will show a blinking Bluetooth® (12) symbol and the device identifier.
2. The Bluetooth® will remain activated to scan for compatible Bluetooth® mobile application to connect to; if no connection is made after 120 seconds, device will show rapidly blinking Bluetooth® (12) symbol to indicate that no connection made, then it will automatically deactivate.
3. Press O/I button (1) while device Bluetooth® is activated to turn off the Bluetooth® function.

Software data transfer via Bluetooth® datalink

1. While device Bluetooth® function is activated, if a compatible mobile application is found over Bluetooth®, device will automatically connect to the software. If a successful connection between device and software is established, the Bluetooth® symbol (12) on the device display will stop blinking and a 4-part circle animation (13) will appear to indicate on-going Bluetooth® connection.
2. When connected, the device will execute data transfer from the device memory to the software via Bluetooth® automatically.
3. After completion of successful data transfer, the device will automatically disconnect from the software and turn off the Bluetooth® function.
4. Refer to the FAQ in the mobile application menu for details.

10. USB datalink function operation

Turn on/off the USB datalink function

1. Connect the device to a computer via USB data cable while the device is turned off.
2. Press O/I button (1) to turn the device on to activate USB datalink function. When USB datalink function is activated, the device will show «USB» blinking on the display and automatically search for compatible software on the computer for data transfer.
3. Disconnect USB data cable to deactivate USB datalink function.

Software data transfer via USB datalink

1. While device USB datalink function is activated, the compatible computer can communicate with the device and execute data transfer per user operation. Refer to software help / user manual for information about how to use the software. The device will show the 4-part cycle animation on the display when communicating to the computer software.
2. After completion of data transfer on the computer software, you can disconnect the device, computer and the USB data cable.

11. Device error and troubleshooting

Device error and troubleshooting

Error Message	Possible cause(s)	Solution
«Er 1»	The result cannot be recorded in memory because device date and time is not set.	Set device date and time.
«NO»	There is no data stored in device memory.	Set device date and time then take measurements.
«Hi»	The PEF value is higher than the PEF range upper limit of 900 ml/min.	Please repeat the measurement.

Error Message	Possible cause(s)	Solution
Bluetooth® symbol (12) blinking rapidly: Bluetooth® connection error.	Bluetooth® connection not established successfully or Bluetooth® connection ends abnormally	Turn off the device and the software, then re-try connection in Bluetooth® datalink function.

Common issues and troubleshooting

Other possible errors and their solutions: If problems occur when using the device, the following points should be checked.

Issues	Possible cause(s)	Solution
The device display remains blank while batteries are installed.	Batteries are installed incorrectly or low in voltage.	1. Check the polarity (+/-) of the battery installation. 2. If the display is erratic or unusual, remove the batteries and re-install new batteries.
The devices fails to provide a reading during measurement, or the measurement reading appears abnormally low or high.	Device does not detect the turbine wind wheel rotation correctly, due to incorrect measurement, incorrect measuring tube installation, or wind wheel being stuck from foreign objects in the turbine.	– Ensure the measurement is taken according to the correct procedure. – Ensure that the measuring tube is correctly installed in the device. – Ensure that the wind wheel of the turbine in the measuring tube is rotating when air is blown through the tube. Any objects, dust, liquids or mucus may interfere with the wind wheel rotation - clean the turbine and measuring tube as needed. – Discuss the abnormal readings with your doctor as needed.
The unit fails to take a measurement, year is blinking.	Device date and time setting is not set.	Ensure the device date and time setting are set.

12. Device maintenance and disposal

Cleaning

Cleaning the mouthpiece and the measuring tube

The peak flow mouthpiece and measuring tube should be cleaned within 30 minutes after each use.

1. Separate the measuring tube (4) and mouthpiece (5) from the device's main body. Disconnect the mouthpiece (5) from the measuring tube (4) by moving it in a forward direction.
2. Immerse the mouthpiece (5) in a prepared soap solution of water and regular dish-washing soap. Thoroughly clean by swirling in the prepared soap solution. Thoroughly rinse the mouthpiece (5) with water (distilled water is recommended). Allow the device to air dry.
NOTE: Do not tap against hard surfaces to avoid damage.
3. Immerse the measuring tube (4) in a prepared soap solution and thoroughly clean by swirling in the prepared soap solution. For cleaning the measuring tube, it is best to use distilled water. Be careful not to damage the turbine assembly. Thoroughly rinse the measuring tube with water then air dry.
4. Visually inspect the measuring tube (4) and mouthpiece (5) after cleaning for residual mucus or debris inside – the measuring tube and mouthpiece should be free of any debris that may interfere with the turbine assembly inside the measuring tube.
NOTE: DO NOT immerse the measuring tube in boiling water!

5. Reconnect the mouthpiece (5) to the measuring tube (4) and reassemble to the main unit as shown. When the tube is pushed completely into the fixed position, it will engage the tab on the body and click into place.
NOTE: If you have «hard water», mineral deposits may be visible on the mouthpiece or measuring tube. Rinse in distilled water, shake and place on a paper towel and allow to air dry before using.

Cleaning the device

Clean the device after use with a clean, damp cloth.

DO NOT immerse the device in water or detergent.

Storage

Store the device in spaces away from direct exposure of sunlight, heat, humidity, and vibrations. For long term storage without use, remove the batteries before storage.

Calibration and support

The device is pre-calibrated based on the measuring tube accompanying the device. The device does not require calibration by the user before use.

DO NOT interchange the measuring tubes between different devices. The device is not designed to be serviced by the user. Contact device distributor for servicing by qualified personnel authorized by the manufacturer.

DO NOT open the device for service or repair.

Disposal

This device is medical electrical equipment. Dispose this device and batteries in accordance with the Waste Electrical and Electronic Equipment (WEEE) directive and applicable local regulations. Do not dispose of the device and batteries with domestic or commercial waste.

13. Specifications and compliance

Technical specifications

 **NOTE:** Technical specifications subjected to change without notice.

Device Type:	Digital Peak Flow Meter
Model number:	PF200B
Reference number	PF 200 BT
Measurement method:	Flow turbine and optical detection
Measurement range:	PEF: 50 to 900 L/min, FEV1: 0.01 to 9.99 L
Measurement resolution:	PEF: 1 L/min, FEV1: 0.0 1L
Measurement specifications:	per ATS Standardization of Spirometry 1994 - Accuracy: ± 5 L/min or $\pm 12\%$ of the reading, whichever is greater FEV1: ± 0.1 L or $\pm 5\%$ of the reading, whichever is greater

Dimensions: 144 x 77 x 48 mm
Weight: 137 g (excluding batteries)

Ingress protection (IP) rating: IP22: Protected against solid foreign objects of 12.5 mm Ø and greater. Protected against vertically falling water drops when enclosure tilted up to 15°

Power source: 2 x 1.5 V AAA batteries

Power rating: DC 3V

Applied part type reference:

Type BF

Electrical protection type: Internally powered ME equipment

Wireless Communication: Bluetooth® Low Energy

Operation mode: Continuous operation

Operating conditions: 10 - 40 °C / 50 - 104 °F
 15 - 90 % relative maximum humidity
 700 hPa ~ 1060 hPa

Storage and transport conditions: -20 - +55 °C / -4 - +131 °F
 15 - 90 % relative maximum humidity

Service life - device: 2 years or 10000 measurements, whichever comes first

Compliance information

This device complies with the requirements of the Medical Device Regulation (EU)2017/745.

Compliant standards

EN 60601-1
 EN 60601-1-2
 EN 60601-1-11
 EN ISO 23747
 ATS standard 1994 update

14. Supplement information for users and patients

Information about measurement and monitoring of peak flow values

A peak flow meter is used to measure a person's «peak expiratory flow», which is the fastest speed a person can blow air out of the lungs after taking in as big a breath as possible. «Peak expiratory flow» is a simple measure of airflow that can tell you how well you are breathing. It tells you how well air is moving through the airways in your lungs. Forced Expiratory Volume (FEV1) is a measure of the volume of air expelled in 1-second. Peak flow meter instructions must be followed carefully to get a correct measure of airflow. If you have a breathing condition such as asthma or COPD, your physician (or other licensed health care professional) may recommend that you use a peak flow meter to watch for changes in your airflow. When the device is used to monitor lung conditions such as asthma and chronic obstructive pulmonary disease (COPD), the user should be under the care of a licensed health care professional. A licensed health care professional's advice is required to understand the meaning and importance of the measurements reported by the device and how to decide on an appropriate treatment plan that defines when to measure.

Your physician (or licensed health care professional) will give you a treatment plan that will tell you what actions to take when you have a change in airflow. In addition, you should record your peak flow measures as recommended by your physician (or other licensed health care professional). Reviewing peak flow measures can help you and your physician (or licensed health care professional) check closely on your asthma or COPD to provide the best treatment for you. The treatment plan given to you by your physician or other licensed health care professional will tell you what action to take when there are changes in your peak flow number. No matter what your peak flow measures are, if you have signs and symptoms such as chest tightness, shortness of breath, coughing or wheezing you should follow your licensed health care professional's advice for contacting him or her.

Which PEF values are normal?

Normal PEF values for males*

Age (years)	Height				
	152 cm	165 cm	180 cm	190 cm	200 cm
20	554	602	649	693	740
25	543	590	636	679	725
30	532	577	622	664	710
35	521	565	609	651	695

Age (years)	Height				
	152 cm	165 cm	180 cm	190 cm	200 cm
40	509	552	596	636	680
45	498	540	583	622	665
50	486	527	569	607	649
55	475	515	556	593	634
60	463	502	542	578	618
65	452	490	529	564	603
70	440	477	515	550	587

Normal PEF values for females*

Age (years)	Height				
	140 cm	152 cm	165 cm	180 cm	190 cm
20	390	423	460	496	529
25	385	418	454	490	523
30	380	413	448	483	516
35	375	408	442	476	509
40	370	402	436	476	502
45	365	397	430	464	495
50	360	391	424	457	488
55	355	386	418	451	482
60	350	380	412	445	475
65	345	375	406	439	468
70	340	369	400	432	461

Normal PEF values for children and adolescents**

Height	Males and females
110 cm	147
112 cm	160
114 cm	173
116 cm	187
120 cm	200
122 cm	214
124 cm	227
127 cm	240
130 cm	254
132 cm	267
134 cm	280
137 cm	293
140 cm	307
142 cm	320
144 cm	334
147 cm	347

Height	Males and females
150 cm	360
152 cm	373
155 cm	387
157 cm	400
160 cm	413
162 cm	427
165 cm	440
167 cm	454

* Leiner GC. et al: Expiratory peak flow rate. *AM Rev Respir Dis* 88:644, 1963

** Polgar G. Promadhat V: *Pulmonary Function Testing in Children: Techniques and Standards*. Philadelphia, W.B. Saunders Company, 1971

Guarantee

This device is covered by a **2 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance.
- Accessories and wearing parts:
 - Batteries
 - Mouthpiece
 - Measuring tube

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website: www.microlife.com/support
 Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

Symbols and definitions

Medical device

CE Marking of Conformity

Importer

Authorized representative
in the European Union

Authorized representative in Switzerland

Manufacturer

Country of manufacture
(Date of manufacture if date printed next to symbol)

Model number

Reference number

Serial number
(YYYY-MM-DD-SSSSS;
year-month-day-serial number)

Lot number (YYYY-MM-DD; year-month-day)

Unique Device Identifier

Caution

General warning sign

Type BF applied part

Direct current

Protected against vertically falling water drops when enclosure tilted up to 15°.

Temperature limitation for operating **and** storage

Humidity limitation for operating **and** storage

Atmospheric pressure limitation

Read instructions for use before operating the device.

Dispose in accordance with waste electrical and electronic equipment (WEEE) directive.

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