Globalcare Medical Technology Co., Ltd
7th Building, 39 Middle Industrial Main Road,
European Industrial Zone, Xiaolan Town,
528415 Zhongshan City, Guangdong Province, PRC
Phone: +86 760 22589901

Donawa Lifescience
Piazza Albania 10, 00153 Rome, Italy
www.donawa.com/contacts

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
- ① Cover of medication container
- ② Medication container
 - a: medication tank
 - b: mesh unit
 - c: connector
- ③ Housing
- ④ ON/OFF button
- ⑤ Operation LED indicator
- ⑥ Status LED indicator
- ⑦ Button for unlocking the medication container
- ⑧ USB-C socket
- ⑨ Battery cover
- ⑩ Child face mask
- ⑪ Adult face mask
- ⑫ Mouthpiece
- ⑬ USB-C cable
- ⑭ Batteries

This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations with regard to the correct disposal of electrical and electronic products. Correct disposal helps to protect the environment and human health.

Read the instructions carefully before using this device.

Type BF applied part

Class II equipment

Serial number

Reference number

Manufacturer

Date of manufacture

ON/OFF button

IP22

Protected against solid foreign objects of 12.5 mm Ø and greater. Protected against vertically falling water drops when enclosure tilted up to 15°.

Authorized representative in the European Union

Authorized representative in Switzerland

Medical device

Importer

Caution

Single patient multiple use (for applied parts only)

Humidity limitation

Temperature limitation

Ambient pressure limitation

Unique Device Identifier

Type number

CE 0123

CE Marking of Conformity

Read these instructions for use carefully. Observe the warnings and safety notes. Keep these instructions for use for future reference.

Make the instructions for use accessible to other users. If the device is passed on, provide the instructions for use to the next user as well.

Intended use:

The device is an aerosol therapy system suitable for domestic use. The device is designed to drive the stainless mesh nebuliser module to produce medical aerosol for lung respiratory disorders.

Intended user:

The user is an adult person without specific knowledge or professional ability. The patient is the intended operator except in the case of children and patients that require special assistance.

Patient population:

The nebuliser can be used on persons over 2 years of age under supervision; it can be used for self-treatment by persons over 12 years of age.

Indications:

Lung acute or chronic diseases of respiratory organ disorders, or inflammation of the upper respiratory system.

Contra-indications:

The device is not indicated to be used with quick-relief medications during life-threatening asthma attacks. No contraindications exist to the administration of aerosols by inhalation. Contraindications related to the medicine used must be checked on the medicine package leaflet. Consult your physician in case of doubts.

Dear Customer,

This nebuliser is an aerosol therapy system suitable for domestic use. This device is used for the nebulisation of liquids and liquid medication (aerosols) and for the treatment of the upper and lower respiratory tract.

If you have any questions, problems or want to order spare parts please contact your local Microlife-Customer Service. Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit the internet at www.microlife.com where you will find a wealth of invaluable information on our products.

Stay healthy – Microlife Corporation!

Table of contents

- 1. Important safety instructions**
 - 2. Preparation before use**
 - Preparing the nebuliser
 - Filling the medication container
 - 3. Using the device for the first time**
 - Inhaling correctly
 - Ending nebulisation
 - 4. LED display color**
 - 5. Cleaning and disinfecting**
 - 6. Malfunctions and actions to take**
 - 7. Maintenance, care and service**
 - 8. Disposal**
 - Battery disposal
 - General disposal
 - Product disposal
 - 9. Guarantee**
 - 10. Technical specifications**
-
- 1. Important safety instructions**

WARNING

- Use the device only as described in this manual and therefore as an aerosol therapy system, following the indications of your doctor. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
- The device is not intended for health care in professional facilities.
- If children or people with physical or mental disabilities use the device, proper supervision by an adult is required.
- Children must be kept under adult supervision to ensure that they do not play with the device.
- Do not use the device on persons who are ventilated and/or unconscious.
- Before use, ensure that there is no visible damage to the device or components. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

- Do not use the device, if you spot damage or if you notice something unusual, if the device does not work properly or if it is making an abnormal sound. If you feel unwell or experience pain, stop using it immediately.
- In case of failure and/or malfunction, read the «Malfunctions and actions to take» section in the user manual. Do not handle or open the device.
- Comply with the safety regulations concerning the electrical devices and in particular:
 - use only original applied parts and components;
 - never submerge the unit in water;
 - do not use the device while taking bath or shower;
 - never wet the device, it is not protected against water penetration;
 - never touch the unit with wet or moist hands;
 - do not leave the unit exposed to the weather elements;
 - the use of this device by children and people with disabilities requires always the close supervision by an adult with full mental faculties;
 - do not pull the power cord or the device itself to unplug it from the power socket;
- Before plugging in the device, make sure that the electrical rating, shown on the rating plate on the bottom of the unit, corresponds to the mains rating.
- Do not leave the unit plugged in when not in use; unplug the device from the wall socket when it is not in operation or fully charged.
- The installation must be carried out according to the instructions of the manufacturer. An improper installation can cause damage to persons, animals or objects, for which the manufacturer cannot be held responsible.
- Before performing any maintenance or cleaning operation, turn off the device and disconnect the plug from the main supply or remove the batteries.
- Some parts of the unit are so small that they may be swallowed by children; keep the equipment out of children's reach.
- If you decide not to use the device any longer, it is recommended to dispose of it according to the current regulations.
- Make sure to:
 - use this device only with medicines prescribed by your doctor;
 - only administer the treatment using the applied part recommended by your doctor depending on the pathology;
- Check the medicine package leaflet for possible contraindications for use with common aerosol therapy systems.
- The device is not a substitute for medical consultation and treatment. Always consult your doctor first in case of pain or illness.
- To avoid strangulation and entanglement, keep cable and air tubes out of reach of young children.
- Do not position the equipment so that it is difficult to operate the disconnection of the device.
- The nebuliser and applied parts are single patient reuse. The device is multi-patient reuse.
- Changes to the device and the components are not permitted.
- Use only applied parts recommended by the manufacturer, use of different components may impair the effectiveness of treatment and damage the device.
- If the USB-C cable of this device is damaged, it must be disposed of. Please contact Customer Services or the retailer.
- Under no circumstances should you open or repair the device yourself. There is a risk of electric shock when opening the device. The device is not disconnected from the power supply unless the adapter is unplugged (and the USB-C cable has no other power connection).
- Never use the device near flammable gases, explosive gas, oxygen or nitrogen oxide.
- Do not operate the device near electromagnetic transmitters.
- Protect the device from impacts.
- Should you require an adapter or extension lead, it must meet the applicable safety requirements. The power limit and the maximum output specified on the adapter must not be exceeded.

2. Preparation before use

Preparing the nebuliser

- If the mesh atomiser has not been assembled, slide it onto the housing (A). See section «Cleaning and disinfecting» for assembly instructions of mesh unit and connector to the medication container.
- **Battery mode:** assemble the batteries into device (B1 – B2) paying attention to polarity
- **USB mode:** connect USB cable ⑬ to device port ⑧ and to USB power supply (not included in delivery. The external power

supply must comply with the European standard EN 60601-1 with protection class 2 and must have the rating specified in section «Technical specifications»)

- For hygiene reasons, it is essential to clean and regularly disinfect the medication container and the atomiser set before using them for the first time and after each treatment. For this purpose, see section «Cleaning and disinfecting».
- If the therapy involves inhaling several different medications one after the other, please be aware that the medication container must be rinsed under warm tap water after each use. See section «Cleaning and disinfecting».

Filling the medication container

- Check that mesh unit is correctly inserted (C2) otherwise liquids may leak from the device. See section «Cleaning and disinfecting» for description of medication container assembly.
- Open the cover of the medication container (D) and fill the medication container (E) as needed. The maximum recommended filling quantity is 10 ml. Avoid overfilling, otherwise liquids may leak from the device.
- Use medication only as instructed by your doctor and ask about the appropriate inhalation period and quantity for you.
- If the prescribed quantity of medication is less than 0.5 ml, it must be topped up with isotonic saline solution. Dilution is also necessary with viscous medications (viscosity < 3). Observe the instructions of your doctor.
- Close the cover of the medication container (F) until the markings on cover and medication container are aligned.
- Securely connect the part of the atomiser set that you would like to use (mouthpiece, adult mask or children's mask) to the connector (2)-c (G). If you are using the adult mask or children's mask, you may attach one of the straps included. Ensure that you securely fix the strap on the mask with a knot.
- Now move the device to your mouth and put your lips firmly around the mouthpiece. When using a mask, place it over nose and mouth, without covering the eyes.
- Ensure that you have added fluid before switching on the device. Start the device by pressing the On/Off button (4).
- If the spray mist is flowing out of the device and the operation LED (5) is lit up the device is operating refer also to «LED display color» section.
- In USB mode, it is possible to change from high-rate therapy to low-rate therapy by pressing the On/Off button (4) for 3 s.

3. Using the device for the first time

Inhaling correctly

- It is important to use the correct breathing technique to ensure that the particles are distributed as widely as possible throughout the airways. To ensure that the particles reach your airways and lungs, you must breathe in slowly and deeply, then hold your breath briefly (5 – 10 seconds) and then breathe out quickly.
- Keep the device as vertical as possible. However, a small angle does not affect usage as the container is leak-proof. Make sure that the medication is in contact with the mesh for proper functioning of the device during use. If the device is tipped gently by 45 degrees with a maximum volume of 10 ml of liquid, the device will work as before.

Ending nebulisation

- If you have finished or want to stop the treatment, you can switch off the device manually using the On/Off button (4).
- The device has an automatic switch-off function, which is activated:
 - After a nebulisation period of 20 minutes.
 - When there is only a very small amount of medication or liquid left or when the substance to be nebulised is no longer in contact with the mesh. The operation LED (5) flashes blue rapidly and the device then switches off.
- Do not use the device if the medication container is empty. The device detects when there is no medication in the medication container and will switch-off automatically.
- For technical reasons, a small amount of medication will remain in the medication container. Do not reuse it.

4. LED display color

LED display color	Meaning
Operation LED Green	• Static ON: device in USB mode and high-rate therapy running • Slow blinking: device in USB mode and low-rate therapy running
Operation LED Blue	• Static ON: device in battery mode and high-rate therapy running • Fast blinking for 3 seconds: absence of liquid inside the liquid container

Operation LED Orange	• Fast blinking: device is running in auto-cleaning mode • Static ON: low charge battery
Operation LED Green-Blue-Orange	• Blinking for 3 seconds: medication container is missing
Status LED Orange	• Static ON: reminder signal that the mesh unit must be changed • Fast blinking: informative signal that auto-cleaning is needed

5. Cleaning and disinfecting

Device(s): Medication container and atomiser set

	Before the first use and after each treatment Follow carefully the cleaning and disinfection instructions of the medication container and atomiser set as they are very important for the device performances, the therapy success and to avoid health risks. • Do not clean the mesh mechanically using a brush or similar device, this could cause irreparable damage. • Please consult your doctor about the additional requirements for hygiene preparation (hand care, handling of medication/inhalation solutions) for high-risk groups (e.g. patients with cystic fibrosis). • The water mentioned later in the cleaning instructions should be as soft as possible.
---	--

Limitations on processing	The maximum times of cleaning and disinfection of medication tank (2)-a and connector (2)-c, mouthpiece (12) and mask(s) (10 – 11) are 366 times. To maintain proper device functioning: • we recommend to replace the mesh unit (2)-b after 1000 operating minutes. The device provides a reminder signal for the mesh unit replacement. To replace the mesh unit, follow the disassembling and reassembling instructions provided in the «Preparation before cleaning» section below. • Before reassembling the medication container, the timer can be reset by pressing the ON/OFF button for 5 s while the device is switched off. • auto-cleaning is recommended periodically. The device provides a reminder for auto-cleaning after 300 operating minutes. The timer is reset when the auto-cleaning function is activated.
----------------------------------	---

Instructions

Preparation before cleaning	Disassembling the mesh atomiser • Press the release button on the rear of the device (7) and vertically pull off the medication container (H). • Remove the mouthpiece or the mask from the medication container: – Open the cover of the medication container (D) – Separate the connector (2)-c, mesh unit (2)-b and medication tank (2)-a (C1) • The reassembly is carried out in reverse order. Note: Insert the mesh unit with electrical contacts towards the medication tank (C1). Note: To avoid microbial growth and residual drug drying, clean and disinfect the medication container and atomiser set immediately after each use.
------------------------------------	---

Cleaning	After each use, the disassembled medication container and the parts of the atomiser set that were used, including mouthpiece or mask, must be cleaned by immersion in a clean container with soft warm water at a max. temperature of 40 °C for approx. 5 minutes. Scrub the entire surface of all components (except the mesh) at least 8 times with a clean, small brush. • Afterwards rinse all parts (except the mesh) thoroughly in running tap water at least 30 seconds. Do not clean the mesh mechanically using a brush or similar device. • Do not wash the mesh under running water. • Do not remove silicone O-ring from cover ①.
Auto-cleaning function for medication container	After manual cleaning, proceed with self-cleaning for the medication container: • Open the cover of the medication container (D) and pour out the remaining fluid. • Pour a sufficient amount (3 – 6 ml) of tap water (soft water) into the medication container (E). Then close the cover (F). Shake the container slightly so that the remaining fluid mixes properly with the cleaning fluid. Let the device run at high nebulisation rate for one minute and then switch it off. • Open the cover of the medication container (D) and pour out all the liquid. • Pour 1.5 – 2.0 ml of tap water or distilled water mixed with white vinegar (1/4 vinegar and 3/4 distilled water) into the connector ②-c of the medication container (I). To avoid leakage through the lateral openings, please do not overfill. Press and hold the On/Off button for 3 seconds when the device switched off, the operation LED ⑤ then flashes blinking red and switched to auto-cleaning mode. The cleaning liquid is now drawn out from the connector ②-c into the medication tank ②-a in reverse mode. During this process, it is possible for droplets of this liquid to leak out of the medication container. Therefore, place some protection under the nebuliser if necessary. Auto-cleaning will automatically run for 3 minutes.
Disinfection	Place the disassembled medication container, the cover, mouthpiece and masks in a clean container. Cover everything with ethyl alcohol (70 – 75%) and leave it for 10minutes. Finally, clean all parts again with sterile water (soft water) and allow to air dry in a clean place.

Disinfection (alternative for cover and mouthpiece only) 	Alternatively, medication tank ②-a, connector ②-c, cover ① and the mouthpiece ⑫ can be placed in boiling water for 5 minutes for disinfection. Note: The mesh ②-b and the mask (⑩ – ⑪) must not be placed in boiling water.
Drying 	- Gently shake the parts a few times so that any water is removed from the tiny holes in the mesh. - Place the individual parts on a clean, dry and absorbent surface and leave them to dry completely (for at least 4 hours). Note: The medication container ② must not be placed in the microwave.
Inspection	Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discolored parts.
Packaging	Pack dry parts in a clean and sealed container when not in use. Do NOT pack wet or damp parts.
Storage 	Storage conditions refer to «Technical specifications». Note: re-clean and re-disinfect the parts if they are stored more than one day.
Transportation	After cleaning and disinfection, always transport the parts within a clean and sealed container. Contamination can be avoided by good hand washing and not touching the inside sections of the parts when taking out and re-assembling the parts for use.

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

6. Malfunctions and actions to take

Use following table for solutions if any problem occur.

Description	Solutions
The nebuliser produces no or too little aerosol.	- Too little medication in the medication container. - The nebuliser is not being held vertically. - Medication liquid has been poured in that is not suitable for nebulisation (e.g. too viscous. The viscosity must not be greater than 3. Ask your pharmacist about the viscosity.) The medication liquid should have been prescribed by a doctor. - The battery is flat, change the batteries or use with USB cable.
The output is too low.	- Particles on the mesh are hindering output. In this case, you should pause the inhalation process and remove the mask or mouthpiece. Then clean and disinfect the medication container (see «Cleaning and disinfecting» section). - The mesh is worn. Replace the mesh part.
What medications are suitable for inhaling?	Naturally, only your doctor can decide which medication should be used for treating your condition. With the NEB 950, you can nebulise medications with a viscosity of less than 3. Medications containing oil (particularly essential oils) must not be used.
There is inhalation solution residue in the medication container	This is normal and due to technical reasons. Stop the inhalation as soon as the atomiser makes a significantly different noise or the device switches off automatically due to a lack of inhalant.

7. Maintenance, care and service

We recommend replacing the nebulising head, face masks and mouthpiece every 6 months (depending on the usage). Order all spare parts from your dealer or pharmacist, or contact Microlife-Service (see foreword).

8. Disposal

Battery disposal

Empty, completely flat batteries must be disposed of in specially designated collection boxes, at recycling points or at electronics retailers. You are legally required to dispose of the batteries.

General disposal

For environmental reasons, do not dispose of the device in household waste at the end of its service life. Dispose of the device at a suitable local collection or recycling point in your country. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Product disposal

Item	Component	Disposal Description
1	Device	The component mainly includes plastic and electronic components. All comply with RoHS and REACH, and all could be safely disposed. This product is subject to European Directive 2012/19/EU on waste electrical and electronic equipment and is marked accordingly. Never dispose of electronic devices with household waste. Please seek out information about the local regulations regarding the correct disposal of electrical and electronic products.
2	Applied parts: medication container (2), masks (11) / (19), mouthpiece (12), replaceable mesh (2)-b	The components consist mainly of plastic and metal. All comply with ROHS and REACH. Before disposal, all the applied parts should be cleaned according to manual instructions and then disinfected by boiling for 5 minutes.

3	USB-C cable (13)	The USB-C cable mainly consists of plastic and metal, fully complies with RoHS and REACH guidelines and can be disposed of safely.
4	2 x 1.5 V alkaline batteries; type AA	The batteries comply with European Directive 2023/1542/EU. Empty, fully discharged batteries must be disposed of via designated collection boxes, recycling points, or electronics retailers, and retailers should be legally obligated to dispose of the batteries.
5	Pouch	The pouch is made of fiber, complies with RoHS and REACH guidelines and can be disposed of safely.

9. Guarantee

This device is covered by a **2 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge. Opening or altering the device invalidates the guarantee. The following items are excluded from the guarantee:

- Transport costs and risks of transport.
 - Damage caused by incorrect application or non-compliance with the instructions for use.
 - Damage caused by accident or misuse.
 - Packaging/storage material and instructions for use.
 - Regular checks and maintenance.
 - Accessories and wearing parts: Batteries, masks, mouthpiece.
- Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website: www.microlife.com/support
- Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

10. Technical specifications

Model:	NEB 950
Type:	GCEC02
Nebulisation rate	Low nebulisation rate: ~ 0.2 ml/min
in USB mode:	High nebulisation rate: ~ 0.4 ml/min
AEROSOL PERFORMANCES ACCORDING TO EN ISO 27427 based on adult ventilatory pattern with Sodium Fluoride (NaF):	
Aerosol output:	0.28 ml
Aerosol output rate:	0.05 ml/min.
Residual Volume:	0.63 ml
GSD (geometric standard deviation):	11.81 ± 0.05
RF (respirable fraction < 5 µm):	63.0 ± 5.86 %
Large particle range (> 5 µm):	37.0 ± 5.86 %
Middle particle range (2-5 µm):	48.3 ± 1.87 %
Small particle range (< 2 µm):	14.7 ± 4.07 %
Particle size:	3.89 ± 0.35 µm (MMAD)
Vibration frequency	108 kHz ± 10%
Sound noise level (average level):	≈ 42 dBA
Nebuliser capacity:	min. 0.5 ml; max. 10 ml
Power input:	2 x 1.5 V alkaline batteries; type AA Input: 5V DC(USB), >500 mA
Operating limits:	20 min. ON/40 min. OFF

Operating conditions:	10 ~ 40 °C / 50 ~ 104 °F 15 – 90 % relative maximum humidity 700 – 1060 hPa Atmospheric pressure
Storage and transport conditions:	-20 ~ +60 °C / -4 ~ 140 °F < 90 % relative maximum humidity 700 – 1060 hPa Atmospheric pressure
Weight:	140 g (with batteries)
Dimensions:	45 x 45 x 116 mm
Automatic switch-off:	Automatic shut down, when the liquid has been used up.
IP Class:	IP22
Expected service life:	Device: 500 hours Mesh unit: 1350 min

Technical alterations reserved.

Class II device as regards protection against electric shocks.

Nebuliser, mouthpiece and masks are type BF applied parts.

Important for the safety of the device and its user: use only power supplies as specified (5V > 500 mA output with USB-C connector) and ensure their certification in accordance with EN 60601-1(with protection class 2) and EN 60601-1-2.

The technical specifications may change without prior notice.

Please report any serious incident that has occurred in relation to the device, injury or adverse event to the local competent authority and to the manufacturer or to the european authorised representative (EU REP).

Vigilance contact point:

<https://ec.europa.eu/growth/sectors/medical-devices/contacts>
AEROSOL CHARACTERISTICS IN ACCORDANCE WITH
REGULATION EN ISO 27427

Standards applied:

Electric safety standards CEI EN 60601-1 electromagnetic compatibility according to CEI EN 60601-1-2.

The device is a Class IIa medical device. **Device complies with European Regulation on medical devices EU MDR 2017/745. Important information regarding Electromagnetic Compatibility (EMC):**

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be restricted if electromagnetic disturbance occurs. This could result in issues such as the failure of the device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of an atomiser set other than that specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Keep portable RF communication devices (including peripheral equipment, such as antenna cables or external antennas) at least 30 cm away from all device parts, including all cables included in delivery.
- Failure to comply with the above could impair the performance of the device.
- Further documentation in compliance with this EMC standard is available from Microlife on www.microlife.com/electromagnetic-compatibility.

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