microlife®

WS300


Diagnostic Scale

EN → 1	RU → 26	RO → 51	TR → 76
FR → 6	UA → 31	DE → 56	GR → 81
IT → 11	PL → 36	NL → 61	AR → 86
ES → 16	HU → 41	SV → 66	FA → 90
PT → 21	CZ → 46	FI → 71	

Preparation

Microlife WS 300



Easy Short Instructions – Simple Scale

1. In order to get accurate measurements, place the scale on a firm level floor.
2. Step on the scale to begin measurement.

Easy Short Instructions – Body Fat Scale

1. In order to get accurate measurements, place the scale on a firm level floor.
2. Power On – Press «**SET**».
Display shows last used User Number.
3. User Number – Press «**SET**», then use ▲/▼ to select **P0** – **P9**. Press «**SET**» to confirm.
4. Gender – Use ▲/▼ to choose female/male. Press «**SET**» to confirm.
5. Height – Adjust 100 – 220 cm with ▲/▼. Press «**SET**» to confirm.
6. Age – Adjust 10 – 80 years with ▲/▼. Press «**SET**» to confirm.
7. Measurement – Step on the scale.
Results display automatically.

Name of Purchaser / Nom de l'acheteur / Nome del rivenditore /
Nombre del comprador / Nome do comprador / Ф.И.О. покупателя
/ П.І.Б. покупця / Imię i nazwisko nabywcy / Vásárló neve /
Jméno kupujícího / Numele cumpărătorului / Name des Käufers
/ Naam koper / Inköparens namn / Ostajan nimi/ Alicinin adi /
Όνοματεπώνυμο αγοραστή / اسم المشتري / نام خریدار

Serial Number / Numéro de série / Numero di serie / Número de
serie / Número de série / Серийный номер / Серійний номер
/ Numer seryjny / Sorozatszám / Výrobní číslo / Număr de serie
/ Serienummer / Serienummer / Sarjanumero / Seri numarası /
Αριθμός σειράς / رقم التسلسل / مدل

Date of Purchase / Date d'achat / Data d'acquisto / Fecha de
compra / Data da compra / Дата покупки / Дата покупки / Data
zakupu / Vásárlás dátuma / Datum năkupu / Data cumpărării /
Kaufdatum / Datum van aankoop / Inköpsdatum / Ostapäivämäärä
/ Satin alma tarihi / Ημερομηνία αγοράς / تاريخ الشراء
/ شماره سریال

Specialist Dealer / Revendeur / Categoria rivenditore / Distribuidor
especializado / Revendedor autorizado / Специализированный
дилер / Специализований дилер / Przedstawiciel / Forgalmazó /
Specializovaný dealer / Distributor de specialitate / Fachhändler
/ Specialist Dealer / Återförsäljare / Alan kauppias / Uzman
satıcı / Εξειδικευμένος αντιπρόσωπος / التاجر المختص
/ تاریخ خرید

- ① SET Button
- ② Up («▲») Button
- ③ Down («▼») Button
- ④ Battery compartment

Display

- ⑤ Low battery indication
- ⑥ Zero Setting
- ⑦ Instant Weight Reading
- ⑧ Error (Err)
- ⑨ User Indication
- ⑩ Gender
- ⑪ Height
- ⑫ Age
- ⑬ User-Linked Weight Reading
- ⑭ Calculation of Body Composition
- ⑮ Body Fat Reading
- ⑯ Body Water Reading
- ⑰ Bone Mass Reading
- ⑱ Muscle Mass Reading
- ⑲ Basal Metabolic Rate (BMR) Reading
- ⑳ Body Mass Index (BMI) Reading
- ㉑ Error (Err2)

Read the instructions carefully before using this device.

Dear Customer,

This diagnostic scale allows you to evaluate your body composition. The knowledge of your body composition allows you to combine a correct nutrition and physical activity as essential conditions to improve the quality of life, promote health, and prevent pathologies.

Detailed user information about our products and about our services can be found at www.microlife.com.

Stay healthy – Microlife Corporation!

Table of contents

1. Important Facts about Body Composition

Body Fat
Body Water
Muscle Mass
Bone Mass
Body Mass Index (BMI)
Basal Metabolic Rate (BMR)

2. Suggestions for Use

3. Using the device for the first time

How to insert the battery
KG / LB / ST Switch
Use as simple scale

4. Entering Personal Data

5. Calculation of Body Composition

6. Important safety instructions

Device care
Disposal

7. Guarantee

8. Technical specifications

1. Important Facts about Body Composition

It is important to know the body composition of an individual in order to be able to estimate his/her health condition. It is not enough to know one's weight but it is necessary to define its quality, which means dividing the human body in its two major component masses:

- Lean body mass: bones, muscles, viscera, water, organs, blood
- Fat body mass: adipose cells

Under the same conditions of weight, these two masses have different volumes.

Body Fat

The total body fat mass is calculated by totalling the essential fat (or primary fat) and the spare fat (or storage fat). The essential fat is necessary for the metabolism and some spare fat, which is theoretically surplus, is necessary for those who carry out physical activities.

For men a body fat percentage between 17 – 26% is considered as normal and for women between 22 – 31%. The body fat % increases with the person's age.

Body Water

Body water % is an important indication of the overall «well being» of a human being as water is the single most important component of body weight. Research shows that 57 % body water is optimum, but normally there is no problem if the reading is 10 % higher or lower.

The body water % is greater in childhood and decreases with age. Women, usually having more adipose tissue than men, have a water ratio of about 55 – 58% while men have about 60 – 62%. Normally below 40% would be seen as too low, and above 70% as too high.

Body weight variations that take place in a few hours are completely normal and are due to water retention, since the extra-cellular water is the only component subject to such rapid variations.

Muscle Mass

Muscle mass is important in determining a healthy body composition. A person with a higher % of muscle mass finds it easier to move, but needs more energy to do it. Exercise is very important in maintaining a healthy body and the muscle mass % is a useful indicator to control it.

The normal muscle mass % on the body weight lies between 38 – 54% for men and between 28 – 39% for women depending on age and physical activity level.

Bone Mass

Bone mass increases rapidly in childhood and reaches its maximum between 30 – 40 years. It decreases slightly afterwards with increasing age.

Healthy nutrition, regular exercise including some level of muscle building can help to reduce bone degeneration. However, it is otherwise difficult to significantly influence bone mass.

Factors like age, gender, weight, and height have a small influence on the bone mass. Average adult bone mass is 15 % for male and 12 % for female. There are no recognized guidelines or recommendations relating to bone mass.

Body Mass Index (BMI)

The BMI is a simple index of weight-for-height that is commonly used to classify underweight, overweight and obesity in adults. It is calculated by dividing the body weight in kg by the square of the height in m of a person.

The WHO BMI Classification is as follows:

- Normal: 18.50 – 24.99 kg/m²
- Underweight: under 18.50 kg/m²
- Overweight: over 24.99 kg/m²
- Obese: over 29.99 kg/m²

Basal Metabolic Rate (BMR)

The BMR in kcal is the amount of energy required by the body when in a state of complete rest to maintain its basic functions. About 70 % of a human's total energy expenditure is due to the basal life processes within the organs of the body. About 20 % of one's energy expenditure comes from physical activity and another 10 % from the digestion of food.

An estimation can be acquired through an equation using age, sex, height, and weight. There are no recognized guidelines or recommendations relating to BMR.

2. Suggestions for Use

- Place the scale on a flat, hard surface. Soft, inaccurate flooring (e.g. rugs, carpets, linoleum) can cause the scale to give uneven reading of your weight.
- Step onto the scale with your feet parallel and weight equally distributed.
- Stand still while the scale measures your weight.
- It is important to weigh yourself with dry, bare feet, correctly positioned on the electrodes. Moreover, for the scale to work properly, avoid your legs touching.
- The scale is safe to use as a simple scale by everyone. The following people should not perform a diagnostic measurement:
 - pregnant women
 - people on cardiovascular medication
 - those fitted with a pace-maker or other forms of electronic, medical equipment
 - those with a temperature
 - children under the age of 10
 - those in dialysis, with oedemas on their legs, people with dimorphism or suffering from osteoporosis
- Always weigh yourself on the same scale each day at the same time, preferably undressed and before breakfast. To get the best results from your scale, weigh yourself twice, and if the two weights are different from each other, your weight is between the two readings. Wait for about 15 minutes after getting up, so that the water can distribute itself around your body.

- The surface of the scale can be slippery, if wet.
- If an error occurs during the analysis of your body composition, the scale will turn off itself automatically.
- The precision of the measurement results can be altered in the following cases:
 - a high alcohol / drug / caffeine level
 - after an intense sports activity
 - during an illness or pregnancy
 - after heavy meals or in presence of dehydration
 - your personal data (height, age, sex) have not been input correctly
 - if you have wet or dirty feet

NOTE: The body composition measured by this scale is to be considered only as approximate information. For more information, consult your physician or dietician.

NOTE: Only people between 10 – 80 years and those who are between 100 – 220 cm tall may use the body composition measurement function of this scale.

3. Using the device for the first time

How to insert the battery

This scale operates with two AAA 1.5 V batteries. When the battery is dead the display shows «LO» ⑤ or the battery symbol. Insert two new AAA batteries into the compartment ④, ensuring the poles match the (+) and (-) symbols as indicated.

KG / LB / ST Switch

1. Turning on the power
Press the SET button to turn on the unit. The display will show the personal settings of the last used user number.
2. Selecting the unit of measurement
Press and hold the SET button. The kg symbol on the LCD will begin to flash. Use the «▲» or «▼» buttons to select lb or st, then press the SET button to confirm.

Use as simple scale

1. In order to get accurate measurements, place the scale on a firm level floor.
2. Step onto the scale and wait without moving; after a few seconds, the display screen will show your weight ⑦.
Do not lean on anything while being weighed as this could alter the reading.

3. Step off the scale, the display screen will continue to show your weight for a few seconds then it will turn off automatically.
4. If the «Err» ⑧ symbol is displayed, it means the weight limit has been exceeded. Please step off the scale immediately to prevent damage to the product. The maximum load capacity is 3 – 180 kg.

4. Entering Personal Data

This scale can memorise the data of 10 different persons.

1. Turning on the power
Press the «SET» button ① to turn on the unit. The display will show the personal settings of the last used user number.
2. Selecting a user number
Press the «SET» button again. The user number will begin to flash. Use the «▲» or «▼» buttons to select a user number from P0 to P9. When the desired user number is displayed, press the «SET» button to confirm.
3. Setting the gender
Use the «▲» or «▼» buttons to choose between female and male, then press the «SET» button to confirm.
4. Setting the height
Use the «▲» or «▼» buttons to adjust your height within the range of 100 cm to 220 cm. When the desired height is displayed, press the SET button to confirm.
5. Setting the age
Use the «▲» or «▼» buttons to adjust your age within the range of 10 to 80 years. When the desired age is displayed, press the «SET» button to confirm.
6. Starting the measurement
Step on the scale to start the measurement. The unit will automatically detect and display your results.

5. Calculation of Body Composition

1. Turning on the power
Press the SET button to turn on the unit. The display will show the personal settings of the last used user number.
2. Selecting a user
Use the «▲» or «▼» buttons to select a user number from P0 to P9. When the desired user number is displayed, press the «SET» button to confirm.

3. Step on the scale making sure you place your bare feet over the two electrodes. The display will show the following data in sequence:
- your weight ⑬ flashing 3 times
 - four half-8 digits «0000» ⑭ appear on the display
 - your body fat percentage ⑮
 - your body water percentage ⑯
 - your bone mass in kg ⑰
 - your muscle mass percentage ⑱
 - your BMR in kcal ⑲
 - your BMI in kg/m² ⑳

NOTE: The alternating display of your measurement results will repeat twice, after which the unit will shut off automatically.

NOTE: Four half-8 digits «0000» will appear on the display to indicate the start of the diagnostic measurement. You can step off the scale, after they have disappeared and your body fat percentage has been displayed.

NOTE: If the «Err 2» ㉑ symbol is displayed, it means the body-fat and body-water percentages cannot be calculated. Please restart the entire procedure from the beginning. If the body-fat percentage is less than 5 % or greater than 60 %, the reading cannot be accurate.

Compare the results obtained with the parameters listed in the section «1. Important Facts about Body Composition»: If your results fall outside the normal values, we suggest you consult your physician.

6. Important safety instructions

- This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical specifications» section.
- Protect it from:
 - water, moisture, and extreme temperatures
 - impact and dropping
 - contamination and dust
 - direct sunlight, heat, and cold
- Do not use this device if you think it is damaged or notice anything unusual.
- Never open this device.

- If the device will not be used for a prolonged period, remove the batteries and ensure nothing is left on the scale.
- Do not mix old and new batteries.
- Do not mix different types of batteries (e.g., alkaline, standard carbon-zinc, or rechargeable nickel-cadmium).
- Read the additional safety information provided within the individual sections of this instruction manual.

Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

Device care

- Clean the device only with a soft, dry cloth.
- Never immerse the device in water or rinse it under running water.

Disposal

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

7. Guarantee

This device is covered by a **2 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge. Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Batteries.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website: www.microlife.com/support

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original

invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

8. Technical specifications

- Equipped with high-precision strain-gauge sensors
- Capacity: 180 kg / 396 lb / 28.8 st
- Resolution: 0.1 kg / 0.2 lb / 1/4 st
- Battery: 3 V (1.5 V AAA x 2)
- Display: LCD with backlight
- Body composition analysis:
 - Body fat in %
 - Body water in %
 - Bone mass in kg
 - Muscle mass in %
 - Basal metabolic rate in kcal (BMR)
 - Body mass index in kg/m^2 (BMI)
- Functions:
 - Automatic zero resetting
 - Automatic switch-off
 - Overload indication
 - Low battery indication
- Technical alterations reserved.

- ① Bouton SET
- ② Bouton («▲») haut
- ③ Bouton («▼») bas
- ④ Compartiment à piles

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