microlife[®]

NC200

Non Contact Thermometer

EN → 1	LV → 45
SV → 10	LT → 54
FI → 18	EE → 63
DA → 27	RU → 72
NO → 36	IS → 82


 Microlife Corporation
 9F, No. 431, RuiGuang Road, NeiHu
 Taipei, 114, Taiwan, China
www.microlife.com




 Microlife UAB
 P. Lukšio g. 32
 08222 Vilnius, Lithuania
eurep@microlife.lt

CE 1639

IB NC 200 N-V10 2026-08-17

Microlife NC 200

Name of Purchaser / Inköparens namn / Ostajan nimi / Forhandlers navn / Kjøpers navn / Pircēja vārds un uzvārds / Pirkėjo pavardė / Ostja nimi / Ф.И.О. покупателя / Nafn kaupanda

Serial Number / Sarjanumero / Serienummer / Serienummer / Sērijas numurs / Serijos numeris / Seerianumber / Серийный номер / Lotunúmer

Date of Purchase / Inköpsdatum / Ostpäivämäärä / Købsdato / Kjøpsdato / Iegādes datums / Pardavimo data / Ostukuupäev / Дата покупки / Kaupdagur

Specialist Dealer / Återförsäljare / Alan kauppias / Special-forhandler / Spesialist forhandler / Specializētais pārstāvis / Pardavusi įstaiga / Ametlik müügiesindaja / Специализированный дилер / Sõluaðili

- ① Measuring sensor
- ② Tracking light
- ③ Self-indicator light
- ④ Display
- ⑤ M-button (memory)
- ⑥ MODE button
- ⑦ START/IO button
- ⑧ Battery compartment cover
- ⑨ All segments displayed
- ⑩ Ready for measuring
- ⑪ Measurement complete
- ⑫ Body mode
- ⑬ Object mode
- ⑭ Low battery indicator
- ⑮ Changing between Celsius and Fahrenheit
- ⑯ Recall mode
- ⑰ Recall the last 30 readings
- ⑱ Measured temperature too high
- ⑲ Measured temperature too low
- ⑳ Ambient temperature too high
- ㉑ Ambient temperature too low
- ㉒ Error function display
- ㉓ Blank display
- ㉔ Flat battery
- ㉕ Date/Time
- ㉖ Beeper function setting
- ㉗ Replacing the battery
- ㉘ Measuring countdown
- ㉙ Protective cap

This Microlife thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, this device can provide a stable, heat-interference-free reading with each measurement. The device performs a self-test every time it is switched on to always guarantee the specified accuracy of any measurement.

This thermometer has been clinically tested and proven to be safe and accurate when used in accordance to the operating instruction manual.

Please read through these instructions carefully in order for you to understand all functions and safety information.

Microlife® is a registered trademark of Microlife Corporation.

Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.

Table of contents

1. Important information

- Device description
- Indications
- Contraindications
- Warning
- Caution
- Electromagnetic compatibility
- Adverse events and reporting

2. Device information

- Package contents

3. Measurement preparation

4. The advantages of this thermometer

5. Control displays and symbols

6. Setting date, time and beeper functions

7. Changing between body and object mode

8. Directions for use

9. Changing between Celsius and Fahrenheit

10. How to recall 30 readings in memory mode

11. Battery installation and replacement

12. Error messages

13. Device maintenance and disposal

Cleaning and disinfecting
Storage
Calibration and support
Disposal

14. Guarantee

15. Technical specifications

16. Compliance information

17. Symbols and definitions

18. www.microlife.com

1. Important information

Device description

A digital infrared thermometer is a medical device that utilizes the principles of human body thermal radiation emission and digital signal processing to compute and provide a body temperature measurement.

Intended purposes:

The device is intended to measure body temperature.

Intended user:

The device is intended to be operated by adults and adolescents with adequate vision, motor functions, and education, capable of understanding the instructions for use and operating general household electrical appliances.

This device is also intended to be operated by healthcare professionals and personnel.

Intended patient:

The intended patients are the general population of all ages.

Intended use environment and conditions:

The device is intended for use in a home healthcare environment (e.g. general household without medically trained personnel) by patients (e.g. for self-measurement) or by a care giver, and in a professional healthcare facility (e.g. physician office).

Indications

This device measures body temperature for indications of:

- Screening for illness and conditions related to high or low body temperature.

Contraindications

- The device is not suitable for use on sites with any breached or compromised skin.

Warning

NOTE: Warning items indicate potentially hazardous situations, if not avoided, may result in death, critical or serious injury to the user or patient.

- This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.
- The measurement result of this device is not a medical diagnosis and not intended to substitute consultation and diagnosis by a qualified professional healthcare provider (e.g., physician, pharmacist, or other licensed health-care professionals). Do not use this device for self-diagnosis or for self-treatment of a medical condition. Seek advice from a health-care professional immediately if the patient is clearly unwell and/or having physiological or medical symptoms.
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed.
- Do not attempt to measure temperature of the body sites that are not specified in this instruction.
- Professional users must follow the country specific applicable medical devices ordinance.

Caution

NOTE: Caution items indicate potentially hazardous situations, if not avoided, may result in minor or negligible injury to the user or patient, or damage to the property or environment.

- Cleaning and/or disinfecting the thermometer is recommended before and after use if the thermometer is dirty or contaminated from use or after storage. Refer to «Cleaning and disinfecting» section for details.
- Do not use this device if you think it is damaged or notice anything unusual.
- Do not attempt to repair, modify, or open this device.
- Vasoconstriction during the initial stages of fever (resulting from infection or other factors) or reduced peripheral perfusion (such as in cases of shock) may result in cooler skin and inaccurate surface temperature measurements. It is recommended to verify core body temperature using alternative methods and consult healthcare professionals for patient assessment if patient's conditions persist.

- If the patient temperature is influenced (e.g. incubator, heating blanket) the measurement result provided by this device shall not be used to determine the presence/absence of high body temperature.

Electromagnetic compatibility

This device is compliant with EN 60601-1-2 Electromagnetic Disturbances standard.

Further documentation in compliance with EN 60601-1-2 EMC standard is available from Microlife on www.microlife.com/electro-magnetic-compatibility.

Do not use this device close to strong electromagnetic fields and portable radio frequency communication devices (for example microwave oven and mobile devices). Keep a minimum distance of 0.3 m from such devices when using this device.

This device is not certified to be used in the vicinity of medical equipment including high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) and computerized tomography (CT) instruments.

Adverse events and reporting

Please report any serious incident that has occurred in relation to the device, injury or adverse event to the local competent authority and to the manufacturer or to the European authorised representative (EU REP).

2. Device information

Package contents

- 1 x Microlife non contact digital thermometer
- 1 x Instruction manual
- 2 x 1.5 V alkaline batteries; type AAA (pre-installed)

3. Measurement preparation

- The patients and thermometer should stay in the same room condition for at least 30 minutes before using.
- Patients should not drink, eat or exercise before/while taking the measurement.
- When taking multiple measurements for comparison, allow a 30 seconds interval between each measurement for measurement reliability.
- Doctors recommend rectal measurement for newborn infants within the first 6 months, as all other measuring methods might lead to ambiguous results.

Caution: In the following situations it is recommended that three temperatures are taken with the highest one taken as the reading:

1. Children under three years of age with a compromised immune system and for whom the presence or absence of high body temperature is critical.
2. When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the device and obtains consistent readings.
3. If the measurement is surprisingly low.

4. The advantages of this thermometer

Measures in a matter of seconds

The innovative infrared technology allows the measurement without even touching the object. This guarantees safe and sanitary measurements within seconds.

Auto measurement with distance control

The device can take a measurement automatically when the device detects the distance is appropriate within 5 cm.

Multiple uses (wide range of measurement)

This thermometer offers a wide range of measurement from 0.1 – 99.9 °C / 32.2 – 211.8 °F, meaning the unit can be used to measure body temperature or it also has a feature allowing it to be used to measure surface temperature of the following examples:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

Accurate and reliable

The unique probe assembly construction incorporates an advanced infrared sensor, ensuring that each measurement is accurate and reliable.

Gentle and easy to use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can even be used on a sleeping child without causing any interruption.
- This thermometer is quick, therefore child-friendly.

Multiple readings recall

Users will be able to recall the last 30 readings with a record of both time and date when entering the recall mode, enabling efficient tracking of temperature variations.

Safe and hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.

High body temperature reminder

10 short beeps and a red LCD backlight reminds the patient that they may have a temperature equal to or higher than 37.5 °C.

Guidance system for self-measurement

A green light on the back shows the user that the device is at the right distance and a measurement will be performed.

5. Control displays and symbols

- **All segments displayed** ⑨: Press the START/IO button ⑦ to turn on the unit; all segments will be shown for 1 second.
- **Ready for measuring** ⑩: When the unit is ready for measuring, the «°C» or «°F» icon will keep flashing while the mode icon (body or object) will be displayed.
- **Measuring countdown** ②⑧: A 3 second countdown will be shown on the display (3, 2, 1), before each measurement.
- **Measurement complete** ⑪: The reading will be shown on the display ④ with the «°C» or «°F» icon and the mode icon steady. The unit is ready for the next measurement as soon as the «°C» or «°F» icon is flashing again.
- **Low battery indicator** ⑭: When the unit is turned on, the  icon will keep flashing to remind the user to replace the batteries.

6. Setting date, time and beeper functions

Setting the date and time

1. After the new batteries are fitted, the year number flashes in the display ②⑤. You can set the year by pressing the M-button ⑤. To confirm and then set the month, press the MODE button ⑥.
2. Press the M-button ⑤ to set the month. Press the MODE button ⑥ to confirm and then set the day.
3. Follow the previously mentioned instructions to set the day, hours and minutes.

4. Once you have set the minutes and pressed the START/IO button ⑦, the date and time are set and the time is displayed.

If no button is pressed for 20 seconds, the device automatically switches to ready for measuring ⑩.

Cancel time setup: Press the START/IO button ⑦ during time setup. The LCD will show Date/Time icons with «--:--». After that press the START/IO button ⑦ to start the measurement. If no further action is taken within 60 seconds, the device will automatically turn off.

Change current date and time: Press and hold the MODE button ⑥ for approx. 8 seconds until the year number starts to flash ②⑤. Now you can enter the new values as described above.

Setting the beeper

1. Press and hold the MODE button ⑥ for 5 seconds to set the beeper ②⑥.
 2. Press the M-button ⑤ to either turn the beeper on or off. The beeper is activated when the beeper icon ②⑥ is shown without a cross.
-  When the beeper setting has been chosen, press the START/IO button ⑦ to enter the «ready for measuring» mode ⑩; otherwise the device automatically switches to ready for measuring after 5 seconds.

7. Changing between body and object mode

1. Press the START/IO button ⑦. The display ④ is activated to show all segments for 1 second.
2. The default mode is body mode ⑫. Press the MODE button ⑥ to switch to object mode ⑬. For switching back to body mode, press the MODE button again.

8. Directions for use

Always remove the protective cap ②⑨ before use.

Measuring in body mode with auto measurement and distance control

1. Press the START/IO button ⑦. The display ④ is activated to show all segments for 1 second.
2. A flashing «°C»/«°F» icon, the blinking blue tracking light ② and a beep indicate that the device is ready for measurement ⑩.

3. Remove any hair, sweat or dirt from the forehead before measuring to ensure the accuracy of the readings.
4. **Aim the thermometer at the center of the forehead with a distance of no more than 5 cm.**
5. **The device will start the measurement automatically, when the measuring sensor ① detects the distance is appropriate within 5 cm.** The display will show a count-down (3, 2, 1); after 3 seconds a long beep will verify the completion of measurement.
6. **Self-measurement:** The green self-indicator light ③ will help to distinguish the start of the measurement. Keep the thermometer pointed towards the forehead until the light switches off.
7. Read the recorded temperature from the LCD display.
8. For the next measurement remove the thermometer from the forehead and wait until the «°C»/«°F» icon is flashing. Follow steps 4 - 5 above.
9. Press and hold the START/IO button ⑦ for 3 seconds to turn off the device; otherwise the device will automatically switch off after approx. 60 seconds.

Measuring in object mode without auto measurement

1. Press the START/IO button ⑦. The display ④ is activated to show all segments for 1 second.
2. Press the MODE button ⑥ to switch to object mode.
3. A flashing «°C»/«°F» icon, the blinking blue tracking light ② and a beep indicate that the device is ready for measurement ⑩.
4. Aim the thermometer at the center of the object you want to measure with a distance of no more than 5 cm. **Press the START/IO button ⑦.** After 3 seconds a long beep will verify the completion of measurement.
5. Read the recorded temperature from the LCD display.
6. For the next measurement wait until the «°C»/«°F» icon is flashing and follow steps 5 - 7 above.
7. Press and hold the START/IO button ⑦ for 3 seconds to turn off the device; otherwise the device will automatically switch off after approx. 60 seconds.

NOTE:

- **Readings from different measuring sites should not be compared as the normal body temperature varies by measuring site and time of day,** being highest in the evening and lowest about one hour before waking up.

- Normal body temperature ranges:
 - Axillar: 34.7 - 37.3 °C / 94.5 - 99.1 °F
 - Oral: 35.5 - 37.5 °C / 95.9 - 99.5 °F
 - Rectal: 36.6 - 38.0 °C / 97.9 - 100.4 °F
 - Microlife NC 200: 35.4 - 37.4 °C / 95.7 - 99.3 °F

9. Changing between Celsius and Fahrenheit

This thermometer can display temperature readings in either Fahrenheit or Celsius. To switch the display between °C and °F, **press and hold** the MODE button ⑥ for 3 seconds; the beeper icon is shown on the display. Press the MODE button again; the current measurement scale («°C» or «°F» icon) will be shown on the display ⑩. Change the measurement scale between °C and °F by pressing the M-button ⑤. When the measurement scale has been chosen, press the START/IO button ⑦ to enter the «ready for measuring» mode; otherwise the device automatically switches to ready for measuring after 10 seconds ⑩.

10. How to recall 30 readings in memory mode

This thermometer can recall the last 30 readings with a record of both time and date.

- **Recall mode ⑩:** Press the M-button ⑤ to enter recall mode when the power is off. The memory icon «M» will flash.
- **Reading 1 - the last reading ⑪:** Press and release the M-button ⑤ to recall the last reading. Number «1» and a flashing «M» are displayed.

Pressing and releasing the M-button ⑤ after the last 30 readings have been recalled will resume the above sequence from reading 1.

11. Battery installation and replacement

This device is supplied with 2 new, long-life 1.5V, size AAA (pre-installed) batteries. Batteries need replacing while the  24 is shown on the display.

Remove the battery cover 27 by sliding it in the direction shown.

Replace the batteries – ensure correct polarity as shown by the symbols in the compartment.

 **Caution:** Do not use expired batteries or mix new and used batteries together.

Caution: Remove batteries if the device is not going to be used for a prolonged period.

12. Error messages

Error	Description	Solutions
Measured temperature too high ⑱ Measured temperature too low ⑲	«H» = too high, measured temperature is higher than 43 °C / 109.4 °F in body mode or 99.9 °C / 211.8 °F in object mode. «L» = too low, measured temperature is lower than 34 °C / 93.2 °F in body mode or 0.1 °C / 32.2 °F in object mode.	Make sure the measuring sensor is clean. Then, take a measurement again.
Ambient temperature too high ⑳ Ambient temperature too low ㉑	«AH» = too high, ambient temperature is higher than 40.0 °C / 104.0 °F. «AL» = too low, ambient temperature is lower than 15.0 °C / 59.0 °F in body mode or lower than 5.0 °C / 41.0 °F in object mode.	Allow the thermometer to rest in a room for at least 30 minutes at ambient temperature (15.0 °C / 59.0 °F in body mode or 5.0 °C / 41.0 °F – 40.0 °C / 104.0 °F in object mode). Then, take a measurement again.
Blank display ㉒	The battery / battery's polarity («+» and «-») has not been loaded correctly.	Check if the battery / battery's polarity («+» and «-») has been loaded correctly.
Flat battery indication ㉓	The battery is low.	The battery should be replaced immediately.

Error	Description	Solutions
Error function display ㉔	The device has a malfunction. «Er 0» / «Er 6»: The system has a malfunction. «Er 2»: Device is directly placed on the forehead / object. Keep the measuring distance of 1 - 5 cm. Do not touch the bottom side (sensing area) of the measuring sensor.	Remove and reload the battery and reboot the device. If the message still appears, please contact the Microlife local retail or distributor.

13. Device maintenance and disposal

Cleaning and disinfecting

Cleaning and/or disinfecting the thermometer is recommended before and after use if the thermometer is dirty or contaminated from use or after storage. Check whether the device is still functional after cleaning and disinfecting. Other cleaning agents or cleaning methods may cause malfunction or damage to the device and will void the warranty.

To clean and disinfect the thermometer:

1. Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring sensor. Ensure that no liquid enters the interior of the device.
2. Let the disinfectant dry off for 10 minutes before next use.

Caution: Do not to scratch the surface of the sensor lens and the display.

Caution: Do not use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the device in water or other cleaning liquids. Especially the display part do not immerse or wipe.

Storage

When not in use: Keep the device and accessories in a dry, cool place away from sunlight, with ambient conditions within the temperature and humidity ranges described in the «Technical specifications». Remove the batteries from the device if the device will not be used for an extended period.

- Protect the device and accessories from the following to avoid damaging the device:
 - water, other liquids, and moisture
 - extreme temperatures
 - impacts and vibrations
 - direct sunlight
 - contamination and dust

Caution: Storing the disused device for an extended period without removing batteries increases the chance of battery fluid leakage, which may lead to device damage and skin irritation when in contact. If your eye or skin is exposed to battery fluid, wash the exposed part immediately with ample clean water. Consult a doctor if irritation or discomfort persists.

Calibration and support

The device is calibrated during manufacturing. In general, it is recommended to have the device verified by your local designated Microlife device distributor every two years, or after mechanical impact, liquid ingress, and/or device malfunctions. For questions related to device measurement accuracy, please contact your local designated Microlife device distributor.

Caution: Do not attempt to service or calibrate the device and accessories yourself.

Disposal

This device is medical electrical equipment. Dispose this device and batteries in accordance with the Waste Electrical and Electronic Equipment (WEEE) directive and applicable local regulations. Do not dispose of the device and batteries with domestic or commercial waste.

14. Guarantee

This device is covered by a **2 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Battery.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website: www.microlife.com/support

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

15. Technical specifications

Device name:	Non contact digital thermometer
Model number:	FR1DG1
Catalogue number:	NC 200
Operating mode:	Adjusted mode
Measuring site:	Forehead
Reference body site:	Sublingual
Measurement range:	Body mode: 34.0 – 43.0 °C / 93.2 – 109.4 °F Object mode: 0.1 – 99.9 °C / Z32.2 – 211.8 °F
Resolution:	0.1 °C / °F

Measurement accuracy (Laboratory):	Body mode: $\pm 0.2^{\circ}\text{C}$, $35.0 \sim 42.0^{\circ}\text{C}$ / $\pm 0.4^{\circ}\text{F}$, $95.0 \sim 107.6^{\circ}\text{F}$ $\pm 0.3^{\circ}\text{C}$, $34.0 \sim 34.9^{\circ}\text{C}$ and $42.1 \sim 43.0^{\circ}\text{C}$ / $\pm 0.5^{\circ}\text{F}$, $93.2 \sim 94.8^{\circ}\text{F}$ and $107.8 \sim 109.4^{\circ}\text{F}$ Object mode: $\pm 1.0^{\circ}\text{C}$, $0.1 \sim 99.9^{\circ}\text{C}$ / $\pm 2^{\circ}\text{F}$, $32.2 \sim 211.8^{\circ}\text{F}$
Clinical results:	Repeatability: 0.24°C Bias: 0.07°C Limits of agreement: 0.90°C
Display:	Liquid Crystal Display, 4 digits plus special icons
Acoustic:	The unit is turned ON and ready for the measurement: 1 short beep. Complete the measurement: 1 long beep. System error or malfunction: 3 short beeps. High body temperature: 10 short beeps.
Memory:	30 readings recall in the memory mode with a record of both time and date.
Backlight:	The display light will be GREEN for 1 second, when the unit is turned ON. The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5°C / 99.5°F . The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5°C / 99.5°F .
Operating conditions:	Body mode: $15 \sim 40.0^{\circ}\text{C}$ / $59 \sim 104.0^{\circ}\text{F}$ Object mode: $5 \sim 40.0^{\circ}\text{C}$ / $41 \sim 104.0^{\circ}\text{F}$ $15 \sim 90\%$ relative maximum humidity $700 \sim 1060$ hPa Atmospheric pressure

Storage and transport conditions:	$-20 \sim +55^{\circ}\text{C}$ / $-4 \sim +131^{\circ}\text{F}$ $15 \sim 90\%$ relative maximum humidity $700 \sim 1060$ hPa Atmospheric pressure
Automatic Switch-off:	Approx. 1 minute after last measurement has been taken.
Battery:	2 x 1.5 V alkaline batteries; size AAA (pre-installed)
Battery lifetime:	approx. 2000 measurements (using new batteries)
Dimensions:	$156.7 \times 43 \times 47$ mm
Weight:	91.5 g (with battery), 68.5 g (w/o battery)
Ingress protection (IP) rating:	IP22: Protected against solid foreign objects of 12.5 mm $\varnothing$ and greater. Protected against vertically falling water drops when enclosure tilted up to 15°
Expected service life:	5 years or 12000 measurements

NOTE: The technical specifications may change without prior notice.

16. Compliance information

This device complies with the requirements of the Medical Device Regulation (EU)2017/745.

Compliant standards:

EN 60601-1, EN 60601-1-2, EN 60601-1-11, EN ISO 80601-2-56

17. Symbols and definitions

Read the instructions carefully before using this device.

Type BF applied part

Caution

Medical device

Reference number

18. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com.

Model number

Serial number (YYYY-MM-DDXXXXX;
YYYY-MM-DD = date of manufacture,
XXXXX = sequence number)

Lot number (YYYY-MM-DD;
YYYY-MM-DD = Date of manufacture)

Manufacturer

Country of manufacture
(Date of manufacture if date printed next to
symbol)

CE Marking of Conformity

Authorized representative in the European Union

Authorized representative in Switzerland

Temperature limitation for operating **or** storage

Humidity limitation for operating **and** storage

Dispose in accordance with waste electrical and
electronic equipment (WEEE) directive.

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