

Guarantee Card


Hypertension

Human Fever

Asthma

Flexible Heating

NC 120

Europe / Middle-East / Africa
Microlife AG
Espenstrasse 139
9443 Widnau / Switzerland
Tel. +41 / 71 727 70 30
Fax +41 / 71 727 70 39
Email admin@microlife.ch
www.microlife.com

Asia
Microlife Corporation.
9F, 431, RuiGang Road, NeiHu
Taipei, 11492, Taiwan, R.O.C.
Tel. 886 2 8797-1288
Fax 886 2 8797-1283
Email service@microlife.com.tw
www.microlife.com

North / Central / South America
Microlife USA, Inc.
1617 Gulf to Bay Blvd., 2nd Floor Ste A
Clearwater, FL 33755 / USA
Tel. +1 727 442 5353
Fax +1 727 442 5377
Email msa@microlifeusa.com
www.microlife.com


CE0044

IB NC 120 V3 0812

Microlife NC 120

- EN** → 1
- RU** → 7
- UA** → 15


microlife

microlife

Microlife NC 120

microlife

Microlife Non Contact Thermometer NC 120

EN

- ① Measuring Sensor
- ② Protection Cap
- ③ START Button
- ④ Display
- ⑤ ON/OFF Button
- ⑥ Mode Switch
- ⑦ Battery Compartment Cover
- ⑧ All Segments displayed
- ⑨ Memory
- ⑩ Ready for Measurement
- ⑪ Measurement complete
- ⑫ Body Mode
- ⑬ Object Mode
- ⑭ Low Battery Indication
- ⑮ Changing between Celsius and Fahrenheit
- ⑯ Recall Mode
- ⑰ Recall the last 30 Readings
- ⑱ Measured Temperature too high

- ⑲ Measured Temperature too low
- ⑳ Ambient Temperature too high
- ㉑ Ambient Temperature too low
- ㉒ Error Function Display
- ㉓ Blank Display
- ㉔ Flat Battery
- ㉕ Replacing the Battery

Read the instructions carefully before using this device.

Type BF applied part

This Microlife Thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, this thermometer can provide a stable, heat-interference-free reading with each measurement. The instrument performs a self-test every time it is switched on to always guarantee the specified accuracy of measurements. This Microlife Thermometer is intended for the periodic measurement and monitoring of human body temperature. **This thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.** Please read through these instructions carefully so that you understand all functions and safety information.

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Table of Contents

1. The Advantages of this Thermometer
2. Important Safety Instructions
3. How this Thermometer measures Temperature
4. Control Displays and Symbols
5. Changing between Body and Object Mode
6. Directions for Use
7. Changing between Celsius and Fahrenheit
8. How to recall 30 readings in Memory Mode
9. Error Messages
10. Cleaning and Disinfecting
11. Battery Replacement
12. Guarantee
13. Technical Specifications
14. www.microlife.com
Guarantee Card (see Back Cover)

1. The Advantages of this Thermometer

Measurement in a matter of seconds

The innovative infrared technology allows measurement without even touching the object. This guarantees safe and sanitary measurements within seconds.

Multiple Uses (Wide Range Measurement)

This thermometer offers a wide measurement range feature from 0 °C to 100.0 °C (32.0 °F to 212.0 °F), meaning the unit can be used to measure body temperature or it can be used to measure surface temperature of the following:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

Accurate and reliable

The unique probe assembly construction, incorporating an advanced infrared sensor, ensures that each measurement is accurate and reliable.

Gentle and Easy to Use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can even be used on a sleeping child, thereby causing no disruption.
- This thermometer is quick, therefore child-friendly.

Auto-Display Memory

The last reading is automatically displayed for 2 seconds when the unit is switched ON.

Multiple Reading Recall

Users will be able to recall the last 30 readings when entering the recall mode, enabling efficient tracking of temperature variations.

Safe and Hygienic

- No direct skin contact.
- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Cleaning of the probe can be done with an alcohol-moistened cotton tissue, making this thermometer completely hygienic for use by the whole family.

Fever Alarm

10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature equal to or higher than 37.5 °C.

2. Important Safety Instructions

- This instrument may only be used for the purposes described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.
- **Never immerse this instrument in water or other liquids. For cleaning please follow the instructions in the «Cleaning and Disinfecting» section.**

- Do not use the instrument if you think it is damaged or notice anything unusual.
- Never open the instrument.
- A basic physiological effect called vasoconstriction can occur in the early stages of fever, resulting in a cool skin effect. The recorded temperature using this thermometer can, therefore, be unusually low.
- If the measurement result is not consistent with the patient's finding or unusually low, repeat the measurement every 15 minutes or double check the result by another core body temperature measurement.
- This instrument comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section!
- Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.
- Protect it from:
 - extreme temperatures
 - impact and dropping
 - contamination and dust
 - direct sunlight
 - heat and cold
- If the instrument is not going to be used for a prolonged period the battery should be removed.

WARNING: Use of this instrument is not intended as a substitute for consultation with your physician. This instrument is NOT waterproof, so never immerse in liquids.

3. How this Thermometer measures Temperature

This thermometer measures infrared energy radiated from the forehead as well as objects. This energy is collected through the lens and converted to a temperature value. Temperature readings obtained by scanning above the eyebrow area will provide the greatest accuracy.

4. Control Displays and Symbols

- **All segments displayed** ⑧: Press the ON/OFF button ⑤ to turn on the unit; all segments will be shown for 2 seconds.
- **Memory** ⑨: The last reading will automatically be shown on the display for 2 seconds.
- **Ready for measurement** ⑩: When the unit is ready for measurement, the «°C» or «°F» icon will keep flashing while the mode icon (body or object) will be displayed.
- **Measurement complete** ⑪: The reading will be shown on the display ④ with the «°C» or «°F» icon and the

mode icon steady. The unit is ready for the next measurement as soon as the «°C» or «°F» icon is flashing again.

- **Low battery indication** ⑭: When the unit is turned on, the battery icon will keep flashing to remind the user to replace the battery.

5. Changing between Body and Object Mode

For changing from body to object mode, slide the mode switch ⑥ at the side of the thermometer downwards. For switching back to body mode, slide the switch up again.

6. Directions for Use

Measuring in Body Mode

1. Press the ON/OFF button ⑤. The display ④ is activated to show all segments for 2 seconds.
2. The last measurement reading will automatically be shown on the display for 2 seconds with the «M» icon ⑨.
3. When the «°C» or «°F» icon is flashing, a beep sound is heard and the thermometer is ready for measurement ⑩.
4. **Aim the thermometer at the center of the forehead with a distance of no more than 5 cm.** If the forehead is covered with hair, sweat or dirt, please remove the obstacle beforehand to improve the reading accuracy.

5. **Press the START button ③ and steadily move the thermometer** from the middle of the forehead to the temple area (about 1 cm above the eyebrow). After 3 seconds a long beep will verify the completion of measurement. If the temple area hasn't been reached before the long beep is heard, repeat the measurement as described above, but move the thermometer a little faster.
6. Read the recorded temperature from the LCD display.

Measuring in Object Mode

1. Follow steps 1-3 above, then aim the thermometer at the center of the object you want to measure with a distance of no more than 5 cm. Press the START button ③. After 3 seconds a long beep will verify the completion of measurement.
2. Read the recorded temperature from the LCD display.

NOTE:

- **Patients and thermometer should stay in steady-state room condition for at least 30 minutes.**
- Don't take a measurement while or immediately after nursing a baby.
- Don't use the thermometer in high humidity environments.
- Patients should not drink, eat, or exercise before/while taking the measurement.

- Don't move the measurement device from the measuring area before hearing the termination beep.
- Use an alcohol swab to carefully clean the probe and wait for 15 minutes before taking a measurement on another patient.
- 10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature equal to or higher than 37.5 °C.
- Always take the temperature in the same location, since temperature readings may vary according to locations.
- Doctors recommend rectal measurement for newborn infants within the first 6 months, as all other measuring methods might lead to ambiguous results. If using a non contact thermometer on those infants, we recommend to always verify the readings with a rectal measurement.
- In the following situations it is recommended that three temperatures are taken and the highest one taken as the reading:
 1. Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 2. When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.

3. If the measurement is surprisingly low.

- **Readings from different measuring sites should not be compared as the normal body temperature varies by measuring site and time of day**, being highest in the evening and lowest about one hour before waking up. Normal body temperature ranges:
 - Axillar: 34.7 - 37.7 °C / 94.5 - 99.1 °F
 - Oral: 35.5 - 37.5 °C / 95.9 - 99.5 °F
 - Rectal: 36.6 - 38.0 °C / 97.9 - 100.4 °F
 - Microlife NC 120: 35.4 °C - 37.4 °C / 95.7 - 99.3 °F

7. Changing between Celsius and Fahrenheit

This thermometer can display temperature measurements in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn OFF the unit, **press and hold** the START button ③ for 5 seconds; after 5 seconds, the current measurement scale («°C» or «°F» icon) will be flashing on the display ⑮. Change the measurement scale between °C and °F by pressing the START button ③. When the measurement scale has been chosen, wait for 5 seconds and the unit will automatically enter the «ready for measuring» mode.

8. How to recall 30 readings in Memory Mode

This thermometer can recall the last 30 readings.

- **Recall mode ⑮:** Press the START button ③ to enter Recall mode when power is off. The memory icon «M» will flash.
- **Reading 1 - the last reading ⑰:** Press and release the START button ③ to recall the last reading. Display 1 alone with memory icon.
- **Reading 30 - readings in succession:** Press and release the START button ③ consecutively to recall the last 30 readings in succession.

Pressing and releasing the START button ③ after the last 30 readings have been recalled will resume the above sequence from reading 1.

9. Error Messages

- **Measured temperature too high ⑱:** Displays «H» when measured temperature is higher than 42.2 °C / 108.0 °F in body mode or 100 °C / 212 °F in object mode.
- **Measured temperature too low ⑲:** Displays «L» when measured temperature is lower than 34.0 °C / 93.2 °F in body mode or 0 °C / 32 °F in object mode.

- **Ambient temperature too high ⑳:** Displays «H» in conjunction with the  when ambient temperature is higher than 40.0 °C / 104.0 °F.
- **Ambient temperature too low ㉑:** Display «L» in conjunction with the  when ambient temperature is lower than 16.0 °C / 60.8 °F in body mode or lower than 5.0 °C / 41.0 °F in object mode.
- **Error function display ㉒:** The system has a malfunction.
- **Blank display ㉓:** Check if the battery has been loaded correctly. Also check polarity (<+> and <->) of the battery.
- **Flat battery indication ㉔:** If the steady battery icon is the only symbol shown on the display, the battery should be replaced immediately.

10. Cleaning and Disinfecting

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the probe lens and the display.

11. Battery Replacement

This instrument is supplied with one lithium battery, type CR2032. Replace with a new CR2032 battery when the flashing battery symbol appears on the display ㉔. Remove the battery cover by sliding it in the direction shown. Remove the battery and replace with a new one ㉕.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

12. Guarantee

This instrument is covered by a **2 year guarantee** from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) confirming date of purchase or the receipt.

- The guarantee covers the instrument. Batteries and packaging are not included.
- Opening or altering the instrument invalidates the guarantee.
- The guarantee does not cover damage caused by improper handling, discharged batteries, accidents or non-compliance with the operating instructions. Please contact Microlife-service.

13. Technical Specifications

Type:	Non Contact Thermometer NC 120
Measurement range:	Body Mode: 34.0-42.2 °C (93.2-108.0 °F) Object Mode: 0-100.0 °C (32-212.0 °F)
Resolution:	0.1 °C / °F
Measurement accuracy:	Laboratory: ±0.2 °C, 36.0 ~ 39.0 °C (±0.4 °F, 96.8 ~ 102.2 °F)
Display:	Liquid Crystal Display, 4 digits plus special icons
Acoustic:	The unit is turned ON and ready for the measurement: 1 short beep Complete the measurement: 1 long beep (1 sec.) if the reading is less than 37.5 °C (99.5 °F), 10 short «beep» sounds, if the reading is equal to or greater than 37.5 °C (99.5 °F). System error or malfunction: 3 short «bi» sounds
Memory:	Auto-Display the last measured temperature 30 readings recall in the Memory Mode

Backlight:	The display light will be GREEN for 4 seconds, when the unit is turned ON. The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5 °C (99.5 °F). The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5 °C (99.5 °F).
Operating temperature:	Body Mode: 16-40.0 °C (60.8-104.0 °F) Object Mode: 5-40.0 °C (41-104.0 °F)
Storage temperature:	-20 °C to +50 °C (-4 °F to 122 °F) 15-95 % relative maximum humidity
Automatic Switch-off:	Approx. 1 minute after last measurement has been taken.
Battery:	CR2032 Battery (X1) 3V - at least 1000 measurements
Dimensions:	120 x 28 x 21 mm
Weight:	39 g (with battery), 36 g (w/o battery)

Reference to standards: ASTM E1965;
IEC 60601-1; IEC 60601-1-2 (EMC)

This device complies with the requirements of the Medical Device Directive 93/42/EEC.
Technical alterations reserved.
According to the Medical Product User Act a biennial technical inspection is recommended for professional users.
Please observe the applicable disposal regulations.

14. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com.

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