
## Guarantee Card


## NC 120


Hypertension

Human Fever

Asthma

Flexible Heating

### Europe / Middle-East / Africa

Microlife AG  
Esenstrasse 139  
9443 Widnau / Switzerland  
Tel. +41 / 71 727 70 30  
Fax +41 / 71 727 70 39  
Email admin@microlife.ch  
www.microlife.com

### Asia

Microlife Corporation.  
9F, 431, RuiGang Road, NeiHu  
Taipei, 11492, Taiwan, R.O.C.  
Tel. +886 2 8797-1288  
Fax +886 2 8797-1283  
Email service@microlife.com.tw  
www.microlife.com

### North / Central / South America

Microlife USA, Inc.  
1617 Gulf to Bay Blvd., 2nd Floor Ste A  
Clearwater, FL 33755 / USA  
Tel. +1 727 442 5353  
Fax +1 727 442 5377  
Email msa@microlifeusa.com  
www.microlife.com

CE0044

IB NC 120 V15-5 1916

## Microlife NC 120

- EN → 1
- FR → 7
- ES → 15
- PT → 23
- DE → 31
- TR → 39
- GR → 45
- RU → 53
- PL → 61
- HU → 67
- BG → 75
- RO → 83
- CZ → 91
- SK → 97
- SL → 103

microlife®


## Microlife NC 120

microlife®

# Microlife NC 120


microlife

# Microlife Non Contact Thermometer NC 120

EN

- ① Measuring sensor
- ② Protection cap
- ③ START button
- ④ Display
- ⑤ ON/OFF button
- ⑥ Mode switch
- ⑦ Battery compartment cover
- ⑧ All segments displayed
- ⑨ Memory
- ⑩ Ready for measuring
- ⑪ Measurement complete
- ⑫ Body mode
- ⑬ Object mode
- ⑭ Low battery indicator
- ⑮ Changing between Celsius and Fahrenheit
- ⑯ Recall mode
- ⑰ Recall the last 30 readings
- ⑱ Measured temperature too high
- ⑲ Measured temperature too low
- ⑳ Ambient temperature too high
- ㉑ Ambient temperature too low
- ㉒ Error function display
- ㉓ Blank display
- ㉔ Flat battery
- ㉕ Replacing the battery

This Microlife thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, this device can provide a stable, heat-interference-free reading with each measurement. The device performs a self-test every time it is switched on to always guarantee the specified accuracy of any measurement.

This Microlife thermometer is intended for the periodic measurement and monitoring of human body temperature. **This thermometer has been clinically tested and proven to be safe and accurate when used in accordance to the operating instruction manual.**

Please read through these instructions carefully in order for you to understand all functions and safety information.

microlife

# Guarantee Card

NC 120

Name of Purchaser / Nom de l'acheteur / Nombre del comprador / Nome do comprador / Name des Käufers / Naam koper / Alıcının Adı / Ονοματεπώνυμο αγοραστή / Ф.И.О. покупателя / Imię i nazwisko nabywcy / Vásárló neve / Име на купувача / Numele cumpărătorului / Jméno kupujícího / Meno zákazníka / Ime in priimek kupca

Serial Number / Numéro de série / Número de serie / Número de série / Serien-Nr. / Seriennummer / Seri Numarası / Αριθμός σειράς / Серийный номер / Numer serijny / Sorozatszám / Серийен номер / Număr de serie / Výrobní číslo / Výrobné číslo / Serijska številka

Date of Purchase / Date d'achat / Fecha de compra / Data da compra / Kaufdatum / Datum van aankoop / Satın Alma Tarihi / Ημερομηνία αγοράς / Дата покупки / Data zakupu / Vásárlás dátuma / Дата на закупуване / Data cumpărării / Datum nákupe / Dátum kúpy / Datum nakupa

Specialist Dealer / Revendeur / Vendedor especializado / Revendedor autorizado / Fachhändler / Specialist Dealer / Uzman Satıcı / Εξειδικευμένος αντιπρόσωπος / Специализированный дилер / Przedstawiciel / Forgalmazó / Специалист дистрибутор / Distribuitor de specialitate / Specializovaný dealer / Špecializovaný predajca / Specializirani trgovec

## Table of Contents

1. The Advantages of this Thermometer
2. Important Safety Instructions
3. How this Thermometer measures Temperature
4. Control Displays and Symbols
5. Changing between Body and Object Mode
6. Directions for Use
7. Changing between Celsius and Fahrenheit
8. How to recall 30 readings in Memory Mode
9. Error Messages
10. Cleaning and Disinfecting
11. Battery Replacement
12. Guarantee
13. Technical Specifications
14. www.microlife.com  
Guarantee Card (see Back Cover)

## 1. The Advantages of this Thermometer

---

### Measures in a matter of seconds

The innovative infrared technology allows the measurement without even touching the object. This guarantees safe and sanitary measurements within seconds.

### Multiple uses (wide range of measurement)

This thermometer offers a wide range of measurement from 0 - 100.0 °C / 32.0 - 212.0 °F, meaning the unit can be used to measure body temperature or it also has a feature allowing it to be used to measure surface temperature of the following examples:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

### Accurate and reliable

The unique probe assembly construction incorporates an advanced infrared sensor, ensuring that each measurement is accurate and reliable.

### Gentle and easy to use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can even be used on a sleeping child without causing any interruption.
- This thermometer is quick, therefore child-friendly.

### Auto-Display Memory

The last reading is automatically displayed for 2 seconds when the unit is switched ON.

### Multiple readings recall

Users will be able to recall the last 30 readings with a record of both time and date when entering the recall mode, enabling efficient tracking of temperature variations.

### Safe and hygienic

- No direct skin contact.
- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Cleaning of the probe can be done with an alcohol-moistened cotton tissue, making this thermometer completely hygienic for use by the whole family.

### Fever alarm

10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature equal to or higher than 37.5 °C.

## 2. Important Safety Instructions

---

- This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.
- **Never immerse this device in water or other liquids. For cleaning please follow the instructions in the «Cleaning and Disinfecting» section.**

- Do not use this device if you think it is damaged or notice anything unusual.
- Never open this device.
- A basic physiological effect called vasoconstriction can occur in the early stages of fever, resulting in a cool skin effect. The recorded temperature using this thermometer can, therefore, be unusually low.
- If the measurement result is not consistent with the patient's finding or unusually low, repeat the measurement every 15 minutes or double check the result by another core body temperature measurement.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.
- The function of this device may be compromised when used close to strong electromagnetic fields such as mobile phones or radio installations and we recommend a distance of at least 1 m. In cases where you suspect this to be unavoidable, please verify if the device is working properly before use.

- Protect it from:
  - extreme temperatures
  - impact and dropping
  - contamination and dust
  - direct sunlight
  - heat and cold
- If the device is not going to be used for a prolonged period the battery should be removed.

**WARNING:** Use of this device is not intended as a substitute for consultation with your physician. This device is NOT waterproof so do not immerse in liquids.

### 3. How this Thermometer measures Temperature

This thermometer measures infrared energy radiated from the forehead as well as objects. This energy is collected through the lens and converted to a temperature value.

### 4. Control Displays and Symbols

- **All segments displayed** (8): Press the ON/OFF button (5) to turn on the unit; all segments will be shown for 2 seconds.
- **Memory** (9): The last reading will automatically be shown on the display for 2 seconds.

- **Ready for measuring** (10): When the unit is ready for measuring, the «°C» or «°F» icon will keep flashing while the mode icon (body or object) will be displayed.
- **Measurement complete** (11): The reading will be shown on the display (4) with the «°C» or «°F» icon and the mode icon steady. The unit is ready for the next measurement as soon as the «°C» or «°F» icon is flashing again.
- **Low battery indicator** (14): When the unit is turned on, the battery icon will keep flashing to remind the user to replace the battery.

### 5. Changing between Body and Object Mode

For changing from body to object mode, slide the mode switch (6) at the side of the thermometer downwards. For switching back to body mode, slide the switch up again.

### 6. Directions for Use

#### Measuring in body mode

1. Press the ON/OFF button (5). The display (4) is activated to show all segments for 2 seconds.
2. The last measurement reading will automatically be shown on the display for 2 seconds with the «M» icon (9).
3. When the «°C» or «°F» icon is flashing, a beep sound is heard and the thermometer is ready for measuring (10).

4. **Aim the thermometer at the center of the forehead with a distance of no more than 5 cm.** Please remove any hair, sweat or dirt from the forehead before measuring to improve the accuracy of the readings.
5. **Press the START button (3) and ensure that the measuring sensor is aimed at the center of the forehead.** After 3 seconds a long beep will verify the completion of measurement.
6. Read the recorded temperature from the LCD display.

#### Measuring in object mode

1. Follow steps 1-3 above, then aim the thermometer at the center of the object you want to measure with a distance of no more than 5 cm. Press the START button (3). After 3 seconds a long beep will verify the completion of measurement.
2. Read the recorded temperature from the LCD display.

#### NOTE:

- **Patients and thermometer should stay in similar room condition for at least 30 minutes.**
- Don't take a measurement while or immediately after nursing a baby.
- Don't use the thermometer in high humidity environments.

- Patients should not drink, eat or exercise before/while taking the measurement.
- Don't move the measurement device from the measuring area before hearing the termination beep.
- Use an alcohol swab to carefully clean the probe and wait for 15 minutes before taking a measurement on another patient.
- 10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature equal to or higher than 37.5 °C.
- Always take the temperature in the same location, since temperature readings may vary according to locations.
- Doctors recommend rectal measurement for newborn infants within the first 6 months, as all other measuring methods might lead to ambiguous results. If using a non contact thermometer on those infants, we always recommend verifying the readings with a rectal measurement.
- In the following situations it is recommended that three temperatures are taken with the highest one taken as the reading:
  1. Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.

2. When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the device and obtains consistent readings.

3. If the measurement is surprisingly low.

- **Readings from different measuring sites should not be compared as the normal body temperature varies by measuring site and time of day,** being highest in the evening and lowest about one hour before waking up.

Normal body temperature ranges:

- Axillar: 34.7 - 37.3 °C / 94.5 - 99.1 °F
- Oral: 35.5 - 37.5 °C / 95.9 - 99.5 °F
- Rectal: 36.6 - 38.0 °C / 97.9 - 100.4 °F
- Microlife NC 120: 35.4 - 37.4 °C / 95.7 - 99.3 °F

#### 7. Changing between Celsius and Fahrenheit

This thermometer can display temperature readings in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn OFF the unit, **press and hold** the START button (3) for 5 seconds; after 5 seconds, the current measurement scale («°C» or «°F» icon) will flash on the display (15). Change the measurement scale between °C and °F by pressing the START button (3). When the measurement scale has been chosen, wait for 5 seconds and the unit will automatically enter the «ready for measuring» mode.

## 8. How to recall 30 readings in Memory Mode

---

This thermometer can recall the last 30 readings.

- **Recall mode** (16): Press the START button (3) to enter Recall mode when power is off. The memory icon «M» will flash.
- **Reading 1 - the last reading** (17): Press and release the START button (3) to recall the last reading. Display 1 alone with memory icon.
- **Reading 30 - readings in succession:** Press and release the START button (3) consecutively to recall the last 30 readings in succession.

Pressing and releasing the START button (3) after the last 30 readings have been recalled will resume the above sequence from reading 1.

## 9. Error Messages

---

- **Measured temperature too high** (18): Displays «H» when measured temperature is higher than 42.2 °C / 108.0 °F in body mode or 100 °C / 212 °F in object mode.
- **Measured temperature too low** (19): Displays «L» when measured temperature is lower than 34.0 °C / 93.2 °F in body mode or 0 °C / 32 °F in object mode.

- **Ambient temperature too high** (20): Displays «H» and ☰ when ambient temperature is higher than 40.0 °C / 104.0 °F.
- **Ambient temperature too low** (21): Displays «L» and ☰ when ambient temperature is lower than 16.0 °C / 60.8 °F in body mode or lower than 5.0 °C / 41.0 °F in object mode.
- **Error function display** (22): The system has a malfunction.
- **Blank display** (23): Check if the battery has been loaded correctly. Also check polarity (<+> and <->) of the battery.
- **Flat battery indication** (24): If the steady battery icon is the only symbol shown on the display, the battery should be replaced immediately.

## 10. Cleaning and Disinfecting

---

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring sensor. Ensure that no liquid enters the interior of the device. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the device in water or other cleaning liquids. Take care not to scratch the surface of the sensor lens and the display.

## 11. Battery Replacement

---

This device is supplied with one lithium battery, type CR2032. Replace with a new CR2032 battery when the flashing battery symbol appears on the display (24). Remove the battery cover (25) by sliding it in the direction shown. Insert the new battery with the + at the top.

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

## 12. Guarantee

---

This device is covered by a **2 year guarantee** from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) confirming date of purchase or the receipt.

- The guarantee covers the device. The battery and packaging are not included.
- Opening or altering the device invalidates the guarantee.
- The guarantee does not cover damage caused by improper handling, a discharged battery, accidents or non-compliance with the operating instructions. Please contact Microlife-service.

### 13. Technical Specifications

---

<b>Type:</b>	Non Contact Thermometer NC 120
<b>Measurement range:</b>	Body mode: 34.0-42.2 °C / 93.2-108.0 °F Object mode: 0-100.0 °C / 32-212.0 °F
<b>Resolution:</b>	0.1 °C / °F
<b>Measurement accuracy:</b>	Laboratory: ±0.2 °C, 36.0 ~ 39.0 °C ±0.4 °F, 96.8 ~ 102.2 °F
<b>Display:</b>	Liquid Crystal Display, 4 digits plus special icons
<b>Acoustic:</b>	The unit is turned ON and ready for the measurement: 1 short beep. Complete the measurement: 1 long beep (1 sec.) if the reading is less than 37.5 °C / 99.5 °F, 10 short «beep» sounds, if the reading is equal to or greater than 37.5 °C / 99.5 °F. System error or malfunction: 3 short «bi» sounds.
<b>Memory:</b>	Auto-Display the last measured temperature 30 readings recall in the memory mode

<b>Backlight:</b>	The display light will be GREEN for 4 seconds, when the unit is turned ON. The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5 °C / 99.5 °F. The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5 °C / 99.5 °F.
<b>Operating conditions:</b>	Body mode: 16-40.0 °C / 60.8-104.0 °F Object mode: 5-40.0 °C / 41-104.0 °F
<b>Storage conditions:</b>	-20 - +50 °C / -4 - +122 °F 15-95 % relative maximum humidity
<b>Automatic Switch-off:</b>	Approx. 1 minute after last measurement has been taken.
<b>Battery:</b>	1 x CR2032 battery 3V
<b>Battery lifetime:</b>	approx. 2000 measurements (using a new battery)
<b>Dimensions:</b>	120 x 28 x 21 mm
<b>Weight:</b>	39 g (with battery), 36 g (w/o battery)
<b>IP Class:</b>	IP21

**Reference to standards:** ASTM E1965;  
IEC 60601-1; IEC 60601-1-2 (EMC);  
IEC 60601-1-11

**Expected service life:** 5 years or 12000 measurements

This device complies with the requirements of the Medical Device Directive 93/42/EEC.  
Technical alterations reserved.  
According to the Medical Product User Act a biennial technical inspection is recommended for professional users.  
Please observe the applicable disposal regulations.

### 14. [www.microlife.com](http://www.microlife.com)

---

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at [www.microlife.com](http://www.microlife.com).

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