Microlife FR 1DM1

3 sec. Measurement

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Fever Alarm / Silent Glow™ Technology

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Forehead Thermometer

Instruction Manual

Please read these instructions carefully before using the instrument and keep them in a safe place.

The Microlife Forehead Thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, the Microlife Forehead Thermometer can provide a stable, heat-interference-free reading with each measurement. The instrument performs a self-test every time it is switched on to always guarantee the specified accuracy of measurements.

The Microlife Forehead Thermometer allows a simple and quick measurement of the body temperature. Please keep in mind that, like any skin temperature, the temperature at the temporal area might differ from the core body temperature (e.g. obtained by a rectal probe), which is due to a basic physiologic effect called vasoconstriction. Since vasoconstriction can occur in the early stages of fever, resulting in a cool skin effect, the recorded temperature using a forehead thermometer may be suspiciously low.

Skin temperature of babies reacts very quickly on ambient temperature. Therefore do not take a measurement while or after nursing; skin temperature may be lower than core temperature.

If the measurement result is not consistent with the patient's finding, always double check the result by another core body temperature measurement, e.g. use an IR – eardrum thermometer or a rectal pen-type thermometer.

Type BF applied part

Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.

Table of Contents

1. The Advantages of your Thermometer

2. Important Safety Instructions

3. Product Description

4. How the Thermometer Measures Body Temperature

5. Control Displays and Symbols

6. Directions for Use

7. Changing from Fahrenheit to Celsius and vice-versa

8. Error Messages

9. How to Recall 12 Readings in Memory Mode

10. Cleaning and Storage

11. Replacing the Battery

12. Guarantee

13. Technical Specifications

14. www.microlife.com

1. The Advantages of your Thermometer

Measurement in 3 seconds

The innovative infrared technology allows measurement simply by scanning the eyebrow area. After a scanning time of 3 seconds the measurement will be displayed on the LCD.

Accurate and reliable

Due to the unique probe assembly construction, incorporating an advanced infrared sensor this device can offer a temperature measurement comparable to one taken from the oral site with a regular digital thermometer.

Gentle and Easy to Use

- Appealing Penguin style delights the children.
- No need for Probe Covers at all, which makes the thermometer simple and easy to use.
- The thermometer can be used without interference to daily lifestyle. A measurement can be taken even while a child is sleeping.
- The thermometer is less threatening to a child than a rectal thermometer and more pleasant to use than an oral thermometer.

Auto-Display Memory

The last reading is automatically displayed for 3 seconds when the unit is switched ON.

Multiple Reading Recalls

Users will be able to recall the last 12 readings when entering the recall mode, enabling efficient tracking of temperature variations.

Safe and Hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.

High Temperature Indication

10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature above 37.5 °C.

2. Important Safety Instructions


- Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions when using on children.
- **Never immerse the thermometer into water or other liquids (not waterproof). For cleaning, please follow the instructions in the «Cleaning and Storage» section.**
- Keep the instrument away from direct exposure to the sun and keep it in a dust-free, dry area at a temperature between 16 - 40 °C (60.8 - 104 °F).
- Do not use the thermometer if there are signs of damage on the measuring tip or on the instrument itself. If damaged, do not attempt to repair the instrument! Please contact your nearest Microlife customer service bureau.
- This thermometer consists of high-quality precision parts. Do not drop the instrument! Protect it from severe impact and shock. Do not twist the instrument and the measuring probe!

WARNING:

- **Use of this thermometer is not intended as a substitute for consultation with your physician.**
- **Thermometer is not waterproof! Please NEVER immerse into liquids!**

3. Product Description

- (1) **Probe**
- (2) **START button**
- (3) **LCD Display**
- (4) **O/I button**
- (5) **Battery Cover**







4. How the Thermometer Measures Body Temperature

The thermometer measures infrared energy radiated from the skin above of the eyebrow area and the surrounding tissue. This energy is collected through the lens and converted to a temperature value.

Temperature readings obtained by scanning above the eyebrow area will provide the greatest accuracy.

5. Control Displays and Symbols

| LCD Display | Display Meaning | Description |
|--|------------------------|--|
|  | All segments displayed | Press the O/I button to turn on the unit, all segments will be shown for 2 seconds. |
|   | Memory | The last reading will be shown on the display automatically for 3 seconds. |
|   | Ready | The unit is ready for measurement, the °C or °F icon will keep flashing. |
|   | Measurement complete | The reading will be shown on the LCD display with the °C or °F icon steady, the unit is ready again for the next measurement as soon as the °C or °F icon is flashing again. |
|  | Low battery indication | When the unit is turned on, the battery icon will keep flashing to remind the user to replace the battery. |

6. Directions for Use

1. Press the O/I button. The LCD is activated to show all segments for 2 seconds.

2. The last measurement reading will be shown on the display automatically for 3 seconds with the «M» icon.

3. When the °C or °F icon is flashing, a beep sound is heard and the thermometer is ready for measurement.

4. Place the probe gently onto the area approximately 1 cm above the center of the eyebrow. If the eyebrow area is covered with hair, sweat or dirt, please remove the obstacle beforehand to improve the reading accuracy.

5. Press the START button once and release it. Start to move the instrument steadily from 1 cm above the center of the eyebrow to 1 cm above the tail of the eyebrow area. A long beep sound will emit after 3 seconds to acknowledge the completion of measurement. If the probe movement is too slow to reach the tail of the eyebrow before the long beep sound heard, repeat the measurement as above from point 4, but simply move the probe a little quicker this time.

6. Read the recorded temperature from the LCD display.

IMPORTANT:

- Patients should stay in steady - state room condition for at least 30 minutes.
- Don't take a measurement while or immediately after nursing a baby.
- Don't use the forehead thermometer in high humidity environments.
- Patients should not drink, eat, or exercise before/while taking the measurement.
- Before placing the thermometer probe onto the measurement area, remove dirt, hair, or sweat.

- Don't move the measurement device from the measuring area before hearing the termination beep.
- Use an alcohol swab to carefully clean the probe and wait for 15 minutes before taking a measurement on another patient.
- 10 short beeps will sound and a red LCD backlight will illuminate when the temperature is higher than 37.5 °C (99.5 °F) in order to alert the patient that he/she may have a temperature above 37.5 °C.
- Always take the temperature in the same location, since the temperature readings may vary from different locations.
- In the following situations it is recommended that three temperatures in the same location be taken and the highest one taken as the reading:
 - 1) New born infants in the first 100 days.
 - 2) Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 - 3) When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.










The thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.

7. Changing from Fahrenheit to Celsius and vice-versa

The thermometer can display temperature measurements in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn OFF the unit, **press and hold** the START button for 5 seconds until «-- -» and the flashing °C or °F icon is shown on the LCD. Switch the measurement scale between °C and °F by pressing the START button again. When the measurement scale has been chosen, wait for 5 seconds and the unit will enter the ready for measuring mode automatically.







8. Error Messages

| Display / Problem | Display Meaning | Possible cause and fault remedy |
|---|-------------------------------|---|
|   | Measured temperature too high | Displays «H» when measured temperature higher than 42.2 °C (108.0 °F). |
|   | Measured temperature too low | Displays «L» when measured temperature lower than 34.0 °C (93.2 °F). |
|   | Ambient temperature too high | Displays «H» in conjunction with the «▲» when ambient temperature is higher than 40.0 °C or 104.0 °F. |
|   | Ambient temperature too low | Display «L» in conjunction with the «▼» when ambient temperature is lower than 16.0 °C or 60.8 °F. |
|  | Error function display | When system has malfunction. |
|  | Blank display | Please check if the battery has been loaded correctly. Also check polarity (<+> and <->) of battery. |
|  | Dead battery indication | If the steady battery icon is the only symbol shown on the display, the battery should be replaced immediately. |

9. How to Recall 12 Readings in Memory Mode


The thermometer can recall the last 12 readings. To recall please follow sequence.

| Display | Display Meaning | Description |
|--|--|--|
|  | Recall Mode | Press the START button to enter Recall Mode when power is off. The memory icon «M» flashes. |
|   | Reading 1 (the last reading) | Press and release the START button to recall the last reading. Display 1 alone with memory icon. |
|   | Reading 2 (the second latest reading) | Press and release the START button to recall the second latest reading. |
|   | Reading 12 (the last reading) | Press and release the START button consecutively to recall readings in succession, up to the last 12 readings. |

Pressing and releasing the START button after the last 12 readings have been recalled will resume the above sequence from reading 1.

10. Cleaning and Storage

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and **never immerse the instrument in water or other cleaning liquids.**


Take care not to scratch the surface of the LCD.

Remove the battery from the instrument if it is not required for extended periods of time in order to avoid damage to the thermometer resulting from a leaking battery.

11. Replacing the Battery

The thermometer is supplied with one lithium battery, type CR2032. Replace with a new CR2032 battery when the flashing battery symbol appears on the LCD display.

To remove the battery, insert the tip of a screwdriver into the slot as indicated to lever it out of the battery cell.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

12. Guarantee

Subject to the following conditions this high-quality measuring instrument is covered by a **2 year guarantee** from the date of purchase. Warranty claims must be lodged within the guarantee period.

This product was manufactured with the utmost of care according to international quality standards. Should you have reason for complaints despite this, please send the instrument, accompanied by the completed Guarantee Card with dealer's stamp as well as original proof of purchase directly or through your medical supplier to your closest Microlife Distributor.

Damage resulting from incorrect use is not covered by the guarantee. Battery and packaging are excluded from the guarantee. Claims beyond this, including claims for damages, are excluded.

Name and address of responsible dealer:

13. Technical Specifications

| | |
|---------------------------------------|---|
| Type: | Forehead Thermometer FR 1DM1 |
| Measuring Range: | 34.0 °C to 42.2 °C (93.2 °F to 108.0 °F) |
| Resolution: | 0.1 °C / °F |
| Accuracy: | Laboratory: ± 0.2 °C, 36.0 ~ 39.0 °C ± 0.4 °F, 96.8 ~ 102.2 °F (Ambient Temperature: 16 - 40 °C) |
| Display: | Liquid Cristal Display , 4 digits plus special icons |
| Acoustic: | a. The unit is turned ON and ready for the measurement: 1 short «beep» sound. b. Complete the measurement: 1 long beep (1 sec.), if the reading is less than 37.5 °C (99.4 °F), 10 short «beep» sounds, if the reading is equal to or greater than 37.5 °C (99.4 °F). c. System error or malfunction: 3 short beeps. |
| Memory: | a. Auto-Display the last measured temperature when switched on. b. 12 readings recall in the Memory Mode. |
| Backlight: | a. The display light will be GREEN for 4 seconds, when the unit is turned ON. b. The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5 °C (99.4 °F). c. The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5 °C (99.4 °F). |
| Operating temperature: | 16 °C to 40 °C (60.8 °F to 104 °F) |
| Storage/transport temperature: | 15 - 95 % relative maximum humidity -20 °C to +50 °C (-4 °F to 122 °F) |
| Automatic Switch-off: | 15 - 95 % relative maximum humidity Approx. 1 minute after last measurement has been taken. |
| Battery: | CR2032 Battery (X1) 3V at least 1000 measurements |
| Dimensions: | 122.36mm (L) x 42.2mm (W) x 43.11mm (H) |
| Weight: | 48g (with battery), 45g (w/o battery) |
| Reference to standards: | ASTM E1965; IEC 60601-1; IEC 60601-1-2 (EMC) |

This device complies with the requirements of the Medical Device Directive 93/42/EEC. Technical alterations reserved.

According to the Medical Product User Act a biennial technical inspection is recommended for professional users.

Please observe the applicable disposal regulations.

14. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com.

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