Microlife IR 1DE1

1 sec. Measurement (Scan-Peak-method)

Mesure en 1 sec. / Medición en 1 seg. / 1 Sek. Messung /
Medição em 1 seg.

Multiple Use (Wide Range 0 °C - 100 °C)

- Body temperature / Température du corps / Temperatura corporal / Körpertemperatur / Temperatura corporal
- Milk surface temperature / Température de surface du lait / Temperatura superficial de la leche / Milch-Oberflächen-temperatur / Temperatura da superfície do leite
- Surface temperature of Baby's bath / Température de surface du bain de bébé / Temperatura superficial del agua del baño del bebé / Badewasser-Oberflächentemperatur / Temperatura superficial del agua del baño
- Ambient temperature / Température ambiante / Temperatura ambiente / Raumtemperatur / Temperatura ambiente

Fever alarm / Silent Glow™ Technology

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Probe Cover Free Ear Thermometer

Instruction Manual

Please read these instructions carefully before using the instrument and keep them in a safe place.

The Microlife Probe Cover Free Ear Thermometer IR 1DE1 is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, the IR 1DE1 can provide a stable, heat-interference-free reading with each measurement. The instrument performs a self-test every time it is switched on to always guarantee the specified accuracy of measurements.

The Microlife Probe Cover Free Ear Thermometer IR 1DE1 is intended for the periodic measurement and monitoring of human body temperature in the home. It is intended for use on people of all ages.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Table of Contents

- 1. The Advantages of this Thermometer**

- 2. Important Safety Instructions**

- 3. Product Description**

- 4. How this Thermometer Measures Ear Temperature**

- 5. Control Displays and Symbols**

- 6. Directions for Use**

- 7. Changing from Fahrenheit to Celsius and vice-versa**

- 8. How to recall 12 readings in Memory Mode**

- 9. Error Messages**

- 10. Cleaning and Storage**

- 11. Replacing the Battery**

- 12. Guarantee**

- 13. Technical Specifications**

- 14. www.microlife.com**

1. The Advantages of this Thermometer

Multiple Use (Wide Range Measurement)

The IR 1DE1 offers a wide measurement range feature from 0 - 100.0 °C (32.0 - 212.0 °F); meaning the unit can be used as an ear thermometer to measure body temperature, but it can also be used to measure surface temperature of following objects:

- Milk surface temperature in baby's bottle.
- Surface temperature of baby's bath.
- Ambient temperature.

Probe Cover Free

Since a probe cover is not required with this thermometer it makes using the IR 1DE1 more user-friendly and more cost-effective.

Measurement in 1 second

The innovative infrared technology allows measurement of ear temperature in only 1 second.

Accurate and reliable

The unique probe assembly construction, incorporating an advanced infrared sensor ensure that each measurement is accurate and reliable.

Gentle and Easy to Use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can be used even on a sleeping child, thereby causing no interruption.
- This thermometer is fast and therefore pleasant to use for children.

Auto-Display Memory

The last reading is automatically displayed for 2 seconds when the unit is switched ON.

Multiple Reading Recalls

Users will be able to recall the last 12 readings when entering the recall mode, enabling efficient tracking of temperature variations.

Safe and Hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Cleaning of the probe tip can be done with an alcohol-moistened cotton tissue, making the IR 1DE1 completely hygienic for use by the whole family.

High Temperature Indication

10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature above 37.5 °C.

2. Important Safety Instructions

- Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions when using on children.
- **Never immerse this thermometer into water or other liquids (not waterproof). For cleaning please follow the instructions in the «Cleaning and Storage» section.**
- Keep the instrument away from direct exposure to the sun and keep it in a dust-free, dry area at a temperature between 5 - 40 °C (41 - 104 °F).
- Do not use the thermometer if there are signs of damage to the measuring tip or to the instrument itself. If damaged, do not attempt to repair the instrument! Please contact your nearest Microlife customer service department.
- Earwax in the ear canal may cause a lower temperature reading. Therefore it is important to ensure the subject's ear canal is clean.
- This thermometer consists of high-quality precision parts. Do not drop the instrument! Protect it from severe impact and shock. Do not twist the instrument and the measuring probe!

WARNING:

- **Use of this thermometer is not intended as a substitute for consultation with your physician.**
- **Thermometer is not waterproof! Please NEVER immerse into liquids without a probe cover.**

3. Product Description


(1) Probe

(2) START button

(3) LCD Display

(4) O/I button

(5) Battery Cover

4. How this Thermometer Measures Ear Temperature

This thermometer measures infrared energy radiated from the eardrum and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. The measured reading obtained directly from the eardrum (Tympanic Membrane) can ensure the most accurate ear temperature.



Measurements taken from the surrounding tissue of the ear canal generate lower readings and may result in misdiagnosis of a fever.

To avoid an inaccurate measurement:

- Switch on the thermometer by pressing the O/I button.
- After one beep is heard (and the temperature scale icon is flashing), straighten the ear canal by gently pulling the middle of the ear back and up.
- Place the probe firmly into the ear canal, press the START button and keep the probe in the ear until the thermometer beeps to identify the completion of the measurement.

This thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.

5. Control Displays and Symbols

LCD Display	Display Meaning	Description
	All segments displayed	Press the O/I button to turn on the unit, all segments will be shown for 2 seconds.
	Memory	The last reading will be shown on the display automatically for 2 seconds.
	Ready	The unit is ready for measurement, the °C or °F icon will keep flashing.
	Measurement complete	The reading will be shown on the LCD display with the °C or °F icon flashing, the unit is ready again for the next measurement.
	Out-of-ear temperature indication	A crossed-ear-icon shown on the LCD if the reading falls outside the range 32.0 ~ 42.2 °C (89.6 ~ 108.0 °F)
	Low battery indication	When the unit is turned on, the battery icon will keep flashing to remind the user to replace the battery.

6. Directions for Use

1. Press the O/I button. The LCD is activated to show all segments for 2 seconds.

2. The last measurement reading will be shown on the display automatically for 2 seconds with the «M» icon.

3. When the °C or °F icon is flashing, a beep sound is heard and the thermometer is ready for measurement.

4. Straighten the ear canal by pulling the ear up and back to give a clear view of the eardrum.

- For children under 1 year:
Pull the ear straight back.

- Children above 1 year and adults:
Pull the ear up and back.

5. While gently pulling the ear, insert the probe snugly into the ear canal and press the START button. Release the button and wait for the beep sound. This is the indication that confirms the end of measurement.

6. Remove the thermometer from the ear canal. The LCD displays the measured temperature.

NOTE:

- 10 short beeps will sound and a red LCD backlight will illuminate when the temperature is higher than 37.5 °C (99.5 °F) in order to alert the patient that he/she may have a temperature above 37.5 °C.
- **In order to assure the accurate readings, please wait at least 30 seconds after 3-5 continuous measurements.**
- Accumulation of ear wax on the probe tip can result in less accurate temperature readings or cross infection between users. **Therefore, it is essential that a clean probe tip is used before each measurement.** Please follow the instructions in the «Cleaning and Storage» section for cleaning.
- **After cleaning the probe tip with alcohol, it is necessary to wait 5 minutes before taking the next measurement,** in order to allow the thermometer to reach its operating reference temperature.

- For an infant, it is best to have the child laying flat with his head sideways so the ear is facing upwards. For an older child or adult, it is best to stand behind and slightly to the side of the patient.
- Always take the temperature in the same ear, since the temperature readings may be different from the right ear and left ear.
- Wait for a few minutes to take the ear temperature after sleeping.
- In the following situations it is recommended that 3 measurements are done in the same ear and the highest one is taken as the reading:
 - 1) New born infants in the first 100 days.
 - 2) Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 - 3) When the user is learning how to use the IR thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.




7. Changing from Fahrenheit to Celsius and vice-versa

The thermometer can display temperature measurements in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn OFF the unit, **press and hold** the START button for 5 seconds until «-- -» and the flashing °C or °F icon is shown on the LCD. Switch the measurement scale between °C and °F by pressing the START button again. When the measurement scale has been chosen, wait for 5 seconds and the unit will enter the ready for measuring mode automatically.

8. How to recall 12 readings in Memory Mode

This thermometer can recall the last 12 readings. To recall please follow sequence.

Display	Display Meaning	Description
	Recall Mode	Press the START button to enter Recall Mode when power is off. The memory icon «M» flashes.
	Reading 1 (the last reading)	Press and release the START button to recall the last reading. Display 1 alone with memory icon.
	Reading 2 (the second latest reading)	Press and release the START button to recall the second latest reading
	Reading 12 (the last reading)	Press and release the START button consecutively to recall readings in succession, up to the last 12 readings.


Pressing and releasing the START button after the last 12 readings have been recalled will resume the above sequence from reading 1.

9. Error messages

Display / Problem	Display Meaning	Possible cause and fault remedy
	Measured Temperature too high	Displays «H» when measured temperature higher than 100.0 °C or 212.0 °F.
	Measured temperature too low	Displays «L» when measured temperature lower than 0 °C or 32.0 °F.
	Ambient temperature too high	Displays «H» in conjunction with the «▲» when ambient temperature is higher than 40.0 °C or 104.0 °F.
	Ambient temperature too low	Displays «L» in conjunction with the «▼» when ambient temperature is lower than 5.0 °C or 41.0 °F.
	Error function display	When system has malfunction.
	Blank display	Please check if the battery has been loaded correctly. Also check polarity (<+> and <->) of batteries.
	Dead battery indication	If the steady battery icon is the only symbol shown on the display, the batteries should be replaced immediately.



10. Cleaning and Storage

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the probe lens and the LCD. Remove the battery from the instrument if it is not required for extended periods of time in order to avoid damage to the thermometer resulting from a leaking battery.

11. Replacing the battery

The thermometer is supplied with one lithium battery, type CR2032. Replace with a new CR2032 battery when the flashing battery symbol appears on the LCD display. Remove the battery cover by sliding it in the direction shown. Take out the battery and put in a new one.

12. Guarantee

Subject to the following conditions this high-quality measuring instrument is covered by a **2 year guarantee** from the date of purchase. Warranty claims must be lodged within the guarantee period.

This product was manufactured with the utmost of care according to international quality standards. Should you have reason for complaints despite this, please send the instrument, accompanied by the completed Guarantee Card with dealer's stamp as well as original proof of purchase directly or through your medical supplier to your closest Microlife Distributor.

Damage resulting from incorrect use is not covered by the guarantee. Battery and packaging are excluded from the guarantee. Claims beyond this, including claims for damages, are excluded. Name and address of responsible dealer:

13. Technical Specifications

Type:	Probe Cover Free Ear Thermometer IR 1DE1
Measuring Range:	0 °C to 100.0 °C (32.0 °F to 212.0 °F)
Accuracy:	Laboratory: ± 0.2 °C, 32.0 ~ 42.2 °C (± 0.4 °F, 89.6 ~ 108.0 °F)
Display:	Liquid Cristal Display with indicating unit 0.1 °C/°F
Accoustic:	a. The unit is turned ON and ready for the measurement: 1 short «bi» sound. b. Complete the measurement: 1 long beep (1 sec.), if the reading is less than 37.5 °C (99.4 °F), 10 short «beep» sounds, if the reading is equal to or greater than 37.5 °C (99.4 °F). c. System error or malfunction: 3 short «bi» sounds.
Memory:	a. Auto-Display the last measured temperature when switched on. b. 12 readings recall in the Memory Mode.
Backlight:	a. The display will be lighted GREEN for 4 sec., when the unit is turned ON. b. The display will be lighted GREEN for 5 sec., when a measurement has been completed with a reading less than 37.5 °C (99.4 °F). c. The display will be lighted RED for 5 sec., when a measurement is completed with a reading equal to or higher than 37.5 °C (99.4 °F).
Operating temperature:	5 °C to 40 °C (41.0 °F to 104 °F)
Storage/transport temperature:	-25 °C to + 55 °C (-13 °F to 131 °F)
Automatic Switch-off:	Approx. 1 minute after last measurement. Foolproof Switch-off: Auto-off if START button is pressed and held for 15 seconds.
Battery:	CR2032 Battery (X1) - at least 1000 measurements
Dimensions:	153 mm (L) x 31 mm (W) x 40 mm (H)
Weight:	53 g (w/ battery), 50 g (w/o battery)
Standards:	Complies with ASTM E-1965 requirements

According to the Medical Product User Act a biennial technical inspection is recommenden for professional users. Please observe the applicable disposal regulations.

14. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com.

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