

# **Demo APP for MicroLifeDeviceSDK (Body Temperature / Thermometer) - Android**

## **Table of Contents**

- Chapter 1    Development Environment**
- Chapter 2    Entry Point and Bluetooth LE Protocol**
- Chapter 3    Bluetooth LE Protocol & APIs**
- Chapter 4    Thermometer APIs**
- Chapter 5    User Interface and Functionality**
- Chapter 6    Demo App**

## Chapter 1 Development Environment

The supported SDK version is as follow:

```
compileSdkVersion 26
buildToolsVersion '26.0.3'

defaultConfig {
    minSdkVersion 19
    targetSdkVersion 26
    versionCode 1
    versionName "1.3"
}
```

- 1.1. Add the library “sdk-release.arr” into the “libs” directory.
- 1.2. In the “build.gradle”, add the description as bellows:

```
compile(name:'sdk-release', ext:'aar')
compile(name:'scalelesdk-v1.4.0', ext:'aar')
```

## Chapter 2 Entry Point and Bluetooth LE Protocol

The “ChoseActivity” is the entry point of the sample application.

The “BtTestActivity” is dedicated to Thermometer.

```
<activity
    android:name=".BPMTestActivity"
    android:screenOrientation="portrait"
    android:windowSoftInputMode="stateHidden" />
<activity
    android:name=".WeightTestActivity"
    android:screenOrientation="portrait"
    android:windowSoftInputMode="stateHidden" />
<activity
    android:name=".BtTestActivity"
    android:screenOrientation="portrait" />
<activity
    android:name=".WBPTTestActivity"
    android:screenOrientation="portrait" />
<activity
    android:name=".ChoseActivity"
    android:screenOrientation="portrait">
    <intent-filter>
        <action android:name="android.intent.action.MAIN" />
        <category android:name="android.intent.category.LAUNCHER" />
    </intent-filter>
</activity>

<activity android:name=".ConnectionActivity">
```

- 2.1 Initialize the instance “wbpProtocol”. This is to fulfill Bluetooth LE features and connection sequence.

```
//Initialize the connection SDK
Global.wbpProtocol = WBPPProtocol.getInstance
    ( aty: this, isSimulation: false, isPrintLog: true, Global.sdkid_WBP);
Global.wbpProtocol.setOnConnectStateListener(this);
Global.wbpProtocol.setOnDataResponseListener(this);
Global.wbpProtocol.setOnNotifyStateListener(this);
Global.wbpProtocol.setOnWriteStateListener(this);
```

- 2.1.1 The “setOnConnectStateListener()” is to get the connection status of device.
- 2.1.2 The “setOnDataResponseListener()” is to get the response from device.
- 2.1.3 The “setOnNotifyStateListener()” is to get the data which is response from device.

2.1.4 The “setOnWriteStateListener()” is to get the data which is sent to device.

2.2 The “isEnabledBt()” or “ isSupportBluetooth() is to check if the smartphone’s Bluetooth is enabled or not. The “isSupportBluetooth()” will prompt a warning message to inform user to turn on Bluetooth if it is disabled.

## Chapter 3 Bluetooth LE Protocol & APIs

### 3.1. Instance of Bluetooth LE Protocol :

#### 3.1.1. Interface :

	public static * Protocol getInstance(Activity aty, boolean isSimulation, boolean isPrintLog, String sdkid)
Definition	Initialize Bluetooth LE Protocol for WatchBP Home A device
Parameter	Activity aty : name of activity or this boolean isSimulation : is simulator or device boolean isPrintLog : is printing log or not. String sdkid : SDK ID of designated device
	<pre>//Initialize the connection SDK Global.bpmProtocol = BPMProtocol.getInstance ( aty: this, isSimulation: false, isPrintLog: true, Global.sdkid);</pre>

### 3.2. Connection State and Result :

#### 3.2.1. Interface :

	public void setOnConnectStateListener(OnConnectStateListener l)
Definition	The “setOnConnectStateListener()” is to get the connection status of device.

#### 3.2.2. Delegate :

	void onBtStateChanged(boolean isEnabled)
Definition	The “onBtStateChanged()” is to monitor the state of Enabled or Disabled.

	void onScanResult(String mac, String name, int rssi)
Definition	This is to get Bluetooth information of devices which discovered in the vicinity.
Parameter	macAddress: MAC of device name: device name RSSI: RSSI

	void onConnectionState(ConnectState state)
Definition	The “onConnectionState()” is to monitor the status of connection.
Parameter	<pre> public enum ConnectState {     ScanFinish,           //Scan finish     Connected,            //Connect success     Disconnect,           //Disconnect     ConnectTimeout,       //Connection timeout     ScaleWake,            //Scale Wake [EBodyProtocol limited]     ScaleSleep            //Scale Sleep [EBodyProtocol limited] } </pre>

### 3.3. Device scanning or discovery :

#### 3.3.1. Interface :

	public void startScan(int timeout)
Definition	The “startScan()” is for device scanning or discovery. The result will be shown with the “onScanResult”.
Parameter	int timeout

	public void stopScan()
Definition	Terminate the scanning process.

#### 3.3.2. Delegate :

	void onConnectionState(ConnectState state)
Definition	The “onConnectionState()” is to monitor the status of scanning.
Parameter	<pre> public enum ConnectState {     ScanFinish,           //Scan finish     Connected,            //Connect success     Disconnect,           //Disconnect     ConnectTimeout,       //Connection timeout     ScaleWake,            //Scale Wake [EBodyProtocol limited]     ScaleSleep            //Scale Sleep [EBodyProtocol limited] } </pre>

### 3.4. Connection :

#### 3.4.1. Interface :

	public void connect(String macAddress)
--	--

Definition	Connect to device with MAC address.
Parameter	macAddress: MAC of device

#### 3.4.2. Delegate :

	void onConnectionState(ConnectState state)
Definition	The “onConnectionState()” is to monitor the status of connection.
Parameter	<pre> public enum ConnectState {     ScanFinish,           //Scan finish     Connected,            //Connect success     Disconnect,           //Disconnect     ConnectTimeout,       //Connection timeout     ScaleWake,            //Scale Wake [EBodyProtocol limited]     ScaleSleep            //Scale Sleep [EBodyProtocol limited] } </pre>

### 3.5. Bonding :

#### 3.5.1. Interface :

	public void bond(String macAddress)
Definition	Binding specified device by MAC
Parameter	macAddress: MAC of device

#### 3.5.2. Delegate :

	void onConnectionState(ConnectState state)
Definition	The “onConnectionState()” is to monitor the status of connection.
Parameter	<pre> public enum ConnectState {     ScanFinish,           //Scan finish     Connected,            //Connect success     Disconnect,           //Disconnect     ConnectTimeout,       //Connection timeout     ScaleWake,            //Scale Wake [EBodyProtocol limited]     ScaleSleep            //Scale Sleep [EBodyProtocol limited] } </pre>

### 3.6. Disconnection :

#### 3.6.1. Interface :

	public void disconnect()
--	--------------------------

Definition	Disconnect device.
------------	--------------------

### 3.6.2. Delegate :

	void onConnectionState(ConnectState state)
Definition	The “onConnectionState()” is to monitor the status of disconnection.
Parameter	<pre> public enum ConnectState {     ScanFinish,           //Scan finish     Connected,            //Connect success     Disconnect,           //Disconnect     ConnectTimeout,       //Connection timeout     ScaleWake,            //Scale Wake [EBodyProtocol limited]     ScaleSleep            //Scale Sleep [EBodyProtocol limited] } </pre>



## Chapter 4 Thermometer APIs

- 4.1. Received device information : Once Bluetooth connects to the thermometer successfully, the first time the thermometer sends the information.

	void onResponseDeviceInfo(String macAddress, int workMode, float batteryVoltage )
Parameter	macAddress : macAddress description  workMode : workMode description  batteryVoltage : batteryVoltage description

- 4.2. Received the current measurement result / data :

	void onResponseUploadMeasureData(ThermoMeasureData data)
Parameter	data : measurement result

## Chapter 5 User Interface and Functionality

### 5.1 Getting Started :


Start the app and then select the button “BODY TEMPERATURE” / “C” to communicate with the designate device Thermometer.



### 5.2 Operation Sequence :

The scanning (discovery) is automatically run to discover devices in the vicinity. If a device is bonded, it will be connected accordingly.

### 5.3 GUI Layout :



5.3.1 Region A : The log window is used to display information about communication handshake between App and device.

5.4 Refer to "BtTestActivity" from the demo application (sample code) to get more detailed.

## Chapter 6 Demo App

### 6.1. Scanning (Discovery), Connection, Data (Readings) :



1. Once starts the demo App, it will be conducting a scanning / discovery process. The device named NC150 BT is displayed in the vicinity. The connection can be created by the MAC address "18:7A:93:6E:68:A2".

2. The received data (included Device information & Readings) is as **Blue** part, and it can be decode in **Green** part. The following are the

(a) Device information:

THERMO : DeviceInfo ->  
macAddress = 187A936E68A2 ,  
workMode = 0 , batteryVoltage = 3.0.

(b) Readings:

THERMO : UploadMeasureData  
-> ThermoMeasureData =  
ThermoMeasureData{ambientTemperature=25.12,  
measureTemperature=36.18,  
mode=0, minute=38, month=3,  
day=12, hour=17, flagErr=0,  
flagFever=0, errorCode=0,  
year=20}.

3. The **Red** part is ACK to inform device the process is complete.

