microlife®

IR210 Ear Thermometer

EN →	1	PT →	60
ES →	11	NL →	70
FR →	21	GR →	80
IT →	31	AR →	91
DE →	41	FA →	100
TR →	51		

 Microlife Corporation
 9F, No. 431, RuiGuang Road, NeiHu
 Taipei, 114, Taiwan, China
www.microlife.com

 Microlife UAB
 P. Lukšio g. 32
 08222 Vilnius, Lithuania
eurep@microlife.lt


 Microlife AG
 Parkweg 2
 9443 Widnau
 Switzerland











CE 1639

IB IR 210 S-V11 2026-08-24

Microlife IR 210

Name of Purchaser / Nombre del comprador/
Nom de l'acheteur / Nome del rivenditore / Name
des Käufers / Nome do comprador / Naam koper
/ Ονοματεπώνυμο αγοραστή / Alicinin adı /
نام خریدار / اسم المشتري

Serial Number / Número de serie/ Numéro de
série / Numero di serie / Serien-Nr. / Número
de série / Serienummer / Αριθμός σειράς / Seri
numarası / شماره سریال / رقم التسلسل

Date of Purchase / Fecha de compra/ Date
d'achat / Data d'acquisto / Kaufdatum / Data
da compra / Datum van aankoop / Ημερομηνία
αγοράς / Satın alma tarihi / تاريخ شراء / تاريخ
اخرید

Specialist Dealer / Distribuidor especializado/
Revendeur / Categoria rivenditore / Fachhändler
/ Revendedor autorizado / Specialist Dealer /
Εξειδικευμένος αντιπρόσωπος / Uzman satıcı
نام فروشنده / التاجر المختص

- ① Measuring sensor
- ② Display
- ③ M-button (memory)
- ④ MODE button
- ⑤ Battery compartment cover
- ⑥ START/IO button
- ⑦ All segments displayed
- ⑧ Memory
- ⑨ Ready for measuring
- ⑩ Correct position indication
- ⑪ Measurement complete
- ⑫ Body mode
- ⑬ Object mode
- ⑭ Incorrect location indication
- ⑮ Error function display
- ⑯ Low battery indicator
- ⑰ Changing between Celsius and Fahrenheit
- ⑱ Recall mode
- ⑲ Recall the last 30 readings
- ⑳ Measured temperature too high
- ㉑ Measured temperature too low
- ㉒ Ambient temperature too high
- ㉓ Ambient temperature too low
- ㉔ Remove and replace the probe cover
- ㉕ Blank display
- ㉖ Flat battery
- ㉗ Date/Time
- ㉘ Beeper function setting
- ㉙ Replacing the battery
- ㉚ Storage holder
- ㉛ Probe cover
- ㉜ How to reload a new probe cover

This Microlife thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, this device can provide a stable, heat-interference-free reading with each measurement. The device performs a self-test every time it is switched on to always guarantee the specified accuracy of any measurement.

This thermometer has been clinically tested and proven to be safe and accurate when used in accordance to the operating instruction manual.

Please read through these instructions carefully in order for you to understand all functions and safety information.

Microlife® is a registered trademark of Microlife Corporation.

Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.

Table of contents

1. Important information

Device description
Indications
Contraindications
Warning
Caution
Electromagnetic compatibility
Adverse events and reporting

2. Device information

Package contents
Device accessories

3. Measurement preparation

4. The advantages of this thermometer

5. Control displays and symbols

6. Setting date, time and beeper functions

7. Directions for use

8. Changing between Celsius and Fahrenheit

9. How to recall 30 readings in memory mode

10. Battery installation and replacement

11. Error messages

12. Device maintenance and disposal

Cleaning and disinfecting
Storage
Calibration and support
Disposal

13. Guarantee

14. Technical specifications

15. Compliance information

16. Symbols and definitions

17. www.microlife.com

1. Important information

Device description

A digital infrared ear thermometer is a medical device that utilizes the principles of human body thermal radiation emission and digital signal processing to compute and provide a body temperature measurement.

Intended purposes:

The device is intended to measure body temperature.

Intended user:

The device is intended to be operated by adults and adolescents with adequate vision, motor functions, and education, capable of understanding the instructions for use and operating general household electrical appliances.

This device is also intended to be operated by healthcare professionals and personnel.

Intended patient:

The intended patients are the general population of all ages.

Intended use environment and conditions:

The device is intended for use in a home healthcare environment (e.g. general household without medically trained personnel) by patients (e.g. for self-measurement) or by a care giver, and in a professional healthcare facility (e.g. physician office).

Indications

This device measures body temperature for indications of:

- Screening for illness and conditions related to high or low body temperature.

Contraindications

- The device is not suitable for use on sites with any breached or compromised skin.

Warning

NOTE: Warning items indicate potentially hazardous situations, if not avoided, may result in death, critical or serious injury to the user or patient.

- This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.
- The measurement result of this device is not a medical diagnosis and not intended to substitute consultation and diagnosis by a qualified professional healthcare provider (e.g., physician, pharmacist, or other licensed health-care professionals). Do not use this device for self-diagnosis or for self-treatment of a medical condition. Seek advice from a health-care professional immediately if the patient is clearly unwell and/or having physiological or medical symptoms.
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed.
- Do not attempt to measure temperature of the body sites that are not specified in this instruction.
- Professional users must follow the country specific applicable medical devices ordinance.

Caution

NOTE: Caution items indicate potentially hazardous situations, if not avoided, may result in minor or negligible injury to the user or patient, or damage to the property or environment.

- Cleaning and/or disinfecting the thermometer is recommended before and after use if the thermometer is dirty or contaminated from use or after storage. Refer to «Cleaning and disinfecting» section for details.
- Do not use this device if you think it is damaged or notice anything unusual.
- Do not attempt to repair, modify, or open this device.
- A basic physiological effect called vasoconstriction can occur in the early stages of high body temperature, resulting in a cool skin effect. The temperature measured during this time will be lower than core body temperature.
- If the patient temperature is influenced (e.g. incubator, heating blanket) the measurement result provided by this device shall not be used to determine the presence/absence of high body temperature.

Electromagnetic compatibility

This device is compliant with EN 60601-1-2 Electromagnetic Disturbances standard.

Further documentation in compliance with EN 60601-1-2 EMC standard is available from Microlife on www.microlife.com/electro-magnetic-compatibility.

Do not use this device close to strong electromagnetic fields and portable radio frequency communication devices (for example microwave oven and mobile devices). Keep a minimum distance of 0.3 m from such devices when using this device.

This device is not certified to be used in the vicinity of medical equipment including high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) and computerized tomography (CT) instruments.

Adverse events and reporting

Please report any serious incident that has occurred in relation to the device, injury or adverse event to the local competent authority and to the manufacturer or to the european authorised representative (EU REP).

2. Device information

Package contents

- 1 x Microlife digital infrared ear thermometer
- 21 x Probe cover
- 1 x Instruction manual
- 1 x Storage holder
- 2 x 1.5 V alkaline batteries; type AAA

Device accessories

Single-use probe cover for hygienic measurement

Microlife probe covers for digital infrared ear thermometer are provided with the device and available as accessory. Use probe covers compatible with the device for measurements to avoid cross contamination between measurements and patients.

3. Measurement preparation

- The patients and thermometer should stay in the same room condition for at least 30 minutes before using.
- Patients should not drink, eat or exercise before/while taking the measurement.

- When taking multiple measurements for comparison, allow a 30 seconds interval between each measurement for measurement reliability.
- Doctors recommend rectal measurement for newborn infants within the first 6 months, as all other measuring methods might lead to ambiguous results.

To avoid an inaccurate measurement

- Fit a new Microlife branded and undamaged probe cover ③ on the measuring sensor ①.
- Switch on the thermometer by pressing the START/IO button ⑥.
- After one beep is heard (and the temperature scale icon is flashing), straighten the ear canal by gently pulling the middle of the ear back and up.
- Place the probe ① firmly into the ear canal. «Good» will be displayed with a short beep to confirm the device has detected the correct position. Press the START/IO button ⑥ and keep the probe in the ear until the thermometer beeps to signal the completion of the measurement.

Caution: DO NOT re-use the probe cover. The probe cover is intended only for single-use to prevent cross-contamination of the device between measurements. Remove used probe cover after measurement and dispose of it in accordance with applicable local regulations.

Caution: Use only Microlife probe covers compatible with the device. Using non-compatible accessories with the device may affect measurement accuracy.

Caution: In the following situations it is recommended that three temperatures are taken with the highest one taken as the reading:

- Children under three years of age with a compromised immune system and for whom the presence or absence of high body temperature is critical.
- When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the device and obtains consistent readings.
- If the measurement is surprisingly low.

4. The advantages of this thermometer

Measurement in 1 second

The innovative infrared technology allows measurement of ear temperature in only 1 second.

Multiple uses (wide range of measurement)

This thermometer offers a wide range of measurement from 0.1 - 99.9 °C / 32.2 - 211.8 °F, meaning the unit can be used to measure body temperature or it also has a feature allowing it to be used to measure surface temperature of the following examples:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

Correct position indication

ACCUsens guidance system confirms the right position in the ear with «good» displayed on the LCD display and a beep.

Probe cover

This thermometer is user-friendly and more hygienic with the usage of a probe cover.

Probe LED

This thermometer includes a probe LED light which enables the user to find the correct ear position in the dark.

Accurate and reliable

The unique probe assembly construction incorporates an advanced infrared sensor, ensuring that each measurement is accurate and reliable.

Gentle and easy to use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can even be used on a sleeping child without causing any interruption.
- This thermometer is quick, therefore child-friendly.

Multiple readings recall

Users will be able to recall the last 30 readings with a record of both time and date when entering the recall mode, enabling efficient tracking of temperature variations.

Safe and hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Using a new probe cover each time, ensures this thermometer is completely hygienic for use by the whole family.

High body temperature reminder

10 short beeps and a red LCD backlight reminds the patient that they may have a temperature equal to or higher than 37.5 °C.

5. Control displays and symbols

- **All segments displayed** ⑦: Press the START/IO button ⑥ to turn on the unit; all segments will be shown for 1 second.
- **Ready for measurement** ⑨: When the unit is ready for measurement, the «°C» or «°F» icon will keep flashing.
- The probe LED light is activated and will keep flashing.
- **Correct position indication** ⑩: The probe LED light will stop flashing (stays illuminated) and «good» will be displayed on the LCD, when the measuring sensor detects an appropriate position.
- **Measurement complete** ⑪: The reading will be shown on the display ② with the «°C» or «°F» icon; the unit is ready for the next measurement, when the «°C» or «°F» icon is flashing again.
- **Low battery indicator** ⑫: When the unit is turned on, the  icon will keep flashing to remind the user to replace the battery.

6. Setting date, time and beeper functions

Setting the date and time

1. After the new battery is fitted, the year number flashes in the display ②. You can set the year by pressing the M-button ③. To confirm and then set the month, press the MODE button ④.
2. Press the M-button ③ to set the month. Press the MODE button ④ to confirm and then set the day.
3. Follow the previously mentioned instructions to set the day, hours and minutes.
4. Once you have set the minutes and pressed the START/IO button ⑥, the date and time are set and the time is displayed.

- ☞ If no button is pressed for 20 seconds, the device automatically switches to ready for measuring ⑨.
- ☞ **Cancel time setup:** Press the START/IO button ⑥ during time setup. The LCD will show Date/Time icons with «--:--». After that press the START/IO button ⑥ to start the measurement. If no further action is taken within 60 seconds, the device will automatically turn off.
- ☞ **Change current date and time:** Press and hold the MODE button ④ for approx. 8 seconds until the year number starts to flash ②7. Now you can enter the new values as described above.

Setting the beeper

1. Press and hold the MODE button ④ for 3 seconds to set the beeper ②8.
2. Press the M-button ③ to either turn the beeper on or off. The beeper is activated when the beeper icon ②8 is shown without a cross.
- ☞ When the beeper setting has been chosen, press the START/IO button ⑥ to enter the «ready for measuring» mode ⑨; otherwise the device automatically switches to ready for measuring after 10 seconds.

7. Directions for use

Measuring in body mode

Important: Before each measurement, fit a new undamaged probe cover ③1 on the measuring sensor ①. Failure to do so will result in incorrect temperature measurements. How to correctly fit a new probe cover ③2 is pictured at the beginning of this instruction.

1. Press the START/IO button ⑥. The display ② is activated to show all segments for 1 second.
2. When the «°C» or «°F» icon is flashing, a beep sound is heard and the thermometer is ready for measuring ⑨.
3. The probe LED light is activated and will keep flashing.
4. Straighten the ear canal by pulling the ear up and back to give a clear view of the eardrum.
 - For children under 1 year; pull the ear straight back.
 - For children of 1 year to adult; pull the ear up and back. Also refer to the short instructions at the front.
5. While gently pulling the ear, insert the probe snugly into the ear canal.

6. The probe LED light will stop flashing (stays illuminated) and «good» will be displayed on the LCD, when the measuring sensor detects an appropriate position.
7. **Immediately** press the START/IO button ⑥. Release the button and wait for the beep sound. This indicates the end of the measurement.
8. Remove the thermometer from the ear canal. The display shows the measured temperature ①1.
9. Replace the probe cover ③1 before starting a new measurement.
10. For the next measurement wait until the «°C»/«°F» icon is flashing and follow steps 5 - 7 above.
11. Press and hold the START/IO button ⑥ for 3 seconds to turn off the device; otherwise the device will automatically switch off after approx. 60 seconds.

Measuring in object mode

Important: Remove the probe cover before each measurement in object mode. Failure to do so will result in incorrect temperature measurements.

1. Press the START/IO button ⑥. The display ② is activated to show all segments for 1 second.
2. Press the MODE button ④ to switch to object mode.
3. Aim the thermometer at the center of the object you want to measure with a distance of no more than 5 cm. **Press the START/IO button ⑥.** After 1 second a long beep will verify the completion of measurement.
4. Read the recorded temperature from the LCD display.
5. For the next measurement wait until the «°C»/«°F» icon is flashing and follow steps 3 - 4 above.

☞ NOTE:

- **Patients and thermometer should stay in similar room condition for at least 30 minutes.**
- **To ensure accurate readings, wait at least 30 sec. after 3-5 continuous measurements.**
- **It is essential that a new undamaged probe cover ③1 is used for each measurement.** Therefore, this device reminds the user to take off the used probe cover when turning off the device. The «probe cover» icon ②4 is displayed and the probe LED light will flash for 3 seconds. For cleaning, follow the instructions in the «Cleaning and disinfecting» section.

- **After cleaning the measuring sensor ① with alcohol, wait 5 minutes before taking the next measurement**, in order to allow the thermometer to reach its operating reference temperature.
- 10 short beeps and a red LCD backlight reminds the patient that they may have a temperature equal to or higher than 37.5 °C.
- For an infant, it is best to have the child lying flat with his/her head sideways so the ear is facing upwards. For an older child or adult, it is best to stand behind and slightly to the side of the patient.
- Always take the temperature in the same ear, since the temperature readings may be different from ear to ear.
- In the following situations it is recommended that three temperatures in the same ear be taken and the highest one taken as the reading:
 - New born infants in the first 100 days.
 - Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 - When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the device and obtains consistent readings.
 - If the measurement is surprisingly low.
- Don't take a measurement while or immediately after nursing a baby.
- Don't use the thermometer in high humidity environments.
- Patients should not drink, eat or exercise before/while taking the measurement.
- Doctors recommend rectal measurement for newborn infants within the first 6 months, as all other measuring methods might lead to ambiguous results.
- **Readings from different measuring sites should not be compared as the normal body temperature varies by measuring site and time of day**, being highest in the evening and lowest about one hour before waking up.
- Normal body temperature ranges:
 - Axillar: 34.7 - 37.3 °C / 94.5 - 99.1 °F
 - Oral: 35.5 - 37.5 °C / 95.9 - 99.5 °F
 - Rectal: 36.6 - 38.0 °C / 97.9 - 100.4 °F
 - Microlife IR 210: 35.4 - 37.4 °C / 95.7 - 99.3 °F

8. Changing between Celsius and Fahrenheit

This thermometer can display temperature readings in either Fahrenheit or Celsius. To switch the display between °C and °F, **press and hold** the MODE button ④ for 3 seconds; the beeper icon is shown on the display. Press the MODE button again; the current measurement scale («°C» or «°F» icon) will be shown on the display ①⑦. Change the measurement scale between °C and °F by pressing the M-button ③. When the measurement scale has been chosen, press the START/IO button ⑥ to enter the «ready for measuring» mode; otherwise the device automatically switches to ready for measuring after 10 seconds ⑨.

9. How to recall 30 readings in memory mode

This thermometer can recall the last 30 readings with a record of both time and date.

- **Recall mode ⑱:** Press the M-button ③ to enter recall mode when the power is off. The memory icon «M» will show.
- **Reading 1 - the last reading ⑲:** Press and release the M-button ③ to recall the last reading. Number «1» and a flashing «M» are displayed.

Pressing and releasing the M-button ③ after the last 30 readings have been recalled will resume the above sequence from reading 1.

10. Battery installation and replacement

This device is supplied with a new 3V CR2032 battery. Battery needs replacing while the  ⑳ is shown on the display. Remove the battery cover ㉑ by sliding it in the direction shown.

Replace the batteries – ensure correct polarity as shown by the symbols in the compartment.

Caution: Remove battery if the device is not going to be used for a prolonged period.

11. Error messages

Error	Description	Solutions
Measured temperature too high (20)	«H» = too high, measured temperature is higher than 43 °C / 109.4 °F in body mode or 99.9 °C / 211.8 °F in object mode.	Make sure the measuring sensor is clean. Make sure the thermometer is properly inserted. Then, take a measurement again.
Measured temperature too low (21)	«L» = too low, measured temperature is lower than 32 °C / 89.6 °F in body mode or 0.1 °C / 32.2 °F in object mode.	
Ambient temperature too high (22)	«AH» = too high, ambient temperature is higher than 40.0 °C / 104.0 °F.	Allow the thermometer to rest in a room for at least 30 minutes at ambient temperature (10,0 °C / 50,0 °F in body mode or 5.0 °C / 41.0 °F – 40.0 °C / 104.0 °F in object mode). Then, take a measurement again.
Ambient temperature too low (23)	«AL» = too low, ambient temperature is lower than 10.0 °C / 50.0 °F in body mode or lower than 5.0 °C / 41.0 °F in object mode.	
Blank display (25)	The battery / battery's polarity («+» and «-») has not been loaded correctly.	Check if the battery / battery's polarity («+» and «-») has been loaded correctly.
Flat battery indication (26)	The battery is low.	The battery should be replaced immediately.
Incorrect location indication (14)	The probe is not correctly inserted in the ear canal.	Please insert the probe as described in this manual.

Error	Description	Solutions
Error function display (15)	The device has a malfunction.	Remove and reload the battery and reboot the device. If the message still appears, please contact the Microlife local retail or distributor.

12. Device maintenance and disposal

Cleaning and disinfecting

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring sensor. Ensure that no liquid enters the interior of the device. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the device in water or other cleaning liquids. Take care not to scratch the surface of the sensor lens and the display.

Storage

When not in use: Keep the device and accessories in a dry, cool place away from sunlight, with ambient conditions within the temperature and humidity ranges described in the «Technical specifications». Remove the batteries from the device if the device will not be used for an extended period.

- Protect the device and accessories from the following to avoid damaging the device:
 - water, other liquids, and moisture
 - extreme temperatures
 - impacts and vibrations
 - direct sunlight
 - contamination and dust

Caution: Storing the disused device for an extended period without removing batteries increases the chance of battery fluid leakage, which may lead to device damage and skin irritation when in contact. If your eye or skin is exposed to battery fluid, wash the exposed part immediately with ample clean water. Consult a doctor if irritation or discomfort persists.

Calibration and support

The device is calibrated during manufacturing. In general, it is recommended to have the device verified by your local designated Microlife device distributor every two years, or after mechanical impact, liquid ingress, and/or device malfunctions. For questions related to device measurement accuracy, please contact your local designated Microlife device distributor.

Caution: Do not attempt to service or calibrate the device and accessories yourself.

Disposal

This device is medical electrical equipment. Dispose this device and batteries in accordance with the Waste Electrical and Electronic Equipment (WEEE) directive and applicable local regulations. Do not dispose of the device and batteries with domestic or commercial waste.

13. Guarantee

This device is covered by a **5 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Battery.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website: www.microlife.com/support
Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

14. Technical specifications

Device name:	Digital infrared ear thermometer
Model number:	IR1DN1
Catalogue number:	IR 210
Operating mode:	Adjusted mode
Measuring site:	Ear
Reference body site:	Sublingual
Measurement range:	Body mode: 32.0 - 43.0 °C / 89.6 - 109.4 °F
Resolution:	0.1 °C / °F
Measurement accuracy (Laboratory):	Body mode: ±0.2 °C, 35.0 ~ 42.0 °C / ±0.4 °F, 95.0 ~ 107.6 °F ±0.3 °C, 32.0 ~ 34.9 °C and 42.1 ~ 43.0 °C / ±0.5 °F, 89.6 ~ 94.8 °F and 107.8 ~ 109.4 °F
Clinical results:	Repeatability: 0.19 °C Bias: 0.03 °C Limits of agreement: 1.33 °C
Display:	Liquid Crystal Display , 4 digits plus special icons
Acoustic:	The unit is turned ON and ready for the measurement: 1 short beep. Complete the measurement: 1 long beep. System error or malfunction: 3 short beeps. High body temperature: 10 short beeps.
Memory:	30 readings recall in the memory mode with a record of both time and date.

Backlight:

The display light will be GREEN for 1 second, when the unit is turned ON. The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5 °C / 99.5 °F.

The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5 °C / 99.5 °F.

Operating conditions:

10 – 40 °C / 50 – 104 °F
15 - 90 % relative maximum humidity
700 – 1060 hPa Atmospheric pressure

Storage and transport conditions:

-20 – +55 °C / -4 – +131 °F
15 - 90 % relative maximum humidity
700 – 1060 hPa Atmospheric pressure

Automatic Switch-off:

Approx. 1 minute after last measurement has been taken.

Battery:

1x CR2032 battery 3V

Battery lifetime:

approx. 800 measurements (using new batteries)

Dimensions:

150 x 43 x 60 mm

Weight:

60 g (with battery), 58 g (w/o battery)

Ingress protection (IP) rating:

IP22: Protected against solid foreign objects of 12.5 mm Ø and greater. Protected against vertically falling water drops when enclosure tilted up to 15°

Expected service life:

5 years or 12000 measurements

NOTE: The technical specifications may change without prior notice.

15.Compliance information

This device complies with the requirements of the Medical Device Regulation (EU)2017/745.

Compliant standards:

EN 60601-1, EN 60601-1-2, EN 60601-1-11, EN ISO 80601-2-56

16.Symbols and definitions

Read the instructions carefully before using this device.

Type BF applied part

Caution

Medical device

Reference number

Model number

Serial number (YYYY-MM-DDXXXXX;
YYYY-MM-DD = date of manufacture,
XXXXX = sequence number)

Lot number (YYYY-MM-DD;
YYYY-MM-DD = Date of manufacture)

Manufacturer

Country of manufacture
(Date of manufacture if date printed next to symbol)

CE Marking of Conformity

Authorized representative in the European Union

Authorized representative in Switzerland

Temperature limitation for operating **or** storage

Humidity limitation for operating **and** storage

Dispose in accordance with waste electrical and electronic equipment (WEEE) directive.

Importer

Unique Device Identifier

Microlife ACCUsens Position Guidance Technology

Microlife Silent Glow™ display technology

Single use only

Keep out of reach of children

Keep away from sunlight

Patient information website

Reminder/Note

Not made with natural rubber latex

Keep dry

17. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com.

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