Microlife Corporation
9F, No. 431, RuiGuang Road,
NeiHu, Taipei, 114, Taiwan, China
www.microlife.com

 Microlife UAB
P. Lukšio g. 32,
08222 Vilnius, Lithuania
eurep@microlife.lt

CE 1639


IB BP W3 Comfort E-V11 2026-04-06

BP W3 Comfort

Wrist Blood Pressure Monitor

EN → 1	CZ → 70	
RU → 17	SK → 86	HU → 134
BG → 36	SL → 102	HR → 151
RO → 54	SR → 118	PL → 167

microlife®

Name of Purchaser / Ф.И.О. покупателя / Име на купувача /
 Numele cumpărătorului / Imię i nazwisko nabywcy /
 Jméno kupujícího / Meno zákazníka / Ime in priimek kupca /
 Ime i prezime kupca / Vásárló neve / Ime i prezime kupca

Serial Number / Серийный номер / Сериен номер /
 Număr de serie / Numer seryjny / Výrobní číslo / Výrobné číslo /
 Serijska številka / Serijski broj / Sorozatszám / Serijski broj

Date of Purchase / Дата покупки / Дата на закупуване /
 Data cumpărării / Data zakupu / Datum nákupu / Dátum kúpy /
 Datum nakupa / Datum kupovine / Vásárlás dátuma /
 Datum kupovine

Specialist Dealer / Специализированный дилер / Специалист
 дистрибутор / Distributor de specialitate / Przedstawiciel /
 Specializovaný dealer / Špecializovaný predajca / Specializirani
 trgovec / Ovlašćeni diler / Forgalmazó / Ovlašteni prodavač

- ① ON/OFF button
- ② Display
- ③ M-button (memory)
- ④ Time button
- ⑤ Battery compartment
- ⑥ Cuff

Display

- ⑦ Systolic value
- ⑧ Diastolic value
- ⑨ Pulse rate
- ⑩ Date/Time
- ⑪ Stored value
- ⑫ Traffic light indicator
- ⑬ Irregular heartbeat (IHB) symbol
- ⑭ Cuff signal indicator «ERR 1»
- ⑮ Cuff fit check
 - A: Suboptimal cuff fit
 - B: Arm movement indicator «ERR 2»
 - C: Cuff pressure check «ERR 3»

- ⑯ Pulse indicator

- ⑰ Battery display

Dear Customer,

This device was developed in collaboration with physicians and clinical tests carried out prove its measurement accuracy to be of a very high standard.*

If you have any questions, problems or want to order spare parts please contact your local Microlife-Customer Service. Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit the internet at www.microlife.com where you will find a wealth of invaluable information on our products.

Stay healthy – Microlife Corporation!

** This device is tested according to the ESH protocol and ISO81060-2:2013.*

Table of contents

1. Introduction

- Document scope
- Disclaimers

2. Important information

Device description
Intended use
Intended user
Intended patient
Intended use environment and conditions
Indications
Contraindications
Side effects
Warning
Caution
Electromagnetic compatibility
Adverse events and reporting

3. Device information

Device accessories
Low battery
Flat battery - replacement

4. Device installation and setup

Inserting the batteries
Setting the date and time

5. Measurement preparation

Before taking a measurement

6. Measurement operation

Starting measurement
How not to store a reading

7. Measurement interpretation

How do I evaluate my blood pressure
Appearance of the irregular heartbeat (IHB) symbol

8. Data memory function

Viewing the stored values
Clearing all values

9. Device error and troubleshooting

10. Device maintenance and disposal

Storage
Calibration and support
Disposal

11. Specifications and compliance

Technical specifications
Compliance information

12. Supplement information for users and patients

Guarantee
Symbols and definitions

1. Introduction

Document scope

Read the instructions carefully before using this device.

This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.

Disclaimers

Microlife® is a registered trademark of Microlife Corporation. Trademarks and trade names are those of their respective owners.

2. Important information

Device description

A digital home-use blood pressure monitor is a medical device that utilizes the principles of cuff-based oscillometric method and digital signal process to compute and provide a blood pressure measurement.

Intended use

This device is intended to measure radial (wrist) blood pressure (systole and diastole) and pulse rate.

Intended user

The device is intended to be operated by adults and adolescents with adequate vision, motor functions, and education, capable of understanding the instructions for use and operating general household electrical appliances.

Intended patient

The intended patients are normotensive and hypertensive adults and adolescents (aged 12 years or older) of the general population.

Intended use environment and conditions

The device is intended for use in a home healthcare environment (e.g. general household without medically trained personnel) by patients (e.g. for self-measurement) or by a care giver.

Indications

This device measures blood pressures for indications of:

- Diagnosis of white-coat hypertension and masked hypertension and identifying white-coat effect and masked uncontrolled hypertension.
- Evaluating blood pressure in response to treatment.
- Confirming the diagnosis of resistant hypertension.
- Detecting morning hypertension.

Contraindications

- The device is not intended for measuring blood pressure in pediatric patients younger than 12 years of age (children, infant, or neonates).
- The device measures blood pressure using a pressured cuff. If the measuring limb suffers from injuries (for example open wounds) or under conditions or treatments (for example intravenous drip) making it unsuitable for surface contact or pressurization, do not use the device, to avoid worsening of the injuries or conditions.

- Avoid taking measurements of patients with conditions, diseases, and susceptible to environment conditions that lead to uncontrollable motions (e.g. trembling or shivering) and inability to communicate clearly (for example children and unconscious patients).
- The device uses oscillometric method to determine blood pressure and requires the measured limb with normal perfusion. The device is not intended to be used on a limb with restricted or impaired blood circulation. Consult with your doctor if you have severe perfusion or blood disorders before using the device.

Side effects

In rare cases, slight bruising may result after measurement due to pressurization of the arm.

Warning

NOTE: Warning items indicate potentially hazardous situations, if not avoided, may result in death, critical or serious injury to the user or patient.

- Avoid taking measurement on the arm on the side of a mastectomy or lymph node clearance.
- Avoid taking measurements on the arm with intravascular access or therapy or an arterio-venous (A-V) shunt. Cuff and pressuriza-

tion may temporarily interfere with blood flow and could result in injury.

- Presence of significant cardiac arrhythmia during measurement may interfere with blood pressure measurement and affect the reliability of blood pressure readings. Consult with your doctor about whether the device is suitable for use in this case.
- DO NOT use this device in a moving vehicle (for example in a car or on an aircraft).
- DO NOT use this device for purposes beyond described in these instructions for use. The manufacturer cannot be held liable for damage caused by incorrect application.
- The measurement result of this device is not a medical diagnosis and not intended to substitute consultation and diagnosis by a qualified professional healthcare provider (e.g., physician, pharmacist, or other licensed health-care professionals).
- DO NOT use this device for self-diagnosis or for self-treatment of a medical condition. Seek advice from a health-care professional immediately if the patient is clearly unwell and/or having physiological or medical symptoms.
- Inspect the device, cuff and other parts for damage. DO NOT USE the device, cuff or parts if they appear damaged or operating abnormally.

- Blood flow of the arm is temporarily interrupted during measurement from cuff pressurization. Extended periods of cuff pressurization reduces peripheral circulation. Beware of signs (e.g tissue discoloration) of impeded peripheral circulation when taking prolonged or multiple measurements. It is recommended to rest between measurements. Abort measurement, loosen the cuff (or disconnect the cuff and device) and rest to restore perfusion.
- DO NOT use this device in oxygen rich environment or near flammable gas.
- DO NOT use this device with other medical electrical (ME) equipment simultaneously. This may cause device malfunction or measurement inaccuracies.
- Use and store the device, cuff and parts in temperature and humidity conditions specified in the «Specifications and compliance». Usage and storage of the device, cuff and parts in conditions outside ranges given in the «Specifications and compliance» may results in device malfunction and the safety of usage.
- Keep the device away from children and people incapable of operating the device. Beware of the risks of accidental ingestion of small parts and of strangulation with the cables and tubes of this device and accessories.
DO NOT let children operate the device alone.

Caution

NOTE: Caution items indicate potentially hazardous situations, if not avoided, may result in minor or negligible injury to the user or patient, or damage to the property or environment.

- The device is not intended to measure pulse rate to check the frequency of a pacemaker.
- DO NOT disassemble or attempt to service the device, accessory, and parts, during use or in storage. Access to the device internal hardware and software is prohibited. Unauthorized access and servicing of the device, during use or in storage, may compromise the safety and performance of the device.
- The device is intended only for measuring blood pressure at wrist. Do not measure other sites because the reading does not reflect your blood pressure accurately.
- After a measurement is completed, loosen the cuff, and rest the wrist to restore limb perfusion, before taking another measurement.
After a measurement is completed, loosen the cuff and rest the wrist for at least 1 minute to restore limb perfusion, before taking another measurement.

- Protect the device and accessories from the following to avoid damaging the device:
 - water, other liquids, and moisture
 - extreme temperatures
 - impacts and vibrations
 - direct sunlight
 - contamination and dust
- This device is reusable. It is recommended to clean and disinfect the device and the accessory prior to use if contamination or cross contamination are possible.
- Stop using this device and cuff and consult with your doctor if you experience skin irritation or discomfort.
- DO NOT use this device, cuff, or parts after the expiration of its stated service life.
- If this monitor is stored at the maximum or minimum storage and transport temperature and is moved to an environment with a temperature of 20 °C, we recommend waiting for approximately 4 hours before using the monitor.
- Storing the device **unused** for an extended period without removing batteries increases the chance of battery fluid leakage, which may lead to device damage and skin irritation when in contact. If your eye or skin is exposed to battery fluid, wash the

exposed part immediately with ample clean water. Consult a doctor if irritation or discomfort persists.

Electromagnetic compatibility

- This device is compliant with electromagnetic disturbances standard.

Further documentation in compliance with EN 60601-1-2 standard is available from Microlife on www.microlife.com/electro-magnetic-compatibility.

- DO NOT use this device in proximity of equipment that may cause electromagnetic disturbance (EMD), such as high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, and computerized tomography (CT) scanners. This device is not certified for operation near these equipments, which may cause device malfunction and measurement inaccuracies.
- Do not use this device close to strong electromagnetic fields and portable radio frequency communication devices (for example microwave oven and mobile devices). Keep a minimum distance of 0.3 m from such devices when using this device.

Caution: The use of non-Microlife or non-compatible accessories may result in increased emissions or decreased immunity of the equipment or system.

Adverse events and reporting

Please report any serious incident, injury or adverse event that has occurred in relation to the device to the manufacturer/ European authorized representative (EU REP), and to the competent authority.

3. Device information

Package Contents

- 1 x Microlife BP W3 Comfort
- 1 x Instruction manual
- 1 x BP diary
- 1 x storage case
- 2 x 1.5 V alkaline batteries; type AAA

CAUTION: Inspect the device, cuff, and other parts for damage. DO NOT USE the device, cuff or parts if they appear damaged or operating abnormally.

Device accessories

Batteries

Use 2 new 1.5 V, size AAA alkaline batteries.

Caution: Do not use expired batteries or mix new and used batteries together.

Caution: Remove batteries if the device is not going to be used for a prolonged period.

You can also operate this device using rechargeable batteries.

Only use «NiMH» type reusable batteries.

Batteries must be removed and recharged when the flat battery symbol appears. They should not remain inside the device as they may become damaged (total discharge as a result of low use of the device, even when switched off).

Batteries cannot be charged in the blood pressure monitor. Recharge batteries in an external charger and observe the information regarding charging, care and durability.

Low battery

When the batteries are approximately $\frac{3}{4}$ empty the battery symbol ⑰ will flash as soon as the device is switched on (partly filled battery displayed). Although the device will continue to measure reliably, you should obtain replacement batteries.

Flat battery - replacement

When the batteries are flat, the battery symbol ⑰ will flash as soon as the device is switched on (flat battery displayed). You cannot take any further measurements and must replace the batteries.

1. Open the battery compartment ⑤ by pulling off the cap.
2. Replace the batteries – ensure correct polarity as shown by the symbols in the compartment.

3. To set date and time, follow the procedure described in Section «Setting the date and time».

 The memory retains all values although date and time must be reset – the year number therefore flashes automatically after the batteries are replaced.

 The memory retains all values although date and time must be reset – the year number therefore flashes automatically after the batteries are replaced.

4. Device installation and setup

Inserting the batteries

After you have unpacked your device, first insert the batteries. The battery compartment ⑤ is on the bottom of the device. Insert the batteries (2 x 1.5 V, size AAA), thereby observing the indicated polarity.

 **Caution:** Inserting the batteries in incorrect polarity orientations may lead to short circuiting and damage the device!

Setting the date and time

1. After the new batteries are fitted, the year number flashes in the display. You can set the year by pressing the M-button ③. To confirm and then set the month, press the time button ④.

2. Press the M-button to set the month. Press the time button to confirm and then set the day.
3. Follow the instructions above to set the day, hour and minutes.
4. Once you have set the minutes and pressed the time button, the date and time are set and the time is displayed.
5. If you want to change the date and time, press and hold the time button for approx. 3 seconds until the year number starts to flash. Now you can enter the new values as described above.

Caution: Make sure date and time settings are correct on the device. Incorrect settings results in misleading data and time records of the measurements.

5. Measurement preparation

Before taking a measurement

- ▶ Avoid heavy activity, eating or smoking immediately before the measurement.
- ▶ Empty your bladder prior to measurement.
- ▶ Sit down on a back-supported chair and relax for 5 minutes. Keep your feet flat on the floor and do not cross your legs.
- ▶ Always measure in a sitting position and on the same wrist. Use the wrist which usually shows higher blood pressure values.
- ▶ Remove any items of clothing and your watch, for example, so that your wrist is free.

- ▶ Always ensure that the cuff is positioned correctly, as shown in the pictures illustrated on the short instruction card.
- ▶ Fit the cuff comfortably but not too tight. The cuff will cover a wrist circumference according to the information in the «Specifications and compliance».

6. Measurement operation

Starting measurement

1. Press the ON/OFF button ① to start the measurement.
2. The cuff fit check ⑮ on the display indicates that the cuff is perfectly placed. If the icon ⑮-A appears, the cuff is fitted suboptimally, but it is still ok to measure.
3. The cuff will now pump up automatically. Relax, do not move and do not tense your arm muscles until the measurement result is displayed. Breathe normally and do not talk.
4. The measurement is performed during the inflation. The inflation speed may vary, this is a normal occurrence.
5. During the measurement, the pulse indicator ⑯ flashes in the display.
6. The result, comprising the systolic ⑦ and the diastolic ⑧ blood pressure and the pulse rate ⑨ are displayed. Note also the explanations on further display symbols in this booklet.

7. Remove and switch off the monitor and enter the result in the enclosed blood pressure pass. (The monitor does switch off automatically after approx. 1 min.).

You can stop the measurement at any time by pressing the ON/OFF button or open the cuff (e.g. if you feel uneasy or an unpleasant pressure sensation).

Caution: Remain still and do not move or talk during measurement. Motions caused by talking, moving, trembling and other vibrations may interfere with the measurement and affect the measurement accuracy!

Caution: You can stop the measurement at any time by pressing the ON/OFF button or open the cuff (e.g. if you feel uneasy or an unpleasant pressure sensation).

How not to store a reading

As soon as the reading is displayed press and hold the ON/OFF button ① until «M» ⑪ is flashing. Confirm to delete the reading by pressing the time button ④.

«CL» is displayed when the reading is deleted from the memory successfully.

7. Measurement interpretation

How do I evaluate my blood pressure

The triangle on the left-hand edge of the display ⑫ points at the range within which the measured blood pressure value lies. The value is either within the optimum (white), elevated (hatched gray) or high (black) range. The classification corresponds to the following ranges defined by international guidelines (ESH, ESC, JSH). Data in mmHg.

NOTE: The blood pressure classification is a general guideline of blood pressure level at home, but diagnosis of hypertension should be made by a healthcare professional based on specific conditions of the patient. Consult with your doctor for questions about the interpretation and classification of your blood pressure values.

Range		Systolic	Diastolic	Recommendation
1.	blood pressure normal	< 120	< 74	Self-check
2.	blood pressure optimum	120 - 129	74 - 79	Self-check
3.	blood pressure elevated	130 - 134	80 - 84	Self-check

Range		Systolic	Diastolic	Recommendation
4.	blood pressure too high	135 - 159	85 - 99	Seek medical advice
5.	blood pressure dangerously high	≥ 160	≥ 100	Urgently seek medical advice!

The higher value is the one that determines the evaluation. Example: a blood pressure value of **140/80** mmHg or a value of **130/90** mmHg indicates «blood pressure too high».

Appearance of the irregular heartbeat (IHB) symbol

This symbol ⑬ indicates that an irregular heartbeat was detected. In this case, the measured blood pressure may deviate from your actual blood pressure values. It is recommended to repeat the measurement.

Information for the doctor in case of repeated appearance of the IHB symbol:

This device is an oscillometric blood pressure monitor that also measures the pulse during blood pressure measurement and indicates when the heart rate is irregular.

8. Data memory function

This device automatically stores the last 30 measurement values.

Viewing the stored values

Press the M-button ③ briefly, when the device is switched off. The display first shows «M» ⑪, and «A» which stands for the average of all stored values.

Pressing the M-button again displays the previous value. Pressing the M-button repeatedly enables you to move from one stored value to another.

 Blood pressure readings with suboptimal cuff fit ⑮-A are not considered in the average value.

 Pay attention that the maximum memory capacity of 30 memories is not exceeded. **When the 30 memory is full, the oldest value is automatically overwritten with the 31 th value.** Values should be evaluated by a doctor before the memory capacity is reached – otherwise data will be lost.

Clearing all values

If you are sure that you want to permanently remove all stored values, hold down the M-button (the device must have been switched off beforehand) until «CL» appears and then release the

button. To permanently clear the memory, press the time button while «CL» is flashing. **Individual values cannot be cleared.**

 **Cancel deletion:** press ON/OFF button ① while «CL ALL» is flashing.

9. Device error and troubleshooting

If an error occurs during the measurement, the measurement is interrupted and an error message, e.g. «ERR 3», is displayed.

Error	Description	Potential cause and remedy
«ERR 1» ⑭	Signal too weak	The pulse signals on the cuff are too weak. Re-position the cuff and repeat the measurement.*
«ERR 2» ⑮-B	Error signal	During the measurement, error signals were detected by the cuff, caused for instance by movement or muscle tension. Repeat the measurement, keeping your arm still.

Error	Description	Potential cause and remedy
«ERR 3» ⑮-C	No pressure in the cuff	An adequate pressure cannot be generated in the cuff. A leak may have occurred. Check that the cuff is correctly connected and is not too loose. Replace the batteries if necessary. Repeat the measurement.
«ERR 5»	Abnormal result	The measuring signals are inaccurate and no result can therefore be displayed. Read through the checklist for taking a reliable measurement and then repeat the measurement.*
«HI»	Pulse or cuff pressure too high	The pressure in the cuff is too high (over 299 mmHg) OR the pulse is too high (over 40 beats per minute). Relax for 5 minutes and repeat the measurement.*
«LO»	Pulse too low	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

** Please immediately consult your doctor, if this or any other problem occurs repeatedly.*

10. Device maintenance and disposal

Cleaning the device

The device can be cleaned when necessary (e.g., between uses by different patients).

Use a soft cloth, dry or wet with detergent, to gently wipe the exterior of the device to remove dust or stains.

Cleaning the cuff

Use a soft cloth, dry or wet with mild detergent, to carefully wipe the cuff to remove dust or stains.

Storage

When not in use:

- Keep the device and accessories in a dry, cool place away from sunlight, with ambient conditions within the temperature and humidity ranges described in the «Specifications and compliance» section.
- Remove the batteries from the device if the device will not be used for an extended period.

Warning: Storing the device **unused** for an extended period without removing batteries increases the chance of battery fluid leakage, which may lead to device damage and skin irritation when in contact. If your eye or skin is exposed to battery fluid, wash the exposed part immediately with ample clean water. Consult a doctor if irritation or discomfort persists.

Calibration and support

The device is calibrated during manufacturing. In general, it is recommended to have the device verified by your local designated Microlife device distributor every two years, or after mechanical impact, liquid ingress, and/or device malfunctions. For questions related to device measurement accuracy, please contact your local designated Microlife device distributor.

Caution: Do not attempt to service or calibrate the device and accessories yourself.

Disposal

This device is medical electrical equipment. Dispose this device and batteries in accordance with the Waste Electrical and Electronic Equipment (WEEE) directive and applicable local regulations. Do not dispose of the device and batteries with domestic or commercial waste.

11. Specifications and compliance

Technical specifications

NOTE: Technical specifications subjected to change without notice.

Device Type:	Digital non-invasive blood pressure monitor
Model number:	BP3KB1-3E
Reference number	BP W3 Comfort
Operating conditions:	10 – 40 °C / 50 – 104 °F 15 – 90 % relative maximum humidity 700 hPa – 1060 hPa
Storage and transport conditions:	-20 – +55 °C 15 – 90 % relative maximum humidity
Weight:	112 g (including batteries)
Dimensions:	81 x 65 x 21 mm
Cuff size:	13.5 – 23 cm (5.31 – 9.05 inches)
Measuring procedure:	Oscillometric, corresponding to Korotkoff method: Phase I systolic, Phase V diastolic
Pressure resolution:	1 mmHg
Cuff pressure display range:	0 – 299 mmHg

Measurement range:	SYS: 60 ~ 255 mmHg DIA: 40 ~ 200 mmHg Pulse: 40 – 199 beats per minute
Static accuracy:	± 3 mmHg
Pulse accuracy:	± 5 % of the readout value
Power source – internal:	2 x 1.5 V AAA batteries
Ingress protection (IP) rating:	IP22: Protected against solid foreign objects of 12.5 mm Ø and greater. Protected against vertically falling water drops when enclosure tilted up to 15°
Applied part type reference:	 Type BF
Service life – device:	5 years or 10000 measurements, whichever comes first
Battery lifetime:	approx. 320 measurements (1.5 V alkaline batteries; size AAA)
Compliance information	This device complies with the requirements of the Medical Device Regulation (EU)2017/745.
Compliant standards:	

EN 60601-1
EN 60601-1-2
EN 60601-1-11
EN IEC 80601-2-30
EN ISO 81060-2

12. Supplement information for users and patients

Guarantee

This device is covered by a **5 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge. Opening or altering the device invalidates the guarantee.

The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.
- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Batteries.

The cuff is covered by a functional guarantee (bladder tightness) for 2 years.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website: www.microlife.com/support

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

Symbols and definitions

Medical device

CE Marking of Conformity

Importer

Authorized representative in the European Union

Manufacturer

Country of manufacture
(Date of manufacture if date printed next to symbol)

Model number

Reference number

Serial number (YYYY-MM-DDXXXXX; YYYY-MM-DD = date of manufacture, XXXXX = sequence number)

Lot number (YYYY-MM-DD; YYYY-MM-DD = Date of manufacture)

Unique Device Identifier

Caution

General warning sign

Type BF applied part

Direct current

Protected against vertically falling water drops when enclosure tilted up to 15°.

Keep dry

Temperature limitation for operating **or** storage

Humidity limitation for operating **and** storage

Atmospheric pressure limitation

Read instructions for use before operating the device.

Dispose in accordance with waste electrical and electronic equipment (WEEE) directive.

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