# Safe, easy and COMFORTABLE

Natural and gentle feeling with vacuum technology

Baby Care Highest comfort with Microlife Breast Pumps.

microlife

A partner for people. For life.



www.microlife.com



# BC 100 Soft

### Manual breast pump

**BPA** free



Light, easy and and ultra-quiet Allows comfortable expressing of milk anywhere



0+

months

### The bottle is completely free of BPA

Ideal for newborns Milk bottle and silicone teat are specially designed for newborns

## BPA O+

### Features

- Vacuum technology simulates natural drinking for gentle and comfortable use
- Extra comfort with 2 expressing settings
- Hygienic and safe accessories are washable and disinfectable
- Fully equipped with 180ml milk bottle (completely BPA free) with cap, bottle stand, two matching screw plugs for bottle and a naturally shaped silicone teat for easy feeding
- Easy storage and transportation with carrying bag
- www.microlife.com/bc100-soft



# BC 200 Comfy

### Electric breast pump



### Individual usage

Highest comfort with 10 different stimulation settings



### 10 different pumping modes

Helps mothers to select the most comfortable pump setting for their milk flow



#### Easy to use

Large illuminated display with timer display easily tracks stimulation and expressing time

#### Features

• Designed to mimic baby's natural sucking rhythm with 2-phase vacuum technology

- Lightweight and compact design for easy expressing at home or on-the-go
- Memory function for personalised usage
- Battery and mains power operated (with battery level indicator)
- Equipped with 180ml milk bottle (completely BPA free) with cap, bottle stand, two matching screw plugs for bottle, naturally shaped silicone teat for easy feeding, medical power adapter and carrying bag
- www.microlife.com/bc200-comfy



## BC 300 Maxi 2 in 1

Dual electric breast pump





Innovative dual mode technology

Suitable for one and two-sided expressing

### Gentle, convenient and time-saving

Double expressing system for more milk in less time



### Highly comfortable

Memory function for personalised usage

#### Features

- Time-saving with dual mode vacuum technology
- Highest comfort with 10 different stimulation and expressing settings
- Large illuminated display with timer display which easily tracks stimulation and expressing time
- Battery and mains power operated (with battery level indicator)
- Equipped with 2 x 180ml milk bottles with caps, bottle stands, two matching screw plugs for bottle, naturally shaped silicone teat for easy feeding, medical power adapter and carrying bag

www.microlife.com/bc300-maxi-2in1

### Overview

Which breast pump is best for you?

	Market Control of Cont	entre entre entre entre	
	BC 100 Soft	BC 200 Comfy	BC 300 Maxi 2 in 1
Manual breast pump	٠		
Single electric breast pump		٠	
Dual electric breast pump			٠
2-phase vacuum technology	•	٠	٠
Usage	Occasionally	Often	Frequent
Weight	135g	330 g	330 g
Size	200×160×66 mm	130 x 120 x 57 mm	130 x 120 x 57 mm
Guarantee	2 years	2 years	2 years
Stimulation mode	2	10	10
Pumping mode	2	10	10
Memory function		٠	٠
-🍟- Display backlight		٠	٠
Battery and mains power operation		٠	٠
Light, compact and quiet	•	٠	٠
O% BPA free	•	٠	٠
Washable accessories	٠	٠	٠
Accessories are ideal for newborns	٠	٠	٠

### Effective and gentle expression with Microlife breast pumps.

Microlife breast pumps allow nursing mothers to express their breast milk easily and in the comfort of their own homes. If not enough breast milk is taken from mothers who have recently given birth, this can lead to the painful build-up of milk, which in the worst-case scenario could cause inflammation of the nipples (mastitis). This build-up of milk can be prevented through targeted expressing of breast milk using the breast pump. In addition, breast pumps help to alleviate cracked or sore nipples.

You can use the integrated stimulation function to draw out flat or inverted nipples, making it easier to then express the breast milk. The breast pump also allows you to store breast milk in the bottle, for example if you will be travelling, or away from your baby for an extended period of time. Babies who cannot be breast-fed directly for health reasons (due to a premature birth or cleft palate) can still be given breast milk by using the breast pump and feeding the baby using a bottle and teat.



Individual solutions for the best possible support for a great and long lasting breastfeeding experience.

#### Manual breast pump

Microlife manual breast pumps are an excellent choice for mothers who express their breast milk occasionally. The innovative design is easy to assemble, lightweight and ultra-silent. No battery or mains power is required as they are operated individually by hand. They are ideal for using on-the-go and are almost silent during use.

#### Single electric breast pump

Microlife electric breast pumps are designed to mimic your baby's natural feeding action, for more comfort and more milk. Ideal for mothers who want to comfortably express their breast milk every day. Powered by batteries or mains adapter to allow comfortable expressing anywhere.

#### Dual electric breast pump

Microlife double pumping systems are the perfect time saver for mothers: The dual breast pump vacuum technology can express more milk in a shorter period of time than conventional breast pumps. Simultaneous expressing of milk at both breasts is proven to be more efficient and may even boost your ability to produce breast milk.

> Microlife AG Espenstrasse 139 9443 Widnau / Switzerland www.microlife.com